



EnzoLyte™ 390 Renin Assay Kit

Fluorimetric

Catalog #	72039
Unit Size	1 Kit
Kit Size	100 Assays (96-well) or 300 assays (384-well)

This kit is optimized to detect the activity of the enzyme renin using a Mc-Ala/Dnp FRET peptide substrate, with the signal monitored at Ex/Em=330/390 nm upon proteolytic cleavage. Ample materials are provided to perform 100 assays in 96-well format or 300 assays in 384-well plate. The kit has the following features:

- **Convenient Format:** Complete kit including all the assay components.
- **Optimized Performance:** Optimal conditions for screening of renin inhibitors.
- **Enhanced Value:** Less expensive than the sum of individual components.
- **High Speed:** Minimal hands-on time.
- **Assured Reliability:** Detailed protocol and references are provided

USA and Canada Ordering Information

AnaSpec Corporate Headquarter

2149 O'Toole Ave.
San Jose, CA 95131
Toll-Free: 800-452-5530
Tel: 408-452-5055
Fax: 408-452-5059
E-mail: service@anaspec.com
Internet: www.anaspec.com

International Ordering Information


A list of international distributors is available at www.anaspec.com.

INTRODUCTION

The renin–angiotensin system (RAS) plays a central role in the regulation of blood pressure and electrolyte homoeostasis.¹ At the first and rate-limiting step of the RAS cascade, renin (EC 3.4.23.15), a highly specific aspartyl protease, cleaves angiotensinogen, produced in the liver, to yield angiotensin I, which is further converted into angiotensin II by ACE (Angiotensin Converting Enzyme). Angiotensin II constricts blood vessels leading to increased blood pressure. It increases the secretion of ADH and aldosterone, and stimulates the hypothalamus to activate the thirst reflex. Since an overactive renin–angiotensin system leads to hypertension, renin is an attractive target for the treatment of this disease.²⁻⁴

The EnzoLyte™ 390 Renin Assay Kit provides a convenient assay for high throughput screening of renin inhibitors and for continuous assay of renin activity using a Mc-Ala/Dnp fluorescence resonance energy transfer (FRET) peptide. The sequence of this FRET peptide is derived from the cleavage site of renin.⁵ In the FRET peptide the fluorescence of Mc-Ala is quenched by Dnp. Upon cleavage into two separate fragments by renin (**Scheme 1**), the fluorescence of Mc-Ala is recovered, and can be monitored at excitation/emission = 330/390 nm. Compared to an EDANS/DABCYL FRET substrate, the Mc-Ala/Dnp FRET substrate has higher sensitivity and can detect the activity of 25 pM renin.

The assays are performed in a convenient 96-well or 384-well microplate format.



Scheme 1. Proteolytic cleavage of Mc-Ala/Dnp FRET peptide by renin.

KIT COMPONENTS, STORAGE AND HANDLING

Note: Store all kit components at -20 °C, except for Component C, which should be stored at -80 °C. Keep Components A and B from light. If used frequently, Component D can be stored at room temperature for convenience.

- Component A:** Renin substrate (2 mM DMSO solution, 50 µL)
Mc-Ala/Dnp FRET peptide
Ex/Em=330 nm/390 nm upon cleavage
- Component B:** Mc-Ala fluorescence reference standard (3 mM, 5 µL)
Ex/Em=330 nm/390 nm
- Component C:** Human recombinant renin (100 µL)
- Component D:** Assay buffer (25 mL)
- Component E:** Renin Inhibitor Ac-HPFV- (Sta)-LF-NH₂ (1mM DMSO solution, 5 µL)

OTHER MATERIALS REQUIRED (BUT NOT PROVIDED)

96-well or 384-well microplate: Black microplates can provide better signal to noise ratios.

Fluorescence microplate reader: Capable of detecting emission at 390 nm with excitation at 330 nm.

STANDARD OPERATION PROTOCOL

Note 1: For fluorometer calibration, please refer to Appendix II - recommended for first-time users.

Screening renin inhibitors using recombinant enzyme

Note: Warm all kit components until thawed to room temperature before starting the experiments.

1. Prepare working solutions.

- Renin substrate solution: Dilute renin substrate (Component A) 1:100 in assay buffer (Component D). For each experiment prepare fresh substrate solution.
- Table 1. Renin substrate solution for one 96-well plate (100 assays).

Components	Volume
Renin substrate (100X, Component A)	50 µL
Assay buffer	4.95 mL
Total volume	5 mL

- Renin diluent: Dilute renin (Component C) 1:100 in assay buffer (Component D). This amount of enzyme is sufficient for a full 96-well plate. If not using the entire plate, adjust the amount of enzyme to be diluted, according to the experimental design.

Note: Prepare enzyme diluents immediately before use. Do not vortex the enzyme solutions. Prolonged storage or vigorous agitation of the diluted enzyme will cause denaturation. Store the enzyme solution on ice.

- Renin inhibitor:⁶ Dilute the 1 mM inhibitor solution (Component E) to 10 µM in assay buffer (Component D). Add 15 µl (96-well plate) or 5 µl (384-well plate) of the 10 µM inhibitor solution into each well.

2. Set up enzymatic reaction.

- Add test compounds and renin solution into the microplate wells. The suggested volume of renin solution for one well of a 96-well plate is 85 µL and test compound is 15 µL, or 28 µL renin solution and 5 µL of test compound for a 384-well plate.

- Set up the following controls at the same time as deemed necessary:
 - Positive control contains diluted renin without test compound.
 - Inhibitor control contains renin and renin inhibitor (provided with the kit).
 - Vehicle control contains renin and vehicle used to deliver test compound (e.g. DMSO).
 - Test compound control contains assay buffer and test compound. Some test compounds have strong autofluorescence and may give false results.
 - Substrate control contains assay buffer.

Note: Using the assay buffer, bring the total volume of all controls to 100 μ L (96-well plate) or 33 μ L (384-well plate)

3. Pre-incubation.

- Incubate the plate at 37°C for 30 min. At the same time, incubate the renin substrate solution at the same temperature.

4. Initiate the enzymatic reaction.

- Add 50 μ L (96-well plate) or 17 μ L (384-well plate) of renin substrate solution into each well. Mix the reagents completely by shaking the plate gently for no more than 30 seconds.

- Measure fluorescence signal:

For kinetic reading: Immediately start measuring fluorescence intensity at Ex/Em=330 nm/390 nm continuously and record data every 1 min. for 15 min.

For end-point reading: Incubate the reaction at 37°C for 15 min. Keep plate from direct light. Then measure fluorescence intensity at Ex/Em=330 nm/390 nm.

- Data analysis: Refer to Appendix I.




Figure 1. The inhibitory curve of renin inhibitor, Ac-HPFV- (Sta)-LF-NH₂ (IC₅₀=0.97 nM)

Appendix I: Data Analysis

- The fluorescence reading from the substrate control well is the background fluorescence. This background reading should be subtracted from the readings of the other wells. The fluorescence readings are expressed in relative fluorescence units (RFU).
- For kinetic reading:
 - Plot data as RFU versus time for each sample. If you want to convert the RFU to the concentration of the product of enzymatic reaction, please refer to [Appendix II](#) for setting up a fluorescence reference standard.
 - Determine the range of initial time points during which the reaction is linear. Typically, 10-15% conversion will be the optimal range.
 - Obtain the initial reaction velocity (V_o) in RFU/min. Determine the slope of the linear portion of the data plot.
 - A variety of data analyses can be done, e.g., determining inhibition %, EC_{50} , IC_{50} , K_m , K_i , etc.
- For endpoint reading:
 - Plot data as RFU versus concentration of test compounds.
 - A variety of data analyses can be done, e.g., determining inhibition %, EC_{50} , IC_{50} , etc.

Appendix II: Instrument Calibration

- **Mc-Ala fluorescence reference standard:** Dilute 3mM Mc-Ala (Component B) to 37.5 μM in deionized water. Do 1:2 serial dilutions to get concentrations of 18.7, 9.4, 4.7, 2.3, 1.15, 0.58 μM , include a water blank. Add 100 μL /well (96-well plate) or 33 μl (384-well plate) of these serially diluted Mc-Ala reference solutions.
- Add 50 μL /well (96-well plate) or 17 μL /well (384-well plate) of the diluted renin substrate solution (refer to SOP step 1 for preparation).

Note: The renin substrate solution is added to the Mc-Ala reference standard to correct for the absorptive quenching by the FRET peptide. If multiple concentrations of substrate are used, this step must be performed for each concentration.

- Plot the Mc-Ala fluorescent reference standard curve as RFU (relative fluorescent units) versus concentration as shown in **Figure 2**.

Note: The final concentration of Mc-Ala reference standard is 25, 12.5, 6.25, 3.125, 1.5, 0.78, 0.39, and 0 μM . This reference standard is used to calibrate the variation of different instruments and different experiments. It is also an indicator of the amount of final product of the renin enzymatic reaction.




Figure 2. Mc-Ala reference standard.

Mc-Ala was serially diluted in assay buffer containing substrate, and the fluorescence recorded at Ex/Em=330 \pm 20 nm/ 390 \pm 20 nm. (Flexstation 384II, Molecular Devices)

References

1. He, FJ. and GA. MacGregor, *J. Renin Angiotensin Aldosterone Syst.* **4**, 11 (2003).
2. Wood, JM. et al., *Hypertension*, **7**, 797 (1985).
3. Shibasaki, M. et al., *Am. J. Hypertens.* **4**, 932 (1991).
4. Wood, JM. et al., *Biochem. Biophys. Res. Comm.* **308**, 698 (2003).
5. Paschalidou, K. et al., *Biochem. J.* **382**, 1031 (2004).
6. Hui, KY. et al., *J. Med. Chem.* **31**, 1679 (1988).