



***E. coli* Ribonuclease H (RNase H), Recombinant**

BACKGROUND

Ribonuclease H (RNase H) is an endoribonuclease which specifically degrades the RNA strand of an RNA/DNA hybrid, leaving the DNA strand and unhybridized RNA intact. ***E. coli* RNase H (=RNaseHI)** was over-expressed in *E. coli* as a recombinant protein and highly purified. MW is 17.6 kDa.

Applications:	1) Removal of mRNA in DNA/RNA hybrid prior to the synthesis of the second strand of cDNA (1, 2) 2) Removal of poly (A) tails from mRNA after hybridization with oligo (dT) (3) 3) Oligodeoxyribonucleotide-directed site-specific cleavage of RNA
Size:	1,000 units
Form:	50 units/ul in 20mM Tris-HCl (pH 7.5), 100mM KCl, 1mM DTT, 50% Glycerol
Specific Activity:	100,000 units/mg protein
Unit Definition:	1 unit is defined as the amount of the enzyme that hydrolyzes 1 nmol of the RNA in 3H labeled M13 DNA/RNA hybrid to acid-soluble ribonucleotides in 20 min at 37°C.
Quality Assurance:	Greater than 95% protein determined by SDS-PAGE (CBB staining) (Fig.1). Endo- and exo-DNase activities and RNase activity were not detected with 100 U/ml RNaseH in 50 ul reaction at 37°C.
Reagents Supplied with Enzyme:	RNaseH Reaction Buffer (10 X): 100 mM Tris-HCl (pH 8.0), 100 mM MgCl ₂ , 500 mM NaCl, 10 mM DTT, 500 ug/ml BSA (Bovine Serum Albumin)
Data Link:	UniProtKB/Swiss-Prot P0A7Y4 (RNH_ECOLI)
Storage:	Store at -20°C
References:	1) Gubler U (1987) "Second-strand cDNA synthesis: mRNA fragments as primers." <i>Method Enzymol</i> 152 : 330-335 PMID: 3309563 2) Sambrook J & Russell DW (2001) <i>Molecular Cloning</i> , Chapter 11 "Preparation of cDNA Libraries and Gene Identification". CSHL Press 3) Vournakis JN <i>et al</i> (1975) "Electrophoretic patterns of deadenylylated chorion and globin mRNAs." <i>Proc.Natl.Acad.Sci.USA</i> 72 : 2959-2963 PMID: 1059086 4) Donis-Keller H (1979) "Site specific enzymatic cleavage of RNA." <i>Nucleic Acids Res.</i> 7 : 179-192 PMID: 386279

***Caution;** To avoid contamination of trace amounts of nucleic acids in BSA, use reaction buffer that does not contain BSA and use RNaseH at higher concentrations.




Fig.1 SDS-PAGE of *E. coli* RNaseHE.

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia, Inc.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

http://www.cosmobio.co.jp/index_e.asp

Phone : +81-3-5632-9617

E-mail: export@cosmobio.co.jp

FAX : +81-3-5632-9618