1	Intended Use
2	Summary and Explanation
3	Principle of the Assay
4	Reagents
5	Caution
6	Storage and Stability
7	Specimen Collection and Preparation
8	Methodology
9	Ring Measurement and Result Processing
10	Limitations of Procedure
11	Expected Values
12	Performance Characteristics
13	Bibliography
14	Summary of Procedure
15	RID Reference Table

Deutsch

Français

## HUMAN IgG, IgA & IgM 'NL' BINDARID™ RADIAL IMMUNODIFFUSION KITS

**For *In vitro* Diagnostic Use Only**

**Product Code: RN004.3, RN010.3, RN012.3  
& RK002**

BINDARID™ is a trademark of The Binding Site Ltd., Birmingham, UK

*Product manufactured by:*  
The Binding Site Ltd, PO Box 11712, Birmingham, B14 4ZB, UK  
[www.bindingsite.co.uk](http://www.bindingsite.co.uk)  
Telephone: +44 (0) 121 436 1000  
Fax: +44 (0) 121 430 7061  
e-mail: [info@bindingsite.co.uk](mailto:info@bindingsite.co.uk)

*FDA (USA) Information*  
Analyte ID Code: 2806 (IgG); 2803 (IgA); 2808 (IgM)  
Test System ID Code: 61067  
Complexity Cat: High

### 1 INTENDED USE

This kit is intended for measuring human IgG, IgA or IgM in serum and other biological fluids.

### 2 SUMMARY AND EXPLANATION

#### **Immunoglobulin G**

In normal adults, IgG constitutes approximately 75% of the total serum immunoglobulins. It has a Mol. Wt. of 150kD, and is composed of two gamma heavy chains and two light chains. The majority of IgG is synthesised in plasma cells, and its ability to cross the placenta gives it a crucial protective role in the first weeks of life when the neonatal immune system is not fully developed. IgG is divided into four subclasses, IgG1,2,3 and 4 which show considerable differences in their properties, including the ability to fix complement, to bind to macrophages and to pass through the placenta (refs. 1,2,3). Normal serum levels of IgG vary with age (see Section 11) and also race. Raised serum levels are associated with chronic liver disease and myeloma; reduced levels may be associated with malignancy and other severe protein-losing conditions (ref. 4).

#### **Immunoglobulin A**

IgA is the major immunoglobulin class of sero-mucous secretions, part of the defence system for external body surfaces. The monomeric form is composed of two alpha heavy chains and two light chains, but a dimeric form, secretory IgA (sIgA) also occurs in secretions and colostrum. sIgA contains a secretory component and a J chain linking the two IgA molecules. Two subclasses of IgA have been identified in humans; IgA1, which accounts for 80-90% of total serum IgA, and IgA2, which accounts for 30-50% of total IgA in secretions such as milk. The two subclasses appear to be regulated independently. Normal serum levels of IgA vary with age (see Section 11). Raised IgA serum levels are associated with breast feeding, as well as chronic infections, liver disease and myeloma; reduced levels may be associated with certain protein losing conditions (ref. 4).

#### **Immunoglobulin M**

IgM is the first class of immunoglobulin synthesised in response to particulate antigens. Each unit consists of two mu heavy chains and two light chains, five units together with a J-chain comprising an IgM molecule. IgM is therefore multivalent, and deals most efficiently with polyvalent antigens such as bacteria and viruses; it also activates complement. On active immunisation IgM rapidly appears in the serum, but levels normally drop after a week, usually in parallel with the increase in IgG. Normal serum levels are age dependent (see Section 11). Raised IgM levels are associated with hepatitis, myeloma, Waldenstrom's macroglobulinaemia and other infections; reduced levels can occur in antibody deficiency syndrome (ref. 4).

Radial immunodiffusion (RID) is a technique that is routinely used for measuring the concentration of various soluble antigens in biological fluids. It is principally derived from the work of Fahey & McKelvey (ref. 5) and Mancini et al. (refs. 6 & 7).

### 3 PRINCIPLE OF THE ASSAY

The method involves antigen diffusing radially from a cylindrical well through an agarose gel containing an appropriate monospecific antibody. Antigen-antibody complexes are formed which, under the right conditions, will form a precipitin ring. The ring size will increase until equilibrium is reached between the formation and breakdown of these complexes, this point being termed 'completion'. At this stage, a linear relationship exists between the square of the ring diameter and the antigen concentration. By measuring the ring diameters produced by a number of samples of known concentration, a calibration curve may be constructed. The concentration of the

antigen in an unknown sample may then be determined by measuring the ring diameter produced by that sample and reading off the calibration curve.

There are three different procedures that may be used with this kit (see Section 8.4). Procedures ONE and TWO require that the rings are measured at completion. A linear calibration curve is constructed for Procedure TWO, whereas for Procedure ONE a reference table (based upon the ideal linear calibration curve) is provided, which converts ring diameters directly to protein concentrations. Using Procedure THREE, ring diameters are measured before completion; the calibration curve produced will be non-linear.

## 4 REAGENTS

- 4.1 **RID plates** (supplied in foil pouches). These contain monospecific antibody to IgG, IgA or IgM in agarose gel. Up to fourteen samples can be run per plate (including calibrators). Preservatives: 0.099% sodium azide, 0.1% E-amino-n-caproic acid (EACA), 0.01% thiomersal (sodium ethylmercurithiosalicylate), 0.01% benzamide.
- 4.2 **Calibrator(s)**. These are supplied in stabilised liquid form as a set of three containing high, medium and low concentrations of the immunoglobulin. (**Note:** the combi kit contains only one calibrator for each of IgG, IgA and IgM). The concentrations given on the vial labels have been obtained by comparison with the CRM470 International Reference Material. Preservatives: 0.099% sodium azide, 0.1% EACA, 0.01% benzamide.
- 4.3 **7% Bovine Serum Albumin (BSA) solution**. This is supplied in stabilised liquid form and is included for use as a diluent. Preservatives: 0.099% sodium azide, 0.1% EACA, 0.01% benzamide.
- 4.4 **Control Serum**. This is supplied in stabilised liquid form. The expected concentrations for IgG, IgA and IgM are marked on the vial label. Preservatives: 0.099% sodium azide, 0.1% EACA, 0.01% benzamide.

## 5 CAUTION

All donors of human serum supplied in this kit have been tested and found negative for Hepatitis B surface antigen (HbsAg) and antibodies to human Immunodeficiency virus (HIV1 & 2) and Hepatitis C virus. However, these tests can not guarantee the absence of infective agents. Proper handling and disposal methods should be established as for all potentially infective material and only personnel adequately trained in such methods should be permitted to perform the procedures.

The plates, calibrators, controls and BSA contain 0.099% sodium azide and plates also contain 0.01% thiomersal as a preservative. Handle with caution - do not ingest or allow contact with skin and mucous membranes. If contact does occur wash with a large volume of water and seek medical advice. Explosive metal azides may be formed with the lead and copper plumbing; on disposal of reagent flush with large volumes of water to prevent azide build up.

Adherence to the given procedure is recommended for the purposes stated and validity or results obtained using methods other than those stated cannot be guaranteed.

Reagents from different batch numbers of kits are NOT interchangeable. If large numbers of test are performed care should be taken to ensure that all reagents are from the same batch.

## 6 STORAGE AND STABILITY

The unopened kits should be stored at 2-8°C and can be used until the expiry date given on the kit box label. DO NOT FREEZE. The expiry dates of the individual components are given on the component labels. RID plates should be stored at 2-8°C and are damaged by temperature extremes. Freezing will destroy the gel, therefore RID plates should be kept away from cooling elements in refrigerators. High temperatures should also be avoided as this will result in moisture loss from the gel, affecting performance. Unopened plates should be stored flat and upside down (pouch label uppermost) to prevent condensation accumulating in the wells. Handle plates with care to prevent gel damage.

Unopened calibrators and controls should be stored at 2-8°C. Once opened they are stable for at least one week at 2-8°C, but for longer storage they should be aliquoted and frozen. All other reagents should be stored at 2-8°C.

## 7 SPECIMEN COLLECTION AND PREPARATION

Use fresh or deep frozen serum samples. Microbially contaminated, haemolysed and very lipaemic serum samples or those containing particulate matter should not be used. Blood samples should be collected by venepuncture, allowed to clot naturally and the serum separated as soon as possible to prevent haemolysis. The serum may be stored at 2-8°C for up to 48 hours prior to assay or for prolonged storage aliquoted and kept at -20°C or below. Repeated freezing and thawing should be avoided.

The BSA included in the kit should be used as diluent when required, as this will maintain the viscosity of the sample. Results can therefore be accurately compared with the calibrators which have a similar viscosity to normal serum.

## 8 METHODOLOGY

(A summary of the entire procedure is given at the end of this instruction leaflet)

### 8.1 Materials provided:

Individual Kits: (Codes: **RN004.3, RN010.3, RN012.3**)

1. 3 x *Human IgG/A/M NL Bindarid* (radial immunodiffusion plates in foil pouches)
2. 8 x *Gel Dividers*
3. 3 x *Human IgG/A/M NL Calibrator* (liquid calibrator set)
4. 1 x 5mL *7% BSA Solution*
5. 1 x *Human IgG, A and M NL Control Serum* (liquid)
6. 1 x Instruction leaflet, including RID Reference table

**Human Immunoglobulins Combi Kit (RK002)**

1. 3 x *Human IgG/A/M NL Bindarid* (radial immunodiffusion plates in foil pouches - one plate for each class)
2. 8 x *Gel Dividers*
3. 3 x *Human IgG/A/M NL Calibrator* (liquid calibrators - 1 for each class)
4. 1 x 5mL *7% BSA Solution*
5. 1 x *Human IgG, A and M NL Control Serum* (liquid)
6. 1 x Instruction leaflet, including RID Reference table

### 8.2 Materials required but not provided:

- 8.2.1 Equipment for collection and preparation of test samples, eg sample tubes, centrifuge etc.
- 8.2.2 Pipettes for accurate dilution of samples, when required.
- 8.2.3 Micropipettes for sample application. These should be capable of accurately delivering 5µL volumes. Binding Site Micropipettes (code AD041) or 'Hamilton' syringes are recommended.
- 8.2.4 Jeweller's Eyepiece (Code AD040) or RID plate reader for magnifying and accurately measuring the precipitin ring diameters to 0.1mm.
- 8.2.5 Graph paper.

### 8.3 Reagent Preparation

#### 8.3.1 RID Plate(s)

To avoid contamination of the gel, plates should be used in a dust-free environment. Take the plate from the foil pouch and remove the lid. If condensation is visible the plate should be kept upside down until the lid has been removed to prevent droplets falling onto the gel. Check the plate to ensure that no damage has occurred in storage or transit, eg splits in the gel. Leave the plate open for 10-15 minutes (or longer if necessary) at room temperature to allow any condensation in the wells or on the gel surface to evaporate. Samples should never be applied to wells in which moisture is still visible.

**Plate partitioning:** The plates may be partitioned into up to four sections using the gel dividers provided prior to use. Each divider should be positioned carefully on the gel, cutting edge downward, with the stabilising arm resting on the central plate label. Press firmly on the arm to cut the gel and leave in position.

Plate partitioning is recommended if only part of the plate is to be used initially or when measuring suspected high concentration samples which could (by diffusing over a wide area) result in antibody depletion occurring elsewhere on the plate. After initial use, partitioned plates should be resealed in their foil pouches and stored at 2-8°C with the gel divider(s) in place. Store partitioned plates right side up and use within four weeks.

#### 8.3.2 Calibrator(s)

The calibrator is prediluted and must be mixed gently before use. It should be applied to the plates neat. The medium and low calibrators should only be used when a calibration curve is required, as for procedures TWO and THREE (see below). With the combi kit, only a single 'high' calibrator for each immunoglobulin class is provided. Dilutions of these must be made using the 7% BSA provided if procedure TWO or THREE is being followed.

These dilutions should normally be a medium dilution (60%, ie 6 parts in 10) and a low dilution (10%, ie 1 part in 10). It is recommended that 120µL of calibrator is mixed with 80µL of the 7% BSA for a 60% dilution, and 25µL of calibrator is mixed with 225µL of the 7% BSA for a 10% dilution.

#### 8.3.3 Control(s)

The liquid control serum should be applied to the plates undiluted, mixing gently immediately before use.

#### 8.3.4 Samples

Samples should not normally require dilution. Sera from patients with myeloma and other clinical conditions, however, can contain high levels of immunoglobulin and may require dilution prior to application. In such cases it is suggested that to obtain adequate accuracy a minimum volume of 25µL of test sample is mixed with the appropriate volume of BSA. For samples having immunoglobulin concentrations below the detection limit of the plates, one of the following is recommended:

- i) Concentrate the sample.
- ii) Make a double fill of the well (see Section 8.5).
- iii) Use a kit with lower measuring range if available, (see catalogue).

### 8.4 Procedures

#### 8.4.1 Procedure ONE: RID Reference Table.

This method does not require the construction of a calibration curve - sample concentrations corresponding to each ring diameter are read directly off the RID Reference Table. Rings must be allowed to develop to completion which will require a minimum diffusion time of 48 hours (72 hours for IgM). The high calibrator should be run on each plate used to ensure all are performing correctly.

#### 8.4.2 Procedure TWO: Calibration Curve at Completion.

In this method, all three calibrators (or the appropriate high calibrator plus the two dilutions for combi kits) are used to produce a linear calibration curve. Rings must be allowed to develop to completion which will require a minimum diffusion time of 48 hours (72 hours for IgM). To conserve wells, one calibration curve can be used for several plates of the same batch used concurrently. In such cases, the neat calibrator should be run on each plate used to ensure all are performing correctly.

#### 8.4.3 Procedure THREE: Calibration Curve prior to Completion.

In this method, all three calibrators (or the high calibrator plus the two dilutions for combi kits) are used to produce a calibration curve which is non-linear, as the rings are measured before completion. The minimum recommended diffusion time is 18 hours. It is advisable to construct a separate calibration curve for each plate used.

## 8.5 Application of Calibrators and Samples

The calibrators, control and test samples should be gently mixed immediately before use. Fill the required number of wells with 5µL of the high calibrator using a micropipette. If Procedure TWO or THREE is being followed, also fill the required number of wells with the medium and low calibrators (medium and low dilutions for combi kits). The remaining wells should then be filled with 5µL of appropriately diluted test samples and controls. Plates should not be left open for long periods during calibrator / test sample application, as this will cause excessive drying of the gel.

**Note:** For those samples suspected of containing low concentrations of the specific proteins, a 'double fill' of the well may be made. The well is initially filled with 5µL of the sample and this is allowed to completely diffuse into the gel, which can take up to 30 minutes. The lid should be kept in place during this period. The second fill (again using 5µL) may then be made and the plate incubated as normal. Results obtained must be corrected for the double sample volume and will be less accurate than those obtained by the normal 'single fill' procedure.

## 8.6 Incubation

After sample application, the lid is tightly closed and the plate stored flat with the lid uppermost at room temperature (approximately 20-24°C). It is essential that the gel is not allowed to dry out during incubation. To minimise evaporation, it is suggested that plates should either be resealed in their foil pouches or stored in a moist box (a sealed plastic box containing damp tissue paper) during incubation. The minimum incubation time for Procedure THREE is 18 hours and for complete diffusion (Procedures ONE and TWO) is 48 hours (72 hours for IgM). Final ring diameters may be affected by temperature; the expected ring size for the high calibrator is 8mm ± 0.3mm (7.5mm ± 0.3mm for IgM) when incubated at 20-24°C. Extremes of temperature should be avoided.

## 8.7 Quality Control

The control serum should be treated exactly like test samples. Values obtained for each control should be within ±10% of the concentration stated on the vial label.

## 9 RING MEASUREMENT AND RESULT PROCESSING

After the required diffusion time, ring diameters should be measured to the nearest 0.1mm, using a jeweller's eyepiece or a RID plate reader. When reading with an eyepiece, use bright side lighting and a dark background. If difficulties are experienced, view the plate macroscopically and mark the edges of the rings on the back of the plate using a needle. The distance between these marks may then be more easily measured.

**Note:** For Procedures ONE and TWO ring diameters must have developed to completion. If there is any doubt, rings should be remeasured after a further 24 hours to ensure there has been no increase in their diameters. The high calibrator should give a ring diameter of 8.0mm ± 0.3mm (7.5mm ± 0.3mm for IgM) at completion. If the ring diameter is outside this range, see Trouble Shooting (Section 10.3).

### Procedure ONE

The concentration of immunoglobulin in each test sample can be read directly from the RID Reference Table, **providing it has been applied neat as recommended.**

Concentrations obtained for samples giving ring diameters greater than the high calibrator should be regarded as approximate, due to the possibility of incomplete diffusion; they may also cause local antibody depletion thereby affecting adjacent ring sizes. Such samples should preferably be diluted appropriately and re-tested. Samples giving ring diameters below the lower limit on the RID Reference Table should be retested in a more concentrated form (see Section 8.3.4). Any change from the recommended sample dilution (ie neat) must be taken into account when calculating the results.

### Example

Test Sample	Dilution	Ring Diameter (mm)	Table Value (mg/L)	Original Sample Conc. (mg/L)
IgA Serum A	Neat	6.4	3120	3120
IgA Serum B	Neat	>10.0	>9080	>9080
IgA Serum B	1/2	7.8	5130	10260*

\* Calculated as follows: Table value x Recommended Diln./Actual Diln., ie 5130mg/L x (1)/(1/2).

### Procedure TWO


Plot the square of the diameters of the precipitin rings formed by the three calibrators (or the high calibrator plus the two dilutions for the combi kits) versus their immunoglobulin concentration (given on the calibrator vial labels). Immunoglobulin concentrations should be along the horizontal (x) axis, ring diameters squared along the vertical (y) axis. A line of best fit is drawn through the three points; the y-intercept should be in the range 10-12mm<sup>2</sup>. The immunoglobulin concentration is determined from the calibration curve; remember to adjust the sample concentration obtained by any dilution factor used.

### Sample Calculation:

IgA calibrators gave the following ring diameters on an IgA test plate at completion:

Calibrator	Conc. (mg/L)	Diameter (D) of ring (mm)	D squared (mm <sup>2</sup> )
High	5450	8.0	64.0
Medium	3270	6.5	42.3
Low	545	3.9	15.2

A calibration curve was plotted using these results:

An unknown sample, applied neat as recommended, gave a 6.1mm diameter ring on this plate. From the above curve, this corresponds to an IgA concentration of 2680mg/L.

### Procedure THREE


Plot the calibration curve as for Procedure TWO. The graph will not be a straight line but a curve, the gradient of which decreases with increasing protein concentration. The y-intercept should be as indicated for Procedure TWO. Test sample protein concentrations are read off the calibration curve; remember to adjust the sample concentration obtained by any dilution factor used.

### Sample Calculation:

IgA calibrators gave the following diameters on an IgA test plate after 18 hours.

Calibrator	Conc. (mg/L)	Diameter (D) of ring (mm)	D squared (mm <sup>2</sup> )
High	5450	7.5	56.3
Medium	3270	6.5	42.3
Low	545	3.9	15.2

A calibration curve was plotted using these results:

An unknown sample, applied neat as recommended, gave a 5.9mm ring on this plate. From the above curve, this corresponds to an IgA concentration of 2450mg/L.

## 10 LIMITATIONS OF PROCEDURE

- 10.1** For Procedure ONE, results generated from ring diameters greater than the high calibrator ring diameter (ie 8mm, or 7.5 for IgM) should be regarded as approximate (see Section 9). For Procedures TWO and THREE, accurate results are limited to the calibration curve between the high and low calibrator (low dilution) values – extrapolation beyond these points is not valid. Samples giving results outside these ranges must be diluted or concentrated as appropriate and retested (see Section 8.3.4).
- 10.2** **FDA (USA) Information – see front page**

### 10.3 TROUBLE SHOOTING

Problem	Possible Causes(s)	Suggested Action(s)
<b>A. No ring for:</b>		
1. Calibrator(s)	Calibrator omitted.	Repeat assay.
2. Test sample	i) Sample omitted.	Repeat assay.
	ii) Concentration too high/low	Dilute/concentrate and reassay.

		Suggested Action(s)	
3. Calibrator(s) and test samples	Plate deterioration	a) Storage damage. Repeat assay using new plate. b) Product expired. Repeat assay using new plate/kit.	
<b>B. Oversize rings for:</b>			
1. High calibrator (more than 8.3mm or 7.8 mm for IgM)	i) Inaccurate ring measurement.	Remeasure using eyepiece or RID plate reader.	
	ii) Incorrect volume applied.	Check 5µL volume applied.	
	iii) Inaccurate volume applied.	a) Micropipette malfunction – check operation and repeat assay.	
		b) Poor technique – repeat assay.	
	iv) Partial evaporation of calibrator on storage.	Repeat assay using new calibrator/kit.	
	v) Plate deterioration.	a) Storage damage. Repeat assay using new plate.	
		b) Product expired. Repeat assay using new kit.	
vi) Local antibody depletion due to adjacent high concentration test samples.	Dilute the sample(s) responsible and repeat assay using new plate.		
2. Test samples (above acceptable range – see Section 10.1)	ii) Incorrect volumes applied.	Check 5µL volume applied.	
<b>C. Undersized rings for:</b>			
1. High calibrator (less than 7.7mm or 7.2mm for IgM)	i) Inaccurate ring measurement.	As for B1 above	
	ii) Incorrect volume applied.		
	iii) Inaccurate volume applied.		
	iv) Calibrator deterioration.	a) Storage damage. Repeat assay using new calibrator. b) Product expired. Repeat assay using new kit.	
	v) Incubation temperature too low (see Section 8.6).	Repeat assay, incubating at 20-24°C.	
2. Test samples (below acceptable range – see Section 10.1)	i) Concentration too low.	See section 8.3.4 and repeat assay.	
	ii) Incorrect volume applied.	Check 5µL volume applied.	
<b>D. Double/Multiple rings:</b>			
i) Non-specific precipitation close to well (due to PEG in gel).		Read outer ring.	
	ii) Poor sample application.	Repeat assay.	
	iii) Calibrator deterioration.	a) Storage damage. Repeat assay using new calibrator.	
		b) Product expired. Repeat assay using new kit.	
iv) Sample deterioration.	Reassay using fresh sample.		
<b>E. Non-circular rings:</b>			
i) Poor sample application.		Repeat assay.	
	ii) Gel dried out before use.	a) Storage damage. Repeat assay using new plate.	
		b) Product expired. Repeat assay using new plate/kit.	
iii) Gel dried out during sample application or incubation.	Repeat assay minimising the time the plate is left open. Incubate with lid on tight in a moist box or sealed foil pouch.		
iv) Local antibody depletion (due to high concentration samples on the plate).		Dilute samples and repeat assay.	
	i) Plate has been frozen.	Repeat assay using new plates. Review storage.	
		ii) Gel dried out before use.	As for E(ii) above.
		iii) Gel dried out during sample application or incubation.	As for E(iii) above.
<b>G. Weak, pitted gel:</b>	Plate has been frozen.	Repeat using new plate. Review storage.	

Problem	Possible Causes(s)	Suggested Action(s)
<b>H. Poor calibration curve:</b>		
1. Curve non-linear (Procedure TWO)	i) Incomplete diffusion.	Incubate for further 24 hours and remeasure the rings.
	ii) Calibrator rings under/oversize.	As for B1 or C1 above. (Similar explanations apply to the medium and low calibrators/dilutions).
	iii) Calibration curve constructed incorrectly.	Check calibration curve construction.
2. y-intercept out-of-range (Section 9)	i) Calibrator rings under/oversize.	As for B1 or C1 above. (Similar explanations apply to the medium and low calibrators/dilutions).
	ii) Calibration curve constructed incorrectly	Check calibration curve construction.

**10.4** Diagnosis cannot be made and treatment must not be initiated on the basis of IgG, IgA and IgM measurements alone. Clinical history and other laboratory findings must be taken into account.

**10.5** If an unexpected result is obtained, the assay should be repeated, preferably with a fresh sample.

If a problem cannot be resolved, please refer to supplier.

## 11 EXPECTED RESULTS

		Mean Conc. (g/L)	Normal Range (g/L)	Number of samples
IgG	Male	9.96	5.66 - 14.25	62
	Female	9.47	6.29 - 12.65	62
IgA	Male	2.05	0.45 - 3.64	63
	Female	1.70	0.49 - 2.91	60
IgM	Male	1.06	0.03 - 2.09	64
	Female	1.11	0.35 - 1.89	60

The above normal ranges are based on the mean concentration  $\pm$  2 standard deviations. The data provided has been obtained from limited numbers of British blood donors and is intended for guidance purposes only. It is strongly recommended that each user should generate his/her own reference ranges for appropriate clinical conditions.

Serum immunoglobulin concentrations are age-dependent - the results below were obtained by a different radial immunodiffusion method (ref. 8).

### Serum Concentrations of IgG, IgA & IgM at Different Ages (g/L)

AGE	IgG		IgA		IgM	
	Mean	Range	Mean	Range	Mean	Range
Cord blood	10.93	7.45 - 16.02	0.01	0.00 - 0.08	0.11	0.04 - 0.26
1/2 - 3 Months	4.86	2.93 - 8.06	0.13	0.03 - 0.57	0.51	0.16 - 1.64
3 - 6 Months	3.60	1.39 - 9.34	0.18	0.04 - 0.78	0.49	0.20 - 1.20
6 - 12 Months	6.66	4.10 - 10.81	0.32	0.13 - 0.82	1.04	0.48 - 2.29
1 - 2 Years	6.30	3.49 - 11.39	0.36	0.13 - 1.02	1.03	0.40 - 2.29
2 - 3 Years	7.61	4.82 - 12.00	0.51	0.22 - 1.18	1.06	0.54 - 2.09
3 - 6 Years	8.50	5.53 - 13.07	0.77	0.33 - 1.80	1.11	0.56 - 2.18
6 - 9 Years	9.68	6.46 - 14.51	1.01	0.28 - 2.22	1.13	0.55 - 2.32
9 - 12 Years	9.62	6.13 - 15.12	1.21	0.57 - 2.56	1.41	0.70 - 2.84
12 - 16 Years	9.89	6.67 - 14.64	1.28	0.77 - 2.19	1.13	0.49 - 2.61
Adults	10.99	6.58 - 18.37	1.61	0.71 - 3.60	1.33	0.40 - 2.63

Concentrations in excess of 100g/L can occur in the serum of myeloma patients.

## 12 PERFORMANCE CHARACTERISTICS

### 12.1 Precision

The precision (repeatability) of this kit is expressed as the mean and the percentage coefficient of variation (CV) which had been determined using serum preparations containing high, medium and low concentrations of IgG, IgA and IgM. All analyses were performed in our laboratory. Each value was calculated from 10 measurements (duplicate determinations on five separate plates from a typical batch) unless otherwise stated. For Procedures ONE and TWO, rings were measured after 48 hours (72 hours for IgM). For Procedure THREE rings were measured after 18 hours.

Sample Pool		Procedure 1 Mean		Procedure 2 Mean		Procedure 3 Mean	
		Conc (g/L)	CV	Conc (g/L)	CV	Conc (g/L)	CV
IgG	High	21.7	2.3%	20.9	2.2%	20.1	3.7%
	Medium	13.5	2.0%	12.5	1.5%	13.5	1.8%
	Low	5.0	3.9%	4.1	4.0%	4.3	9.7%
IgA	High	4.84	3.0%	4.9	2.9%	4.56	4.7%
	Medium	2.99	1.9%	2.95	1.6%	3.05	1.9%
	Low	1.16	3.9%	1.00	6.6%	1.00	6.6%
IgM	High	2.45	1.9%	2.46	1.7%	2.48	6.9%
	Medium	1.52	2.7%	1.51	2.3%	1.57	1.3%
	Low	0.53	9.9%	0.54	2.5%	0.55	5.0%

### 12.2 Within plate and inter-batch variation:

The within plate variation is expressed as the mean  $\pm$  standard deviation of determinations of CV made using 3 plates from separate batches. Six measurements were made per plate, using a human serum pool as the sample.

The interbatch variation is expressed as the CV of mean diameter values obtained from recent batches of plates. The mean diameter for each batch was calculated using

the ring diameter at completion obtained using a human serum pool as the sample, applied to two (or more) plates from each batch (six ring measurements per plate).

**15 RID REFERENCE TABLE**

	Within-plate variation	Interbatch variation
	Mean CV % ± SD (N=3)	CV (%) (N=3)
IgG	0.72 ± 0.27	0.12
IgA	0.83 ± 0.56	1.07
IgM	0.72 ± 0.23	1.12

**13 BIBLIOGRAPHY**

- Hamilton, R H (1987). Human IgG subclass measurements in the clinical laboratory. Clin. Chem. **33**, 1707 - 1725.
- Schur, P H (1987), IgG subclasses - a review. Annals of Allergy, **58**, 89-99.
- Vlug, A & Van Remortel, P (1989). The structure and function of human IgG subclasses. Eur. Clin. Lab. **8**(5), 26-35.
- Zilva, J F & Pannall, P R (1984). Clinical chemistry in diagnosis and treatment. Publ. Lloyd-Luke (Medical Books) Ltd., London. pp 348-352.
- Fahey, J L & McKelvey, E M (1965). Quantitative determination of serum immunoglobulins in antibody-agar plates. J. Immunol., **94**, 84-90.
- Mancini, G, Vaerman, J P et al. (1963). Peters H. (ed), Amsterdam, Elsevier Publishing Co., p370-373. Protides of the biological fluids. XI Colloquium.
- Mancini, G, Carbonara, A O et al. (1965). Immunochemical quantitation of antigens by single radial immunodiffusion. Immunochem. **2**, 235-254.
- Lentner, C (ed.) (1984) Geigy Scientific Tables, 8<sup>th</sup> edition, vol 3. Publ. Ciba-Geigy Ltd., Basle, Switzerland. Table 36, p154.

**14 SUMMARY OF PROCEDURE**

- Select Procedure ONE, TWO or THREE. Procedure THREE must be used if results are required quickly.
- Combi Kit only:** Prepare calibrator dilutions if using Procedure TWO or THREE.
- Prepare sample dilutions; this is only required for samples with known high immunoglobulin concentrations.
- Allow condensation to evaporate from RID plate(s).
- Apply calibrator(s), control and samples to RID plate(s) in 5µL volumes.
- Replace lid(s) and incubate at room temperature (20-24°C) for fixed time period - minimum 18 hours (Procedure THREE) or until rings are complete (minimum 48 hours for IgG, IgA; 72 hours for IgM).
- Measure the ring diameters.
- Read results off RID Reference Table (Procedure ONE) or plot calibration curve and read off results (Procedures TWO and THREE).

RID Reference Table for Human Immunoglobulin Concentration (mg/L)

Diameter of ring (mm)	IgG	IgA	IgM
4.0	2500	605	344
4.1	2840	687	390
4.2	3190	773	440
4.3	3540	859	488
4.4	3900	945	537
4.5	4260	1030	586
4.6	4640	1120	639
4.7	5030	1220	692
4.8	5440	1320	749
4.9	5820	1410	802
5.0	6240	1510	859
5.1	6650	1610	916
5.2	7100	1720	977
5.3	7540	1830	1040
5.4	7980	1930	1100
5.5	8430	2040	1160
5.6	8900	2160	1230
5.7	9370	2270	1290
5.8	9850	2380	1360
5.9	10300	2490	1420
6.0	10800	2620	1490
6.1	11300	2740	1560
6.2	11900	2870	1630
6.3	12400	2990	1700
6.4	12900	3120	1770
6.5	13400	3250	1850
6.6	14000	3390	1920
6.7	14500	3520	2000
6.8	15100	3650	2080
6.9	15700	3800	2160
7.0	16200	3930	2230
7.1	16800	4070	2310
7.2	17400	4220	2400
7.3	18000	4370	2480
7.4	18600	4510	2560
7.5	19200	4660	2650
7.6	19900	4820	2740
7.7	20500	4970	2830
7.8	21200	5130	2910
7.9	21800	5290	3000
8.0	22500	5450	3100
8.1	23200	5610	3190
8.2	23800	5770	3280
8.3	24500	5940	3370
8.4	25200	6110	3470
8.5	25900	6280	3570
8.6	26600	6450	3670
8.7	27300	6620	3770
8.8	28100	6800	3870
8.9	28800	6980	3970
9.0	29600	7160	4070
9.1	30300	7350	4180
9.2	31100	7530	4280
9.3	31900	7720	4390
9.4	32600	7910	4490
9.5	33400	8090	4600
9.6	34200	8290	4710
9.7	35000	8480	4820
9.8	35800	8680	4930
9.9	36700	8880	5050
10.0	37500	9080	5160

**Note:** The above values assume that test samples are applied undiluted in 5µL volumes. The high calibrator should give a ring diameter of 8.0 ± 0.3mm (IgG, IgA) or 7.5 ± 0.3mm (IgM) at completion when incubated at 20-24°C.

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