



**MONOCLONAL ANTIBODY**

*For research use only, Not for diagnostic use.*

**Catalog No. RIT-M001**

**Anti KS [Keratan Sulfate] (R-10G)**  
**- lacking oversulfated structures -**

**BACKGROUND**

Monoclonal antibody R-10G recognises Keratan Sulfate lacking oversulfated structures in oligosaccharide segments of Keratan Sulfate glycosaminoglycan chains on hiPS cells. Digestion with Keratanase II, Keratanase or Endo-β-galactosidase removes these epitopes from Keratan Sulfate proteoglycans.

R-10G is useful as a potent tool for the evaluation and standardization of hiPS cells in regenerative medicine.


<b>Product type</b>	Primary antibodies
<b>Immunogen</b>	Human iPS cell (Tic)
<b>Host Species</b>	Mouse (C57BL/6)
<b>Fusion Partner</b>	P3U1
<b>Clone Designation</b>	R-10G
<b>Isotype</b>	IgG1
<b>Source</b>	Ascites
<b>Purification</b>	Affinity purified by Protein G
<b>Form</b>	Liquid
<b>Formulation Buffer</b>	Phosphate buffered saline. *NOTE:PBS doesn't contain preservative.Preservative is added based on the research purpose
<b>Concentration</b>	1 mg / mL
<b>Volume</b>	100 ul
<b>Label</b>	Unlabeled
<b>Specificity</b>	R-10G epitopes are oligosaccharide segments of Keratan Sulfate glycosaminoglycan chains consisting of disaccharide-repeating units with galactose and 6-sulfated <i>N</i> -acetyl-glucosamine, lacking oversulfated structures. Digestion with Keratanase II, Keratanase or Endo-β-galactosidase removes these epitopes from Keratan Sulfate glycosaminoglycan chains. <b>R-10G epitopes are expressed on human iPS/ES cells but not on human EC cells.</b>
<b>Cross species reactivity</b>	Human, Other species have not been tested.
<b>Storage Conditions</b>	Store at -70°C. Aliquot to avoid cycles of freeze / thaw.

<b>Application notes</b>	• <b>Western blotting</b> : 3 µg / mL
Recommended dilutions	• <b>Immunofluorescence</b> : 10 µg / mL

Other applications have not been tested.  
 Optimal dilutions/concentrations should be determined by the end user.


<b>References</b>	1) Kawabe K, Tateyama D, Toyoda H, Kawasaki N, Hashii N, Nakao H, Matsumoto S, Nonaka M, Matsumura H, Hirose Y, Morita A, Katayama M, Sakuma M, Kawasaki N, Furue MK, and Kawasaki T (2013). A novel antibody for human-induced pluripotent stem cells and embryonic stem cells recognizes a type of keratan sulfate lacking oversulfated structures. <i>Glycobiology</i> 23 (3), 322–336. PubMed: <a href="https://pubmed.ncbi.nlm.nih.gov/23154990/">23154990</a>
-------------------	---

[ ANTIBODY CHARACTERIZATION ]



**Fig.1 Screening of R-10G mAb by western blotting.** Tic cell lysates in the complete RIPA buffer were desolved by SDS-PAGE on a 4-15% gradient gel under nonreducing, followed by immunoblot detection with R-10G.

- A :** GelCode\* Blue staining of SDS-PAGE of the Tic cell lysates (10ug protein).  
**B :** Tic cell lysates were analyzed by western blotting with R-10G.  
 The molecular mass markers are shown on the left.



**Fig. 2. Localization of the R-10G and TRA-1-81 epitopes on cultured Tic cells visualized on laser confocal microscopy.** Tic cells cultured on Millicell\* EZ slides were double-stained first with R-10G and Alexa Fluor 488-conjugated secondary (anti-mouse IgG1) antibody, followed by with TRA-1-81 and Alexa Fluor 555-conjugated secondary (anti-mouse IgM) antibody. Cells were observed at two different magnifications: ×100 (upper panel) and ×400 (lower panel).

- (I)** Nomarski imaging. **(II)** Nuclear counterstaining with TO-PRO\*-3. **(III)** Antigens for R-10G (green). **(IV)** Antigens for TRA-1-81 (conventional hiPS marker antibody) (red). **(V)** Merged image of (III) and (IV). **(VI)** Close-up view of V (×400).

GelCode\* is Thermo Fisher Scientific Co. and the subsidiary company's trademark.  
 Millicell\* is a Merck Millipore ®  
 TO-PRO\* is a Molecular Probes ®

*For research use only, Not for diagnostic use.*



**COSMO BIO CO., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9613

FAX : +81-3-5632-9619