

---

Product Manual


# OxiSelect™ Nitrotyrosine ELISA Kit

## Catalog Number

STA-305	96 assays
STA-305-5	5 x 96 assays

**FOR RESEARCH USE ONLY**  
Not for use in diagnostic procedures


---



**CELL BIOLABS, INC.**  
*Creating Solutions for Life Science Research*

## **Introduction**

The modification of tyrosine residues in proteins to 3-nitrotyrosine by peroxynitrite (Figure 1) or other potential nitrating agents has been detected in biological systems that are subject to oxidative stress. Detection of nitrotyrosine-containing proteins has been reported in many human and animal diseases or cellular models of disease. While all tyrosine residues in proteins may theoretically be targets for nitration, presumably the efficiency of tyrosine nitration is dependent on various biological conditions such as the local production and concentration of the reactive species, the existence and availability of antioxidants and scavengers, the accumulation of inflammatory cell and the presence of pro-inflammatory cytokines, as well as the proximity and compartmentation of these components.



**Figure 1. 3-Nitrotyrosine Formation**

Cell Biolabs' Nitrotyrosine ELISA Kit is a competitive enzyme immunoassay developed for rapid detection and quantitation of 3-nitrotyrosine in protein sample. The quantity of 3-nitrotyrosine in protein sample is determined by comparing its absorbance with that of a known nitrated BSA standard curve. The kit has a nitrotyrosine detection sensitivity range of 20 nM to 8.0  $\mu$ M. Each kit provides sufficient reagents to perform up to 96 assays, including standard curve and unknown protein samples.

## **Assay Principle**

The nitrotyrosine quantitation kit is a competitive ELISA. The unknown protein nitrotyrosine sample or nitrated BSA standards are first added to a nitrated BSA preabsorbed EIA plate. After a brief incubation, an anti-nitrotyrosine antibody is added, followed by an HRP conjugated secondary antibody. The protein nitrotyrosine content in unknown sample is determined by comparing with a standard curve that is prepared from predetermined nitrated BSA standards.

## **Related Products**

1. STA-303: OxiSelect™ Nitrotyrosine Immunoblot Kit
2. STA-310: OxiSelect™ Protein Carbonyl ELISA Kit

3. STA-318: OxiSelect™ AOPP Assay Kit
4. STA-332: OxiSelect™ MDA ELISA Kit
5. STA-334: OxiSelect™ HNE Adduct ELISA Kit
6. STA-320: OxiSelect™ Oxidative DNA Damage ELISA Kit (8-OHdG)
7. STA-325: OxiSelect™ Oxidative RNA Damage ELISA Kit (8-OHG)

### **Kit Components**

1. Nitrotyrosine Coated EIA Plate (Part No. 230501): one strip well 96-well plate.
2. Anti-Nitrotyrosine Antibody (Part No. 230502): One 20  $\mu$ L vial of anti-nitrotyrosine Rabbit IgG.
3. Secondary Antibody, HRP Conjugate (Part No. 231003): One 20  $\mu$ L vial.
4. Assay Diluent (Part No. 310804): One 50 mL bottle.
5. 10X Wash Buffer (Part No. 310806): One 100 mL bottle.
6. Substrate Solution (Part No. 310807): One 12 mL amber bottle.
7. Stop Solution (Part. No. 310808): One 12 mL bottle.
8. Nitrated BSA Standard (Part No. 230503): One 500  $\mu$ L vial of 1 mg/mL Nitrated BSA in PBS with a nitrotyrosine content of 40  $\mu$ M (2.7 mole of nitrotyrosine per mole of BSA). The protein nitrotyrosine level is predetermined by a spectrophotometric method as described by Ischiropoulos et al (See Ref. 3).

### **Materials Not Supplied**

1. Protein samples such as purified protein, plasma, serum, cell lysate
2. 10  $\mu$ L to 1000  $\mu$ L adjustable single channel micropipettes with disposable tips
3. 50  $\mu$ L to 300  $\mu$ L adjustable multichannel micropipette with disposable tips
4. Multichannel micropipette reservoir
5. Microplate reader capable of reading at 450 nm (620 nm as optional reference wave length)

### **Storage**

Upon receipt, aliquot and store the Nitrated BSA Standard at -20°C to avoid multiple freeze/thaw cycles. Store all other kit components at 4°C until their expiration dates.

### **Preparation of Reagents**

- 1X Wash Buffer: Dilute the 10X Wash Buffer Concentrate to 1X with deionized water. Stir to homogeneity.

- Anti-Nitrotyrosine Antibody and Secondary Antibody: Immediately before use dilute the Anti-Nitrotyrosine Antibody 1:1000 and Secondary Antibody 1:1000 with Assay Diluent. Do not store diluted solutions.

## **Preparation of Standard Curve**

Prepare a dilution series of nitrated BSA standards in the nitrotyrosine concentration range of 0 nM – 8000 nM by diluting the Nitrated BSA stock solution in Assay Diluent (Table 1).

<b>Standard Tubes</b>	<b>Nitrated BSA Standard (μL)</b>	<b>Assay Diluent (μL)</b>	<b>Nitrated BSA (μg/mL)</b>	<b>Nitrotyrosine (nM)</b>
1	60	240	200	8000
2	100 of Tube #1	300	50	2000
3	100 of Tube #2	300	12.5	500
4	100 of Tube #3	300	3.125	125
5	100 of Tube #4	300	0.78	31.25
6	100 of Tube #5	300	0.195	7.81
7	100 of Tube #6	300	0.049	1.95
8	0	300	0	0

**Table 1. Preparation of Nitrated BSA Standards**

## **Assay Protocol**

1. Prepare and mix all reagents thoroughly before use. Each protein sample including nitrated BSA and blank should be assayed in duplicate.
2. Add 50 μL of unknown protein sample or nitrated BSA standard to the wells of the EIA plate. Incubate at room temperature for 10 minutes on an orbital shaker.
3. Add 50 μL of the diluted anti-nitrotyrosine antibody to each well, incubate at room temperature for 1 hour on an orbital shaker.
4. Wash microwell strips 3 times with 250 μL 1X Wash Buffer per well with thorough aspiration between each wash. After the last wash, empty wells and tap microwell strips on absorbent pad or paper towel to remove excess 1X Wash Buffer.
5. Add 100 μL of the diluted Secondary Antibody-Enzyme Conjugate to all wells.
6. Incubate at room temperature for 1 hour on an orbital shaker.
7. Wash microwell strips 3 times according to step 4 above. Proceed immediately to the next step.
8. Warm Substrate Solution to room temperature. Add 100 μL of Substrate Solution to each well, including the blank wells. Incubate at room temperature on an orbital shaker. Actual incubation time may vary from 2-30 minutes.

*Note: Watch plate carefully; if color changes rapidly, the reaction may need to be stopped sooner to prevent saturation.*

9. Stop the enzyme reaction by adding 100  $\mu$ L of Stop Solution into each well, including the blank wells. Results should be read immediately (color will fade over time).
10. Read absorbance of each microwell on a spectrophotometer using 450 nm as the primary wave length.

### Example of Results

The following figures demonstrate typical Nitrotyrosine ELISA results. One should use the data below for reference only. This data should not be used to interpret actual results.




Figure 2: Nitrotyrosine ELISA Standard Curve.




Figure 3: Protein Nitration by tetranitromethane. STO (MEF) cells were lysed in 25mM HEPES, pH 7.5, 150 mM NaCl, 1% NP-40, 10 mM MgCl<sub>2</sub>, 1 mM EDTA, 2% Glycerol. Cell Lysate was nitrated with tetranitromethane (TNM). The protein 3-nitrotyrsonne levels were determined as described in the assay instructions.

## **References**

1. Gow A. J., Farkouh C. R., Munson D. A., Posencheg M. A and Ischiropoulos H. (2004) *Am J Physiol Lung Cell Mol Physiol.* **287**, L262-8.
2. Bian, K., Ke, Y., Kamisaki, Y. and Murad, F. (2006) *J Pharmacol Sci Vol.* **101**, 271-279.
3. Ischiropoulos, H., Zhu, L., Chen, J., Tsai, M., Martin, J. C., Smith, C. D. and Beckman, J. S. (1992) *Archiv. Biochem. Biophys.* **298**, 431-437.

## **Recent Product Citations**

1. Cheah, F.-C. et al. (2008). Airway inflammatory cell responses to intra-amniotic lipopolysaccharide in a sheep model of chorioamnionitis. *Am. J. Physiol. Lung Cell Mol. Physiol.* **296**:L384-L393.
2. Li, X. et al. (2008). Lipoamide protects retinal pigment epithelial cells from oxidative stress and mitochondrial dysfunction. *Free Radic. Biol. Med.* **44(7)**:1465-1474.

## **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

## **Contact Information**

Cell Biolabs, Inc.  
7758 Arjons Drive  
San Diego, CA 92126  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2007-2009: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.