## **SAFETRANS ( $\alpha$ -CDE)**

### **BACKGROUND**

Gene and oligonucleotide therapies are emerging as a potential strategy for the treatment of genetic diseases, cancers, cardiovascular diseases and infectious diseases. Recently, polyamidoamine (PAMAM) dendrimer (generation (G) 3) conjugate with  $\alpha$ -cyclodextrin ( $\alpha$ -CDE) has been known as novel DNA, shRNA and siRNA carriers.  $\alpha$ -CDE works as a superior nucleic acids carrier without severe cytotoxicity *in vitro* and *in vivo*.

**Storage**                      Stored at -20°C

**Package**                     2 mg

**pDNA Transfection**      【In the case of charge ratio 100 (carrier/pDNA)】

**Protocol**

1. Seed the cells ( $1 \times 10^4 \sim 2 \times 10^5$  cells/well) on a 24-well-plate, and cultured for 1 day at 5% CO<sub>2</sub> incubator. (40~60 % confluent should be set for transfection.)
2. Add the serum-free medium 200 uL/well into 1.5 mL tube, and add the pDNA (2 ug/well) dissolved in TE buffer (pH 7-8)
3. Add the  $\alpha$ -CDE Transfection reagent (1 mg/mL) 189.3 uL into pDNA solution of 2 and mixing by pipetting (5~6 times) or slight vortexing (10 s).
4. Incubation for 10-15 min at room temperature to form the complex with pDNA.
5. During the complexation in 4, cells were washed with PBS or serum free medium.
6. Add the serum-free medium 200 uL/well and pDNA complex solution (1 ug/uL) prepared in 3, then incubated for 1-3 h at 5% CO<sub>2</sub> incubator<sup>1)</sup>.
7. After removal of the medium containing pDNA complex, washed with PBS or serum-free medium for 1 time.
8. Add the culture medium (containing serum) and further incubated for 24 h, and determine the gene expression.

**1) Incubation time can be changed for customer's use. Negligible cytotoxicity should be observed until 24 h incubation.**

**Note: We recommend for the use of adhesion cells cultured in 24 well-plates in this transfection protocol. In transfection with other charge ratios, please see the content of pDNA and the transfection reagent indicated the table 1.**

**Table 1. Examples of Transfection**

**【In the case of transfection with pDNA (2 ug)/well】**

Charge ratio	$\alpha$ -CDE Transfection reagent (ug/well)
20	37.9 ug
50	94.6 ug
100	189.3 ug

#### **siRNA Transfection**

**【In the case of a charge ratio of 100 (carrier/siRNA)】**

#### **Protocol**

1. Seed the cells ( $1 \times 10^4 \sim 2 \times 10^5$  cells/well) on a 24-well-plate, and cultured for 1 day at 5% CO<sub>2</sub> incubator. (40~60 % confluent should be set for transfection.)
2. Add the serum-free medium 50 uL/well into tube, and add the siRNA (0.4 ug/well) dissolved in RNase-free water.
3. Add the  $\alpha$ -CDE Transfection reagent (1 mg/mL) 37.9 uL into siRNA solution of 2 and mixing by pipetting (2~3 times).
4. Incubation for 10-15 min at room temperature to form complex with siRNA.
5. During the complexation in 4, cells were washed with PBS or serum free medium.
6. Add the serum-free medium 220 uL/well and siRNA complex solution prepared in 3, then incubated for 1-3 h at 5% CO<sub>2</sub> incubator<sup>1)</sup>.
7. Add the serum (30 uL/well) and further incubated for 24-48 h, and determine the RNAi effects.

**1) Incubation time can be changed for customer's use. Negligible cytotoxicity should be observed until 24 h incubation.**

**Note: We recommend for the use of adhesion cells cultured in 24 well-plates in this transfection protocol. In transfection with other charge ratios, please see the content of pDNA and transfection reagents indicated the table 2.**

**Table 2. Example of Transfection**

【In the case of transfection with siRNA (0.4 ug )/well】

Charge ratios	$\alpha$ -CDE Transfection reagent (ug/well)
20	7.6 ug
50	18.9 ug
100	37.9 ug

**shRNA Transfection  
Protocol (example for  
pSilencer)**

【In the case of Charge ratio 100 (carrier/pSilencer (shRNA expressing vector))】

1. Seed the cells ( $1 \times 10^4 \sim 2 \times 10^5$  cells/well) on a 24-well-plate, and cultured for 1 day at 5% CO<sub>2</sub> incubator. (40~60 % confluent should be set for transfection.)
2. Add the serum-free medium 200 uL/well into 1.5 mL tube, and add the pSilencer (2 ug/well) dissolved in TE buffer (pH 7-8)
3. Add the  $\alpha$ -CDE Transfection reagent (1 mg/mL) 189.3 uL into pSilencer solution of 2 and mixing by pipetting (5~6 times) or slight vortexing (10 s).
4. Incubation for 10-15 min at room temperature to form complex with pSilencer.
5. During the complexation in 4, cells were washed with PBS or serum free medium.
6. Add the serum-free medium 200 uL/well and  $\alpha$ -CDE/pSilencer complex solution prepared in 3, then incubated for 1-3 h at 5% CO<sub>2</sub> incubator<sup>1)</sup>.
7. Add the serum (22.2 uL/well) and further incubated for 24-48 h, and evaluate the RNAi effects.

**1) Incubation time can be changed for customer's use. Negligible cytotoxicity should be observed until 24 h incubation.**

**Note: We recommend for the use of adhesion cells cultured in 24 well-plates in this transfection protocol. In transfection with other charge ratios, please see the content of pSilencer and transfection reagents indicated the table 3.**

**Table 3. Examples of Transfection**

**【In the case of transfection with pSilencer (2 ug)/well】**

Charge ratio	$\alpha$ -CDE Transfection reagent (ug/well)
20	37.9 ug
50	94.6 ug
100	189.3 ug

**Reverse Transfection Protocol**

1. Add the serum-free medium 200  $\mu$ L/well into tube, and add the siRNA or miRNA (0.4 ug/well) dissolved in RNase-free water.
2. Add the  $\alpha$ -CDE Transfection reagent 37.9  $\mu$ L (stock conc. 1 mg/mL) into siRNA or miRNA solution of 1 and mixing by pipetting (2~3 times).
3. Incubation for 10-15 min at room temperature to form complex with siRNA or miRNA.
4. Add  $\alpha$ -CDE Transfection reagent/siRNA or miRNA complexes (7.47 ug) 200  $\mu$ L/well.
5. Add the cell suspension ( $5 \times 10^4$  cells/mL) 200  $\mu$ L/well, and incubate for 1-3 h at 5% CO<sub>2</sub> incubator.<sup>1)</sup>
6. Exchange the culture medium (containing 10% serum) 500  $\mu$ L, and further incubate for 21 h at 5% CO<sub>2</sub> incubator.
7. Assay the transfected cells.

**1) Incubation time can be changed for customer's use. Negligible cytotoxicity should be observed until 24 h incubation.**

**Note: We recommend for the use of adhesion cells cultured in this transfection protocol.**

**Trouble Shooting**

**【Low transfection efficient】**

- ① Charge ratio between  $\alpha$ -CDE Transfection reagent and pDNA may not be optimized.

When a charge ratio between  $\alpha$ -CDE Transfection reagent and pDNA is not optimized, the net charge of the complex provides negative or far positive, resulting in the low transfection efficiency. In that case, please optimize a charge ratio between  $\alpha$ -CDE Transfection reagent and pDNA.

- ② The amount of the complex of  $\alpha$ -CDE/pDNA may be not enough.

In the case of low transfection efficiency and no cytotoxicity of the complex, please increase the amount of the complex of  $\alpha$ -CDE/pDNA.

- ③ The medium containing serum was used during the mixing of  $\alpha$ -CDE Transfection reagent and pDNA.

Please use serum-free medium during the complex preparation.

- ④ It took time to use the complex of  $\alpha$ -CDE Transfection reagent and pDNA.

Please perform transfection after complex formation as soon as possible.

- ⑤ Incubation time for transfection may not be optimized.

Optimized incubation time is dependent on the kind of cells. Please optimize incubation time for each experiment.

- ⑥ Problem in vectors.

The factors such as promoter, origin of replication and size of pDNA may sometimes affect gene expression.

- ⑦ Cell density may not be proper. 40~60 % confluent should be beset for transfection.

- ⑧ Purity of pDNA may not be high.

- ⑨ Problem in reporter assay. Please perform the reporter assay with a positive control.

### **【Very high cytotoxicity】**

- ① Incubation time with the complex of  $\alpha$ -CDE Transfection reagent and pDNA may be too long.

- ② The amount of the complex of  $\alpha$ -CDE Transfection reagent and pDNA may be too much.

- ③ Stress for cells. Washing process should be done quickly and gently.

- ④ Long time culture with serum-free medium.

After transfection for 2-3 h, further incubation should be done with medium containing 10% serum.

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- ① Same cell density should be used in each experiment.
- ② Contamination of mycoplasma.
- ③ Cell passage may be too high.
- ④ Different qualification of serum.

**【Precipitation was observed when mixing with pDNA】**

- ① Low purification of pDNA
- ② Mixing ratio of  $\alpha$ -CDE Transfection reagent and pDNA may not be proper.
- ③ The composition of the solution dissolving pDNA may not be proper.
- ④  $\alpha$ -CDE Transfection reagent may be expired.

**Note**

This reagent is for research use only. Do not use in humans or diagnostic procedures.

**Reference**

1. Kihara F, Arima H, Tsutsumi T, Hirayama F, Uekama K. In Vitro and In Vivo Gene Transfer by an Optimized  $\alpha$ -Cyclodextrin Conjugate with Polyamidoamine Dendrimer. *Bioconjug. Chem.*, 14, 342-350, 2003.
2. Kihara F, Arima H, Tsutsumi T, Hirayama F, Uekama K. Effects of Structure of Polyamidoamine Dendrimer on Gene Transfer Efficiency of the Dendrimer Conjugate with  $\alpha$ -Cyclodextrin. *Bioconjug. Chem.*, 13, 211-219, 2002.
3. Arima H, Kihara F, Hirayama F, Uekama K. Enhancement of Gene Expression by Polyamidoamine Dendrimer Conjugates with  $\alpha$ -,  $\beta$ -,  $\gamma$ -Cyclodextrins. *Bioconjug. Chem.*, 12, 476-484, 2001.
4. Tsutsumi T, Hirayama F, Uekama K, Arima H. Potential use of Polyamidoamine Dendrimer/ $\alpha$ -Cyclodextrin Conjugate (generation 3, G3) as a Novel Carrier for Short Hairpin RNA-Expressing Plasmid DNA. *J. Pharm. Sci.*, 97, 3022-3034, 2008.
5. Tsutsumi T, Hirayama F, Uekama K, Arima H. Evaluation of polyamidoamine dendrimer/ $\alpha$ -cyclodextrin conjugate (generation 3, G3) as a novel carrier for small interfering RNA (siRNA). *J. Control. Release*, 119, 349-359, 2007.
6. Tsutsumi T, Hirayama F, Uekama K, Arima H. Potential use of polyamidoamine dendrimer/ $\alpha$ -cyclodextrin conjugate (generation 3, G3) as a novel carrier for short hairpin RNA-expressing plasmid DNA, *J. Pharm. Sci.*, 97, 3022-3034, 2008.

*For research use only. Not for clinical diagnosis.*

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1. Kihara F, Arima H, Tsutsumi T, Hirayama F, Uekama K. In Vitro and In Vivo Gene Transfer by an Optimized  $\alpha$ -Cyclodextrin Conjugate with Polyamidoamine Dendrimer. Bioconjug. Chem., 14, 342-350, 2003.
2. Kihara F, Arima H, Tsutsumi T, Hirayama F, Uekama K. Effects of Structure of Polyamidoamine Dendrimer on Gene Transfer Efficiency of the Dendrimer Conjugate with  $\alpha$ -Cyclodextrin. Bioconjug. Chem., 13, 211-219, 2002.
3. Arima H, Kihara F, Hirayama F, Uekama K. Enhancement of Gene Expression by Polyamidoamine Dendrimer Conjugates with  $\alpha$ -,  $\beta$ -,  $\gamma$ -Cyclodextrins. Bioconjug. Chem., 12, 476-484, 2001.
4. Tsutsumi T, Hirayama F, Uekama K, Arima H. Potential use of Polyamidoamine Dendrimer/ $\alpha$ -Cyclodextrin Conjugate (generation 3, G3) as a Novel Carrier for Short Hairpin RNA-Expressing Plasmid DNA. J. Pharm. Sci., 97, 3022-3034, 2008.
5. Tsutsumi T, Hirayama F, Uekama K, Arima H. Evaluation of polyamidoamine dendrimer/ $\alpha$ -cyclodextrin conjugate (generation 3, G3) as a novel carrier for small interfering RNA (siRNA). J. Control. Release, 119, 349-359, 2007.
6. Tsutsumi T, Hirayama F, Uekama K, Arima H. Potential use of polyamidoamine dendrimer/ $\alpha$ -cyclodextrin conjugate (generation 3, G3) as a novel carrier for short hairpin RNA-expressing plasmid DNA, J. Pharm. Sci., 97, 3022-3034, 2008.

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COSMO BIO CO., LTD.