

The RI Assay of anti Androstenedione-11Alpha-Succ-BSA according to FKA-106

The Cross reactivities of anti Androstenedione-11Alpha-Succ-BSA (FKA-106) by RI Assay

Antigen	Androstenedione-Succ-BSA	
Antibody	FKA-106 Serum	Freeze-drying of 200µl / Vial
Dilution of Serum	2000test for RI assay	250µl/tube use
Binding of Total	60%	
Androstenedione dose of B0/2	200 pg	
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.	

Steroid

Cross Reaction

Androstenedione	100 %
Androstenetrione	200 %
5Alpha- Androstanedione	32.30 %
DHA	5.00 %
Progesterone	3.90 %
Testosterone	2.90 %
Cortisol	0.20 %
Corticosterone	0.20 %
Estradiol	0.01 %

For research use only; not for use as a diagnostic.

COSMO BIO CO.,LTD

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
 TEL: +81-3-5632-9605 FAX: +81-3-5632-9614 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp