
Medaka Vitellogenin (Vg) ELISA Kit


In recent years the influence of Endocrin-Disrupting-Chemicals (EDC) upon living organisms has been one of the hot topics in environmental study. However, determination of the chemical compound species that have EDC-like influence has not yet been completed, and it is absolutely necessary to take actions to complete such a task.

Vitellogenin (Vg), a precursor of yolk protein, is a female-specific protein that is expressed only in the blood of higher vertebrate oviparous organisms. Vg is normally expressed in liver under regulation of estrogen (female sex hormone) during egg yolk formation process, and when it reaches ovum via blood stream it contributes to the egg yolk formation. Since estrogen has been known to have an effect upon male species to induce unusual expression of Vg, an ELISA screening method detecting the Vg in EDC-exposed samples has become the most attractive bioscreening assay system.

[Principle]

This kit is based on the ELISA method which is used with specific anti- medaka polyclonal antibody. Anti- medaka polyclonal antibody is coated on the microplate, so that Vg in the sample and standard is captured. After reacting biotin-conjugated anti-medaka polyclonal antibody and HRP-conjugated streptavidin, the complex of coated-antibody—antigen—antibody is formed. Finally, you could measure the Vg concentration depend on the coloring strength which has developed with enzymatic reaction.

[Kit Contents]

(1) Antibody-coated microtiter plate (8 wells×12)	1 plate
(2) The standard, Medaka Vg, lyophilized	200µl × 2
(3) Dilution solution	30mL × 2
(4) Biotin conjugated anti-medaka Vg antibody (×100)	80µl × 1
(5) HRP-conjugated streptavidin (×500)	20µl × 1
(6) OPD (o-phenylenediamin) tablets	2 tab.
(7) Substrate buffer	15 mL × 1
(8) Stop solution	6 mL × 1
(9) Wash solution (×20)	30 mL × 1

[Required Apparatuses]

- (1) A microplate reader
- (2) A micropipet
- (3) A microplate washer

[Assay Method]

- (1) Preparation of reagents to be used

- ① Wash solution

Make sure that the solution does not contain any crystallized material, upon reaching room temperature. Working solution is prepared by dilution of 30 mL of the stock wash solution with distilled deionized water to 570 mL. This solution is stable for 14 days when stored in a refrigerator.

- ② Antibody-coated microtiter plate

Apply 300µL in each well and wait another 10 to 60 minutes.

- ③ Standard dilution

Dilute the concentrated standard solution (256 ng/mL) to make 256, 64, 16, 4, 1 ng/mL with Dilution solution.

- ④ Biotin conjugated anti-medaka Vg polyclonal antibody (×100)

The stock Biotin conjugated anti-medaka Vg polyclonal antibody can be prepared for 96 well reaction upon dilution of 60µL stock solution with 6mL Dilution solution just before use.

- ⑤ HRP-conjugated Streptavidin (×500)

The stock HRP-conjugated Streptavidin can be prepared for 96 well reaction upon dilution of 12µL stock solution with 6mL Dilution solution just before use.

- ⑥ Coloring solution

Dissolve two OPD tablets with 13 mL of Substrate buffer.

- (2) Preparation of samples to be used

The homogenate of Medaka liver should be diluted more than 10 times with Dilution solution.

- (3) Assay procedure

- ① Discard the wash solution from the wells completely, but do not keep wells in dry condition for reliable data.



- ② Apply 50µL of each concentration of standard solution and samples into the wells and incubate for 2 hours.

- ③ After incubation, discard the reaction solution and wash with 300µL wash solution. Repeat this washing procedure another 2 times.

- ④ Apply 50 μ l Biotin conjugated anti-medaka Vg polyclonal antibody and incubate for 1 hour.
- ⑤ After the incubation, wash the well 3 times with wash solution (See ③).
- ⑥ Apply 50 μ l HRP-conjugated streptavidin and incubate for 1 hour.
Substrate buffer should be restored to room temperature from refrigerator at this time.
- ⑦ After the incubation, wash the well 3 times with wash solution (See ③).
- ⑧ Apply 100 μ l Coloring solution into each well for 5 minutes in room temperature.
- ⑨ Apply 50 μ l of Stop solution for stopping the enzymatic reaction.
- ⑩ Read absorbance with a microplate reader at 490nm or 492nm.
- ⑪ Measure the medaka Vg concentration using standard curve.

[Standard curve]

[Reproducibilities]

Inter assay (n=20)

Sample	A	B
Mean (ng/mL)	30.1	14.0
SD	1.4	0.4
CV (%)	4.7	3.1

Intra assay (n=20)

Sample	A	B
Mean (ng/mL)	32.5	14.8
SD	1.9	1.2
CV (%)	6.0	8.4

[Specificity: cross reactivity]

	Medaka	Carp	Mummichog
Immunoblotting	(+)	(-)	(+)

[Recovery test]

In the recovery study, the recovery percentages 106 was obtained for 10 times dilutions of the liver homogenate sample.

[Usage notes]

- ① The Reagents should be stored at recommended temperature.
The standard, medaka Vg: -80°C , Other reagents: $2\sim 8^{\circ}\text{C}$
- ② Do not use the reagents which is expired the date of usage.
- ③ Do not leave the standard, medaka Vg, for long time under room temperature.
- ④ Biotin conjugated anti-medaka Vg polyclonal antibody and HRP-conjugated streptavidin is not stable under sun light, so that do not leave it for long time under sun light.
- ⑤ The apparatus for making coloring solution should be clean one.
- ⑥ OPD (o-phenyldiamin) is harmful, so that take good care to it.
- ⑦ Stop solution, 2N H_2SO_4 , is strong acid, so that take good care to it.
- ⑧ The kit is constructed with well-adjusted combination for each lot. Replace combination among different lots may cause unexpected results.
- ⑨ This kit is only for research use. Do not use for medicine or any other things.
- ⑩ When using the reagents, take care to avoid them from touching to skin, mucous membrane, clothes, and getting into eye.
- ⑪ If the reagents happen to get into eye or mouth, wash out them and go to see a doctor if you need.
- ⑫ After using the kit, wash your hand very carefully.
- ⑬ If you find that the packages of the reagents are broken or something trouble, do not use them.
- ⑭ When you store the reagents, make sure to avoid them from vaporizing, falling down.
- ⑮ After using the reagents, the packages should be discarded under established rule.
- ⑯ We can not grantee for any other usage of the packages and the accompaniments.

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