

## **PURIFIED RABBIT ANTI-PORCINE TYPE I COLLAGEN**

(For research use only)

Catalog n° 20191

**Description** Affinity chromatography purified rabbit antibody to type I collagen

extracted, purified from porcine tendon.

Format Purified, freeze-dried antibody in 0.5ml vials. Reconstitute with

0.5ml diH2O and store aliquots (≥100µl) at -20°C.

Stability Unreconstituted: 24 months at -20°C; avoid freeze-thaw cycles

by aliquoting before use and storage at -20°C.

Reconstituted: 6 months (mini) at -20°C.

**Applications** IFA, IHC(p), ELISA, (SP)RIA (tested).

SpecificityPorcine collagen, type I:100(%; at a 1: 200Porcine collagen, types II, IV, V:< 0.1</th>RIA dilution)Porcine collagen, type III:< 1.0</th>Porcine elastin, fibronectin:< 0.1</th>Human, mouse collagens, type I:< 0.1</th>Bovine, rat collagens, type I:< 1.0</th>

**Use** For the immunostaining of extra or intracellular components in light

microscopy and the quantitation of porcine type I collagen in

biological fluids or culture supernatants.

Working dilutions . IFA\* (with fluorescein anti-rabbit IgG conjugate); use diluted at

≥ 1:80 on frozen porcine tissues (skin, liver ...).

. IHC(p)\* (immunoperoxidase technique with the Envision / Dako kit\*); use diluted at ≥ 1:500 on fixed, paraffin-embedded porcine

tissues (skin, liver ...).

**. ELISA**; use diluted at  $\geq$  1:200 (OD  $\geq$  500).

**References** (177) Rousseau, Gagnieu, Biomaterials, 2002, 23, 1503-1510

<sup>\*:</sup> Our recommended protocols [IHC(f) / IFA(2), 4-06/04; IHC(p) / Envision, 4-01/04) are available on request.