FASTKIT™ Immunochromato PEANUT <Instruction Manual>

※ Please read this instruction before use.

Background

Amendments to Food Sanitation Law in Japan have rendered it mandatory to label 5 specific raw food materials, (Eggs, Milk, Wheat, Buckwheat and Peanuts), identified to have a high possibility of triggering food allergies, as of April 2002. Verification of the Manufacturer's Data Report and verification of contamination of specified ingredients by testing is critical in order to carry out labeling correctly.

The FASTKIT IMMNOCHOMATO is a detection kit for peanut proteins in food using immunochromatography, with simple operation and results obtainable in a short period of time.

[Characteristics]

- 1) The FASTKIT Immunochromato can simultaneously identify the presence of multiple peanut proteins.
- 2) The FASTKIT Immunochromato is widely applicable to samples ranging from raw ingredients to processed foods.
- 3) One step simple operation; no background expertise necessary.
- 4) Results can be quickly obtained (15 mins).

Components of the Kit

| | |
|---|-----------------|
| A : Test plate | 20tests |
| B : Dilution buffer | 50mL× 1 bottle |
| C : Extraction/concentration buffer (×10 concentrate) | 100mL× 1 bottle |

Purpose and Performance of the Kit


- Detection of peanut proteins contained in food samples.
- The Minimum detection sensitivity is several ppm, expressed as concentration in food (calculated from measured values that have been obtained using "FASTKIT Elisa PEANUT").*

Screening Principles


When the sample solution is dropped onto the sample dropper area of the test plate, colloidal gold labeled rabbit anti-peanut antibodies (labeled (2) in diagram) are dissolved, creating a complex with the peanut proteins (labeled (1) in diagram) contained in the sample solution. These complexes move to the development area by capillary action, and are trapped by the anti-peanut antibodies (labeled (3) in diagram) fixed at interpretation area (T), with the colloidal gold creating a red-purple colored line at interpretation area (T). Interpretation in this Kit is carried out by observation, which interprets the presence of peanut proteins contained in the sample solution.

On the other hand, regardless of the presence of peanut proteins contained in the sample, the excess colloidal gold labeled antibodies move to the development area and are trapped by the anti-rabbit immunoglobulin antibodies (labeled (4) in diagram) fixed at interpretation area (C), with the colloidal gold creating a red-purple colored line at interpretation area (T). Normal movement of the sample solution to the development area is verified by a red-purple colored line at interpretation area (C).

<Separate parts of the test plate>

<Principle diagram>

Procedure 1 (Preparation of Reagents)

[Equipment]

- A homogenizer (a food cutter, a grinder or a mill)
- A centrifuge (one which can be operated at higher than 3,000 x g and 4°C is preferable), centrifuge tubes, filter paper, funnels, graduated cylinders and beakers
- Micropipettes (those which can dispense 10 to 1000-μL of liquid are preferable)

[Preparation of Reagents]

- 1) Return dilution buffer and extraction/concentration buffer to room temperature before use.
- 2) Dilute extraction/concentration buffer × 10 with purified water before use.

Note 1) When sedimentation has occurred in the extraction/concentration buffer, return to room temperature and dissolve sedimentation before diluting with purified water.

Note 2) The dilution buffer and extraction/concentration buffer can be used irrespective of the type of Kit as they are of identical composition throughout the FASTKIT Immunochromato Series.

[An example of sample-extract preparation]

Conditions for the optimal extraction of peanut proteins may differ depending on the food sample, therefore it is necessary to analyze the composition of each sample.

Following is a general example of an extraction procedure.

- 1) Each of the food samples should be ground into a homogenous state or paste.
- 2) Each 2-g sample should be homogenized with 38 mL of the diluted extraction buffer for 30 sec. This procedure should be repeated three times for each sample.
- 3) The homogenate should be centrifuged for 20 min at higher than 3,000 x g and 4°C.
- 4) The supernatant should be filtered through filter paper.
- 5) The filtrate should be diluted 10 parts to 1 with the dilution buffer to provide a sample extract.

Note 1) It is absolutely necessary to prevent the peanut components from being contaminated throughout the process. Therefore, equipment such as the homogenizer should be washed after each procedure. (Wash the equipment with a neutral detergent and then soak them in an alkaline detergent overnight, or wash them in an ultrasonic washer for 30 min.)

Note 2) The sample extracts should be adjusted to neutral pH (6.0 to 8.0). As it is difficult to extract the peanut proteins from certain foods, the samples should be well processed and the extraction process carried out with great care.

Note 3) As the presence of insoluble substances such as fat are the cause of non-specific reactions it may be necessary to alter the centrifuging variables. (e.g. centrifuge for 30 min at 4°C at higher than 10,000 x g)

Note 4) Store the sample extracts at 4°C and subject them to analysis as soon as possible. Avoid freezing.

<Procedure of sample extraction>

Sample

| ←Grind into a homogeneous state

2 g of the sample

| ←Add 38 mL of extraction buffer and homogenize

Sample Extract

| ←Centrifuge for 20 min at 3,000 x g and 4°C

Supernatant

| ←Filter through filter paper

Filtrate

| ←Dilute 10 parts to 1 with the dilution buffer

Extract for analysis

Procedure 2 (IMMUNOCHROMATO Procedure)

Warning: Correct results can become unobtainable due to the effect of moisture absorption on the immunochromatography test plate, therefore, only remove from the aluminum pouch immediately before use and only after returning it to room temperature.

[Preparation of Reagents]

Return the test plate to room temperature while still inside the aluminum pouch, then remove from the aluminum pouch directly before use. Make sure not to touch the test plate's sample dropper or the interpretation area with your finger etc.


After sufficiently mixing the sample solution prepared in accordance with operation method 1 (preparation of sample solution), drop 100μL on the test plate's sample dropper area using a micropipette, and leave sitting on a level surface for 15 minutes and make a visual interpretation.

Note 1) The micropipette tip must be replaced after each test.


[Interpretation of results]

Warning: The interpretation result is obtained 15 minutes after the commencement of the test. Also, irrespective of the shade, test are interpreted as positive when the red-purple colored line is observed.


Positive: Positive test results are indicated by a red-purple line at interpretation area (T) and interpretation area (C). A weak positive test result is designated as + w, positive as +, and a strong positive as ++.

Negative: A negative test result is indicated by a red-purple line at interpretation area (C).

Interpretation deferral (re-test): When no red-purple colored line is observed in both interpretation area (T) and interpretation area (C), or when the red-purple colored line is only observed in interpretation area (T), then subject to a retest.

False Positive (Cross-Reaction or Non-Specific Reaction) and False Negative Reactions

- 1) With the exception of the peanut no cross-reaction occurs with the other specific ingredients, (eggs, cow milk, wheat and buckwheat).
- 2) Cross-reactions may occur with the protein of almond, cashew nut, macadamia nut, hazelnut, soybean* and red kidney bean “Taisyo Kintoki”, resulting in false-positive results*.
- 3) Nonspecific reactions may occur, in the presence of highly viscous foods* and extremely concentrated proteins. In such cases, carry out analysis with optimum concentrations of dilutions.
- 4) The following foods may yield false-negative results: foods meat products sterilized by high-temperature and high-pressure treatment, canned whale products, retort-pouched foods, and baked foods

Correlation between FASTKIT ELISA and FASTKIT Immunochromato*

| Overall agreement ratio 244/249specimens (96.8%) | | FASTKIT ELISA peanut | | | | Total |
|--|----------|----------------------|---------|--------|----------|--------------|
| | | More than 10ppm | 5~10ppm | 1~5ppm | Negative | |
| FASTKIT Immunochromato Peanut | Positive | 9 | 0 | 3 | 3 | 15 |
| | Negative | 1 | 0 | 1 | 232 | 234 |
| Total | | 10 | 0 | 4 | 235 | 249specimens |

- 1) Specimens indicating 1ppm or greater using FASTKIT ELISA Peanut are positive.
- 2) Detection sensitivity depends on the components of food.

Precautions for the Use

[General Warnings]

- 1) Read the instruction leaflet carefully. Use the kit according to the procedures described in the instruction leaflet.
- 2) Do not use any kit that has passed its expiry date. The expiration date is printed on a label on the outer case and test plate's aluminum bag.
- 3) FASTKIT has been developed to measure the allergic substances in selected foods. Therefore, it is not a diagnostic reagent for food allergies. There has been no proven correlation between the test results of this kit and the allergy development and or symptoms.
- 4) Do not determine the presence or absence of allergic substances according not only to the results of the FASTKIT but also on other information obtained by reviewing the records of the raw ingredients, manufacturing processes and so on.
- 5) This instruction is prepared as a guideline for inspection personnel. It is the responsibility of the user to determine the usefulness of the kit and or the suitability to each of the selected foods.
- 6) Kit specifications may be altered without prior notice.

[Safety Precautions]

- 1) Do not smear any reagent from the kit on your skin, mucosa, clothing, etc.
- 2) If any of the reagents are swallowed or come into contact with the eyes, rinse thoroughly with tap water and seek medical attention.

Storage and Expiration Date

[Storage]

The kits should be stored in a refrigerator (2 to 8°C) in the absence of light. Avoid freezing.

[Expiration Date]

Valid for 12 months from the production date in an unopened aluminum bag. The expiration date is printed on the outer case and test plate's aluminum bag.

[References]

- 1) Ministry of Health, Labour and Welfare, "Test Methods of Food Containing Allergy Substances" (SHOKUHATSU No. 1106001)
- 2) Eiji Marui: Shokuhin Eisei Kenkyu (Food Sanitation Research) Vol.52 (5), p25-31, 2002
- 3) Hiroshi Akiyama, Masatake Toyoda: Shokuhin Eisei Kenkyu (Food Sanitation Research) Vol.52 (6), p65-73, 2002

※ The FASTKIT Immunochromato Series applies the basic technology of immunochromatograph held by Becton Dickinson and Company, Ltd. and is made with the antibodies utilized in Nippon Meat Packers, Inc.**

[Distribution Source and Contact] **

Indicate side surface of outer case

[Manufacturer] **

Research and Development Center
Nippon Meat Packers, Inc.
3-3 Midorigahara, Tsukuba, Ibaraki, 300-2646, Japan
Phone : +81-29-847-7825, Fax : +81-29-847-7824
URL : <http://www.rdc.nipponham.co.jp>

E2-C410520C71

COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN

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