

# Labiase

Cat#OZK-OZ30-EX

from *Streptomyces fulvissimus* TU-6.

**Labiase**, produced by a submerged culture of *Streptomyces fulvissimus* TU-6, is a new enzyme preparation that lyses effectively cell walls of numerous lactic acid bacterium. This preparation is used as a tool for studies of cell walls structure of lactic acid bacterium and the preparation of plasmid DNA, intracellular enzyme and protoplast from lactic acid bacterium and so on.

## SPECIFICATIONS:

Appearance.....Lyophilized powder (containing lactose)  
 Package Size 500mg  
 Activity..... $\geq 10$  units/g ( $\beta$ -N-Acetyl-D-glucosaminidase)  
 Other activities contained Lysozyme  $\geq 2 \times 10^4$  units/g  
 Optimum pH and temperature.....See Fig. 1,2  
 pH and thermal stability .....See Fig. 1,2




Fig.1 pH-stability and activity  
 Stability:18 hours treatment at 25°C  
 Activity:37°C




Fig.2 Thermal-stability and activity  
 Stability:10 minutes treatment  
 Activity:pH4.0

## ASSAY FOR $\beta$ -N-ACETYL-D-GLUCOSAMINIDASE ACTIVITY:

### Unit Definition

One unit of enzyme catalyzes the release of 1  $\mu$  mole of *p*-nitrophenol from *p*-nitrophenyl N-acetyl- $\beta$ -D-glucosaminide per minute at 37°C, pH4.0.

### Method

Reaction mixture

Substrate.....5mM *p*-nitrophenyl N-acetyl- $\beta$ -D-glucosaminide in 0.2M

Citrate-phosphate buffer, pH 4.0 0.95ml

Enzyme.....suitably diluted enzyme 0.05ml

Total volume 1.00ml

### Procedure

After incubation for 10 minutes at 37°C, terminate by 2ml of 1M Na<sub>2</sub>CO<sub>3</sub>. Read A<sub>405</sub>.

### Calculation

$$\text{Enzyme unit} = \frac{A_{405} \times 3}{17.7 \times 10^3 \times 10} \times 10^3$$

The molar extinction coefficient of free *p*-nitrophenol under these conditions is 17,700.

## ASSAY FOR LYTIC ACTIVITY TOWARD LACTIC ACID BACTERIUM

### Method

Reaction mixture

Buffer.....50mM Citrate -sodium citrate buffer, pH 4.0	1ml
Substrate...Lactic acid bacterium suspension( $A_{660} \doteq 5.0$ )	1ml
Enzyme.....500 $\mu$ g-protein/ml solution	1ml
Distilled water	2ml
Total volume	5ml

### Procedure

Before adding the enzyme solution preincubation is carried out for 5 minutes at 37°C. After incubation for 2 hours at 37°C with gentle shaking,  $A_{660}$  of the mixture is determined. As a reference, 1 ml of distilled water is used instead of enzyme solution.

### Calculation

$$\text{Percentage decrease in } A_{660} = \frac{(d_0 - dt) - (D_0 - Dt)}{d_0} \times 100$$

d(0 or t):  $A_{660}$  of reaction mixture after 0 or 2 hours

D(0 or t):  $A_{660}$  of reference mixture after 0 or 2 hours

### Result

See Table 1.

The extent of lysis of lactic acid bacterium cells by **Labiase** varies with lactic acid bacterium strain, growth stage of lactic acid bacterium, or cultural condition.

### STORAGE:

Lyophilized preparation is stable for at least 1 year when stored at 4°C.

Table 1. Lysis of lactic acid bacterium cells by **Labiase**.

Strains	Lysis(%)
<i>Leuconostoc mesenteroides</i> IFO 3832	92.8
<i>Leuconostoc lactis</i> IFO 12455	99.6
<i>Streptococcus faecalis</i> IFO 3971	96.5
<i>Lactobacillus acidophilus</i> IFO 13951	97.9
<i>Lactobacillus plantarum</i> IAM 1041	96.8
<i>Lactobacillus casei subsp. casei</i> IFO 3533	59.4
<i>Lactobacillus casei subsp. ramosus</i> IAM 1118	91.8
<i>Lactobacillus brevis</i> IFO 3345	95.6
<i>Lactobacillus fermentum</i> IFO 3071	70.4
<i>Lactobacillus sp. (buchneri)</i> IFO 3961	89.4
<i>Lactobacillus fructivorans</i> IFO 13954	96.9
<i>Lactobacillus delbrueckii subsp. bulgaricus</i> JCM 1002	97.7
<i>Lactobacillus delbrueckii subsp. delbrueckii</i> JCM 1012	92.0
<i>Lactobacillus helveticus</i> JCM 1120	98.7
<i>Lactobacillus helveticus</i> JCM 1005	84.5
<i>Lactobacillus kefir</i> JCM 5818	93.0
<i>Lactobacillus sanfrancisco</i> JCM 5818	82.1
<i>Bifidobacterium bifidum</i> JCM 1255	62.9
<i>Streptococcus salivarius subsp. thermophilus</i> IFO 13957	94.0
<i>Streptococcus mutans</i> IFO 5705	73.9

### REFERENCES:

- 1) Ohbuchi, K., Hasegawa, K., Hamachi, M., Ozeki, K., and Kumagai, C. (2001) J. Biosci. Bioeng. 91, 487.

#### Distributor



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME

KOTO-KU, TOKYO 135-0016, JAPAN

TEL : +81-3-5632-9617

FAX : +81-3-5632-9618

URL : <http://www.cosmobio.co.jp/>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

#### Manufacturer

**ÖZEKI** OZEKI CORPORATION

4-9, Imazu Dezaike-Cho, Nishinomiya-shi

Hyogo 663-8227 JAPAN

TEL. : +81-798-32-2171

FAX. : +81-798-34-7475

URL : <http://www.ozeki.co.jp>