




White Adipocyte Culture kit, V-1 (Rat)



(Rat Cryopreserved Cell and culture medium)

Principle


White adipocytes, a type of subcutaneous adipocytes, play an important role in energy storage. White Adipocyte Culture Kit, V-1 (Rat) contains preadipocytes isolated from rat subcutaneous adipose tissues and culture medium that induces differentiation of precursor cells into mature adipocytes. The kit provides a convenient system for studying the mechanism of adipogenesis as well as for analyzing effects of drugs on metabolic syndrome such as obesity, diabetes and hypertension.



White adipocytes
Day 9



Visceral adipocytes
Day 9



Epididymal adipocytes
Day 9

Components

Components	Size	Quantity	Storage Conditions	Shelf Life
Subcutaneous White Preadipocytes, rat	3.0 x 10 ⁶ cells	1	Liquid Nitrogen	1 year
			-80°C Freezer	1 year
White Adipocyte Culture Medium (PMC-SACMR-COS)	250 mL	1	-80°C Freezer	1 year
			-20°C Freezer	6 months

Shipping : dry ice

Materials required but not provided

- Variable volume pipettes
- Culture plate



COSMO BIO CO., LTD.

Inspiration for Life Science


Precautions

- Read the instructions carefully before beginning the culture.
- This kit is for research use only, not for human or diagnostic use.
- Always wear gloves and lab coat when handling the cell culture.


Protocols <Cultured with the 24well-plate>

1. Thaw Adipocyte Culture Medium in a 37°C water bath with gentle shaking .
2. Quickly place Common Precursor Cell vial in a 37°C water bath until the contents are thawed.
3. Transfer thawed cells into a 15 ml centrifuge tube containing 10 ml of Adipocyte Culture Medium. Mix gently and centrifuge for 5 minutes at 4°C at 170 x g for 5 minutes.
4. After removing the supernatant, resuspend cells in 10 ml of Adypocyte Culture medium and centrifuge for at 4°C at 170 x g for 5 minutes.
5. After removing the supernatant, resuspend cells in 12.5 ml of Adypocyte Culture medium.
6. Dispense 0.5 ml of cell suspension to each well of 24-well plate.
7. Incubate the plate at 37°C under 5 % CO₂, 100 % humidity.
8. After 1 Day culture, add 0.5 ml of Adipocyte Culture Medium gently to each well of 24-well plate.
9. Change the medium every 2 days. Be gentle not to disturb the cell layer.
 - i. Approximately 3 days of culture, preadipocyte culture becomes confluent.
 - ii. Approximately 4 -5 days of culture, cells turn to mature adipocytes.
 - iii. Approximately 7 days of culture, cells become hypertrophic.
 - iv. Appoximately 8 days of culture, start detaching from the bottom of the well.


Cellular morphology




1 Day



2 Days



3 Days



8 Days

References

- (1) Hashimoto, T., Igarashi, J., Kosaka, H. Sphingosine Kinase is Induced in Mouse 3T3-L1 cells and Promotes Adipogenesis. J. Lipid Res. 50, 602-610 (2009)
- (2) Takahashi, K., Yoshina, S., Masashi, M., Ito, W., Inoue, T., Shiwaku, H., Arai, H., Mitani, S., Okazawa, H. Nematode Homologue of PQBP1, a Mental Retardation Causative Gene, Is Involved in Lipid Metabolism. PLoS One. 4, e4104 (2009)
- (3) Oguri, A., Tanaka, T., Iida, H., Meguro, K., Takano, H., Oonuma, H., Nishimura, S., Morita, T., Yamasoba, T., Nagai, R., Nakajima, T. Involvement of Cav3.1 T-type Calcium Channels in Cell Proliferation in Mouse Preadipocytes. Am. J. Physiol. Cell Physiol. 298, C1414-C1423 (2010)

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619