COSMO BIO Co., LTD.  
Inspiration for Life Science

CF001 Ver.1.0

Jul. 2010

Cat No. PRX-CF001

# Remarkable Yield Translation System Trial Kit 【RYTS Trial Kit】

## Instruction Manual (Version 1.0)

For research purposes only.

Before using this product, please read carefully  
this instruction manual.

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FAX : +81-3-5632-9619

## Table of Contents

1. Product overview.....	2
2. Product features.....	2
3. Kit contents .....	3
4. Before using this kit .....	3
4.1 Template considerations .....	3
4.2 Requirement of nuclease free environment.....	3
5. Construction and preparation of expression template .....	4
5.1 DNA Template.....	4
5.2 RNA Template.....	5
5.3 Related products.....	5
6. Protein synthesis .....	6
6.1 Additional material and equipment required .....	6
6.2 Standard procedure.....	6
6.3 Evaluation of produced protein activity.....	7
7. Combined application of RYTS and Clover <i>Direct</i> <sup>™</sup> .....	8
7.1 About Clover <i>Direct</i> <sup>™</sup> .....	8
7.2 Standard procedure using with Clover <i>Direct</i> <sup>™</sup> .....	8
7.3 Expression of site-directly labeled proteins .....	9
7.4 Related products to Site-Directed Protein Functionalization .....	10
8. Troubleshooting.....	11
9. References.....	12
10. Contact and support .....	14

## 1. Product overview

The Remarkable Yield Translation System (RYTS) Kit is a cell-free protein synthesis system. This kit includes an *E.coli* extract and all essential components for coupled transcription/translation reaction.

The RYTS Kit contains the *E.coli* Lysate, which is prepared according to unique method developed for highly efficient production of extensive protein by RIKEN (1-4). ProteinExpress Co., Ltd. has been licensed to distribute this method by RIKEN under JP patent No. 4477923, US No. 7615344 and EP No.1655375.

## 2. Product features

Synthesis of recombinant protein using RYTS Kit is highly efficient, convenient and quick. Maximum yield of protein can reach 150 µg per the 300 µl reaction.

Furthermore, this kit is available for the expression of protein with unnatural amino acid by exploitation of CloverDirect™ tRNA Reagents for Site-Directed Protein Functionalization. CloverDirect™ allows the incorporation of unnatural amino acids at defined positions of proteins using *in vitro* translation (see page 8 for details). If you need more information on CloverDirect™, please visit our web site.

([http://www.cosmobio.co.jp/export\\_e/products/proteins\\_and\\_peptides/products\\_prx\\_20091030.asp?entry\\_id=2769](http://www.cosmobio.co.jp/export_e/products/proteins_and_peptides/products_prx_20091030.asp?entry_id=2769)).

### 3. Kit contents

No.	cap color	Contents	Q'ty	Storage
1	white	<i>E.coli</i> Lysate	1 tube	−80 °C
2	pink	2× Reaction Mix	1 tube	
3	blue	Methionine	1 tube	
4	orange	Enzyme Mix	1 tube	−20 °C
5	green	CAT Control Vector 250ng/μl	1 tube	
6	clear	Nuclease Free Water	1 tube	

#### Storage condition

Store at appropriate temperature indicated in the Kit contents. Avoid multiple freeze-thawing cycles and vortex, which may cause for decrease of protein productivity.

### 4. Before using this kit

#### 4.1 Template considerations

It is able to use circular DNA, linear DNA and mRNA for the protein synthesis using RYTS Kit. Use a DNA template and mRNA with high purity (Abs260 / Abs280 > 1.8).

#### 4.2 Requirement of nuclease free environment

To eliminate nuclease contamination, you should wear gloves when preparing reaction mixture, reaction tubes and pipette tips should be nuclease-free. Avoid a template DNA or mRNA solution from contamination of nucleases.

## 5. Construction and preparation of expression template

### 5.1 DNA template

The DNA template (circular DNA or linear DNA) for protein synthesis should be contained that a protein-coding sequence be under the control of a T7 promoter and located downstream of a ribosomal binding site (RBS) sequence. Additionally, T7 terminator sequence is located at downstream of the coding sequence. Distance between T7 promoter and start ATG should not exceed 100 base pairs, and also distance between the RBS sequence and start ATG should not be more than 5-8 base pairs. Please refer to Figure 1 for an example of the DNA template.

Figure 1. An example of template

In case of very low protein yield, optimization of nucleotide sequence (codon usage, addition of N-terminal tags, etc.) is required to improve protein productivity. We recommend the use of ProX tag, which is original N-terminal peptide tag developed for efficient expression in *E.coli* cell-free protein synthesis system. ProX tag sequence is indicated at Figure 2.

(A)

5' -	AUG	UCU	AAA	CAA	AUC	GAA	GUA	AAC	UUU	UCU	AAU	GAG	- 3'
	Met	Ser	Lys	Gln	Ile	Glu	Val	Asn	Phe	Ser	Asn	Glu	

(B)

Figure 2. (A) ProX tag sequence, and (B) an example of template

### 5.2 mRNA template

Messenger RNA, which is generally synthesized from circular or linear DNA template using commercial transcription kit or RNA polymerase, is also available for RYTS reaction. In such case, the DNA template for transcription of mRNA should be also constructed according to section [5.1 DNA template](#). But, the T7 promoter, which is recognized by T7 RNA polymerase, can be replaced to other promoter recognized by other RNA polymerase (e.g., SP6 promoter recognized by SP6 RNA polymerase).

### 5.3 Related products

Cloning vectors for the construction of circular template DNA with ProX tag

- pROX-FL92.1ttt / Product Code: # TS013

Linear template generation kits for the construction of DNA template with or without ProX tag and His-tag

- RYTS Linear Template Set for *E.coli* (His-tag) / Product Code: #TS001

- RYTS Linear Template Set for *E.coli* (ProX-tag) / Product Code: #TS002

Transcription kit for preparation of mRNA

- T7 RNA polymerase / Product Code: #RP6T7

## 6. Protein synthesis

### 6.1 Additional material and equipment required

- DNA template (circular or linear) or mRNA encoding the protein of interest
- Nuclease-free microcentrifuge tubes and pipette tips
- Incubator, water bath or shaker with thermal controller

### 6.2 Standard procedure

1. Remove the reagents of this kit from freezer, and allow them to thaw on ice.
2. Set up the following reaction on ice in a nuclease-free microcentrifuge tube.

Cap No. / Cap color / Regents	volume
1 / white / <i>E.coli</i> Lysate	100 $\mu$ l
2 / pink / 2 $\times$ Reaction Mix	150 $\mu$ l
3 / blue / Methionine	3 $\mu$ l
4 / orange / Enzyme Mix	15 $\mu$ l
a) circular DNA or linear DNA b) mRNA	a) 0.4 – 1.0 $\mu$ g b) 50 – 100 $\mu$ g
6 / clear / Nuclease Free Water	up to 300 $\mu$ l


3. Mix thoroughly by pipetting several times, then incubate the reaction at 30°C for two hours.
4. After reaction, store the reaction mixture at 4°C for short term storage.

#### Note;

We recommend performing additional reaction using "CAT control vector" in this kit for positive control of the reaction. CAT assay, an enzymatic assay for quantification of chloramphenicol acetyl transferase (CAT) activity, is available for measuring typical productivity of this kit. If you wish to use CAT assay, please refer to the previous released article (5).

### 6.3 Evaluation of produced protein activity

To evaluate performance of RYTS Kit in comparison to a competitor's cell-free protein synthesis kit (Kit A), activities of Green Fluorescent Protein (GFP) resulting from protein syntheses by these two kits were determined. After these reactions, the fluorescent intensities of each reaction mixtures were measured by a fluorescence scanner. It was shown that the obtained fluorescent intensity from RYTS reaction mixture was 3-fold greater than from Kit A reaction mixture (see Figure 3).

**Figure 3. Comparison of produced GFP activity**

Applied volume : 50  $\mu$ l of translational reaction mix

Fluorescence detection (Ex : 488nm / Em : 520 nm)

## 7. Combined application of RYTS and CloverDirect™

### 7.1 About CloverDirect™

CloverDirect™ tRNA Reagent for Site-Directed Protein Functionalization allows the incorporation of unnatural amino acids at defined positions of proteins using a cell-free protein synthesis system (6-10). The protein with unnatural amino acid will be obtained within a few hours just by adding CloverDirect™ reagents and DNA template having an amber stop codon (UAG) or a four-base codon (CGGG) to a cell-free protein synthesis system. In our product lineup, transfer RNAs charged with unnatural amino acids containing fluorescent groups, biotin, PEG, photo-cross-linker, etc. are available (see page 10). For more information about CloverDirect™, please visit our web site.

([http://www.cosmobio.co.jp/export\\_e/products/proteins\\_and\\_peptides/products\\_prx\\_20091030.asp?entry\\_id=2769](http://www.cosmobio.co.jp/export_e/products/proteins_and_peptides/products_prx_20091030.asp?entry_id=2769))

### 7.2 Standard procedure using with CloverDirect™

CloverDirect™ reagent contains lyophilized unnatural aminoacyl-tRNA and tRNA buffer. Dissolve unnatural aminoacyl-tRNA with tRNA buffer, and mix with template DNA (or mRNA) and reagents of RYTS Kit. Typical reaction using with CloverDirect™ is indicated below.

Cap No. / Cap color / Reagents	Volume
1 / white / <i>E.coli</i> Lysate	100 µl
2 / pink / 2× Reaction Mix	150 µl
3 / blue / Methionine	3 µl
4 / orange / Enzyme Mix	15 µl
a) circular DNA or linear DNA b) mRNA	a) 0.4 – 1.0 µg b) 50 – 100 µg
Dissolved unnatural aminoacyl-tRNA	10 µl
6 / clear / Nuclease Free Water	up to 300µl

### 7.3 Expression of site-directly labeled proteins

A typical result of the protein labeling using RYTS and CloverDirect™ TAMRA (Product Code, #CLD02 and #CLD06) is shown in Figure 4.


Figure 4. In-gel detection of produced TAMRA labeled proteins

2UAG: UAG codon is inserted after initiator AUG codon    ProX tag\*: ProX tag is fused to the N-terminus

Applied volume: 0.25 µl of RYTS translational reaction mix    White arrow: TAMRA labeled proteins

\*ProX tag is original N-terminal peptide tag developed for incorporation of unnatural amino acids (11).

## 7.4 Related products to Site-Directed Protein Functionalization

Cloning vectors for the construction of circular template DNA with ProX tag

- pROX-FL92.1amber / Product Code: # TS011
- pROX-FL92.1cggg / Product Code: # TS012

CloverDirect™ tRNA Reagents for Site-Directed Protein Functionalization

### Site-Directed Fluorescence Labeling

CloverDirect™ CR110-X-AF-tRNA	[5-CR110-X : Abs/Em = 498/521nm]
CloverDirect™ HiLyte Fluor™ 488-AF-tRNA	[HiLyte Fluor™ 488 : Abs/Em = 497/525nm]
CloverDirect™ TAMRA-X-AF-tRNA	[5(6)-TAMRA-X : Abs/Em = 546/575nm]
CloverDirect™ ATTO 633-AF-tRNA	[ATTO633 : Abs/Em = 629/657nm]
CloverDirect™ ATTO 655-X-AF-tRNA	[ATTO655-X : Abs/Em = 633/684nm]

### Site-Directed Biotin Labeling

CloverDirect™ Biotin-AF-tRNA	[Biotin]
CloverDirect™ Biotin-X-AF-tRNA	[Biotin-X]
CloverDirect™ Biotin-XX-AF-tRNA	[Biotin-XX]

### Site-Directed Post-translational Modification

CloverDirect™ Tyr(PO <sub>3</sub> H <sub>2</sub> )-tRNA	[O-phospho-Tyr]
CloverDirect™ Lys(Me)-tRNA	[ε-methyl-Lys]
CloverDirect™ Lys(Me <sub>2</sub> )-tRNA	[ε-dimethyl-Lys]
CloverDirect™ Lys(Ac)-tRNA	[ε-acetyl-Lys]

### Site-Directed Unnatural Mutagenesis

PEGylated amino acids

CloverDirect™ PEG4-AF-tRNA	[Methyl-PEG4]
CloverDirect™ PEG8-AF-tRNA	[Methyl-PEG8]
CloverDirect™ PEG12-AF-tRNA	[Methyl-PEG12]

Cross-linking amino acids

CloverDirect™ BPA-tRNA	[p-benzoyl-phenylalanine]
CloverDirect™ AcPhe-tRNA	[p-acetyl-phenylalanine]

Photo-isomerizable amino acid

CloverDirect™ azoAla-tRNA	[p-phenylazophenyl-alanine]
---------------------------	-----------------------------

For more information about CloverDirect™, please visit our web site.

([http://www.cosmobio.co.jp/export\\_e/products/proteins\\_and\\_peptides/products\\_prx\\_20091030.asp?entry\\_id=2769](http://www.cosmobio.co.jp/export_e/products/proteins_and_peptides/products_prx_20091030.asp?entry_id=2769))

## 8. Troubleshooting

Low or no protein production	
1) Expired Kit	Use before expiration date.
2) Inappropriate storage condition	Store this product at appropriate temperature, and avoid multiple freeze-thawing.
3) Nuclease contamination	Nucleases markedly influence on productivity of protein (See section 4.2, page 3).
4) Difficult to produce the target protein using <i>E.coli</i> cell-free system	Addition of ProX tag to the protein N-terminus may be effective (See page 4).

Low or no activity of protein	
1) Disulfide bond is essential for the protein activity	Addition of oxidized glutathione may be effective. The adequate concentration of oxidized glutathione is ten millimolar.
2) Post-translational modification is required for the protein activity	In <i>E.coli</i> cell-free system, post-translational modification activity such as glycosylation and phosphorylation is lacked.
3) Protein folding is not correct	Lower reaction temperature may be effective for correct folding. Otherwise, please try to investigate refolding condition.

If your trouble is not improved according to the troubleshooting, please contact us (See page 14, Contact and support).

## 9. References

1) Cell-free protein synthesis system from *Escherichia coli* cells cultured at decreased temperatures improves productivity by decreasing DNA template degradation.

Eiko Seki, Natuko Matsuda, Shigeyuki Yokoyama and Takanori Kigawa.

*Analytical Biochemistry*, 377, 156-161 (2008)

2) A highly efficient cell-free protein synthesis system from *Escherichia coli*.

Dong-Myung Kim, Takanori Kigawa, Cha-Yong Choi and Shigeyuki Yokoyama.

*Eur. J. Biochem.*, 239, 881 -886 (1996)

3) Preparation of *Escherichia coli* cell extract for highly productive cell-free protein expression.

Takanori Kigawa, Takashi Yabuki, Natsuko Matsuda, Takayoshi Matsuda, Rie Nakajima, Akiko Tanaka and Shigeyuki Yokoyama.

*Journal of Structural and Functional Genomics*, 5, 63–68 (2004)

4) Automated system for high-throughput protein production using the dialysis cell-free method.

Masaaki Aoki, Takayoshi Matsuda, Yasuko Tomo, Yukako Miyata, Makoto Inoue, Takanori Kigawa and Shigeyuki Yokoyama.

*Protein Expression and Purification*, 68, 128-136 (2009)

5) Chloramphenicol Acetyltransferase from Chloramphenicol-Resistant Bacteria.

W. V. Shaw.

*Methods Enzymol.*, 43, 737-755 (1975)

6) FRET analysis of protein conformational change through position-specific incorporation of fluorescent amino acids

Daisuke Kajihara, Ryoji Abe, Issei Iijima, Chie Komiyama, Masahiko Sisido and Takahiro Hohsaka

*Nature Methods.*, 3, 923-929 (2006).

7) Position-specific incorporation of fluorescent non-natural amino acids into maltose-binding protein for detection of ligand binding by FRET and fluorescence quenching

Issei Iijima and Takahiro Hohsaka

*ChemBioChem.*, 2009, 10, 999-1006.

8) Position-specific incorporation of biotinylated non-natural amino acids into a protein in a cell-free translation system

Takayoshi Watanabe, Norihito Muranaka, Issei Iijima and Takahiro Hohsaka

*Biochem. Biophys. Res. Commun.*, 361, 794-799 (2007)

9) Comprehensive screening of amber suppressor tRNAs suitable for incorporation of non-natural amino acids in a cellfree translation system

Hikaru Taira, Yosuke Matsushita, Kenji Kojima, Kaori Shiraga and Takahiro Hohsaka

*Biochem. Biophys. Res. Commun.*, 374, 304-308 (2008).

10) Efficient Incorporation of Nonnatural Amino Acids with Large Aromatic Groups into Streptavidin in In Vitro Protein Synthesizing Systems

Takahiro Hohsaka, Daisuke Kajihara, Yuki Ashizuka, Hiroshi Murakami and Masahiko Sisido

*J. Am. Chem. Soc.*, 121, 34-40 (1999).

11) Incorporation of fluorescent non-natural amino acids into N-terminal tag of proteins in cell-free translation and its dependence on position and neighboring codons

Ryoji Abe, Kaori Shiraga, Shogo Ebisu, Hiroaki Takagi and Takahiro Hohsaka

*J. Biosci. Bioeng.*, 110, 32-38 (2010)

Manufactured by  ProteinExpress

COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619

# Remarkable Yield Translation System Trial Kit 【RYTS Trial Kit】

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