



## Instruction Manual for usage of Hybri-Bags

### Description

Hybri-bags are easy-to-use clear polyethylene bags available 2 types of hardness (soft and hard). They can be sealed with a standard heat-sealer and are available in 200 x 300 (mm) size. Gridded hybri-bags are also available to avoid any hassle of having to use a ruler to cut the bags straight. They are heat-sealable, resistant to temperatures up to at least 80°C(soft) and 120°C(hard), and do not curl when cut.

| Catalog Number | Product Name                                     | Size      |
|----------------|--|-----------|
| SE-S-1021-EX   | Hybridization Bags<br>Hybri-Bag Soft             | 50 Sheets |
| SE-S-1001-EX   | Hybridization Bags<br>Hybri-Bag Hard             | 50 Sheets |
| SE-S-1002-EX   | Hybridization Bags<br>Hybri-Bag Hard (1 cm grid) | 50 Sheets |


All hybrid-bags are coated with Nylon outside and Polyethylene and Polypropylene on the inside for soft and hard types respectively.

### Application

Hybridization bags are useful clear bags used for hybridization and incubation steps. They can be used for membrane incubation steps in Southern, Northern, and Western blotting protocols for applications such as hybridizations and washes, antibody and conjugate incubations. These tough, durable bags are boilable and freezable.

### How to use this product

1. Cut the Hybri-bag to adjust the size to the membrane.
2. Put membrane and buffer including antibody into the Hybri-bag, and seal with a heat sealer. (Note that you need to remove the bubbles from the hybri-bag when you seal to perform even blocking.)
3. Perform blocking on the shaker.



### Storage

No expiry date, however store at room temperature and avoid storing in a humid environment which may cause Hybri-bags to curl.

For research use only; not for diagnostic purposes



COSMO BIO Co., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.com>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618