## Anti-LexA antibody, rabbit serum

Immunized Animal: Rabbit  
Polyclonal antiserum

Product No.: BAM-61-001 50ul  
BAM-61-002 250ul

*E. coli* LexA protein binds specifically to the SOS-box sequence) and represses the genes belonging to the SOS regulon. In response to DNA damage, RecA protein is activated by ss-DNA accumulated in the damaged cells and promotes autocleavage of LexA repressor by its coprotease activity. As the results, DNA repair genes and error prone polymerases are induced, and DNA damage is repaired and mutation is induced (1).

The *lexA* gene is used for yeast two-hybrid experiments as a bait to identify the protein-protein interaction in vivo (2).

This product was prepared by immunizing rabbit with full-size highly-purified recombinant LexA protein. Using this antibody, 23 kD LexA protein was identified in the *E. coli* whole-cell lysate (Fig 1) and the expression of bait constructs was identified in yeast extracts by western blotting.

### Specifications


Form: antiserum added with 0.05% sodium azide  
Storage: 4 C (longer period, -70 C)

### Usages

- 1) Studies on the SOS regulation in *E. coli* (3). For western blotting; 1000~3000 fold dilution.
- 2) Construction and expression of a bait protein fused to LexA protein can be examined by western blotting of the yeast extracts, using the antiserum.

Purified LexA protein is available from BioAcademia (#01-002) to be used as a positive control for western blotting.

Fig.1 Detection of LexA repressor in the whole cell lysate of *E. coli* by this antiserum

### References

1. Friedberg EC, et al. DNA Repair and Mutagenesis 2<sup>nd</sup> Ed., ASM Press (2005)
2. Sambrook J & Russell DW, Molecular Cloning 3<sup>rd</sup> Ed. Cold Spring Harbor Press (2001)
3. Hishida T, et al., Genes Dev. 18, 1886-1897 (2004)

*\*For research use only; not for use as a diagnostic.*

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