Reagent for electrophoresis  
**2D-SILVER STAIN-II**

● **Introduction**

For detection of proteins and nucleic acids migration by polyacrylamide gel electrophoresis, silver staining is now being highlighted as a sensitive staining method. However the original silver staining method is not necessarily easy to work and takes long-time. **2D-SILVER STAIN-II** is an improved reagent kit developed to provide simple and fast staining.

● **Features**

- 1) Rapid staining results in a short time after electrophoresis.
- 2) Staining with higher sensitivity  
For proteins: 50 to 100 times more sensitive than CBB staining  
For nucleic acids: 50 to 100 times more sensitive than EB staining
- 3) Simple preparation and simple operation.
- 4) No use of such materials as heavy metals controlled by regulations.
- 5) Capable of restaining CBB-stained gels. (Double-staining)
- 6) Staining process can be stopped any time during development for desired chromatic figures.

● **Applications**

**2D-SILVER STAIN-II** is applicable to the detection of proteins and nucleic acids subjected to polyacrylamide and SDS-polyacrylamide slab gel electrophoresis.

● **Components**

Reagents	Main ingredients	Volume
① Fixing Reagent	Thiourea	100mL×1
② Pretreatment reagent	Dithiothreitol, Glutar-aldehyde, Thiourea	100mL×1
③ Staining Solution A	Silver nitrate	100mL×1
④ Staining Solution B	Ammonium hydroxide, Sodium hydroxide	100mL×1
⑤ Concentrated developer	Citric acid, Formaldehyde, Sodium thiosulfate	100mL×1
⑥ Stopper	Citric acid	100mL×1

● **Procedure**

For protein staining

In addition to the reagents supplied with the kit, the following materials are required.

- 1) Methanol
- 2) Acetic acid
- 3) Deionized water (less than  $10^{-6}$  mho conductivity)
- 4) Measuring cylinder
- 5) Flat-bottomed tray for developing gels
- 6) Measuring pipette

1. **Preparation of reagents**

Prepare the following solutions just before staining the gel. The volume of reagents used for preparation shown below is based on the size of 140 mm × 140 mm × 1.0 mm slab gel. For any different sizes of slab gels, the volume of reagents should be prepared in proportion to the size of the gel to be used.

1) **Fixing solution I**

Mix 100mL of methanol and 20 mL of acetic acid with 80 mL of deionized water. Stir and use the mixture as fixing solution I.

2) **Fixing solution II**

Mix 60 mL of methanol, 20 mL of acetic acid and 10 mL of ① Fixing Reagent with 110 mL of deionized water. Stir and use the mixture as fixing solution II.

3) **Pretreatment solution**

Mix 100 mL of methanol, 10 mL of ② Pretreatment Reagent and 90 mL of deionized water. Stir and use the mixture as pretreatment solution.


4) **Silver staining solution**

Mix 10 mL of ③ Staining Solution A and 10 mL of ④ Staining Solution B and stir. Then add 180 mL of deionized water and stir. Use the mixture as silver staining solution.


5) **Developer**

Add 10 mL of ⑤ Concentrated Developer to 190 mL of deionized water and stir. Use the mixture as developer.

2. **Staining Procedure**

- 1) Soak the gel in fixing solution I in a clean flat-bottomed tray, and shake it for 10 min-utes.
- 2) Remove the solution, and then pour 200 mL of fixing solution II to the tray and shake it for 15 minutes.
- 3) Remove the solution, and then pour 200 mL of pretreatment solution to the tray and shake it for 10 minutes.
- 4) Remove the solution, and then pour 200 mL of the deionized water to the tray and shake it for 5 minutes.

- 5) Remove the solution, and then pour 200 mL of silver staining solution to the tray and shake it for 15 minutes.
- 6) Remove the silver staining solution into a container. Immediately add 2 ~ 3 mL of concentrated hydrochloric acid to the waste solution, to convert the silver into silver chloride. (This procedure must be followed.) Pour 200 mL of deionized water to the tray and shake it for 2 minutes, and then remove the water. Repeat this step three times.
- 7) Pour 200 mL of developer to the tray and shake it well.
- 8) When appropriate chromatic figures are obtained, pour 10 mL of ⑥ Stopper Solution to the same tray and shake it for a moment.
- 9) Leave it about 10 minutes and then wash the gel with deionized water two or three times, and keep it.

In case of nucleic acids staining, Please use 10% Trichloroacetic acid and 50% Methanol in stead of fixing solution and pretreatment solution supplied.

Reagents	Volume (For 140mm×140 mm×1.0 mm slab gel)	Duration of shaking**		
		Gel ≤ 1mm	Gel ≤ 2 mm**	IEF
1. Fixing solution I	200 mL	10 min.	20 min.	60 min.**
2. Fixing solution II	200 mL	15 min.	30 min.	30 min.
3. Pretreatment solution	200 mL	10 min.	20 min.	30 min.
4. Deionized water	200 mL	5 min.	10 min.	5 min.
5. Silver staining solution	200 mL	15 min.	25 min.	15 min.
6. Deionized water	200 mL×3 times	2 min. every time	5 min. every time	2 min. every time
7. Developer	200mL	5 to 10 min.	5 to 10 min.	5 to 10 min.
8. Stopper	10 mL			

Cautions: \*a The duration of shaking is based on the temperature of 20 to 23°C.

\*b The quantity of the each reagent be twice that of gel ≤ 1 mm.

\*c In case of staining after IEF, use 200 mL of 3.5% Sulfosalicylic acid/11.5% TCA solution as fixing solution I.

● **Notes**

- 1) As the silver staining solution produces explosive silver amide when left to stand, the solution must be treated with hydrochloric acid or sodium chloride, to convert the silver into silver chloride. (See Procedure 6.)
- 2) Since Staining Solutions B among the components of 2D-Silver Stain-II, contains ammonium hydroxide and sodium hydroxide as an alkali substance respectively, take care to prevent them from contacting the eye and the skin. In case of eye contact, rapidly wash the eye with flowing water, and consult a doctor. In case of contact with the skin and clothes, rapidly wash away with water, and consult a doctor if inflammation is caused.
- 3) Staining Solution A and other components may also cause a stimulus to the mucous membrane and the skin. In case of contact with the eye, skin, and clothes, rapidly wash away with water, and consult a doctor if inflammation should be caused.
- 4) Addition of the silver staining solution to the sample gel may color or make turbid the solution. Such coloration and turbidity do not affect the stain. (See Procedure 5.)
- 5) Deionized water or distilled water of less than  $10^{-6}$  mho conductivity is recommended.
- 6) The tray to be used for the staining should be clean and have a smooth surface.
- 7) When a sample gel stained by the CBB staining is to be subjected to the silver staining, the sample gel should be washed with deionized water thoroughly, (3 hours or longer while changing deionized water several times prior to the silver staining.)
- 8) Care should be taken not to damage the gel during operation.
- 9) Care should be taken to see that the gel may not come out of the surface of solution.
- 10) Use a pair of disposable gloves in staining operation not to touch the gel.
- 11) If a proper shaker is not available, shake the tray by hand.
- 12) The sensitivity of a gel thicker than 1 mm declines somewhat as compared with a thinner gel.
- 13) Keep the reagents in a tightly closed container and refrigerated, and use within 6 months after unsealing.
- 14) Use the reagents (methanol, acetic acid) of superior quality.

● **Package size**

Enough to stain 10 slab gels (140 mm×140 mm×1.0mm)

● **Storage**

The kit should be stored in a dark place at 2-10°C.

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