



Anti Ecdysone-6-CMO-BSA polyclonal antibody

Product No. FKA-612E

The Cross reactivities of anti Ecdysonel-6-carboxymethyl-oxime-BSA IgG (FKA-612E) by ELISA with Ecdysone-6-CMO-HRP (FKA-611)

Antigen	Ecdysone-6-CMO-BSA		
Antibody(FKA-612E)	Purified by Protein A	lyophilized	0.2µl / Vial
Dilution of IgG	2000test for ELISA	50µl/well use	
Label (FKA-611)	Ecdysone-6-CMO-HRP		1ml / Vial
	1,000test for ELISA	50µl/well use	
O.D. (492/600nm) of B0	about 1.0		
Ecdysone dose of B0/2	5 ng/ml (250 pg/well)		
Strage/Handling	Maintain at 4°C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.		

Steroid	Cross Reaction
Ecdysone	100 %
20beta-OH-ecdysone	1.50 %
Cholesterol	0.00 %
Testosterone	0.00 %
Pregnenolone	0.00 %
Progesterone	0.00 %
Cortisol	0.00 %
The other Steroids	0.00 %

Note: The ELISA of Ecdysone carries out Antibody 50 µl and HRP conjugate 50 µl and Ecdysone sample 50 µl on microplate coated with an anti rabbit IgG (H+L) serum and then incubate at room temperature overnight.

For research use only; not for use as a diagnostic.



**FKA
#612E**

Cosmo FKA 611と612E による EcdysoneのELISAの特性

The Cross reactivities of anti Ecdysone-6-CMO-BSA serum by EIA with Ecdysone-6-CMO-HRP

抗体(FKA 612) anti Ecdysone-6-CMO-BSA serum
 rabbit No.N-22-3 0.2 μ l/vial
 Dilution of antiserum 1000, 000倍希釈液 50 μ l/well使用

標識(FKA 611) Ecdysone-6-CMO-HRP
 10 μ g/1 μ l/1ml/Vial 1000倍希釈液
 50, 000倍希釈液 50 μ l/well使用

O.D. 492nm of B0 1.0前後
 Ecdysone dose of B₀/2 5 ng/ml前後 (250 pg/well)

Cross reaction

Ecdysone	100 %
20 β -OH-ecdysone	1.5 %
Cholesterol	0 %
Testosterone	0 %
Pregnenolone	0 %
Progesterone	0 %
Cortisol	0 %