Code No.KAL-KM018

For research use only


## Thermosensitive TRP Channel Anti Rat TRPV1 (VR-1) Polyclonal Antibody

Capsaicin, a pungent ingredient of hot peppers, receptor has six transmembrane domains. It's a non-selective channel with high permeability of  $Ca^{2+}$ . Capsaicin, fat-soluble pain stimulus substance, has vanillyl group and is classified into the family of vanilloids. This receptor was named vanilloid receptor subtype 1 (VR-1) first, and is now named TRPV1 (transient receptor potential vanilloid subfamily member 1).


TRPV1 is activated not only by capsaicin but also by heat (over 43°C) or proton, and is found in the study of pain reception or transmission of stimuli.

This polyclonal antibody is specific for TRPV1 of rat, and has been proved to be useful for the immunohistochemistry.

Package Size	5 $\mu$ g (50 $\mu$ L/vial)
Format	Rabbit polyclonal antibody purified by antigen G affinity chromatography.
Buffer	1%BSA as a stabilizer, containing 0.1% Proclin as a bacteriostat
Storage	Store below $-20^{\circ}C$ Once thawed, store at $4^{\circ}C$ . Repeated freeze-thaw cycles should be avoided
Purification method	This antibody was purified from rabbit serum by Protein G affinity chromatography.
Working dilution for immunohistochemistry:	0.1 $\mu$ g/mL;

dorsal root ganglion (DRG) of lumbar region  
(normal rat), 30  $\mu$ m of thickness  
Fukuoka, T.  
Second Department of Anatomy, Hyogo college  
of medicine, Hyogo, Japan

dorsal root ganglion (DRG) of lumbar region  
(normal rat), 30  $\mu$ m of thickness  
Fukuoka, T.  
Second Department of Anatomy, Hyogo college  
of medicine, Hyogo, Japan

Preparation of antibodies and instruction  
Tominaga, M.  
Department of Physiology, Faculty of Medicine,  
Mie University, Japan


## Thermosensitive TRP Channel Anti Rat TRPV1 (VR-1) Polyclonal Antibody

### 【References】

1. Tominaga M, & Julius D (2000): capsaicin receptor in the pain pathway. *Jpn J Pharmacol* 83(1): 20-24
2. Tominaga M, Caterina MJ, Malmberg AB, Rosan TA, Gilbert H, Skinner K, Raumann BE, Basbaum AI, & Julius D (1998): The cloned capsaicin receptor integrates multiple pain-producing stimuli. *Neuron* 21(3): 531-543
3. Davis JB, Gray J, Gunthorpe MJ, Hatcher JP, Davey PT, Overend P, Harries MH, Latcham J, Clapham C, Atkinson K, Hughes SA, Rance K, Grau E, Harper AJ, Pugh PL, Rogers DC, Bingham S, Randall A, & Sheardown SA (2000): Vanilloid receptor-1 is essential for inflammatory thermal hyperalgesia. *Nature* 405(6783): 183-187
4. Hwang SW, Cho H, Kwak J, Lee SY, Kang CJ, Jung J, Cho S, Min KH, Suh YG, Kim D, & Oh U (2000): Direct activation of capsaicin receptors by products of lipoxygenases: endogenous capsaicin-like substances. *Proc Natl Acad Sci USA* 97(11): 6155-6160
5. Mandadi S. et al. (2006): Increased sensitivity of desensitized TRPV1 by PMA occurs through PKCepsilon-mediated phosphorylation at S800. *Pain*. Jul; 123(1-2): 106-16.\*

\* Application Reference

Chemical structure of capsaicin

vanillyl group

Additional: Anti Thermosensitive TRP Channel antibodies available from TRANSGENIC INC.

TRPV1	KM018	Anti Rat TRPV1 (VR-1) Polyclonal Antibody
TRPV2	KM019	Anti Rat TRPV2 (VRL-1) Polyclonal Antibody
TRPM8	KM060	Anti Rat TRPM8 (CMR1) Polyclonal Antibody
phospho-TRPV1	KM112	Anti Rat phospho TRPV1 (VR-1) Polyclonal Antibody
TRPV4	KM119	Anti Mouse TRPV4 Polyclonal Antibody
TRPA1	KM120	Anti Mouse TRPA1 Polyclonal Antibody

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

1. 富永 真琴:多刺激痛み受容体としてのカプサイシン受容体の構造と生理機能実験医学 Vol. 18 No.17 : 2325-2330, 2000
2. Tominaga M, & Julius D (2000): capsaicin receptor in the pain pathway. *Jpn J Pharmacol* 83(1): 20-24
3. Tominaga M, Caterina MJ, Malmberg AB, Rosan TA, Gilbert H, Skinner K, Raumann BE, Basbaum AI, & Julius D (1998): The cloned capsaicin receptor integrates multiple pain-producing stimuli. *Neuron* 21(3): 531-543
4. Davis JB, Gray J, Gunthorpe MJ, Hatcher JP, Davey PT, Overend P, Harries MH, Latcham J, Clapham C, Atkinson K, Hughes SA, Rance K, Grau E, Harper AJ, Pugh PL, Rogers DC, Bingham S, Randall A, & Sheardown SA (2000): Vanilloid receptor-1 is essential for inflammatory thermal hyperalgesia. *Nature* 405(6783): 183-187
5. Hwang SW, Cho H, Kwak J, Lee SY, Kang CJ, Jung J, Cho S, Min KH, Suh YG, Kim D, & Oh U (2000): Direct activation of capsaicin receptors by products of lipoxygenases: endogenous capsaicin-like substances. *Proc Natl Acad Sci USA* 97(11): 6155-6160
6. Mandadi S. et al. (2006): Increased sensitivity of desensitized TRPV1 by PMA occurs through PKCepsilon-mediated phosphorylation at S800. *Pain*. Jul; 123(1-2): 106-16.\*

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