MONOCLONAL ANTIBODIES PRODUCED  
BY THE RAT LYMPH NODE METHOD

For research use only

## Anti Collagen IV $\alpha$ 5(IV) Chain, Human (Mono)

Cat#: SGE-C452 Quantity: 500  $\mu$ l

**This product can be used for immunostaining of cryostat sections and for western blotting**


### How to stain human cryostat sections with C-452-EX (H52)

#### **Materials**

1. Acid solution : KCl-HCl buffer, pH 1.5, Ion strength=0.1 (Mix equal volumes of 0.1 M KCl and 0.1 M HCl.)
2. Blocking solution: 10% goat serum (normal goat serum diluted with PBS)
3. First antibody: H52 monoclonal Ab is usually diluted to 1/10 - 1/100 with PBS containing 1 mg/ml of bovine serum albumin before use.
4. Secondary Ab: Affinity purified Ab to rat IgG (Cappel #55760) is usually diluted to 1/50- 1/100 with PBS containing 5% human normal serum before use.

#### **Procedure**

1. Cryostat sections: 4  $\mu$ m
2. Treatment with acetone for 10 min
3. Acid-solution treatment with 0.1 M KCl- HCl (pH1.5) for 10 min
4. Washing with PBS
5. Blocking with 10% normal goat serum for 10 min
6. Washing
7. First Ab for 60 min
8. Washing
9. Secondary Ab for 30 min
10. Washing
11. Mount with Perma Fluor (Tanner #11611)

Human cryostat section stained with monoclonal antibody H52.

Glomerular basement membrane (BM), part of tubular BM, and Bowman's capsular BM are stained.

COSMO BIO CO., LTD.  
Inspiration for Life Science

Engl. 060523

MONOCLONAL ANTIBODIES PRODUCED  
BY THE RAT LYMPH NODE METHOD

Anti Collagen IV  $\alpha$ 5(IV) Chain, Human (Mono)

For research use only

- Product name: Anti Collagen IV  $\alpha$ 5(IV) Chain, Human (Mono)  
Code: SGE-C452  
Quantity: 500  $\mu$ l  
Appearance: Solution. Diluted ascitic solution obtained from nude mice; diluted with culture medium containing 10 % fetal bovine serum and 0.1%  $\text{NaN}_3$ .  
Clone names: H52 (rat IgG2b/kappa)  
Specificity of antibodies: H52 is specific to NC1 domain of  $\alpha$ 5(IV) (Reference 1).  
Reactive with NC1 domain of  $\alpha$ 5(IV) of rat, rabbit, guinea pig, dog, and bovine.  
Preparation of antibodies : Monoclonal antibodies were prepared by the rat lymph node method developed by Shigei Medical Research Institute with a synthetic peptide of type IV collagen as an immunogen.  
use : Staining of human cryostat sections by indirect immunofluorescence (Acid treatment is necessary.); Western blotting (dilution x100)  
Storage : Store at below  $-30^\circ\text{C}$ . Stable at these conditions for several years  
Reference : 1) Sado et al. (1995) Establishment by the rat lymph node method of epitope-defined monoclonal antibodies recognizing the six different a chains of human type IV collagen. Histochem. Cell Biol. 104: 267-275.

**Distributor**

COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME  
KOTO-KU, TOKYO 135-0016, JAPAN  
TEL : +81-3-5632-9617  
FAX : +81-3-5632-9618  
URL : <http://www.cosmobio.co.jp/>  
e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

**Manufacturer**

**Shigei Medical Research Institute**

2117 Yamada , Okayama  
701-0202, Japan  
TEL: +81-86-282-3113  
FAX:+81-86-282-3115  
E-mail: [inst@shigei.or.jp](mailto:inst@shigei.or.jp)

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

MONOCLONAL ANTIBODIES PRODUCED  
BY THE RAT LYMPH NODE METHOD

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