



EnzoLyte™ 490 HIV-1 Protease Assay Kit

Fluorimetric

Catalog #	71127
Unit Size	1 Kit
Kit Size	500 Assays (96-well) or 1250 Assays (384-well)

This kit is optimized to detect Human Immunodeficiency Virus-1 (HIV-1) protease activity using an EDANS/DABCYL FRET peptide as a substrate. The fluorescence signal can be monitored at Ex/Em=340 nm/490 nm upon proteolytic cleavage. Ample materials are provided to perform 500 assays in a 96-well format or 1250 assays in a 384-well format. The kit has the following features:

- **Convenient Format:** Complete kit including all the assay components.
- **Optimized Performance:** Optimal conditions for the detection of HIV-1 protease activity.
- **Enhanced Value:** Less expensive than the sum of individual components.
- **High Speed:** Minimal hands-on time.
- **Assured Reliability:** Detailed protocol and references are provided.

USA and Canada Ordering Information

AnaSpec Corporate Headquarter

2149 O'Toole Ave.
 San Jose, CA 95131
 Toll-Free: 800-452-5530
 Tel: 408-452-5055
 Fax: 408-452-5059
 E-mail: service@anaspec.com
 Internet: www.anaspec.com

Technical Support

Tel: 408-452-5055
 Fax: 408-434-9266
 E-mail: assay@anaspec.com

International Ordering Information


A list of international distributors is available at www.anaspec.com.

INTRODUCTION

The 10~12 kD aspartic protease of Human Immunodeficiency Virus-1 (HIV-1) is required for the post-translational cleavage of the precursor polypeptides, Pr^{gag} and Pr^{gag-pol}¹. Since these cleavages are essential for the maturation of HIV infectious particles, this protease has become one of the key targets for developing anti-AIDS drugs.

The EnzoLyte™ 490 HIV-1 Protease Assay Kit provides a convenient assay for high throughput screening of HIV-1 protease inhibitors and for continuous quantification of HIV-1 protease activity using a fluorescence resonance energy transfer (FRET) peptide. The sequence of this FRET peptide is derived from the native p17/p24 cleavage site on Pr^{gag} for HIV-1 protease. In this FRET peptide, the fluorescence of EDANS is quenched by DABCYL. Upon cleavage into two separate fragments by the HIV-1 protease at the Tyr-Pro bond (see **Scheme 1**), the fluorescence of EDANS is recovered, and can be monitored at excitation/emission = 340 nm/490 nm.

The assays are performed in a convenient 96 or 384-well microplate format.



Scheme 1. Proteolytic cleavage of EDANS/DABCYL FRET peptide by HIV-1 protease.

KIT COMPONENTS, STORAGE AND HANDLING

Note: Store all kit components at -20 °C, and keep components A and B away from light.

- Component A:** HIV-1 protease substrate (600 µL)
EDANS/DABCYL FRET peptide, Ex/Em=340nm/490nm upon cleavage
- Component B:** EDANS, fluorescence reference standard (100 µM DMSO solution, 20 µL)
Ex/Em=340nm/490nm
- Component C:** Pepstatin A, a characterized HIV-1 protease inhibitor (27.4 µg powder)
- Component D:** 2X Assay buffer (50 mL)
- Component E:** Stop solution (30 mL)
- Component F:** DMSO (100 µL)
- Component G:** DTT (1 M, 200 µL)

OTHER MATERIALS REQUIRED (BUT NOT PROVIDED)

HIV-1 protease: HIV-1 protease can be produced either from *E. coli*² or by chemical synthesis³.

96-well microplate: Black microplate provides better signal to noise ratio.

Fluorescence microplate reader: Capable of excitation at 340±30 nm and emission at 490±30 nm.

STANDARD OPERATION PROTOCOL

Note 1: For fluorescence instrument calibration, please refer to Appendix II (recommended for first-time users).

Note 2: Warm all kit components until thawed to room temperature before starting the experiments.

Please choose protocol A or B based on your needs.

Protocol A Screening HIV protease inhibitors using purified HIV-1 protease

1. Prepare working solutions

- Assay buffer: Dilute 2X assay buffer (component D) to 1X with deionized water. Add 1 M DTT (component G) 1:1000 into 1X assay buffer. Use this DTT-containing 1X assay buffer in ***all*** the following steps.

Note: Prepare fresh DTT-containing 1X assay buffer for each experiment.

- HIV-1 protease substrate solution: Dilute HIV-1 protease substrate (component A) 1:50 in assay buffer. Mix the reagents well.

Note: Prepare fresh substrate solution for each experiment.

- HIV-1 protease diluent: Dilute HIV protease to an appropriate concentration in assay buffer.

Note: Prepare enzyme diluent immediately before use. Do not vortex enzyme. Prolonged storage of diluted enzyme or vigorously vortexing will denature the enzyme. Keep the enzyme on ice.

- Test compound: Dilute test compounds with deionized water or an appropriate vehicle.

- Pepstatin A⁴ (control inhibitor): Add 20 μ L of DMSO (component F) into one vial of Pepstatin A (component C) to get a concentration of 2 mM. Dissolve completely by vortexing. Dilute 2 mM Pepstatin A to 2 μ M in assay buffer.

Note: 2 mM Pepstatin A should be stored at -20 °C for further use. It is not stable, and should be prepared fresh for each experiment.

2. Set up enzymatic reaction

- Add test compounds and HIV-1 protease diluent into a microplate. The suggested total volume of HIV-1 protease diluent and test compound for a 96-well plate is 50 μ L. The suggested total volume of HIV-1 protease and test compound for a 384-well plate is 20 μ L.
- Set up the following controls at the same time:
 - Positive control contains HIV-1 protease diluent without test compound.
 - Inhibitor control contains HIV-1 protease diluent and known inhibitor, Pepstatin A (10 μ L for 96-well plate or 4 μ L for 384-well plate).
 - Vehicle control contains HIV-1 protease diluent and vehicle used to deliver test compound (e.g. DMSO).
 - Test compound control contains assay buffer and test compound. Some test compounds have strong autofluorescence and may give false results.
 - Substrate control contains assay buffer.

Note: Bring the total volume of all the controls to 50 μ L/well for a 96-well plate or 20 μ L/well for a 384-well plate with assay buffer.

3. Pre-incubation

- Incubate the plate at the desired temperature for enzymatic reaction (e.g. 25°C or 37°C) for 10-15 min. In the meantime, also incubate the HIV-1 protease substrate solution at the same temperature.

4. Initiate the enzymatic reaction

- Add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of HIV-1 protease substrate. Mix the reagents completely by shaking the plate gently for 30-60 seconds.
- Measure fluorescence signal:

For kinetics reading: Immediately start measuring fluorescence intensity at Ex/Em=340 \pm 30 nm/490 \pm 30 nm continuously and record data every 5 minutes for 30 to 60 minutes.

For end-point reading: Incubate the reaction at room temperature for 30 to 60 minutes. Keep the plate away from direct light. (Optional) add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of stop solution (component E). Mix the reagents, then measure fluorescence intensity at Ex/Em=340 \pm 30 nm/490 \pm 30 nm.

- Data analysis: Refer to Appendix I.

Protocol B Measuring HIV-1 protease activity in biological samples

1. Prepare working solutions

- 2X Assay buffer: Add 1 M DTT (component G) 1:500 into 2X assay buffer (component D). Use this DTT-containing 2X assay buffer in ***all*** the following steps.
Note: Prepare fresh DTT -containing assay buffer for each experiment.
- HIV-1 protease substrate solution: Dilute HIV-1 protease substrate (component A) 1:50 in 2X assay buffer. Mix the reagents well.
Note: Prepare fresh substrate solution for each experiment.

2. Set up enzymatic reaction

- Add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of HIV-1 protease containing biological sample.
- Set up the following controls at the same time:
 - Positive control contains HIV-1 protease standard.
 - Negative control contains biological sample without HIV-1 protease.
 - Substrate control contains deionized water.

Note: Bring the total volume of all the controls to 50 μL /well for a 96-well plate or 20 μL /well for a 384-well plate.

3. Initiate the enzymatic reaction

- Add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of HIV-1 protease substrate. Mix the reagents completely by shaking the plate gently for 30-60 seconds.
- Measure fluorescence signal:
 - For kinetics reading: Immediately start measuring fluorescence intensity at $\text{Ex/Em}=340\pm30$ nm/ 490 ± 30 nm continuously and record data every 5 minutes for 30 to 60 minutes.
 - For end-point reading: Incubate the reaction at room temperature for 30 to 60 minutes. Keep the plate away from direct light. (Optional) add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of stop solution (component E). Mix the reagents, then measure fluorescence intensity at $\text{Ex/Em}=340\pm30$ nm/ 490 ± 30 nm.
- Data analysis: Refer to Appendix I.

Appendix I: Data Analysis

- The fluorescence reading from the substrate control well is the background fluorescence. Subtract this background reading from the readings of the other wells to get the relative fluorescence unit (RFU).
- For kinetics reading:
 - Plot data as RFU versus time for each sample. To convert RFU to concentration of the product of enzymatic reaction, please refer to [Appendix II](#) on how to set up the fluorescence reference standard.
 - Determine the range of initial time points during which the reaction is linear. 10-15% conversion appears to be the optimal range.
 - Obtain the initial reaction velocity (V_o) in RFU/min. Determine the slope of the linear portion of the data plot.
 - A variety of data analyses can be done, e.g., determining inhibition %, EC_{50} , IC_{50} , K_m , K_i , etc.
- For endpoint reading:
 - Plot data as RFU versus concentration of test compounds.
 - A variety of data analyses can be done, e.g., determining inhibition %, EC_{50} , IC_{50} , etc.

Appendix II: Instrument Calibration

- **EDANS fluorescence reference standard:** Dilute 100 μM EDANS (component B) to 1 μM in deionized water. Perform 1:2 serial dilutions to get concentrations of 500, 250, 125, 62.5, 31.25, 15.63 and 0 nM. Add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of the serially diluted EDANS from 1 μM to 0 nM.
- Add 50 μL /well (96-well plate) or 20 μL /well (384-well plate) of HIV-1 protease substrate solution (refer to protocol B step 1 for preparation).

Note: The HIV-1 protease substrate solution should be added to the EDANS reference standard to correct for the fluorescence inner filter effect.

- (Optional) If the stop solution (component E) was added into the enzymatic reaction before taking the end-point reading, the same volume of stop solution should be added to the reference standard wells to obtain a better comparison.
- Plot EDANS fluorescence reference standard as RFU (relative fluorescence unit) versus concentration as in **Figure 1**.

Note: The final concentrations of the EDANS reference standard solutions are 500, 250, 125, 62.5, 31.25, 15.63, 7.81, and 0 nM. This reference standard curve is used to calibrate for the variation of different instruments and the different batches of experiments. It is also an indicator of the amount of HIV-1 protease enzymatic reaction final product.




Figure 1. The EDANS reference standard calibration curve.

EDANS was diluted in assay buffer containing HIV-1 protease substrate. 100 μL of EDANS at each concentration was added into a black 96-well microplate. The fluorescence signal was measured by a fluorescence microplate reader (FLx800, Bio-Tek Instruments) with a filter set of Ex/Em=360 \pm 40 nm/460 \pm 40 nm. (Samples were done in duplicates).

REFERENCES

1. S. Seelmeier, H. Schmidt, V. Turk, H. K. von der, *Proc.Natl.Acad.Sci.U.S.A* 85, 6612-6616 (1988).
2. H. Gehring et al., *J.Virol.Methods* 109, 143-152 (2003).
3. J. Schneider and S. B. Kent, *Cell* 54, 363-368 (1988).
4. S. Seelmeier, H. Schmidt, V. Turk, H. K. von der, *Proc.Natl.Acad.Sci.U.S.A* 85, 6612-6616 (1988).