



DHL™ Cell Cytotoxicity Assay Kit

Fluorimetric

Catalog #	71302
Unit Size	1 kit
Kit Size	500 assays (96-well) or 1000 assays (384-well)

This kit provides a convenient fluorescent method to detect lactate dehydrogenase (LDH) released from damaged cells. Released cytoplasmic LDH is a well-established indicator of cell membrane damage and cell cytotoxicity. The fluorescence signal can be monitored at excitation/emission=530-560 nm/590 nm. The kit has the following features:

- **Convenient Format:** All essential assay components are included.
- **Optimized Performance:** Optimal conditions for the detection of cell cytotoxicity.
- **Enhanced Value:** Less expensive than the sum of individual components.
- **High Speed:** Minimal hands-on time.
- **Assured Reliability:** Detailed protocol and references are provided

USA and Canada Ordering Information

AnaSpec Corporate Headquarter

2149 O'Toole Ave.
 San Jose, CA 95131
 Toll-Free: 800-452-5530
 Tel: 408-452-5055
 Fax: 408-452-5059
 E-mail: service@anaspec.com
 Internet: www.anaspec.com

Technical Support

Tel: 408-452-5055
 Fax: 408-434-9266
 E-mail: assay@anaspec.com


International Ordering Information

A list of international distributors is available at www.anaspec.com.

INTRODUCTION

The damage of cell membrane leads to release of cytoplasmic enzymes. The measurement of lactate dehydrogenase (LDH) release is a well-accepted assay to estimate cell membrane integrity and quantify cell cytotoxicity¹⁻⁵. LDH release assay has been proven to correlate very well with traditional ⁵¹Cr release assay and trypan blue staining^{2,6}.

The EnzoLyte™ Cell Cytotoxicity Assay Kit uses resazurin as a fluorogenic indicator to measure the activity of LDH released from damaged cells (Scheme 1). The fluorescent signal is proportional to the number of dead cells ($r^2 > 0.95$). Meanwhile, the viable cells produce negligible fluorescent signal under same condition (Figure 1). Therefore, the assay can be performed in a mixture of damaged and viable cells. The kit is suitable for high throughput screening of cytotoxicity of a variety of compounds.



Scheme 1. In a cell population that co-exists viable cells and damaged cells, the DHL™ cell cytotoxicity assay kit only detects the dehydrogenases (e.g. LDH) activity released from damaged cell, not those in live cells. By using an enzyme-coupled reaction, the dehydrogenases in the medium will convert non-fluorescent resazurin to strongly fluorescent resorufin, which can be detected at emission of 590 nm with the excitation at 530-560nm.




Figure 1. The increment of fluorescence signal is correlated with the increment of lysed cells. 3T3 cells were seeded into 96-well plate. 10 μ L of lysis solution was added to the cells to release cytoplasmic LDH (lysed 3T3), while 10 μ L of growth medium was added to another sets of wells (viable 3T3). 50 μ L of LDH assay solution was added. The fluorescence signal was monitored at Ex/Em=530 \pm 30/590 \pm 30 nm 10 minutes later. The assay can detect as little as 97 lysed 3T3 cells ($> \pm 3S.D.$), while the living cells produce little fluorescence signal.

KIT COMPONENTS, STORAGE AND HANDLING

Note: Store all the kit components at -20 °C, and keep component A and B from light.

Component A: Assay mixture (1 bottle)

Component B: Assay buffer (30 mL)

Component C: Lysis solution (10 mL)

Component D: Stop solution (20 mL)

OTHER MATERIALS REQUIRED (BUT NOT PROVIDED)

96 or 384-well microplate: Tissue culture microplate with black wall and clear bottom.

Fluorescence microplate reader: Capable of detecting emission at 590 nm with excitation at 530-560 nm.

STANDARD OPERATION PROTOCOL

Note: Please use protocol A or B based on your needs.

Protocol A. Assay the cytotoxicity effect of test compounds.

1. Prepare cell culture.

- Seed 1×10^4 cells per well in a microplate. Add test compounds and then culture cells in a 37°C incubator for the desired exposure period. The total volume suggested for a 96-well plate is 100 μ L. The total volume suggested for a 384-well plate is 50 μ L. All samples must have at least four parallel wells.
- Set up the following controls at the same time. All controls must have at least four parallel wells.
 - Positive control contains cells and known proliferation or cytotoxicity inducer.
 - Negative control contains cells but no test compounds.
 - Vehicle control contains cells and the vehicle used to deliver test compounds.
 - Non-cell control contains growth medium but no cells.

Note: LDH contained in serum will contribute to background fluorescence.

- Test compound control contains growth medium and test compound. Some test compounds have strong autofluorescence and may give false results.

Note: Match the total volume of all the controls to 100 μ L for a 96-well plate or 50 μ L for a 384-well plate by growth medium.

2. Prepare LDH assay solutions.

Note: Warm up all the kit components until thawed at 22-25 °C before starting the experiments. Warm up the stop solution (component D) at 37 °C to dissolve all the precipitation.

- LDH assay solution: add 25 mL of assay buffer (component B) to one vial of assay mixture (component A), which is sufficient for five 96-well plate. Mix the reagents completely. Store unused portion of LDH assay solution at -20°C.

3. Assay LDH activity.

- Retrieve the cells from 37°C incubator and incubate at 22-25°C for 20-30 minutes.

Note: It is important to equilibrate the temperature of the cells to 22-25 °C before the assay.

- For the 96-well plate, add 10 μ L/well of lysis solution (component C) to half of the parallel wells of samples and controls. Meanwhile, add 10 μ L/well of Hank's balance solution (HBSS) or phosphate-buffered saline (PBS) to the rest half of the parallel wells of samples and controls. For

the 384-well plate, add 5 μL /well of lysis solution (component C) to half of the parallel wells of samples and controls. Meanwhile, add 5 μL /well of Hank's balance solution (HBSS) or phosphate-buffered saline (PBS) to the rest half of the parallel wells of samples and controls.

- Incubate the plate on a microplate shaker for 1-2 min at 50-100 rpm to facilitate cell lysing.
- For 96-well plate, add 50 μL of LDH assay solution to each well. Mix the reagents by shaking gently for 30 seconds. (Add 25 μL of LDH assay solution for a 384-well plate).
- Incubate the reaction at 22-25°C for 10 minutes.
- (Optional) Add 20 μL of stop solution (component D) to each well (5 μL for 384-well plate).
- Measure fluorescence intensity at Ex/Em=530-560 nm/590 nm immediately.

4. Perform data analysis:

- The fluorescence reading from the non-cell control well is the background fluorescence. The readings from other wells need to be subtracted with this background fluorescence.
- Calculate the percentage of cytotoxicity for samples and controls according to the following formula:

$$\% \text{ Cytotoxicity} = 100 \times (LDH_{\text{cytotoxic}}) / (LDH_{\text{maximal}})$$

The fluorescence reading from unlysed cells represents $LDH_{\text{cytotoxic}}$.

The fluorescence reading from lysed cells represents LDH_{maximal} .

Note: Negative control wells still have $LDH_{\text{cytotoxic}}$ reading caused by spontaneous release of LDH^2 or cells autolysis.

Protocol B. Cytotoxicity Assay for Co-culture System.

1. Prepare cells according to Table 1.

- Seed effector cell and target cells at desired testing concentration and ratio in a microplate. Match the total volume to 100 μ L per well for a 96-well plate (50 μ L for a 384-well plate).
- Set up the following controls at the same time.
 - Target cell control contains target cells at the concentration used in test wells. Prepare at least four parallel wells
 - Effector cell control contains effector cells at the concentration used in test wells.
 - Non-cell control contains growth medium but no cells.

Note: LDH contained in serum will contribute to background fluorescence.

Note: Match the total volume of all the controls to 100 μ L for a 96-well plate or 50 μ L for a 384-well plate by growth medium.

- Culture cells in a 37°C incubator for the desired exposure period.

2. Prepare LDH assay solutions.

Note: Warm up all the kit components until thawed at 22-25 °C before starting the experiments.

- LDH assay solution: add 25 mL of assay buffer (component B) to one vial of assay mixture (component A), which is sufficient for five 96-well plate. Mix the reagents completely. Store unused portion of LDH assay solution at -20°C.

3. Initiate enzymatic reaction.

- Retrieve the cells from 37°C incubator and incubate at 22-25°C for 20-30 minute.
 - Note: It is important to equilibrate the temperature of the cells to 22-25 °C before the assay.*
- For the 96-well plate, add 10 μ L/well of lysis solution (component C) to half of the parallel wells of the target cell control. For the 384-well plate, add 5 μ L/well of lysis solution (component C) to half of the parallel wells of the target cell control.
- Incubate the plate on a microplate shaker for 1-2 min at 50-100 rpm to facilitate cell lysing.
- For the 96-well plate, add 50 μ L of LDH assay solution to each well. Mix the reagents by shaking gently for 30 seconds. (Add 25 μ L of LDH assay solution for the 384-well plate).
- Incubate the reaction at 22-25°C for 10 minutes.
- (Optional) Add 20 μ L of stop solution (component D) to each well (5 μ L for 384-well plate). Mix the reagent by shaking 10 second.
- Measure fluorescence intensity at Ex/Em=530-560 nm/590 nm immediately.
- Continue to data analysis.

4. Perform data analysis:

- The fluorescence reading from the non-cell control well is the background fluorescence. The readings from other wells need to be subtracted with this background fluorescence.
- Calculate percent cytotoxicity according to the following formula⁷:

$$\% \text{ Cytotoxicity} = 100 \times$$

$$\frac{(\text{LDH}_{\text{cytotoxic}} - \text{Spontaneous LDH}_{\text{target cell}} - \text{Spontaneous LDH}_{\text{effector cell}})}{(\text{Maximal LDH}_{\text{target cell}} - \text{Spontaneous LDH}_{\text{target cell}})}$$

The fluorescence reading from testing wells containing both target and effector cell represents $\text{LDH}_{\text{cytotoxic}}$.

The fluorescence reading from lysed target cell control represents $\text{Maximal LDH}_{\text{target cell}}$.

The fluorescence reading from unlysed target cell control represents $\text{Spontaneous LDH}_{\text{target cell}}$.

The fluorescence reading from unlysed effector cell control represents $\text{Spontaneous LDH}_{\text{effector cell}}$.

REFERENCES

1. T. Mosmann, *J.Immunol.Methods* 65, 55-63 (1983).
2. C. Korzeniewski and D. M. Callewaert, *J.Immunol.Methods* 64, 313-320 (1983).
3. B. Arechabala, C. Coiffard, P. Rivalland, L. J. Coiffard, Y. Roeck-Holtzhauer, *J.Appl.Toxicol.* 19, 163-165 (1999).
4. N. W. Roehm, G. H. Rodgers, S. M. Hatfield, A. L. Glasebrook, *J.Immunol.Methods* 142, 257-265 (1991).
5. C. Legrand et al., *J.Biotechnol.* 25, 231-243 (1992).
6. D. Marzin, *Cell Biol.Toxicol.* 15, 359-365 (1999).
7. G. Konjevic, V. Jurisic, I. Spuzic, *J.Immunol.Methods* 200, 199-201 (1997).