



Taq DNA Polymerase, Economy

BACKGROUND

Thermus aquaticus DNA polymerase (**Taq DNA polymerase**) gene was expressed in *E.Coli* in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and the MW is 94 kDa, same as that of the natural enzyme.

- This enzyme is suitable for PCR reactions; capable of amplifying DNA with various primers.

Applications:

- 1) High-throughput PCR
- 2) Colony PCR
- 3) Incorporation of dUTP, dITP, and fluorescence-labeled nucleotides
- 4) Primer extension
- 5) Addition of a single nucleotide (adenosine) at the 3'-blunt ends

Size:

200 U (5U/ μ l)

Concentration:

5 units/ μ l, where one unit is defined as the amount of enzyme that can incorporate 10 nmols of total dNTPs into an acid-insoluble material in 30 minutes at 74°C when activated salmon sperm DNA was used as template/primer.

Form:

20mM Tris-HCl (pH 8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.5% Tween20, 0.5% Igepal CA-630

Quality Assurance:

Greater than 95% of protein determined by SDS-PAGE (CBB staining) (Fig.1) The absence of endonucleases and exonucleases was confirmed.

PCR Test:

Good amplification result was obtained in PCR reaction using λ DNA as a template (Fig.2).

Reagents Supplied with Enzyme:

10 x Reaction Buffer (*Taq*): 100mM Tris-HCl (pH 8.3), 500mM KCl, 15mM MgCl₂

Storage:

Store at -20°C

References:

Related Products

BAM-02-001-EX	Taq DNA Polymerase(+dNTPs)
BAM-02-021-EX	Pfu DNA Polymerase(+dNTPs)



General composition of PCR reaction mixture (total 50ul)

Taq DNA polymerase (5 units/ul)	*0.25 ul
10 x Reaction Buffer (Taq)	5 ul
2.5mM (each) dNTPs	4ul
Template	<500ng
Primer 1	0.2 ~ 1.0uM (final conc.)
Primer 2	0.2 ~ 1.0uM (final conc.)
Sterile distilled water	up to 50ul

*Use of excess amount is not recommended.

PCR condition

98°C	10 sec	} 25 cycles
57°C	30 sec	
72°C	8 min	

(2 min in the case of 2 kb DNA)




Fig.1 SDS-PAGE of Taq DNA polymerase




Fig.2 Amplification of λ DNA

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia, Inc.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619