ACTIFLUOR™ ADAMTS13 Activity*

REF 812

EC	REP
----	-----

Obelis s.a
 Bd. Général Wahis 53, 1030 Brussels, BELGIUM

INTENDED USE

The ACTIFLUOR™ ADAMTS13 Activity assay is a fluorescence resonance energy transfer (FRET) assay for the measurement of ADAMTS13 in human plasma. The assay is for *in vitro* diagnostic use.

EXPLANATION OF THE TEST

ADAMTS13, also known as von Willebrand Factor (vWF) cleaving protease, is a zinc metalloproteinase that cleaves ultra large vWF multimers (UL-vWF) at the Tyr(1605)-Met(1606) bond located in the A2 region of vWF.¹ Studies have shown that low levels of ADAMTS13 activity are associated with Thrombotic Thrombocytopenia Purpura (TTP), a life-threatening hematological condition characterized by low platelet count, microvascular thrombi, red cell fragmentation, CNS and renal complications.^{2,3} An ADAMTS13 activity level below 5% of normal leads to an accumulation of UL-vWF multimers in plasma.⁴ The UL-vWF multimers bind to receptors on platelets inducing platelet aggregation and formation of intravascular microthrombi.

Congenital TTP is a rare heritable disorder caused by mutations within the ADAMTS13 gene which result in the production of non-functional ADAMTS13 protein.^{5,6} The acquired form of TTP is an autoimmune-like disorder caused by the development of autoantibodies to ADAMTS13 that inhibit enzyme activity.⁷ Measurement of the ADAMTS13 activity level has been useful in identifying patients with TTP from other thrombocytopenic conditions such as hemolytic uremic syndrome (HUS) and idiopathic thrombocytopenic purpura (ITP).⁸

The existing assays for measuring ADAMTS13 activity including multimer analysis, collagen binding and vWF ELISAs can be long, difficult to perform and highly variable. A rapid and easy to perform FRET assay using a synthetic VWF73 peptidyl substrate was described by Kokame, *et al.*⁹ Recently, Zhang, *et al.*, developed a novel recombinant VWF86 FRET substrate based upon autoquenching of fluorescein.¹⁰

PRINCIPLE OF THE METHOD

ACTIFLUOR ADAMTS13 is a FRET assay that measures the amount of ADAMTS13 activity in human plasma. A citrated plasma sample is assayed for ADAMTS13 protease activity using recombinant VWF86-ALEXA FRET substrate. Proteolytic cleavage of the VWF86-ALEXA FRET substrate between the tyr/met residues by ADAMTS13 uncouples the ALEXA fluorochromes resulting in an increase in fluorescence. The increase in fluorescence over time (V_{max}) due to cleavage of the substrate by ADAMTS13 is monitored at 37°C using a spectrofluorometer ($Ex=485\text{ nm}$; $Em=535\text{ nm}$). A standard curve is constructed using normal plasma with a known concentration of ADAMTS13. The ADAMTS13 activity in the plasma is determined by interpolation of the V_{max} values from the standard curve.

* United States Patent No. 7,270,976

REAGENTS – Sufficient for assaying 16 samples in duplicate

- R1 Fluorescence Microwells**, 48 (white)
R2 ADAMTS13 Standard, 250 µL, 2 vials (lyophilized)
R3 ADAMTS13 Positive Control, 150 µL (lyophilized)
R4 DMSO, 0.5 mL
R5 ALEXA488-VWF86 FRET Substrate, 300 µL (lyophilized)
R6 Assay Buffer, 10 mL, 2 vials (lyophilized)
R7 ADAMTS13 Inactivated Plasma, 0.5 mL, 3 vials (lyophilized)

WARNINGS AND PRECAUTIONS

This product contains human source material that has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using registered methods. As no known test method can provide complete assurance that products derived from human specimens will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

Do not use kit components beyond the expiration date. Do not mix reagents from different kits. Avoid microbial contamination of the kit components. Do not mouth pipette.

FRET substrate	Warning		H315, H319, P264, P280, P302 + P352, P305 + P351 + P338, P332 + P313, P337 + P313	
Assay Buffer	Danger	  	CONT	4-Morpholinopropansulphonic acid, Sodium azide H302, H311, H315, H319, H335, H411, P261, P264, P270, P273, P280, P302 + P352, P305 + P351 + P338, P332 + P313, P337 + P313

Hazard Statements:

H302 Harmful if swallowed.
H311 Toxic in contact with skin.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.
H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements:

P261 Avoid breathing dust.
P264 Wash thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

Precautionary Statements

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P332 + P313 If skin irritation occurs: Get medical advice/attention.
P337 + P313 If eye irritation persists: Get medical advice/attention.

REAGENT PREPARATION AND STORAGE

R2 ADAMTS13 Standard

Reconstitute a vial with 250 µL of filtered deionized or distilled water. Let stand for 15 minutes at room temperature (18°-25°C) before gently mixing. See the vial label for the lot specific standard concentration "C". The ADAMTS13 Standard may be held at room temperature until use. Store reconstituted standard at 2°-8°C if not used within 2 hours. The standard is stable for up to 4 hours when at 2°-8°C and for up to 12 weeks when stored at -70°C or colder. Avoid multiple freeze-thaws.

R3 Positive Control

Reconstitute with 150 µL of filtered deionized or distilled water. Let stand for 15 minutes at room temperature (18°-25°C) before gently mixing. Positive Control may be held at room temperature until use. Store reconstituted control at 2°-8°C if not used within 2 hours. The Positive Control is stable for up to 4 hours when at 2°-8°C and for up to 12 weeks when stored -70°C or colder. Avoid multiple freeze-thaws.

R5 ALEXA488-VWF86 FRET Substrate

Prepare a 25% DMSO solution by diluting the DMSO reagent 1:4 with filtered deionized water. Reconstitute the substrate with 300 µL of the 25% DMSO solution. Let stand for 15 minutes at room temperature (18°-25°C) before gently mixing. DO NOT VORTEX! NOTE: Protect from light. The substrate may be held at room temperature until use. The ALEXA488-VWF86 FRET Substrate is stable for up to 2 weeks when protected from light and stored at 2°-8°C and for up to 12 months when stored at -20°C or colder. Avoid multiple freeze-thaws.

R6 Assay Buffer

Reconstitute a vial of Assay Buffer with 10 mL of filtered deionized or distilled water. Let stand for 15 minutes before briefly vortexing. A vial contains sufficient Assay Buffer to perform 24 assays. Assay Buffer is stable for up to 1 day at room temperature (18°-25°C) and for up to 12 weeks when stored at -20°C.

R7 ADAMTS13 Inactivated Plasma

Reconstitute a vial of ADAMTS13 Inactivated Plasma with 500 µL of filtered deionized or distilled water. Let stand for 15 minutes at room temperature (18°-25°C) before gently mixing. The ADAMTS13 Inactivated Plasma may be held at room temperature until use. Store reconstituted plasma at 2°-8°C if not used within 2 hours. The ADAMTS13 Inactivated Plasma is stable for up to 4 hours when at 2°-8°C and for up to 12 weeks when stored at -70°C or colder. Avoid multiple freeze-thaws.

SPECIMEN COLLECTION AND PREPARATION

Only citrate collected platelet poor plasma may be used for this assay. Do Not Use EDTA. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guidelines-Fifth Edition", CLSI Document H21-A5, Vol. 28, No. 5, January 2008. Plasma collection should be performed as follows:

1. Collect 9 parts of blood into 1 part of 3.2% (0.109 M) trisodium citrate anticoagulant solution.
2. Centrifuge the blood sample at 1,500 x g for 15 minutes.
3. Plasma should be stored at room temperature and assayed within 2 hours. Alternatively, plasma may be stored at -70°C for up to 6 months
4. Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at room temperature and assayed within 2 hours.

PROCEDURE

Materials Provided – See Reagents

Materials Required But Not Provided

filtered deionized or distilled water
8-channel pipette covering 50-200 µL
1-channel pipette covering 10-200 µL
microcentrifuge tubes, 1.5 mL
microwell plate fluorometer

Preparation of Standards and Reagents

A. Preparation of ADAMTS13 Standards

Prepare seven (7) ADAMTS13 standards for generating the standard curve by making serial dilutions of the ADAMTS13 Standard in microfuge tubes as follows:

Tube	Standard Concentration	Volume of	Added to
1	"C" ng/mL	120 µL of Standard	0 µL
2	C/2 ng/mL	60 µL from Tube 1	60 µL of ADAMTS13 Inactivated Plasma
3	C/4 ng/mL	60 µL from Tube 2	60 µL of ADAMTS13 Inactivated Plasma
4	C/8 ng/mL	60 µL from Tube 3	60 µL of ADAMTS13 Inactivated Plasma
5	C/16 ng/mL	60 µL from Tube 4	60 µL of ADAMTS13 Inactivated Plasma
6	C/32 ng/mL	60 µL from Tube 5	60 µL of ADAMTS13 Inactivated Plasma
7	0 ng/mL	0 µL	60 µL of ADAMTS13 Inactivated Plasma

B. Preparation of Sample Plasma and Positive Control

Sample plasma and Positive Control are each diluted 1:2 with ADAMTS13 Inactivated Plasma by adding 50 µL of citrated plasma or Positive Control to 50 µL of Inactivated Plasma and mixing gently. Hold at room temperature.

C. Preparation of ALEXA488-VWF86 FRET Substrate


Pre-warm the Assay Buffer to 37°C. Just prior to use in Step 5 below, dilute the reconstituted ALEXA488-VWF86 FRET Substrate 1:25 with Assay Buffer (e.g. add 0.2 mL of Substrate to 4.8 mL of Assay Buffer). Note: The assay requires 1.6 mL of prepared substrate for every 16 microwells used.

Assay Procedure (Kinetic Method)

1. Add 80 µL of Assay Buffer to each microwell to be used.
2. Add 20 µL of each ADAMTS13 Standard from Step A to wells A1-A2 through G1-G2 respectively. Add 20 µL of 1:2 diluted Positive Control from Step B to wells H1-H2.
3. Add 20 µL of 1:2 diluted plasma samples from Step B to each of two microwells.
4. Place the microwell plate in a fluorescence plate reader set at Ex=485 nm and Em=535 nm, with a cut-off at 530 nm, and the temperature set at 37°C. Use a low setting for photomultiplier tube. For readers using filters, use an FITC filter combination set. Allow the standards and samples to warm for 3 minutes.
5. Add 100 µL of 1:25 diluted ALEXA488-VWF86 FRET Substrate from Step C to each microwell.
6. Measure the increase in fluorescence for 20 minutes, collecting data at 20 second intervals. Calculate the rate of change in fluorescence, RFU/sec or V_{max} .

RESULTS

A standard curve is constructed by plotting the RFU/sec or V_{max} for each standard versus the corresponding concentration of ADAMTS13 in ng/mL. A standard curve should be generated each time the assay is performed. A representative standard curve using a SpectraMax® Spectrofluorometer (Molecular Devices, CA) is shown below.

CALCULATION (Kinetic Method)

Interpolate the ADAMTS13 concentration in the diluted sample and Positive Control directly from the standard curve via the RFU/sec or V_{max} of the sample. Multiply this concentration by 2 to obtain the ADAMTS13 concentration in the neat sample and Positive Control. If the RFU/sec or V_{max} of the sample falls above the range of the standard curve, dilute the plasma further with ADAMTS13 Inactivated Plasma and repeat the assay.


Alternative Assay Procedure (Endpoint Method)

If the fluorescence plate reader is not suitable for measurement in the kinetic mode, the endpoint method can be used.

1. Add 80 μ L of Assay Buffer to each microwell to be used.
2. Add 20 μ L of each ADAMTS13 Standard from Step A to wells A1-A2 through G1-G2 respectively. Add 20 μ L of 1:2 diluted Positive Control from Step B to wells H1-H2.
3. Add 20 μ L of 1:2 diluted plasma samples from Step B to each of two microwells.
4. Place the microwell plate in a fluorescence plate reader set at $Ex=485$ nm and $Em=535$ nm, with a cut-off at 530 nm, and the temperature set at 37°C. Use a low setting for photomultiplier tube. For readers using filters, use an FITC filter combination set. Allow the standards and samples to warm for 3 minutes.
5. Add 100 μ L of 1:25 diluted ALEXA488-VWF86 FRET Substrate from Step C to each microwell.
6. Immediately read the fluorescence in the microwells ($T=0$).
7. Incubate the microwells for 20 minutes at 37°C. Read the fluorescence again ($T=20$). Calculate the net change in fluorescence, ΔRFU ($RFU_{T=20} - RFU_{T=0}$).

RESULTS

A standard curve is constructed by plotting the ΔRFU for each standard versus the corresponding concentration of ADAMTS13 in ng/mL. A standard curve should be generated each time the assay is performed. A representative standard curve using a SpectraMax® Spectrofluorometer (Molecular Devices, CA) is shown below.

CALCULATION (Endpoint method)

Interpolate the ADAMTS13 concentration in the diluted sample and Positive Control directly from the standard curve via the ΔRFU . Multiply this concentration by 2 to obtain the ADAMTS13 concentration in the neat sample and Positive Control. If the ΔRFU of the sample falls above the range of the standard curve, dilute the plasma further with ADAMTS13 Inactivated Plasma and repeat the assay.

EXPECTED VALUES

In a study of normal individuals ($n=234$), the mean (± 1 SD) ADAMTS13 activity level was found to be 666 ± 135 ng/mL.

QUALITY CONTROL

The Positive Control should be included whenever the assay is performed. If the Positive Control fails to fall within the specified range, the test should be repeated. Contact BioMedica Diagnostics if the control repeatedly fails to fall within the specified range.

TRACEABILITY OF CALIBRATORS AND CONTROLS

Information on traceability of calibrators and controls is available upon request.

PERFORMANCE CHARACTERISTICS

Precision

The intra-assay coefficient of variation was determined to be 4.1%

The inter-assay coefficient of variation was determined to be 4.4%.

REFERENCES

1. Furlan, M., Robles, R. and Lämmle, B. L. *Blood* 1996, **87**: 4223-4234.
2. Tsai, H. M. *Blood* 1996, **87**: 4235-4244.
3. Furlan, M., *et al.* *N Engl J Med* 1998, **339**: 1578-1584.
4. Moake, J. L., *et al.* *N. Engl J Med* 1982, **307**: 1432-1435.
5. Zheng X, *et al.* *J Biol Chem* 2001, **276**: 41059-41063.
6. Levy, G. G., Nichols, W. C. and Lian, E. C. *Nature* 2001, **413**: 488-494.
7. Tsai, H. M., Rice, L. and Sarode, R. *Ann Intern Med* 2000, **132**: 794-799.
8. Unpublished data from BioMedica Diagnostics.
9. Kokame K, Nobe Y, Kokubo Y, Okayama, A. and Miyata, T. *Br. J. Haematol.* 2005, **129**: 93.
10. Zhang, L., Lawson, H. L., Harrish, V. C., Huff, J. D., Knovich, M. A. and Owen, J. *Anal Biochem* 2006, **358**: 298-300.

SpectraMax is a trademark of Molecular Devices, Sunnyvale, CA, USA

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