



MONOCLONAL ANTIBODY

*For research use only. Not for clinical diagnosis.***Catalog No. AGE-M03****Anti GA-pyridine****BACKGROUND**

Glycolaldehyde formed as a result of the myeloperoxidase-H₂O₂ (MPO) reaction can react with proteins to yield various AGEs. Recently, a novel specific GA-derived AGE, called GA-pyridine, has been described in foam cells and the extracellular matrix of human atherosclerotic fibrotic lesions, glomerular mesangial and Bruch's membrane and choroid.

Product type	Primary antibodies
Immunogen	GA-LDL
Raised in	Mouse
Myeloma	P3U1
Clone number	2A2
Isotype	IgG1
Host	Mouse
Source	Ascites
Purification	Protein G
Buffer	PBS containing 0.1% proclin as a preservative
Concentration	0.2 mg / ml
Volume	100 ul
Label	Unlabeled
Specificity	GA-pyridine
Cross reactivity	-
Storage	Store below -20°C (below -70°C for prolonged storage) Aliquot to avoid cycles of freeze/thaw.

Application notes	• Western blotting: 1/200 - 1/2000
Recommended dilutions	• Immunohistochemistry: 1/200 - 1/400
	• ELISA: 1/200 - 1/400

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References	1) Nagai R., Hayashi CM., Xia L., Takeya M., Horiuchi S: Identification in human atherosclerotic lesions of GA-pyridine, a novel structure derived from glycolaldehyde-modified proteins. J Biol Chem. 277, 48905-48912 (2002) PMID: 12377783
	2) Glenn JV., Mahaffy H., Wu K., Smith G., Nagai R., Simpson DAC., Boulton ME., Stitt AW. Advanced Glycation End Product (AGE) Accumulation on Bruch's Membrane: Links to Age-Related RPE Dysfunction. Invest. Ophth. Vis. Sci. 50, 441-451 (2009) PMID: 18676633

ANTIBODY CHARACTERIZATION




Fig.1 GA-pyridine structure




Fig.2 Immunoreactivity of the GA-pyridine (2A2) monoclonal antibody to GA-HAS, GA-BSA, MG-HAS, GO-HSA and HSA

RELATED PRODUCTS:

Product Name	Quantity	Maker	Cat#
Anti N ^ε -(carboxymethyl) lysine [CML] (2G12) Monoclonal Antibody	100 ul	CAC	AGE-M01
Anti N ^ε -(carboxyethyl) lysine [CEL] (CEL-SP) Monoclonal Antibody	100 ul	CAC	AGE-M02
Anti GA-pyridine (2A2) Monoclonal Antibody	100 ul	CAC	AGE-M03
Anti N ^ω -(carboxymethyl) arginine [CMA] (2G12) Monoclonal Antibody	100 ul	CAC	AGE-M04
CML-BSA/N ^ε -(carboxymethyl) lysine-BSA	200 ul	CSR	AGE-GP01
CEL-BSA/N ^ε -(carboxyethyl) lysine-BSA	200 ul	CSR	AGE-GP02
GA-BSA/Glycolaldehyde-BSA	200 ul	CSR	AGE-GP03
Ribose-gelatin	500 ul	CSR	AGE-GP04
Mild-AGE-BSA	200 ul	CSR	AGE-GP05

ELISA protocol

Coating

- 1) Distribute 100 ul / well of the sample solution (1 ug/mL in PBS) to 96 well microtiter plates (Thermo, MaxiSorp).
- 2) Incubate the plates 2h at RT or overnight at 4 degrees.
- 3) Discard the supernatant of sample solution.
- 4) Wash the plates three times with washing.buf.(PBS/0.05%Tween 20)

Blocking

- 1) Distribute 200 ul / well of 0.5% gelatin-PBS to 96 well microtiter plates
- 2) Incubate the plates 1h at RT.
- 3) Discard the the supernatant of 0.5% gelatin-PBS
- 4) Wash the plates three times with washing.buf.(PBS/0.05%Tween 20)

Primary antibody

- 1) Distribute 100 ul / well of Primary antibodies diluted with washing buf. as suggested in the APPLICATIONS to each well.
- 2) Incubate the plates 1h at RT.
- 3) Discard the supernatant of Primary antibody solution.
- 4) Wash the plates three times with washing.buf.(PBS/0.05%Tween 20)

Secondary antibody

- 1) Distribute 100 ul / well of secondary antibodies (HRP-anti mouse IgG) diluted with washing buf. as suggested in the APPLICATIONS to each well.
- 2) Incubate the plates 1h at RT.
- 3) Discard the supernatant of secondary antibody.
- 4) Wash the plates three times with washing.buf.(PBS/0.05%Tween 20)

OPD color reaction

- 1) Reaction for 2-10 minutes at RT.
- 2) Distribute 100 uL / well of 2M H₂SO₄ to each well and stop enzyme reaction.
- 3) After gentle mixing, determine the absorbance at 492 nm of each well by a spectrophotometer.

For research use only. Not for clinical diagnosis.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619