



## MONOCLONAL ANTIBODY

*For research use only. Not for clinical diagnosis.*

**Catalog No. KMU-M01**

# Anti-Nitroguanosine monoclonal antibody

## BACKGROUND

8-Nitroguanosine is a nitrated nucleic acid which is formed by peroxynitrite, myeloperoxidase, nitrite, and peroxide. It is known that the nitration of guanine is enhanced in virus infection<sup>1, 2</sup>, bacterial infection<sup>3, 4</sup>, inflammatory disease<sup>5</sup>, cancers<sup>5</sup>, and diseases associated with smoking<sup>6</sup>. 8-nitroguanosine is thought to be one of the makers of DNA damage caused by oxidative stress. Cyclic GMP (cGMP) is one of the important substances for the signal transfer. On the other hand, 8-Nitro-cGMP (nitrated cGMP) has been identified *in vivo*<sup>3</sup>. Therefore, 8-Nitro-cGMP can potentially act as a mediator for reactive oxygen signaling<sup>3, 7-9</sup>. Although the product, Anti-Nitroguanosine monoclonal antibody (NO<sub>2</sub>G52), does not cross-react with normal nucleotide bases, it selectively reacts with nitro functionalized nucleotides such as nitro-cGMP, nitro-GMP, and nitro-GTP. Therefore, Anti-Nitroguanosine monoclonal antibody is universal antibody of nitrated guanine which modified 8<sup>th</sup> position of guanine with nitro group. NO<sub>2</sub>G52 is commonly used for immunostaining. In addition, NO<sub>2</sub>G52 can potentially utilized for affinity purification of nitroguanine derivative<sup>3, 6</sup>.

<b>Product type</b>	Primary antibodies
<b>Host</b>	Mouse (BALB/c)
<b>Form</b>	Liquid 1 mg/ml PBS solution; 0.1% ProClin as a preservative
<b>Volume</b>	200 ug
<b>Clone</b>	NO <sub>2</sub> G52
<b>Isotype</b>	IgG1

---

<b>Application notes</b>	IHC, ELISA Immunohistochemistry, 10 ug/ml ELISA, 1 ug/ml Optimal dilutions/concentrations should be determined by the end user.
<b>Cross reactivity</b>	Strongly reacts (10 umol/l) 8-NO <sub>2</sub> -guanosine 8-NO <sub>2</sub> -guanine  Slightly cross-reacts ( >1 mmol/l) 8-Br-guanosine 8-Br-guanine 8-Cl-guanine



No cross-reaction

guanosine, guanine, 8-OH-guanine, 8-OH-deoxyguanosine, xanthine, adenine, adenosine, thymine, deoxythymidine, uracil, uridine, 3-NO<sub>2</sub>-tyrosine, 2-NO<sub>2</sub>-imidazole, cytosine




Fig. 1. The reactivity of Monoclonal Antibody(NO<sub>2</sub>G52).




Fig. 2. Production of 8-nitroguanine by airway epithelial cell with Idiopathic Fibrosis (IPS).

**A.** 8-nitroguanine (green) and iNOS (red) **B.** Bright-field **C.** 8-nitroguanine (green) and 8-oxoguanine (Red) **D.** 8-nitroguanine (green) and mitochondria (red)




Fig. 3. Endogenous guanine nitration in RAW 264.7 cells, a murine macrophage cell line



## Storage

Store below -20°C (below -70°C for prolonged storage).

Aliquot to avoid cycles of freeze/thaw.

## References

- 1) T. Akaike, S. Okamoto, T. Sawa, J. Yoshitake, F. Tamura, K. Ichimori, K. Miyazaki, K. Sasamoto and H. Maeda, 8-nitroguanosine formation in viral pneumonia and its implication for pathogenesis, *Proc. Natl. Acad. Sci. USA*, **100**, 685-690 (2003).
- 2) J. Yoshitake, T. Akaike, T. Akuta, F. Tamura, T. Ogura, H. Esumi, and H. Maeda, Nitric oxide as an endogenous mutagen for Sendai virus without antiviral activity, *J. Virol.*, **78**, 8709-8719 (2004).
- 3) T. Sawa, M. H. Zaki, T. Okamoto, T. Akuta, Y. Tokutomi, S. Kim-Mitsuyama, H. Ihara, A. Kobayashi, M. Yamamoto, S. Fujii, H. Arimoto, and T. Akaike, Protein S-guanylation by the biological signal 8-nitroguanosine 3',5'-cyclic monophosphate, *Nat. Chem. Biol.*, **3**, 727-735 (2007).
- 4) M. H. Zaki, S. Fujii, T. Okamoto, S. Islam, S. Khan, K. A. Ahmed, T. Sawa, and T. Akaike, Cytoprotective function of heme oxygenase 1 induced by a nitrated cyclic nucleotide formed during murine salmonellosis, *J. Immunol.*, **182**, 3746-3756 (2009).
- 5) Y. Terasaki, T. Akuta, M. Terasaki, T. Sawa, T. Mori, T. Okamoto, M. Ozaki, M. Takeya and T. Akaike, Guanine nitration in idiopathic pulmonary fibrosis and its implication for carcinogenesis, *Am. J. Respir. Crit. Care. Med.*, **174**, 665-673 (2006).
- 6) T. Sawa, M. Tatemichi, T. Akaike, A. Barbin and H. Ohshima, Analysis of urinary 8-nitroguanine, a marker of nitrative nucleic acid damage, by high-performance liquid chromatography-electrochemical detection coupled with immunoaffinity purification: association with cigarette smoking, *Free Radic. Biol. Med.*, **40**, 711-720 (2006).
- 7) M. Feelisch, Nitrated cyclic GMP as a new cellular signal, *Nat. Chem. Biol.*, **3**, 687-688 (2007).
- 8) K. A. Ahmed, T. Sawa, T. Akaike, Protein cysteine S-guanylation and electrophilic signal transduction by endogenous nitro-nucleotides, *Amino Acids*, **41**, 123-130 (2011).
- 9) T. Sawa, H. Arimoto and T. Akaike, Regulation of redox signaling involving chemical conjugation of protein thiols by nitric oxide and electrophiles, *Bioconjugate Chem.*, **21**, 1121-1129 (2010).

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)