Anti-Acetylaminofluorene-DNA adducts (Anti-AAF-DNA adducts)

BACKGROUND

DNA adducts in mammalian cells exposed to *N*-acetoxy-2-acetylaminofluorene* (NA-AAF), an activated derivative of the potent carcinogen 2-acetylaminofluorene (2-AAF), play significant roles in cell killing, chromosome aberration, sister-chromatid exchange, gene mutation and neoplastic transformation (1,2). NA-AAF binds covalently to guanine in the DNA of mammalian cells and produces three different DNA adducts. The C-8 adducts, *N*-(deoxyguanosin-8-yl)-2-acetylaminofluorene (dG-C8-AAF) and deacetylated *N*-(deoxyguanosin-8-yl)-2-aminofluorene (dG-C8-AF), account for the major portion of the DNA-bound products, while the minor N2 adduct, 3-(deoxyguanosin-N2-yl)-2-acetylaminofluorene (dG-N2-AAF), accounts for the rest of them. The relative induction levels of the two major C-8 adducts vary among cell types. These adducts distort the DNA helix as do UV-induced cyclobutane pyrimidine dimers (CPD), and therefore they are repaired by nucleotide excision repair in human cells. Iwamoto *et al* (3) have established monoclonal antibodies against dG-C8-AAF in denatured DNA. These antibodies enable one to detect dG-C8-AAFs in DNA from cultured cells using an enzyme-linked immunosorbent assay (ELISA) and to visualize them in cultured cells or rodent tissues using an immunofluorescence (IF). This technology would contribute to understanding of molecular mechanisms in AAF-related research fields including cancer research, anticancer research or toxicology.

*NA-AAF is available from NCI Chemical Carcinogen Reference Standards Repository (Midwest Research Institute, Kansas City, MO)

Product type	Primary antibody
Clone number	AAF-1
Isotype	IgG1 (kappa)
Host	Mouse
Source	The hybridoma was established by fusion of mouse myeloma cells with Balb/c mouse splenocytes immunized with NA-AAF-modified single-stranded DNA conjugated with methylated BSA. This hybridoma (clone AAF-1) culture supernatant was collected and precipitated with ammonium sulfate. After centrifugation, the pellet dissolved in small volume of double-distilled water was dialysed against PBS. The dialysate was then lyophilized.
Form	This antibody is lyophilized form. Reconstitute with 100 µl of distilled water. No preservative is contained.
Volume	100 uL (Lyophilized Powder)
Specificity	AAF-1 primarily binds to dG-C8-AAF in denatured DNA, although dG-C8-AF in DNA is also recognized with slightly less efficiency.
Storage	Lyophilized form (Before reconstitution) : store at -20°C. Reconstituted form : store at -20°C. After reconstitution, it is stable for at least 1 year when stored at -20°C. It should be divided into small quantity to avoid freezing and thawing many times.

Application notes Recommended dilutions

- **Immunohistochemistry:** 1/25
- **ELISA:** 1/100

Other applications have not been tested.

Optimal dilutions/concentrations should be determined by the end user.

References

- 1) R.H. Heflich and R.E. Neft, Genetic toxicity of 2-acetylaminofluorene, 2-aminofluorene and some of their metabolites and model metabolites. *Mutation Res.* 318 (1994) 73-174.
- 2) E.Kriek, Fifty years of research on *N*-acetyl-2-aminofluorene, one of the most versatile compounds in experimental cancer research. *J. Cancer Res. Clin. Oncol.* 118 (1992) 481-489.
- 3) T. Iwamoto *et al.*, In situ detection of acetylaminofluorene-DNA adducts in human cells using monoclonal antibodies. *DNA Repair* 3 (2004) 1475-1482.

ANTIBODY CHARACTERIZATION


Fig.1 The dose-dependent formation of NA-AAF-induced DNA adducts in human cells.
Cells were exposed to NA-AAF for 0.5 h and the formation of DNA adducts in denatured DNA (500 ng/well) was determined using a sensitive-direct-binding ELISA with AAF-1 (1/100). (Details are shown in Ref. 3.)


Fig.2 The formation of NA-AAF-induced DNA adducts in human cells.
Cells were exposed to 200 μM NA-AAF or solvent for 0.5 h. After permeabilization and fixation, DNA adducts were visualized by sequential treatment of AAF-1 (1/25) and Alexa Fluor 488 goat anti-mouse IgG conjugate. Nuclear DNA was counterstained with DAPI. (Details are shown in Ref. 3.)

RELATED PRODUCTS:

Product Name	Maker	Cat#
Anti cyclobutane pyrimidine dimers (CPDs) Monoclonal Antibody (Clone:TDM-2)	CAC	NM-DND-001
Anti (6-4) photoproducts (6-4PPs) Monoclonal Antibody (Clone:64M-2)	CAC	NM-DND-002
Anti Dewar photoproducts (DewarPPs) Monoclonal Antibody (Clone:DEM-1)	CAC	NM-DND-003
Anti Acetylaminofluorene(AAF)-DNA adducts Monoclonal Antibody (Clone:AAF-1)	CAC	NM-MA-001
Anti XPA Monoclonal Antibody (Clone:A-2)	CAC	KUP-TM-M01
Anti XPA Monoclonal Antibody (Clone:5F12)	BAM	70-032
Anti XPF Monoclonal Antibody (Clone:A-2)	CAC	KUP-TM-M02
Anti XPG Monoclonal Antibody (Clone:A-2)	CAC	KUP-TM-M03
Anti ERCC1 Monoclonal Antibody (Clone:A-2)	CAC	KUP-TM-M04
PROTAMINE SULFATE COATED ELISA PLATE 96	CSR	NM-MA-P001
PROTAMINE SULFATE COATED ELISA PLATE 96 x 5	CSR	NM-MA-P002
PROTAMINE SULFATE COATED ELISA PLATE 96 x 10	CSR	NM-MA-P003

PROTOCOLS:

ELISA

A. The coating of microtiter plates by protamine sulfate

- 1) Prepare 0.003% protamine sulfate solution in distilled water and stir for 1 hour.
- 2) Distribute 50 μ L / well of the solution to 96 well microtiter plates (Polyvinylchloride flat-bottom, Thermo, Cat. No. 2801, Milford, MA).
- 3) Incubate the plates at 37°C overnight and coat protamine sulfate on plates by drying completely.
- 4) Wash the plates three times with 100 μ L / well of distilled water.
- 5) These plates can be stored for long times in dark.

B. Cell culture and NA-AAF exposure

- 6) Plate cells in 10-cm dishes and culture one or two days.
- 7) Prepare 100 mM NA-AAF in DMSO and expose cells to NA-AAF (for example ; 0, 25, 50, 100, 150 μ M) for 30 minutes at 37°C.
- 8) Wash cells twice by 10 mL of Dulbecco's PBS (DPBS), and then cells were harvested by a cell scraper from the dishes and centrifuged at 10,000 x g for 15 seconds at 4°C.
- 9) Cell pellets were stored at -80°C until processing.

C. DNA isolation

- 10) Genomic DNA was purified using a QIAamp Blood Kit (QIAGEN, Cat. No. 51104 or 51106). DNA concentrations were calculated from the absorbance at 260 nm.

D. DNA sample coating to the microtiter plates precoated with protamine sulfate

- 11) Prepare sample DNA solutions in PBS at the concentration of 10 μ g / mL
- 12) To denature DNA, heat DNA solutions in a hot plate at 100°C for 10 minutes and chill rapidly in an ice bath for 15 minutes.
- 13) Distribute 50 μ L / well of each denatured DNA solution to protamine sulfate precoated 96 well microtiter plates (use 4 wells for each sample) and dry completely overnight at 37°C.

E. DNA adduct detection

- 14) Wash the DNA-coated plates 5 times with 150 μ L / well of PBS-T (0.05% Tween-20 in PBS).
- 15) Distribute 150 μ L / well of 2% FBS in PBS to each well to prevent non-specific antibody binding.
- 16) Incubate 30 minutes at 37°C
- 17) Wash the plates 5 times with 150 μ L / well of PBS-T.
- 18) Distribute 100 μ L / well of AAF-1 antibodies diluted with PBS as suggested in the **APPLICATIONS** to each well and incubate 30 minutes at 37°C.
- 19) Wash the plates 5 times with 150 μ L / well of PBS-T
- 20) Distribute 100 μ L / well of 1:2000 Biotin-F(ab')₂ fragment of anti-mouse IgG (H+L) (Rockland, Code 710-106-019) diluted with PBS to each well and incubate 30 minutes at 37°C.
- 21) Wash the plates 5 times with 150 μ L / well of PBS-T.
- 22) Distribute 100 μ L / well of 1:10000 Peroxidase-Streptavidin (Zymed, Cat. No. 43-4323) diluted with PBS to each well and incubate 30 minutes at 37°C.
- 23) Wash the plates 5 times with 150 μ L / well of PBS-T.
- 24) Wash the plates once with 150 μ L / well of Citrate-phosphate buffer (pH5.0) [Citric acid monohydrate 5.10 g, Na₂HPO₄ 7.30 g, Distilled water 1000 ml]. Keep the buffer solution in the plates until the next substrate solution is ready.
- 25) After throwing the buffer away, distribute 100 μ L / well of the substrate solution [o-Phenylene diamine 8 mg, H₂O₂ (35%) 4 μ l, Citrate-phosphate buffer (pH 5.0) 20 ml] to each well and incubate 30 minutes at 37°C]
- 26) Distribute 50 μ L / well of 2M H₂SO₄ to each well and stop enzyme reaction.
- 27) After gentle mixing, determine the absorbance at 492 nm of each well by a spectrophotometer.

PROTOCOLS:

Immunofluorescence microscopy

A. Cell culture and NA-AAF exposure

- 1) Culture the cells in the appropriate condition in 35-mm glass-bottom dishes (MatTek, Ashland, MA). (For example, inoculate 4×10^5 cells per dish, then incubate for one or two days in a CO₂ incubator.)
- 2) Prepare 100 mM NA-AAF in DMSO and expose cells to NA-AAF (for example; 0, 200 μ M) for 30 minutes at 37°C. Wash cells 2 times with 2 mL of DPBS.

B. Cell fixation and permeabilization

- 3) Pour 1 mL of 0.5% Triton X-100 in PBS, and permeabilize the cells for 15 minutes on ice.
- 4) Wash the cells 2 times with 2 mL of DPBS.
- 5) Pour 2 mL of methanol-aceton (1 : 1) into each dish, and fix the cells for 10 minutes at -20°C.
- 6) Dry the dishes at room temperature for 5-10 minutes.

C. Indirect Immunofluorescence

- 7) Add 70 μ L of 2M HCL and denature cellular DNA for 30 minutes at room temperature.
- 8) Wash the cells 5 times with 2 mL of PBS.
- 9) Add 70 μ L of 20% FBS in PBS to prevent non-specific antibody binding.
- 10) Incubate 30 minutes at 37°C with gentle shaking.
- 11) Wash the cells 5 times with 2 mL of PBS.
- 12) Add 70 μ L of AAF-1 antibodies diluted with PBS containing 5% FBS as suggested in the **APPLICATIONS** onto the cells and incubate for 30 minutes at 37°C with shaking (Optimization of antibody concentration or incubation condition is recommended if necessary.)
- 13) Wash the cells 5 times with 2 mL of PBS. (Subsequent steps must be done in the dark.)
- 14) Add 70 μ L of 1:100 Alexa Fluor 488-F(ab')₂ fragment of anti-mouse IgG (H+L) (Molecular Probes, Cat. No. A-11017) diluted with PBS containing 5% FBS and incubate for 30 minutes at 37°C with shaking. Distribute 100 μ L / well of 1:10000 Peroxidase-Streptavidin (Zymed, Cat. No. 43-4323) diluted with PBS to each well and incubate 30 minutes at 37°C.
- 15) Wash the cells 5 times with 2 mL of PBS.
- 16) Add 70 μ L of 0.05 μ g/ mL DAPI in PBS and incubate for 5 minutes at 37°C with shaking.
- 17) Wash the cells 5 times with 2 mL of PBS.
- 18) Promptly add 20 μ L of Vectashield mounting medium (Vector, Cat. No. H-1000) onto the cells, then put a cover slip on them.
- 19) Observe fluorescence images by a fluorescence microscope.

120221-1

For research use only, Not for diagnostic use.

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FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

Monoclonal Antibodies against DNA Damage

Powerful tools for studying DNA damage and its biological effects

Monoclonal antibodies against UV-induced DNA Damage

Anti Cyclobutane Pyrimidine Dimers (CPDs) [Clone : TDM-2]

Anti (6-4) photoproducts (6-4PPs) [Clone : 64M-2]

Anti Dewar photoproducts (DewarPPs) [Clone : DEM-1]

Prolonged exposure to solar UV radiation may result in acute and chronic health effects to the skin, eye, and immune system, including skin cancers. These harmful effects are suggested to be closely related to DNA damage. The major types of DNA damage induced by solar UV radiation are cyclobutane pyrimidine dimers (CPDs), (6-4) photoproducts (6-4PPs), and Dewar photoproducts (DewarPPs), which are formed between adjacent pyrimidine nucleotides on the same strand of DNA. These helix-distorting DNA lesions are repaired exclusively by a nucleotide excision repair system in humans. Mori *et al.* have developed and characterized monoclonal antibodies specific for CPDs and for 6-4PPs (1). Matsunaga *et al.* have established and characterized monoclonal antibodies against DewarPPs (2). These antibodies enable one to quantitate photoproducts in DNA purified from cultured cells or from the skin epidermis using an enzyme-linked immunosorbent assay (ELISA) and to visualize and measure photoproducts in DNA in cultured cells or the skin using indirect immunofluorescence. Thus, this technology will contribute to understanding the molecular mechanisms of cellular responses to UV light and DNA damage in many research fields including cancer research, photobiology, dermatology, ophthalmology, immunology, and cosmetology.


Features

- Highly specific for the target lesion
- Research applications include ELISA, IF and IHC
- Useful for research in DNA damage and repair
- Allows visualization of the DNA repair process
- Applicable to a broad range of research fields including cancer research, photobiology, dermatology, ophthalmology, immunology, and cosmetology


(1) Toshio Mori, Misa Nakane, Tsuyoshi Hattori, Tsukasa Matsunaga, Makoto Ihara, Osamu Nikaïdo, Simultaneous establishment of monoclonal antibodies specific for either cyclobutane pyrimidine dimer or (6-4) photoproduct from the same mouse immunized with ultraviolet-irradiated DNA. *Photochem. Photobiol.*, 54: 225-232 (1991).

(2) Tsukasa Matsunaga, Yuri Hatakeyama, Michi Ohta, Toshio Mori and Osamu Nikaïdo, Establishment and characterization of a monoclonal antibody recognizing the Dewar isomers of (6-4) photoproducts. *Photochem. Photobiol.*, 57: 934-940 (1993).


UV-induced major DNA damage

Cyclobutane pyrimidine dimer (CPD)

(6-4) photoproduct (6-4PP)

Dewar photoproduct (DewarPP)

Description	Host	Clone	Application	Cat. No.	Quantity
Anti CPDs	Mouse	TDM-2	ELISA / IC	CAC-NM-DND-001	1 vial
Anti 6-4PPs	Mouse	64M-2	ELISA / IC	CAC-NM-DND-002	1 vial
Anti DewarPPs	Mouse	DEM-1	ELISA / IC	CAC-NM-DND-003	1 vial

COSMO BIO CO., LTD.

Monoclonal Antibodies against DNA Damage

Monoclonal antibodies against UV-induced DNA Damage


Anti CPDs [Clone : TDM-2]

Anti 6-4PPs [Clone : 64M-2]

Anti DewarPPs [Clone : DEM-1]

Immunocytochemistry


In situ visualization of XPB (TFIIH) and RPA at CPD sites after micropore UV irradiation

Katsumi et. al., *J. Invest. Dermatol.* 117: 1156-1161, 2001

The technique of micropore UV irradiation combined with fluorescent antibody labeling is very powerful for examining whether a protein of interest is recruited to the sites of UV-induced DNA damage. Micropore UV irradiation induces UV-damage at localized areas of nuclei using a polycarbonate isopore membrane filter. The polycarbonate blocks UV radiation, and cells are exposed only through the 5 μ m pores of the filter. 0.5 h after micropore UV irradiation, cells were fixed and immunofluorescent double staining for DNA damage and repair protein were performed.

In situ Visualization of XPB and CPD 30 min after micropore UV irradiation


Nishiwaki et. al., *J. Invest. Dermatol.* 122: 526-532, 2004.

Cells were doubly stained for XPB and for CPD 0.5 h after local UV irradiation. In normal MSU-1 cells, XPB foci overlapped with the corresponding CPD foci, indicating that XPB is quickly recruited to the sites of DNA damage for repair. In contrast, no or less bright XPB foci at the DNA damage sites were observed in repair deficient TTD cell lines.


ELISA

A sensitive ELISA for measuring UV-induced DNA damage

Genomic DNA is purified from UV-damaged cells and denatured DNA is used to coat wells of a 96 well plate. The binding of TDM-2 or 64M-2 to DNA damage is detected by sequential treatment with biotinylated 2nd antibody and streptavidin-peroxidase. Then, the absorbance of colored products derived from OPD is measured at 492 nm.

Quantification of DNA damage repair by ELISA

Nakagawa et. al., *J. Invest. Dermatol.* 110: 143-148, 1998.


Normal human cells repair 90% of the initial 6-4PP within 3 h after UV irradiation, while they remove 50% of the initial CPD within 24 h. Both damage are repaired by the same nucleotide excision repair (NER) pathway, but 6-4PP forms bigger distortion in DNA than CPD does, resulting in much more efficient repair. In contrast, repair deficient XP-C cells can not repair both damage at all.

Anti Acetylaminofluorene-DNA Adducts Monoclonal Antibody


Anti AAF-DNA adducts [Clone : AAF-1]

DNA adducts in mammalian cells exposed to N-acetoxy-2-acetylaminofluorene (NA-AAF), an activated derivative of the potent carcinogen 2-AAF, play significant roles in cell killing, chromosome aberration, gene mutation and neoplastic transformation. NA-AAF binds covalently to guanine in the DNA of mammalian cells and produces three different DNA adducts. The C-8 adducts dG-C8-AAF and deacetylated dG-C8-AF account for the major portion of the DNA-bound products, while the minor N2 adduct dG-N2-AAF accounts for the remainder. The relative induction levels of the two major C-8 adducts vary among cell types. These adducts distort the DNA helix and therefore are repaired by nucleotide excision repair in human cells. Our AAF-1 antibodies bind most efficiently to dG-C8-AAF and less efficiently to dG-C8-AF in denatured DNA. The antibodies enable one to detect AAF-DNA adducts in DNA from cultured cells using an enzyme-linked immunosorbent assay (ELISA) and to visualize them in cultured cells or rodent tissues by immunofluorescence (IF). This technology will contribute to understanding of molecular mechanisms in AAF-related research fields including cancer research, anticancer research and toxicology.

AAF-DNA adducts recognized by AAF-1


The dose-dependent formation of NA-AAF-induced DNA adducts in human cells.

Cells were exposed to NA-AAF for 0.5 h and the formation of DNA adducts in denatured DNA (500 ng/well) was determined using a sensitive-direct-binding ELISA with AAF-1 (1/100).

Visualization of NA-AAF-induced DNA adducts in human cells.

Cells were exposed to 200 μM NA-AAF or solvent for 0.5 h. After permeabilization and fixation, DNA adducts were visualized by sequential treatment of AAF-1 (1/25) and ALEXA FLUOR® 488 goat anti-mouse IgG conjugate. Nuclear DNA was counterstained with DAPI.

Description	Host	Clone	Application	Cat. No.	Quantity
Anti AAF-DNA adducts	Mouse	AAF-1	ELISA / IC	CAC-NM-MA-001	1 vial

Useful for ELISA assays with DNA damage antibodies


PROTAMINE SULFATE COATED ELISA PLATE

New


Protamine sulfate is a small cationic protein that binds to negatively charged DNA. Protamine sulfate coated wells capture sample DNA more efficiently; a critical step in the accurate and reproducible determination of DNA damage detection by ELISA.

- Steady DNA binding
- High signal detection of a small amount (low concentration) sample
- Room temperature preservation

* Plate seal 1 sheet

Protamine coating increases DNA-binding

UV- or mock-irradiated DNA (20 ng) was added to plates either coated, or not coated, with protamine sulfate. CPDs were then detected by ELISA using TDM-2 antibody. Protamine sulfate coated wells produced strong dose-dependent CPD signals whereas non-coated wells produced very poor signals.

Description	Cat. No.	Quantity
PROTAMINE SULFATE COATED ELISA PLATE 96	CSR-NM-MA-P001	1 plate
PROTAMINE SULFATE COATED ELISA PLATE 96×5	CSR-NM-MA-P002	5x1 plate
PROTAMINE SULFATE COATED ELISA PLATE 96×10	CSR-NM-MA-P003	10x1 plate

Monoclonal Antibodies against DNA Damage

Antibodies against Nucleotide excision repair (NER) factors


Anti XPA [Clone : A-2]

Anti XPG [Clone : G-26]


Anti XPF [Clone : 19-16]

Anti ERCC1 [Clone : E1-44]

Nucleotide excision repair (NER) is a major repair system for removing a variety of DNA lesions including UV-induced cyclobutane pyrimidine dimer and (6-4) photoproduct as well as chemical-induced bulky base adducts. Defects in the NER system give rise to xeroderma pigmentosum (XP), an autosomal recessive disease characterized by a predisposition to skin cancer and in some cases neurological abnormalities. The early process of human NER, from damage recognition to dual incision (removal of damage-containing oligonucleotides), is accomplished by six core NER factors, XPC-RAD23B, TFIIH, XPA, RPA, XPF-ERCC1 and XPG *in vitro*.

Current Model for the Dual Incision Process of NER

Description	Host	Clone	Application	Cat. No.	Quantity
Anti XPA XPA has an ability to bind to DNA with some preference to damaged DNA and interacts with most of other NER factors. XPA appears to be involved in a proper assembly of preincision complex and verification of damaged DNA strand.	Mouse	A-2	WB	CAC-KUP-TM-M01	100 µl
	Mouse	5F12	WB / ELISA	BAM-70-031	50 µg
	Mouse	5F12	WB / ELISA	BAM-70-032	250 µg
Anti XPF XPF harbors a nuclease domain and forms a stable complex with ERCC1. The ERCC1-XPF complex has a unique ability to make a nick on the DNA strand which makes the transition from duplex to single-stranded DNA in the 5' to 3' direction. In the NER process, ERCC1-XPF is responsible for 5'-incision at a dual incision step.	Mouse	19-16	WB / IF	CAC-KUP-TM-M02	100 µl
Anti XPG XPG is a structure-specific endonuclease with an opposite polarity to ERCC1-XPF and makes a nick on the DNA strand which makes the transition from single-stranded to duplex DNA in the 5' to 3' direction. In the NER process, XPG is responsible for 3'-incision at a dual incision step.	Mouse	G-26	WB	CAC-KUP-TM-M03	100 µl
Anti ERCC1 ERCC1 forms a stable complex with XPF and the heterodimer has an ability to make a nick on the DNA strand which makes the transition from duplex to single-stranded DNA in the 5' to 3' direction. In the NER process, ERCC1-XPF complex is responsible for 5'-incision at a dual incision step.	Mouse	E1-44	WB	CAC-KUP-TM-M04	100 µl

For research use only, Not for diagnostic use.

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- (1) Toshio Mori, Misa Nakane, Tsuyoshi Hattori, Tsukasa Matsunaga, Makoto Ihara, Osamu Nikaido, Simultaneous establishment of monoclonal antibodies specific for either cyclobutane pyrimidine dimer or (6-4) photoproduct from the same mouse immunized with ultraviolet-irradiated DNA. *Photochem. Photobiol.*, 54: 225-232 (1991).
- (2) Tsukasa Matsunaga, Yuri Hatakeyama, Michi Ohta, Toshio Mori and Osamu Nikaido, Establishment and characterization of a monoclonal antibody recognizing the Dewar isomers of (6-4) photoproducts. *Photochem. Photobiol.*, 57: 934-940 (1993).

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