
Product Manual


OxiSelect™ Oxygen Radical Antioxidant Capacity (ORAC) Activity Assay

Catalog Number

STA-345

192 assays

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures



CELL BIOLABS, INC.
Creating Solutions for Life Science Research


Introduction

Oxidative stress is a physiological condition where there is an imbalance between concentrations of reactive oxygen species (ROS) and antioxidants. However, excessive ROS accumulation will lead to cellular injury, such as damage to DNA, proteins, and lipid membranes. The cellular damage caused by ROS has been implicated in the development of many disease states, such as cancer, diabetes, cardiovascular disease, atherosclerosis, and neurodegenerative diseases. Under normal physiological conditions, cellular ROS generation is counterbalanced by the action of cellular antioxidant enzymes and other redox molecules. Because of their potential harmful effects, excessive ROS must be promptly eliminated from the cells by this variety of antioxidant defense mechanisms. Antioxidants include both hydrophilic and lipophilic molecules for metabolizing ROS.

Although the products of ROS-induced oxidative stress are extensively used to monitor the effects of oxidative stress, it is also important to evaluate the antioxidant capacity of biological fluids, cells, and extracts. The Oxygen Radical Absorbance Capacity (ORAC) Assay is a classic tool for measuring the antioxidant capacity of biomolecules from a variety of samples. The ORAC Activity Assay is based on the oxidation of a fluorescent probe by peroxy radicals by way of a hydrogen atom transfer (HAT) process. Peroxyl radicals are produced by a free radical initiator, which quenches the fluorescent probe over time. Antioxidants present in the assay work to block the peroxy radical oxidation of the fluorescent probe until the antioxidant activity in the sample is depleted. The remaining peroxy radicals destroy the fluorescence of the fluorescent probe. This assay continues until completion, which means both the antioxidant's inhibition time and inhibition percentage of free radical damage is a single value. The sample antioxidant capacity correlates to the fluorescence decay curve, which is usually represented as the area under the curve (AUC). The AUC is used to quantify the total peroxy radical antioxidant activity in a sample and is compared to an antioxidant standard curve of the water soluble vitamin E analog Trolox™ (see Assay Principle below).

Cell Biolabs' OxiSelect™ ORAC Activity Assay is a fast and reliable kit for the direct measurement of ORAC antioxidant capacity from cell lysate, plasma, serum, tissue homogenates, and food extracts. Each kit provides sufficient reagents to perform up to 192 assays, including blanks, antioxidant standards and unknown samples. The assay is designed for use in single plate microplate readers as well as readers with high-throughput capabilities. Please read the complete kit insert prior to performing the assay.

Assay Principle



Related Products

1. STA-346: OxiSelect™ HORAC Activity Assay
2. STA-305: OxiSelect™ Nitrotyrosine Protein ELISA Kit
3. STA-310: OxiSelect™ Protein Carbonyl ELISA Kit
4. STA-350: OxiSelect™ Comet Assay (3-Well Slides)
5. STA-320: OxiSelect™ Oxidative DNA Damage ELISA Kit (8-OHdG Quantitation)
6. STA-324: OxiSelect™ Oxidative DNA Damage Quantitation Kit (AP Sites)
7. STA-330: OxiSelect™ TBARS Assay Kit (MDA Quantitation)
8. STA-334: OxiSelect™ HNE-His Adduct ELISA Kit
9. STA-337: OxiSelect™ 8-iso-Prostaglandin F2α ELISA Kit (96 Assays)
10. STA-340: OxiSelect™ Superoxide Dismutase Activity Assay
11. STA-340: OxiSelect™ Superoxide Dismutase Activity Assay
12. STA-341: OxiSelect™ Catalase Activity Assay

Kit Components

1. 96-well Microtiter Plate (Part No. 234501): Two 96-well clear bottom black plates.
2. Fluorescein Probe (100X) (Part No. 234502): One 0.5 mL vial.
3. Free Radical Initiator (Part No. 234503): One 5.0 mL amber bottle.
4. Antioxidant Standard (Trolox™) (Part No. 234504): One 1.5 mL vial of a 1 mM solution.
5. Assay Diluent (4X) (Part No. 234505): One 50 mL bottle.

Materials Not Supplied

1. Sample extracts for testing
2. Deionized water
3. 37°C incubator
4. Bottles, flasks, and conical or microtubes necessary for reagent preparation
5. Reagents and materials necessary for sample extraction and purification
6. 10 µL to 1000 µL adjustable single channel micropipettes with disposable tips
7. 50 µL to 300 µL adjustable multichannel micropipette with disposable tips
8. Multichannel micropipette reservoir
9. Fluorescence microplate reader equipped with a 485 nm excitation filter and 530 nm emission filter

Storage

Upon receipt, store the Fluorescein Probe (100X), Free Radical Initiator and Antioxidant Standard frozen at -20°C and Assay Diluent at 4°C. Aliquot as necessary to avoid multiple freeze/thaws. Store all remaining kit components at room temperature until their expiration dates.

Preparation of Reagents

- **Fluorescein Probe (100X):** Dilute the Fluorescein Probe 1:100 with Assay Diluent. Mix to homogeneity. Label this as 1X Fluorescein Solution. Use only enough Fluorescein Probe as necessary for immediate applications.

Note: Do not store diluted Fluorescein Probe solutions.

- **Assay Diluent (4X):** Dilute the Assay Diluent 1:4 with deionized water. Mix to homogeneity. Use this for all sample and standard dilutions. Store the 1X Assay Diluent at 4°C.

Sample Preparation

Note: Samples should be stored at -70°C prior to performing the assay. Sample should be prepared at the discretion of the user. The following recommendations are only guidelines and may be altered to optimize or complement the user's experimental design.

- **Deproteinized Fractions:** Samples can be deproteinized and have their non-protein fractions assayed. Mix samples with 0.5 M perchloric acid (1:1, v/v), centrifuge at 10,000 x g for 10 minutes at 4°C. Remove the supernatant for measuring the non-protein fraction in the assay.
- **Cell Culture:** Wash cells 3 times with cold PBS prior to lysis. Lyse cells with sonication or homogenation in cold PBS and centrifuge at 10,000 x g for 10 minutes at 4°C. Aliquot and store the supernatant for use in the assay.
- **Lipophilic Fractions:** Lipophilic samples must be treated in order to be soluble in an aqueous environment. Dissolve samples in 100% acetone and then dilute in a solution of 7% beta-cyclodextrin and 50% acetone. Incubate the mixture for 1 hour at room temperature with mixing. Further dilute samples as necessary prior to testing.

- **Plasma:** Collect blood with heparin and centrifuge at 4°C for 10 minutes. Remove the plasma and aliquot samples for testing. Blood plasma or serum should be diluted 100-fold or more with Assay Diluent prior to performing the assay.
- **Tissue Lysate:** Sonicate or homogenize tissue sample on cold PBS and centrifuge at 10,000 x g for 10 minutes at 4°C. Aliquot and store the supernatant for use in the assay.
- **Plasma or Serum:** Dilute 100-fold with Assay Diluent immediately before use.
- **Urine:** Test neat or diluted with Assay Diluent if appropriate.
- **Nutrition Samples:** Results may vary depending on sample source and purification. Dilution and preparation of these samples is at the discretion of the user, but use the following guidelines:

Solid or High Protein Samples: Weigh solid sample and then homogenize after adding deionized water (1:2, w/v). Centrifuge the homogenate at 10-12,000 x g for 10 minutes at 4°C. Recover the supernatant which is the water-soluble fraction. Separately recover the insoluble fraction (pulp) and wash with deionized water. Combine this wash with the supernatant. The pooled supernatant can be diluted with Assay Diluent and used directly in the assay. The pulp is further extracted by adding pure acetone (1:4, w(solid pulp)/v) and mixing at room temperature for 30-60 minutes. Centrifuge the extract/solid at 12,000 x g for 10 minutes at 4°C. Recover the acetone extract and dilute with Assay Diluent as necessary prior to running the assay. The total ORAC value is calculated by combining the results from the water-soluble fraction and the acetone extract from the pulp fraction.

Aqueous Samples: Centrifuge the sample at 5-10,000 x g for 10 minutes at 4°C to remove any particulates. Dilute the supernatant as necessary prior to running the assay. Certain liquids such as juice extracts may be tested without dilution.

Preparation of Antioxidant Standard Curve

1. Prepare fresh standards by diluting the 1 mM Antioxidant Standard stock solution in Assay Diluent. Use only enough Antioxidant Standard as necessary for immediate applications. (Example: Add 50 µL of Antioxidant Standard stock tube to 950 µL of Assay Diluent)
2. Prepare a series of the remaining antioxidant standards according to Table 1 below.

Tubes	Trolox™ Antioxidant Standard (µL)	Assay Diluent (µL)	Trolox™ Concentration (µM)
1	50	950	50
2	40	960	40
3	30	970	30
4	20	980	20
5	10	990	10
6	5	995	5
7	2.5	997.5	2.5
8	0	1000	0

Table 1. Preparation of Standards

Note: Do not store diluted Antioxidant Standard solutions.

Assay Protocol

Note: Each Antioxidant Standard and samples should be assayed in duplicate or triplicate. A freshly prepared standard curve should be used each time the assay is performed.

1. Add 25 μ L of the diluted Antioxidant Standard or samples to the 96-well Microtiter Plate.
2. Add 150 μ L of the 1X Fluorescein Solution to each well. Mix thoroughly. Incubate the plate for 30 minutes at 37°C.
3. Add 25 μ L of the Free Radical initiator into each well using either a multichannel pipette or a plate reader liquid handling system.
4. Shake the plate immediately after adding the Free Radical Initiator for 15 seconds to ensure homogeneity.
5. Immediately begin reading sample and standard wells with a fluorescent microplate reader with an excitation wavelength of 480nm and an emission wavelength of 530nm. Read the wells in increments between 1 and 5 minutes for a total of 60 minutes. Save values for Calculation of Results below.

Note: The final assay values should be less than 10% of the initial values in order for the assay to be completed. If assay values are not less than 10% of initial values, then the assay should be repeated with further dilution of samples.

Example of Results

The following figure demonstrates typical OxiSelect™ ORAC Activity Assay results. One should use the data below for reference only. This data should not be used to interpret or calculate actual sample results.




Figure 1: ORAC Activity Assay Standard Curve.

Calculation of Results

Note: A spreadsheet application or plate reader software can be used to perform the calculations.



1. Calculate the area under the curve (AUC) for each sample and standard using the final assay values and the linear regression formula below.

The AUC can be calculated from the equation below:

$$\text{AUC} = 1 + \text{RFU}_1/\text{RFU}_0 + \text{RFU}_2/\text{RFU}_0 + \text{RFU}_3/\text{RFU}_0 + \dots + \text{RFU}_{59}/\text{RFU}_0 + \text{RFU}_{60}/\text{RFU}_0$$


RFU₀ = relative fluorescence value of time point zero.

RFU_x = relative fluorescence value of time points (eg. RFU₅ is relative fluorescence value at minute five)



2. Calculate the Net AUC by subtracting the Blank AUC from the AUC of each sample and standard.

$$\text{Net AUC} = \text{AUC (Antioxidant)} - \text{AUC (blank)}$$



3. Graph the Net AUC on the y-axis against the Trolox™ Antioxidant Standard concentration on the x-axis (see Figure 1)
4. Calculate the μMole Trolox™ Equivalents (TE) of unknown sample by comparing the standard curve. Results (ORAC value) may be expressed as TE per L or g of sample.

Calculation Example:

25 μL of 10-fold diluted sample is assayed along with 25 μL of each Trolox™ antioxidant standard including blank as described in Assay Protocol. The average AUC is 4.3 for blank and 9.1 for sample.

$$\text{Net AUC} = \text{AUC (Antioxidant)} - \text{AUC (blank)} = 9.1 - 4.3 = 4.8$$

Based on the Trolox™ antioxidant standard curve, the equivalent Trolox concentration is 20 μM, therefore: ORAC value (Sample) = 20 μM x 10 (dilution factor) = 200 μM TE = 200 μMole TE/L

References

1. Ames, B.N., Shigenaga, M.K., and Hagen, T.M. (1993) *Proc. Natl. Acad. Sci. USA* 90: 7915-7922.
2. Cao, G., and Prior, R. (1999) *Methods Enzymol.* 299: 50-62.
3. Huang, D., Ou, B., Hampsch-Woodill, M., Flanagan, J., and Prior, R. (2002) *J. Agric. Food Chem.* 50: 4437- 4444.
4. Huang, D., Ou, B., & Prior, R. (2005) *J. Agric. Food Chem.* 53: 1841-1856.
5. Ou, B., Hampsch-Woodill, M., and Prior, R. (2001) *J. Agric. Food Chem.* 49: 4619-4626.
6. Rice-Evans, C., and Miller, N.J. (1994) *Methods Enzymol.* 234: 279-293.

Trolox™ is a trademark of Hoffman-LaRoche

Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

Contact Information

Cell Biolabs, Inc.
7758 Arjons Drive
San Diego, CA 92126
Worldwide: +1 858 271-6500
USA Toll-Free: 1-888-CBL-0505
E-mail: tech@cellbiolabs.com
www.cellbiolabs.com

©2008: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.