




MS-4 Agarose

A molecular screening agarose for improved resolution of DNA fragments with 500 bp or less, especially sized- primer fragments. MS-4 is a type of agarose unique to Conda.

At 3 % concentration, MS-4 Agarose gives a resolution of DNA fragments similar to gels made with polyacrylamide at concentrations of 8 %. While MS-4 may be dissolved carefully by microwaving, gels are best prepared by autoclaving.



Features:

- Excellent resolution of DNA fragments lower than 500 bp, specially smaller sized- primer fragments.
- Forms a very clear, transparent gel, even at concentrations of 5% or higher.
- Efficient mechanical handling at all concentrations. The chances of gel breaking or cracking when handled are greatly minimized.

Functional Tests:

- DNA resolution: bands appear sharp and finely resolved.
- DNase/RNase activity: none detected.
- Gel background: very low after EtBr staining.
- DNA binding: very low.

Specifications:

	3%	5%
Moisture	≤ 7%	
Ash	≤ 0.3%	
EEO	≤ 0.12%	
Sulfate	≤ 0.11%	
Clarity (NTU)	≤ 7.5	
Gel Strength (g/cm ²)	≥ 500	≥ 1000
Gelling Temperature (°C)	≤ 31	
Melting Temperature (°C)	≤ 70	

Ranges of separation:

3%	80 - 500 bp
4%	30 - 300 bp
5%	10 - 200 bp

These ranges are approximate and have been calculated in 1XTAE buffer.

To achieve the best resolution of MS - 4 Agarose gels, store at 4° - 8°C for 20 hours before use.