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- [1] Givens MD. A clinical, evidence-based approach to infectious causes of infertility in beef cattle. *Theriogenology* (2006) 66, 648-654.
- [2] Heuer C, Healy A, Zerbini C. Economic effects of exposure to bovine viral diarrhea virus on dairy herds in New Zealand. *J Dairy Sci* (2007) 90, 5428-5438.
- [3] Miles DG. Overview of the North American beef cattle industry and the incidence of bovine respiratory disease (BRD). *Anim Health Res Rev* (2009) 10, 101-103.
- [4] Tsuchiaka S, Masuda T, Sugimura S, Kobayashi S, Komatsu N, Nagai M, Omatsu T, Furuya T, Oba M, Katayama Y, Kanda S, Yokoyama T, Mizutani T. Development of a novel detection system for microbes from bovine diarrhea by real-time PCR. *J Vet Med Sci* (2016) 78, 383-389.
- [5] Hilton WM. BRD in 2014: where have we been, where are we now, and where do we want to go? *Anim Health Res Rev* (2014) 15, 120-122.
- [6] Kishimoto M, Tsuchiaka S, Rahpaya SS, Hasebe A, Otsu K, Sugimura S, Kobayashi S, Komatsu N, Nagai M, Omatsu T, Naoi Y, Sano K, Okazaki-Terashima S, Oba M, Katayama Y, Sato R, Asai T, Mizutani T. Development of a one-run real-time PCR detection system for pathogens associated with bovine respiratory disease complex. *J Vet Med Sci* (2017) 79, 517-523.
- [7] Rahpaya SS, Tsuchiaka S, Kishimoto M, Oba M, Katayama Y, Nunomura Y, Kokawa S, Kimura T, Kobayashi A, Kirino Y, Okabayashi T, Nonaka N, Mekata H, Aoki H, Shiokawa M, Umetsu M, Morita T, Hasebe A, Otsu K, Asai T, Yamaguchi T, Makino S, Murata Y, Jan Abi AJ, Omatsu T, Mizutani T. Dembo polymerase chain reaction technique for detection of bovine abortion, diarrhea, and respiratory disease complex infectious agents in potential vectors and reservoirs. *J Vet Sci* (2018) 19, 350-357.
- [8] Mahlum CE, Haugerud S, Shivers JL, Rossow KD, Goyal SM, Collins JE, Faaborg KS. Detection of bovine viral diarrhea virus by TaqMan reverse transcription polymerase chain reaction. *J Vet Diagn Invest* (2002) 14, 120-125.
- [9] Suo B, He Y, Tu SI, Shi X. A Multiplex Real-Time Polymerase Chain Reaction for Simultaneous Detection of *Salmonella* spp., *Escherichia coli* O157, and *Listeria monocytogenes* in Meat Products. *Foodborne Pathogens and Disease* (2010) 7, 619-628.

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- [10] Logar K, Kopinč R, Bandelj P, Starič J, Lapanje A, Ocepek M. Evaluation of combined high-efficiency DNA extraction and real-time PCR for detection of *Mycobacterium avium* subsp. *paratuberculosis* in subclinically infected dairy cattle: comparison with faecal culture, milk real-time PCR and milk ELISA. *BMC Vet Res* (2012) 8, 1-10.
- [11] Horwood PF, Mahony TJ. Multiplex real-time RT-PCR detection of three viruses associated with the bovine respiratory disease complex. *J Virol Methods* (2011) 171, 360-363.
- [12] Boxus M, Letellier C, Kerkhofs P. Real Time RT-PCR for the detection and quantitation of bovine respiratory syncytial virus. *J Virol Methods* (2005) 125, 125-130.
- [13] Sachse K, Salam HSH, Diller R, Schubert E, Hoffmann B, Hotzel H. Use of a novel real-time PCR technique to monitor and quantitate *Mycoplasma bovis* infection in cattle herds with mastitis and respiratory disease. *Vet J* (2010) 186, 299-303.
- [14] Marques LM, Amorim AT, Martins HB, Rezende IS, Barbosa MS, Lobão TN, Campos GB. Timenetsky J. A quantitative TaqMan PCR assay for the detection of *Ureaplasma diversum*. *Vet Microbiol* (2013) 167, 670-674.
- [15] Shirafuji H, Yazaki R, Shuto Y, Yanase T, Kato T, Ishikura Y, Sakaguchi Z, Suzuki M, Yamakawa M. Broad-range detection of arboviruses belonging to Simbu serogroup lineage 1 and specific detection of Akabane, Aino and Peaton viruses by newly developed multiple TaqMan assays. *J Virol Methods* (2015) 225, 9-15.
- [16] Rodríguez-Lázaro D, Hernández M, Scotti M, Esteve T, Vázquez-Boland JA, Pla M. Quantitative detection of *Listeria monocytogenes* and *Listeria innocua* by real-time PCR: assessment of *hly*, *iap*, and *lin02483* targets and AmpliFluor technology. *Appl Environ Microbiol* (2004) 70, 1366-1377.
- [17] Luong ML, Clancy CJ, Vadnerkar A, Kwak EJ, Silveira FP, Wissel MC, Grantham KJ, Shields RK, Crespo M, Pilewski J, Toyoda Y, Kleiboeker SB, Pakstis D, Reddy SK, Walsh TJ, Nguyen MH. Comparison of an *Aspergillus* real-time polymerase chain reaction assay with galactomannan testing of bronchoalveolar lavage fluid for the diagnosis of invasive pulmonary aspergillosis in lung transplant recipients. *Clin Infect Dis* (2011) 52, 1218-1226.
- [18] Stoddard RA, Gee JE, Wilkins PP, McCaustland K, Hoffmaster AR. Detection of pathogenic *Leptospira* spp. through TaqMan polymerase chain reaction targeting the *LipL32* gene. *Diagn Microbiol Infect Dis* (2009) 64, 247-255.
- [19] Wernike, K., Hoffmann, B., Kalthoff, D., König, P. and Beer, M. 2011. Development and validation of a triplex real-time PCR assay for the rapid detection and differentiation of wild-type and glycoprotein E-deleted vaccine strains of Bovine herpesvirus type 1. *J. Virol. Methods* 174: 77-84.