# Extracellular Vesicles derived from *Saccharomyces cerevisiae* 1611-4 strain

Cat. No. YSEV-R5

Updated on November 4, 2025

[www.cosmobiousa.com](http://www.cosmobiousa.com)

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Exosomes released by cells in multicellular organisms are a type of extracellular vesicles (EVs). They are expected to be applied in pharmaceuticals and cosmetics, and developments are underway in various fields. Meanwhile, similar EVs are also released from unicellular microorganisms, which are spherical structures ranging from 20 to 400 nm. These EVs are involved in inter-microbial and microbe–host cell communication, similar to exosomes.<sup>1</sup> EV production is an essential function for microorganisms, and elucidating the roles of EVs in bacterial interactions or host interactions is expected to lead to applications in vaccine development, gut microbiota research, and drug delivery systems (DDS).<sup>2</sup>

Yeasts, which are unicellular fungi, have also been reported to produce EVs, and it has been suggested that they are involved in the transport of various bioactive substances.<sup>3-5</sup> *Saccharomyces cerevisiae* is widely used as a versatile commercial yeast with high alcohol tolerance, selected for brewing and baking. This product is EVs purified from the culture supernatant of the *Saccharomyces cerevisiae* 1611-4 strain using ultrafiltration.\*

\* The strain *Saccharomyces cerevisiae* 1611-4 was independently isolated, cultured, and identified from plant materials and local food ingredients by Professor Akihiro Yamaguchi, Ms. Kurumi Kameda, and laboratory members at the Laboratory of Applied Microbiology, Rakuno Gakuen University.

## 【 II 】 Product Information

COSMO BIO CO., LTD (CSR)

Cat. No.	Description	Quantity	Particle Concentration	Storage
YSEV-R5	Extracellular Vesicles derived from <i>Saccharomyces cerevisiae</i> 1611-4 strain	200 $\mu$ L	$> 1 \times 10^{10}$ particles/mL in PBS Filter sterilized	4°C

Note: Particle number varies by lot. Please refer to the attached CoA.

## 【 III 】 Usage

- Sterilized with a 0.22  $\mu$ m membrane filter.
- Dilute with appropriate buffer or medium ( $10 \times$  –  $2000 \times$ ) according to your experiment.
- Store at 4°C : use as soon as possible after opening.

## 【IV】 Reference Data


Figure 1. Particle size distribution of EVs derived from *Saccharomyces cerevisiae* 1611-4 strain

## 【V】 References

1. Obana, N., Kurosawa, M., Toyofuku, M. & Nobuhiko, N. Biogenesis and Functions of Membrane Vesicles Actively Produced by Microbes. *KAGAKU TO SEIBUTSU* **54**, 812–819 (2016).
2. Obana, N. & Nomura, N. Functions and biosynthesis of membrane vesicles produced actively by Gram-positive bacteria. *Japanese J. Lact. Acid Bact.* **27**, 10–16 (2016).
3. Oliveira, D. L. et al. Characterization of yeast extracellular vesicles: Evidence for the participation of different pathways of cellular traffic in vesicle biogenesis. *PLoS One* **5**, e11113 (2010).
4. Zhao, K. et al. Extracellular vesicles secreted by *Saccharomyces cerevisiae* are involved in cell wall remodelling. *Commun. Biol.* **2**, (2019).
5. Rizzo, J., Rodrigues, M. L. & Janbon, G. Extracellular Vesicles in Fungi: Past, Present, and Future Perspectives. *Front. Cell. Infect. Microbiol.* **10**, (2020).

COSMO BIO Co., LTD.

【JAPAN】  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [support@cosmobiousa.com](mailto:support@cosmobiousa.com)  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600

13304

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# Extracellular Vesicles derived from *Saccharomyces cerevisiae* 1611-4 strain

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Extracellular Vesicles derived from Yeast *Saccharomyces cerevisiae* 1611-4 strain  
Cat. No. YSEV-R5

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