EV-Capture™ EV Purification Spin Column Kit

Cat. No. EVP01-010

Updated on November 4, 2025

www.cosmobiousa.com

【 I 】 Background and Measurement Principle

Extracellular vesicles (EVs) are lipid bilayer-enclosed particles released from cells, and they play an important role in intercellular communication. They are being actively studied. Ultracentrifugation is the most commonly used method for purifying EVs from biological samples or culture supernatants, but it requires expensive equipment and lengthy procedures. Therefore, various alternative methods have been developed, making it essential to select the most appropriate method depending on the application.

This product separates EVs from other components by taking advantage of the fact that the EV surface is negatively charged. It is a spin column type that can be used with a benchtop centrifuge, allowing for inexpensive and easy recovery of intact EVs. The recovered EVs can be used for particle confirmation, characterization, and functional analysis.

【 II 】 Kit Components

Storage temperature : Room temperature

Component	Quantity
EV-Capture™ EV Purification Spin Column	10
2mL capless tubes	10
Bottom caps	10
Sample Buffer (X10)	10 mL × 1
Elution Buffer	8 mL × 1

Materials required but not provided

- Compact Benchtop Centrifuge
- Rotator
- Tween20

Depending on the sample, recovery yield may be low. If it does not affect subsequent experiments, adding Tween20 to the elution buffer at a final concentration of 0.01% can improve recovery.

- Microtubes

Tubes for collecting the purified product. If necessary, please use low-binding tubes such as those with hydrophilic coatings.

- Dialysis Tubes or Ultrafiltration Columns

The purified product contains a high concentration of salt. If required, perform solvent exchange using dialysis or ultrafiltration.

【III】 Protocol

- ① Dilute Sample Buffer (X10) 10-fold with ultrapure water to prepare the required volume.
※ Note: 4 mL of Sample Buffer is used per sample for equilibration and washing. Calculate the total volume required, including dilution.
- ② Break off the bottom tab of the EV-Capture™ column and set it in a 2 mL capless tube
- ③ Centrifuge at 5000 rpm for ~15 sec to remove the storage solution. Discard the collected storage solution.
- ④ Remove the top cap and reset the column in a 2 mL capless tube
- ⑤ Add 300 µL of Elution Buffer, vortex for 10 sec to soak the resin.
- ⑥ Centrifuge at 5000 rpm for ~15 sec to remove the Elution Buffer. Discard the collected solution.
- ⑦ Add 500 µL of Sample Buffer, centrifuge at 5000 rpm for 15 sec, and discard the buffer. Repeat 3 times to equilibrate the resin. dilute the sample with buffer before application.

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- ⑧ Add sample to the column, close the cap.
 - * Use supernatant obtained by centrifugation (2000–10000 rpm) to remove larger particles.
 - * Filter samples with appropriate pore size membranes.
 - * If the pH of the sample exceeds 9, the recovery yield may decrease. Please neutralize or dilute the sample with buffer before application.
 - * Apply 300–700 µL of sample per column.
 - * Serum samples should be diluted 2–6 fold with Sample Buffer; do not exceed 200 µL net volume
- ⑨ Mix by hand up and down to disperse resin.
- ⑩ Rotate at room temperature for 10–15 min.
- ⑪ Remove bottom cap and reset in a 2 mL capless tube.
- ⑫ Centrifuge at 5000 rpm for ~15 sec to remove flow-through. Discard the flow-through.
- ⑬ Add 500 µL Sample Buffer, vortex for 10 sec, centrifuge at 5000 rpm for 15 sec, discard buffer. Repeat twice.
- ⑭ Wash resin with 500 µL Sample Buffer, centrifuge at 5000 rpm for 15 sec, discard. Repeat 3 times.
 - * If resin adheres to the upper lid or the column wall, wash it down so that the resin falls to the bottom of the column.
- ⑮ Discard the 2 mL cap-less tube and attach a microtube to the column for collecting the purified product.
- ⑯ Add 150 µL of Elution Buffer (optionally add 0.01% Tween20 to improve yield).
- ⑰ Vortex for 10 sec.
- ⑱ Centrifuge at 5000 rpm for ~30 sec to recover purified EVs
- ⑲ Repeat twice and combine eluates as purified EV sample.
 - * Note: Resin fragments may contaminate eluates; use supernatant after centrifugation for nanoparticle tracking analysis.

【IV】 Additional Information

Example applications are available online. Please visit COSMO BIO's website (<https://www.cosmobio.co.jp>) and search for 'EV-Capture'.

COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: support@cosmobiousa.com
URL: www.cosmobiousa.com
Phone/FAX: (+1) 760-431-4600

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EV-Capture™ EV Purification Spin Column Kit

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