Gram-positive bacterial EV ELISA Kit

Cat. No. EVEL01

Updated on November 4, 2025

www.cosmobiousa.com

【 I 】 Background and Principle

Extracellular vesicles (EVs) are particles with lipid bilayers released from cells. They are considered to play important roles in intercellular communication¹ and are actively studied.

It has also been clarified that EVs are responsible for communication not only between multicellular organisms but also between microbes, and even between microbes and host cells.² EV production is an essential function for microorganisms, and elucidation of the functions of EV-mediated bacterial-bacterial and bacterial-host interactions is expected to contribute to vaccine development, host-microbe interaction studies in the gut, drug delivery systems, and other applied fields.³

For quantification of EVs from multicellular organisms, ELISA kits using antibodies against marker proteins have been developed. However, for microbes, marker proteins have not been well identified, and quantification by nanoparticle tracking analysis is mainstream, which requires expensive equipment.

This product captures EVs using a plate coated to carry positive charges on the surface and detects them with an anti-LTA (Lipoteichoic acid) antibody as a probe.

This allows sensitive relative quantification of Gram-positive bacterial EVs without requiring expensive equipment or special techniques.

【 II 】 Product Features

- Enables relative quantification of microbe-derived EVs using the same procedure as conventional ELISA
- Standards are included, allowing immediate start of experiments.
- High sensitivity enables measurement with small sample volumes.

Limit of detection : 4.0×10^5 particles/mL

【III】 Kit Components

Storage temperature : 4°C

No	Component	Volume	Quantity
1	EV Detection Plate	8-well × 12 strips	1 plate
2	Lactobacillus EVs Standard (<i>Lactobacillus paracasei</i> 180913-R1 strain) 1.28 × 10 ⁷ particles/mL	1.0 mL	1 vial*
3	Assay Buffer	35 mL	1 bottle
4	Wash Buffer (10X)	40 mL	1 bottle
5	Anti-LTA Antibody (100X, red cap)	120 µL	1 vial
6	HRP-labeled Secondary Antibody (200X, green cap)	150 µL	1 vial
7	Substrate Solution	12 mL	1 bottle
8	Stop Solution	6 mL	1 bottle
9	Plate Seal		3 sheets

* Standard curve can be prepared twice (n=2)

Materials required but not provided

- Micropipettes (10 – 1000 µL)
- Multichannel pipette
- Reservoir
- Plate shaker
- Plate reader (capable of measuring at 450 nm)
- Plate washer

【IV】 Preparation of Reagents and Samples

【IV – 1】 EV Standards (prepare for 2 wells per plate)

	Concentration (particles/mL)	EVs Standard	Assay Buffer	Dilution
A	1.28 × 10 ⁷			1
B	6.40 × 10 ⁶	250 µL of A	250 µL	2
C	3.20 × 10 ⁶	250 µL of B	250 µL	2
D	1.60 × 10 ⁶	250 µL of C	250 µL	2
E	8.00 × 10 ⁵	250 µL of D	250 µL	2
F	4.00 × 10 ⁵	250 µL of E	250 µL	2

Preparation of Standard Solutions

1. Add 250 µL of Assay Buffer to 250 µL of the EVs Standard Solution included in the kit (A in the table above) to make a 2-fold dilution. Mix well to obtain solution B.
2. Prepare subsequent 2-fold serial dilutions in the same manner.
3. Use 100 µL of each dilution per well.

Note: Since each concentration is measured in duplicate (n = 2), 200 µL of each dilution will be required.

Caution: Prepare the EVs Standard Solution freshly at the time of use, and only in the amount needed.

【IV】 Reagents and Sample Preparation (cont.)

【IV – 2】 Wash Buffer

- Dilute Wash Buffer (10X) 10-fold with purified water.

Example: For one plate, dilute 40 mL Wash Buffer (10X) with 360 mL purified water.

【IV – 3】 Anti-LTA Antibody

- Dilute Anti-LTA Antibody (100X) 100-fold with Assay Buffer.

Example: For one plate, dilute 100 µL Anti-LTA Antibody (100X) with 10 mL Assay Buffer.

Caution: Prepare the antibody solution freshly at the time of use, and only in the required amount.

【IV – 4】 HRP-labeled Secondary Antibody

- Dilute HRP-labeled Secondary Antibody (200X) 200-fold with Assay Buffer.

Example: For one plate, dilute 50 µL HRP-labeled Secondary Antibody (200X) with 10 mL Assay Buffer.

Caution: Prepare the antibody solution freshly at the time of use, and only in the required amount.

【IV – 5】 Samples

- Dilute samples at least 2-fold with Assay Buffer so that the total reaction volume becomes 100 µL
- EVs may not be captured if the solvent contains high salt concentrations. Replace or dilute the solvent with phosphate buffer or an equivalent buffer to reduce the salt concentration to 150 mM or lower.
- For some samples, the relationship between concentration and absorbance may not be linear even if within the calibration curve. Prepare multiple dilutions to confirm linearity.

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- ① Bring EV Detection Plate and reagents to room temperature.
- ② Prepare EV standards as described in [IV-1].
- ③ Add 100 µL of prepared EV standards (4.00×10^5 to 1.28×10^7 particles/mL) or sample solution to each well.
- ④ Seal the plate and shake with a plate shaker (800 rpm, 30 seconds).
- ⑤ Incubate at room temperature for 2 hours or overnight at 4° C.
- ⑥ Remove reaction solution completely, add 300 µL Wash Buffer (see [IV-2]) to each well, and wash. Repeat 3 times.
- ⑦ Add 100 µL of diluted Anti-LTA Antibody (see [IV-3]) to each well.
- ⑧ Seal the plate, shake at 800 rpm for 30 seconds, and incubate at room temperature for 2 hours.
- ⑨ Remove antibody solution completely, wash each well 3 times with 300 µL Wash Buffer.
- ⑩ Add 100 µL of diluted HRP-labeled Secondary Antibody (see [IV-4]) to each well.

【V】 Assay Procedure (cont.)

- ⑪ Seal the plate, shake at 800 rpm for 30 seconds, and incubate at room temperature for 2 hours.
- ⑫ Remove antibody solution completely, wash each well 3 times with 300 μ L Wash Buffer.
- ⑬ Add 100 μ L of Substrate Solution to each well and incubate at room temperature for about 20 minutes, or until sufficient color develops.
- ⑭ After confirming the color intensity, add 50 μ L of Stop Solution to each well.
- ⑮ Measure absorbance of each well with a plate reader at 450 nm.
- ⑯ Draw a calibration curve using the concentrations and absorbance of standards, and calculate the concentrations of the samples.

【VI】 References

1. Y. Naito, Y. Yoshioka, Y. Yamamoto and T. Ochiya: Cell Mol Life Sci., 74, 697 (2017).
2. Obana, N., Kurosawa, M., Toyofuku, M. & Nobuhiko, N. Biogenesis and Functions of Membrane Vesicles Actively Produced by Microbes. KAGAKU TO SEIBUTSU 54, 812–819 (2016).191.
3. Obana, N. & Nomura, N. Functions and biosynthesis of membrane vesicles produced actively by Gram-positive bacteria. Japanese J. Lact. Acid Bact. 27, 10–16 (2016).

【VII】 Notes

Examples are available online.

Please search “Gram-positive bacterial EV ELISA” on our website.

COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: support@cosmobiousa.com
URL: www.cosmobiousa.com
Phone/FAX: (+1) 760-431-4600

13746

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Gram-positive bacterial EV ELISA Kit

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