# Tear Mucin Assay Kit (O-Glycan Assay Method)

For Research Use Only

Cat. No. MUC01

Edition Date: 2024/11/01

[www.cosmobiousa.com](http://www.cosmobiousa.com)

## 【 1 – 1 】 Background and Principle of Measurement

Mucin is a glycoprotein that is a major component of mucus, including tears, saliva, gastric fluid, and intestinal fluid. The basic structure of mucin is a core protein with a repeating structure of 10-80 amino acid residues, to which O-glycans are attached. The molecular weight of mucin ranges from 1 to 10 million, and 50-80% of it is composed of glycans. The sugar chains attached to mucins are composed of N-acetylgalactosamine (GalNAc), which is the reducing end of the sugar chain, attached to serine and threonine of the core protein, forming a branched structure with a wide variety of structures. The negatively charged ends of the mucin sugar chains keep the surface of the cornea hydrophilic (Fig. 1).


Fig. 1 Structure of Mucin

Some of these sugar chains are known to have specific recognition functions for viruses and fungi<sup>1)</sup>. Because of these properties, mucins are also known to prevent viruses, pathogens, and fungus-derived toxins from migrating beyond the corneal epithelial layer into the corneal tissue, i.e., they are positioned as a corneal barrier function substance.

In addition, mucin in tear fluid plays an important role in tear fluid stability, and a decrease in mucin is thought to lead to the onset of dry eye<sup>2),3)</sup> (Figure 2). Mucins in tear fluid can be broadly classified into secretory-type mucins and membrane-type mucins. Two types of secretory-type mucins are known: MUC5AC (secreted from conjunctival goblet cells) and MUC7 (secreted from the tear film), while membrane-type mucins include MUC1, MUC4 and MUC16 (on the cell membrane of the corneal epithelium)<sup>4)</sup>.

This kit can measure mucins in tear fluid under alkaline conditions by measuring the fluorescence intensity obtained by  $\beta$ -elimination of O-glycans from the core protein and fluorescence labeling of the reducing end of the sugar chain<sup>5),6)</sup>.



Fig. 2 Types of mucins expressed on the ocular surface and barrier function of mucins in tear fluid

## 【 1 – 2】 Kit components

Storage temperature: 4 to 10°C

This product can measure 50 samples.

Component	Volume	Quantity	Precautions for handling
Elution Buffer	30mL	1	Wear protective equipment such as glasses and gloves when handling, and take care to avoid contact with the body.
Slurry for spin column	45mL	1	
Standard solution GalNAc 100 µg/mL	1mL	1	
Reagent A	0.3mL	1	
Reagent B	1.5mL	1	
Stop solution	15mL	1	
Empty column	-	50	
Centrifuge tube	-	50	


## What to prepare

- Schirmer test strips
- Micro test tubes (2.0 mL, 1.5 mL, 0.5 mL)
- A fluorescence spectrophotometer (using a microcell) or a fluorescence plate reader is required.
- When using a fluorescence plate reader, please prepare a 96 well black plate

## 【II – 1】 Preparation of standard solution

Dilute the standard solution (N-acetylgalactosamine 100 µg/mL) with the extraction solution as follows.

- ① 100 µg/mL   ② 50 µg/mL   ③ 25 µg/mL   ④ 12.5 µg/mL   ⑤ 6.25 µg/mL   ⑥ 3.125 µg/mL   ⑦ 1.5625 µg/mL  
⑧ Blank(Elution Buffer)

Number the tubes from ② to ⑧ and dispense 500 µL of the extract into a new 1.5 mL micro tube.

Take 500 µL from ① and add it to ②, and mix well. Take 500 µL from ② and add to ③. Repeat the same procedure up to ⑦ to make a standard solution. The blank in ⑧ is the Elution Buffer only.

The prepared standard solution can be stored at 4°C for 3 weeks.

## 【II – 2】 Sample Preparation Method

Collect tear fluid for 5 minutes using Schirmer test strips. When storing the collected tear fluid on Schirmer test strips, place the strips in 2.0 mL microtubes and store at 4°C to avoid freezing. If stored for a long period of time, store at -20°C or below.

### 【III】 Assay method

- 1 Fold a small piece of Schirmer test paper, press it into the bottom of a 2.0 mL micro tube, add 200  $\mu\text{L}$  of the extract solution, and mix it with a test tube mixer for 30 seconds. Extract at room temperature for 1 hour, stirring occasionally, while the Schirmer test strips are fully immersed in the extraction solution.
- 2 During the 1-hour extraction, prepare an empty column for purification. After the slurry gel filtration carrier has been well stirred, 800  $\mu\text{L}$  is placed in the Empty column and set in a centrifuge tube without lid. Centrifuge at  $800 \times g$  for 10 minutes at room temperature using a swinging rotor centrifuge. At this time, remove any liquid in the centrifuge tube (see Figure 3).
- 3 Charge 50  $\mu\text{L}$  of the extract from the Schirmer test paper obtained in step 1 to the top of the gel in the Empty column prepared in step 2. Centrifuge at  $800 \times g$  for 10 minutes at room temperature. Collect the eluted solution in a centrifuge tube. Approximately 50  $\mu\text{L}$  of the eluted solution can be collected, this solution is used as the sample solution (see Figure 3).
- 4 Take 20  $\mu\text{L}$  of Sample Mucin Solution and Standard Solution and add in a 0.5 mL microtest tube. Just before use, mix Reagent A: Reagent B at a ratio of 1:5 (v/v)\* 1, add 24  $\mu\text{L}$  of this solution, mix, and then heat at  $100^\circ\text{C}$  for 30 minutes.
- 5 Allow the mixture to leave at room temperature, and after confirming once it has cooled to room temperature, add 200  $\mu\text{L}$  of the reaction stopper solution and stir.
- 6 Dispense 100  $\mu\text{L}$ /well onto a black plate, and measure at excitation wavelength of 336 nm and fluorescence wavelength of 383 nm with a fluorescence plate reader (if a fluorescence plate reader is not available, measure one sample at a time with a fluorescence spectrophotometer using a micro cell).
- 7 Create a standard curve (4-parameter logistic curve: 4PL) from the fluorescence values of the standard solution and calculate the concentration of mucin in the sample solution

\* Use the mixture of Reagent A and Reagent B immediately.


Fig. 3 Purification method by spin column

## 【IV】 Standard curve



Fig. 4 Standard Curve (4PL)

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Reference data 1 Changes in excitation and fluorescence wavelengths in standard solutions at various concentrations

**A:** Excitation spectrum for fluorescence wavelength (Em) 383 nm

**B:** Fluorescence spectrum for excitation wavelength (Ex) 336 nm

The recommended measurement wavelengths are Ex 336 nm and Em 383 nm, but in the case of fluorescence plate readers that use interference filters, it may not be possible to measure at the recommended measurement wavelengths depending on the half-value bandwidth of the filter.

In this case, please increase the wavelength difference between the excitation wavelength and the fluorescence wavelength and measure.

## 【VI】 References

---

- [1] Gipson IK, Hori Y, Argüeso P. Character of Ocular Surface Mucins and Their Alteration in Dry Eye Disease. *Ocul Surf.* 2004;2(2):131-148.
- [2] Argüeso P, Balaram M, Spurr-Michaud S, Keutmann HT, Dana MR, Gipson IK. Decreased Levels of the Goblet Cell Mucin MUC5AC in Tears of Patients with Sjögren Syndrome. *Invest Ophthalmol Vis Sci.* 2002;43(4):1004-1011.
- [3] Uchino Y, Uchino M, Yokoi N, et al. Alteration of Tear Mucin 5AC in Office Workers Using Visual Display Terminals: The Osaka Study. *JAMA Ophthalmol.* 2014;132(8):985-992.
- [4] Inatomi T, Spurr-Michaud S, Tisdale AS, Zhan Q, Feldman ST, Gipson IK. Expression of secretory mucin genes by human conjunctival epithelia. *Invest Ophthalmol Vis Sci.* 1996;37(8):1684-1692.
- [5] Honda Susumu, Matsuda Yoshikazu, Takahashi Masaye, Kakehi Kazuaki, Ganno Shigetake. Fluorimetric determination of reducing carbohydrates with 2-cyanoacetamide and application to automated analysis of carbohydrates as borate complexes. *Anal Chem.* 1980;52(7):1079-1082.
- [6] Crowther RS, Wetmore RF. Fluorometric assay of O-linked glycoproteins by reaction with 2-cyanoacetamide. *Anal Biochem.* 1987;163(1):170-174.

COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioua.com](mailto:info@cosmobioua.com)  
URL: [www.cosmobioua.com](http://www.cosmobioua.com)  
Phone/FAX: (+1) 760-431-4600

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# Tear Mucin Assay Kit (O-Glycan Assay Method)

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[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

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- [1] Gipson IK, Hori Y, Argüeso P. Character of Ocular Surface Mucins and Their Alteration in Dry Eye Disease. *Ocul Surf.* 2004;2(2):131-148.
- [2] Argüeso P, Baram M, Spurr-Michaud S, Keutmann HT, Dana MR, Gipson IK. Decreased Levels of the Goblet Cell Mucin MUC5AC in Tears of Patients with Sjögren Syndrome. *Invest Ophthalmol Vis Sci.* 2002;43(4):1004-1011.
- [3] Uchino Y, Uchino M, Yokoi N, et al. Alteration of Tear Mucin 5AC in Office Workers Using Visual Display Terminals: The Osaka Study. *JAMA Ophthalmol.* 2014;132(8):985-992.
- [4] Inatomi T, Spurr-Michaud S, Tisdale AS, Zhan Q, Feldman ST, Gipson IK. Expression of secretory mucin genes by human conjunctival epithelia. *Invest Ophthalmol Vis Sci.* 1996;37(8):1684-1692.
- [5] Honda Susumu, Matsuda Yoshikazu, Takahashi Masaye, Kakehi Kazuaki, Ganno Shigetake. Fluorimetric determination of reducing carbohydrates with 2-cyanoacetamide and application to automated analysis of carbohydrates as borate complexes. *Anal Chem.* 1980;52(7):1079-1082.
- [6] Crowther RS, Wetmore RF. Fluorometric assay of O-linked glycoproteins by reaction with 2-cyanoacetamide. *Anal Biochem.* 1987;163(1):170-174.

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