# IMMUNO SHOT -Platinum-

Protein-free immuno-reaction enhancing solution

## Instruction

### ---- Contents ----

(1) Introduction .....	2
(2) Products .....	2
(3) How to use IMMUNO SHOT -Platinum- .....	2
(4) Western Blotting .....	3
(5) ELISA .....	3
(6) Trouble shooting .....	4
(7) Contact Information .....	4

#### Cautions

1. Research use only. Do not use for medical purpose.
2. Do not dilute or add other agents in IMMUNO SHOT -Platinum- solution to get the best result.

#### (1) Introduction

IMMUNO SHOT -Platinum- is, like IMMUNO SHOT, an enhancer of antigen-antibody reaction. In Western blotting and ELISA, researchers often experiences weak signal or high background. IMMUNO SHOT -Platinum- improves these problems by just using it as antibody diluents. Due to the principle of working mechanism, IMMUNO SHOT -Platinum- can be used in many assay systems that use antigen-antibody reaction.

IMMUNO SHOT -Platinum- is, unlike IMMUNO SHOT, consisted of all chemical components and contains no proteins. So this product is useful for assay systems that have to be avoided protein contamination.

#### How IMMUNO SHOT -Platinum- Works

IMMUNO SHOT -Platinum- contains polymers. One polymer, by changing the physicochemical properties of antigen and antibody, enhances the mutual accessibility, and facilitate the specific reaction. The other polymer reduces non-specific binding of antibody. Thus, IMMUNO SHOT -Platinum- enhances the antigen-antibody reaction while reducing background.

#### • Features of IMMUNO SHOT -Platinum- •

##### 1. High signal with low background

IMMUNO SHOT -Platinum- enhances the antigen-antibody reaction. Comparing with the method using detergent-containing buffer, several to over 10-fold stronger signal can be obtained while the background level is low. Thus, you can get much higher S/N ratio than usual method.

##### 2. Effective for saving antibody usage and time of reaction time

Because higher signal can be obtained using IMMUNO SHOT -Platinum-, you can reduce the amount of antibody used and the time required for reactions.

##### 3. Can be used for many reactions

IMMUNO SHOT -Platinum- can be used not only for Western blotting and ELISA, but also for other assay systems using antigen-antibody reactions. In addition, IMMUNO SHOT -Platinum- does not affect activities of HRP or ALP, and can be used for assay systems using these enzymes.

##### 4. Easy to use

IMMUNO SHOT -Platinum- is one solution format and formulated as to Ready to Use. Just exchange your dilution buffer of antibodies to IMMUNO SHOT -Platinum- solution.

#### (2) Products

IMMUNO SHOT family has following products. This manual applies all these products.

Cat number	Product name	Quantity
CSR-IS-P-250	IMMUNO SHOT -Platinum-	250ml
CSR-IS-P-500	IMMUNO SHOT -Platinum-	500ml

#### (3) Related products

Some of related products are listed below. To make 100% protein-free assay system we recommend Blocking Reagent for WB / ELISA-Chemically Defined-.

Cat number	Product name	Quantity
CSR-IS-CD-500W	Blocking Reagent for WB-Chemically Defined-	500ml
CSR-IS-CD-500E	Blocking Reagent for ELISA -Chemically Defined-	500ml

(4) Principle of usage of IMMUNO SHOT -Platinum-

- IMMUNO SHOT -Platinum- is one solution format. Just use the same solution for the dilution of the 1st and 2nd antibody. No change is required for the other assay protocol. For details see the later instructions.
- Do not dilute the IMMUNO SHOT -Platinum- solution. The solution is optimized. Dilution causes loss of IMMUNO SHOT -Platinum-'s ability.
- IMMUNO SHOT -Platinum- can be used in Western blotting, antibody sandwich ELISA with either the 1st or 2nd antibody-labeled type, antigen sandwich ELISA.

(5) Western Blotting (WB)

WB is a method to detect proteins by specific antibodies. Usually, proteins are separated by SDS-PAGE and transferred to membrane made by nitrocellulose or PVDF. The way to use IMMUNO SHOT -Platinum- in WB is described below.

- 1) SDS-PAGE and transfer of protein to PVDF membrane should be done by usual method.
- 2) Blocking and the washing should be done by usual method.
- 3) Dilute the 1st antibody with IMMUNO SHOT -Platinum-. The dilution factor is influenced by many factors, such as antibody species, amount of antigen, etc. Though you can reduce antibody concentration by using IMMUNO SHOT -Platinum-, we recommend performing a pre-test to determine the best antibody concentration.
- 4) Dilute the 2nd antibody with IMMUNO SHOT -Platinum-. The best dilution factor is influenced by many factors, such as antibody species, amount of antigen, etc. Refer to the supplier's instruction to determine the best antibody concentration.
- 5) For visualization, many users use HRP-, or AP-labeled antibody. In both cases, please watch the strength of staining or luminescence and stop the reaction. Longer reaction gives you high background or appearance of extra band.

(6) ELISA

ELISA is a method to determine the amount of antigen or antibody in samples by using labeled antigen or antibody. The sandwich ELISA is most widely used, where antigen sample is applied on solid phase antibody and bound antigen is reacted by HRP-labeled detection antibody. The way to use IMMUNO SHOT -Platinum- in the sandwich ELISAs are described below.

- 1) Antibody attachment (solid phase), blocking and the washing procedure should be done by usual method.
  - 2) Dilute the antigen with PBS-T. Pour appropriate amount of these diluents into each well, mix and incubate for appropriate time. IMMUNO SHOT -Platinum- may be used for antigen dilution in order to facilitate reaction of antigen and solid phase antibody.
  - 3) Dilute the HRP-labeled detection antibody with IMMUNO SHOT -Platinum-. The best dilution factor is influenced by many factors, such as antibody species, amount of antigen, etc. Refer to the supplier's instruction to determine the best antibody concentration.
  - 4) For substrate reaction, follow the ordinary method.
  - 5) In some case 1st antibody and HRP-labeled 2nd antibody system is used. In such a case use IMMUNO SHOT -Platinum- for dilution of both antibodies.
- In other ELISA system, use IMMUNO SHOT -Platinum- for reaction buffer for antigen-antibody.

(7) Trouble shooting

Trouble	Cause and resolutions
Western blotting	
Weak signal	1. Low antigen conc.: Use higher antigen conc.
	2. Low antibody conc.: Survey best antibody conc.
	3. Not enough transfer: Use higher current or longer transfer time.
	4. Blocking too strong: Do not use long time blocking.
	5. Too much transfer: When using nitrocellulose, proteins pass though the membrane by strong transfer manipulation. Check the procedure or exchange membrane to PVDF.
Partial whitening (lumines.)	6. Too much antigen or antibody: Over signaling often suppress the luminescence and cause partial whitening of a band. Control the amount of antigen or antibody concentration.
Too many extra-bands	7. Too much higher antibody conc.: Higher antibody conc. often causes non-specific signaling. Control the antibody conc.
	8. Too much antigen: Higher antigen often causes non-specific signaling. Control the amount of antigen.
	9. Not enough blocking: Some antigen and antibody have preference of blocking agents, change the blocking agents or check the blocking conditions.
	10. Not enough washing: Increase the number and time of washing.
High background	11. High antibody conc. or too long incubation: Reduce the antibody conc. or shorten the incubation time.
ELISA	
Weak signal	1. Too low antigen or antibody conc.: Increase the concentrations.
Too strong signal	2. Too high antigen or antibody conc.: Check the antigen and antibody concentration by performing the titration.
	3. Too long incubation: shorten the incubation time.
High background	4. Too high antigen or antibody conc.: check antigen and antibody concn.
	5. Not enough blocking: Some antigen and antibody preference of blocking agents, change the blocking agents or check the blocking conditions.
	6. Not enough washing or too much washing: check the number and time of washing.

COSMO BIO Co., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU, TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)

— 3 —

— 4 —

# IMMUNO SHOT -Platinum-

Protein-free immuno-reaction enhancing solution

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