



*ELISA kit for measuring  
UV-induced DNA damage*

**High Sensitivity 6-4PP/  
(6-4)Photoproducts  
ELISA kit**  
**(with mAb clone 64M-2)**

Catalog Number : NM-MA-K002(96 tests)


*For research use only, Not for diagnostic use.*

- Please read this manual thoroughly before use -

## INTRODUCTION

Prolonged exposure to solar UV radiation may result in harmful acute and chronic effects to the skin (including skin cancers), eye, and immune system. These harmful effects appear to be closely related to UV-induced DNA damage. The major types of DNA damage induced by solar UV radiation are cyclobutane pyrimidine dimers (CPDs), (6-4) photoproducts (6-4PPs), and Dewar photoproducts (DewarPPs), which are formed between adjacent pyrimidine nucleotides on the same DNA strand. These helix-distorting DNA lesions are repaired exclusively by a nucleotide excision repair system in humans. To better study molecular events surrounding UV-induced DNA damage and repair, Mori *et al.* previously developed and characterized monoclonal antibody (mAb) specific for CPDs and mAb specific for 6-4PPs (1) while Matsunaga *et al.* developed and characterized mAb specific for DewarPPs (2). Three of these antibodies (CPDs: clone TDM-2; 6-4PPs: clone 64M-2; DewarPPs: clone DEM-1) continue to be cited frequently in the literature, often for use in ELISA by a recommended procedure.

This High Sensitivity (6-4)photoproducts (6-4PPs) ELISA Kit is the only commercially available ELISA utilizing anti-6-4PPs clone 64M-2 and has been optimized for high sensitivity detection of 6-4PPs in DNA purified from cultured cells or from skin epidermis. This ELISA detects 6-4PPs from dipyrimidines in all DNA sequence contexts (i.e., TT, TC, CT and CC). Thus, the availability and convenience of this ELISA Kit will contribute to further understanding molecular mechanisms involved in cellular responses to UV radiation and DNA damage with applications across many research fields including cancer research, photobiology, dermatology, ophthalmology, immunology, and cosmetics science.



**Figure 1: Structures of UV-induced DNA damage in thymine-thymine sequence**

## ASSAY PRINCIPLE

The format of this assay is ELISA with colorimetric detection. In brief, genomic DNA purified from UV-damaged cells is heat denatured and applied to microtiter wells pre-coated with protamine sulfate. 6-4PP specific monoclonal antibody clone 64M-2 (Cosmo Bio Cat. No. CAC-NM-DND-002) is then added to each well for thirty minutes and unbound antibody is removed by rinsing. The amount of 64M-2 antibody remaining in each well is then measured by sequential treatment of wells with biotinylated 2nd antibody, streptavidin-peroxidase, and o-phenylenediamine (OPD). The reaction between peroxidase, H<sub>2</sub>O<sub>2</sub> and OPD produces a yellow orange color, the strength of which is generally proportional to the amount of 64M-2 antibody remaining bound to the plate. The color development reaction is stopped and the absorbance of each well at 492 nm is measured with a spectrophotometer.




Figure 2: An ELISA for 6-4PPs

## REACTIVITY

- 1) Anti-6-4PPs monoclonal antibody clone 64M-2 recognizes 6-4PPs on single-stranded DNA.
- 2) 64M-2 binds to 6-4PPs formed each dipyrimidine sequence context (TT, TC, CT and CC).
- 3) 64M-2 stably binds to 6-4PPs in DNA longer than eight bases.
- 4) 64M-2 binds to 6-4PPs in UV-irradiated DNA purified from a wide range of sources the prokaryote and eukaryote irradiated with UV.

## KIT COMPONENTS

Item	Amount
ELISA plate precoated with protamine sulfate (12 x 8 well strips)	1 plate
Positive control Calf thymus DNA, UVC irradiated (10 J/m <sup>2</sup> )	1 vial (20 µg/mL, 500 µL)
Negative control Calf thymus DNA, not irradiated	1 vial (20 µg/mL, 500 µL)
Assay Diluent Concentrate (10X)	1 vial (10 mL)
Wash Buffer Concentrate (20X)	2 x 15 mL vials
Blocking Reagent Concentrate (50X)	1 vial, lyophilized. Reconstitute with 400 µL purified water before use.
Anti-6-4PPs Monoclonal Antibody (clone 64M-2) (100X)	1 vial, lyophilized. Reconstitute with 150 µL purified water for a 100X working solution.
Biotinylated 2nd antibody (100X)	1 vial (150 µL)
Streptavidin-peroxidase (100X)	1 vial (150 µL)
OPD Tablet (5mg)	3 tablets
OPD Diluent Concentrate (10X)	1 vial (1.7 mL)
Stop Solution	1 vial (12 mL)
Plate Cover Film	3 covers Two films are spares; it is not necessary to use all three films.
Instruction Manual	1 manual

## MATERIALS TO BE SUPPLIED BY THE USERS

- DNA samples
- DNA Purification Kit (for sample preparation)  
Recommended: QIAamp Blood Kit (QIAGEN, Cat. No. 51104 or 51106)
- 100 °C Heating Block
- Ice bath (Crush ice)
- Purified water
- 10 µL - 1000 µL adjustable single channel micropipettors and disposable tips
- 50 µL - 150 µL adjustable multichannel micropipettors and disposable tips
- Reservoir for Wash Solution
- 1.5 mL tubes (for diluting samples)
- 15 mL or 50 mL tubes (for dilutions)
- 37°C Incubator (non-humidified)
- Absorbance microplate reader capable of reading 492 nm
- Vortex mixer
- Desktop centrifuge

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

## STORAGE AND EXPIRATION

<b>Unopened kit</b>	: 4 °C
<b>Opened kit</b>	
Reconstituted solutions	: -20 °C
Reconstituted Antibody solution	: avoid repeated thaw and freeze cycles
Positive and Negative Controls	: -20 °C
Other reagents	: 4 °C
ELISA plate precoated with protamine sulfate	: room temperature, protect from light
<b>Expiration date</b>	
6 months from the shipping date.	

## PREPARATION OF REAGENTS

➤ **Bring all reagents to room temperature (18-25 °C) before use.**

### **1. ELISA plate precoated with protamine sulfate**

Bring to room temperature (18-25 °C) before use. Return unused wells to foil pouch.

### **2. Positive and Negative 6-4PP Controls**

The concentration of the control solutions is 20 µg/mL. Prepare 6-4PPs DNA Control solutions to a concentration of 4 µg/mL with 1X Assay Diluent.

### **3. Assay Diluent**

Dilute 10 mL of Assay Diluent concentrate (10X) with 90 mL purified water to make 100 mL of Assay Diluent (1X).

### **4. Wash Buffer**

Dilute 15 mL of Wash Buffer concentrate (20X) with 285 mL purified water to make 300 mL of Wash Buffer (1X).

### **5. Blocking Reagent**

The *Blocking Reagent* is lyophilized. Reconstitute with 400 µL of purified water. Upon reconstitution, the solution is a 50X concentrate. Dilute 1:50 with Assay Diluent to prepare *Blocking Reagent Working Solution*.

### **6. Anti-6-4PPs**

The *Anti-6-4PPs* antibody is lyophilized. Reconstitute with 150 µL of purified water. Upon reconstitution, the *Anti-6-4PPs* solution is a 100X concentrate. Dilute 1:100 with Assay Diluent to prepare *Anti-6-4PPs Working Solution*.

### **7. Biotinylated 2<sup>nd</sup> Antibody**

The *Biotinylated 2<sup>nd</sup> Antibody* solution is a 100X concentrate. Dilute 1:100 with Assay Diluent to prepare *Biotinylated 2<sup>nd</sup> Antibody Working Solution*.

### **8. Streptavidin-peroxidase**

The *Streptavidin-peroxidase* solution is a 100X concentrate. Dilute 1:100 with Assay Diluent to prepare *Streptavidin-Peroxidase Working Solution*.

### **9. OPD Diluent**

*OPD Diluent* is provided as a 10X concentrate. Dilute 500 µL of *OPD Diluent* concentrate (10X) with 4.5 mL of purified water to prepare 5 mL of 1X *OPD Diluent*. Prepare immediately before use.

### **10. OPD Substrate Solution**

Dissolve one OPD tablet in 5 mL 1X OPD Diluent to make Working *OPD Substrate Solution*. Prepare immediately before use.

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

# ASSAY PROTOCOLS

## **A. Cell culture and UV irradiation**

1. Plate cells in 10 cm dishes and culture for one or two days.
2. Wash cells once with Dulbecco's PBS (DPBS) and irradiate with UV (e.g., 0, 2.5, 5, 7.5, 10 J/m<sup>2</sup> at 254 nm). To study DNA repair, irradiate cells with 10 J/m<sup>2</sup> and incubate for various amounts of time before harvesting (e.g., 3, 6, 12 hours) to allow repair.
3. Wash cells with 10 mL DPBS. Harvest by scraping cells from dish. Centrifuge at 10,000xg for 15 seconds at 4 °C.
4. Store cell pellets at –80 °C until ready for DNA isolation.

## **B. DNA isolation**

5. Purify genomic DNA using a QIAamp Blood Kit (QIAGEN, Cat. No. 51104 or 51106) or similar. DNA concentrations are calculated by absorbance at 260 nm.

## **C. DNA sample coating to the ELISA plate precoated with protamine sulfate**

6. Prepare sample DNA or 6-4PPs DNA Control solutions to a concentration of 4 µg/mL with 1X *Assay Diluent*. Denature DNA solutions by heating to 100°C for 10 minutes, then chill rapidly in an ice bath for 15 minutes.
7. Apply 50 µL/well of denatured DNA solution to the ELISA plate wells precoated with protamine sulfate (duplicates recommended) and dry completely overnight by incubation at 37 °C.

## **D. DNA damage detection**


8. Wash the DNA-coated plates 5 times with 150 µL/ well of 1X *Wash Buffer*. Discard *Wash Buffer*, and tap the plate on the paper towel to remove the solution completely.
9. Add 150 µL/well *Blocking Reagent Working Solution* to each well to prevent non-specific binding of antibody.
10. Seal the wells with a Plate Cover Film and incubate 30 minutes at 37 °C.
11. Wash the plates 5 times with 150 µL/well of 1X *Wash Buffer*. Discard *Wash Buffer*, and tap the plate on the paper towel to remove the solution completely.
12. Add 100 µL/well of *Anti-6-4PP Working Solution*, seal the wells with a Plate Cover Film and incubate 30 minutes at 37 °C.
13. Wash the plates 5 times with 150 µL/well of 1X *Wash Buffer*. Discard *Wash Buffer*, and tap the plate on the paper towel to remove the solution completely.
14. Add 100 µL/well *Biotinylated 2<sup>nd</sup> Antibody Working Solution*, seal the wells with a Plate Cover Film and incubate 30 minutes at 37 °C.
15. Wash the plates 5 times with 150 µL/well of 1X *Wash Buffer*. Discard *Wash Buffer*, and tap the plate on the paper towel to remove the solution completely.
16. Add 100 µL/well of *Streptavidin-Peroxidase Working Solution*, seal the wells with a Plate Cover Film and incubate 30 minutes at 37 °C.
17. Wash the plates 5 times with 150 µL/well of 1X *Wash Buffer*. Discard *Wash Buffer*, and tap the plate on the paper towel to remove the solution completely.
18. Add 100 µL/well Working *OPD Substrate Solution* to each well, seal the wells with a Plate Cover Film and incubate 30 minutes at 37 °C.
19. Add 100 µL/well *Stop Solution* to each well to stop enzyme reaction.
20. Mix gently and immediately determine the absorbance at 492 nm of each well using a spectrophotometer.

## NOTES


- Do not mix or substitute reagents with those from other lots or sources.
- If a precipitate appear in *Assay Diluent* concentrate or *Wash Buffer* concentrate, warm it gently to dissolve the precipitate before use.

[www.cosmobio.co.jp](http://www.cosmobio.co.jp)

## EXAMPLE OF RESULTS



**Figure 3: UV-induced 6-4PPs in DNA measured by ELISA**



**Figure 4: Formation and repair of UV-induced 6-4PPs in human cells measured by ELISA**

UVC radiation induces 6-4PPs in DNA of HeLa cells in dose-dependent manner. The initial level of 6-4PPs induced by 10 J/m<sup>2</sup> of UVC decreases over time as 6-4PPs are repaired, indicating the capacity of nucleotide excision repair in HeLa cells.

## SELECTED REFERENCES

- 1) Mori, T., *et al.*, *Photochem. Photobiol.* **54**, 225-232 (1991).
- 2) Matsunaga, T., *et al.*, *Photochem. Photobiol.* **54**, 403-410 (1991).
- 3) Potten, C.S., *et al.*, *Int. J. Radiat. Biol.* **63**, 313-324 (1993)
- 4) Kobayashi, N., *et al.*, *J. Invest. Dermatol.* **101**, 685-689 (1993)
- 5) Todo, T., *et al.*, *Nature* **361**, 371-374 (1993)
- 6) Nakane, H., *et al.*, *Nature* **377**, 165-168 (1995)
- 7) Komatsu, Y., *et al.*, *Nucleic Acids Res.* **25**, 3889-3894 (1997)
- 8) Kobayashi, N., *et al.*, *J. Invest. Dermatol.* **110**, 806-810 (1998)
- 9) Nakagawa, A., *et al.*, *J. Invest. Dermatol.* **110**, 143-148 (1998)
- 10) Otsoshi, E., *et al.*, *Cancer Res.* **60**, 1729-1735 (2000)
- 11) Katsumi, S., *et al.*, *J. Invest. Dermatol.* **117**, 1156-1161 (2001)
- 12) Kobayashi, N., *et al.*, *Pigment Cell Res.* **14**, 94-102 (2001)
- 13) Wakasugi, M., *et al.*, *J. Biol. Chem.*, **277**, 1637-1640 (2002).
- 14) Imoto, K., *et al.*, *J. Invest. Dermatol.* **119**, 1177-7782 (2002)
- 15) Nishiwaki, Y., *et al.*, *J. Invest. Dermatol.* **122**, 526-532 (2004)
- 16) Sugasawa, K., *et al.*, *Cell* **121**, 387-400 (2005).
- 17) Yasuda, G., *et al.*, *Mol. Cell. Biol.*, **27**, 6606-6614 (2007).
- 18) Matsumoto, M., *et al.*, *J. Cell Sci.*, **120**, 1104-1112 (2007).
- 19) Yamamoto, A., *et al.*, *DNA Repair*, **6**, 649-657 (2007).

More than 200 papers using 64M-2 antibodies have been published so far.

## RELATED PRODUCTS

Product Name	Maker	Cat#
Anti cyclobutane pyrimidine dimers (CPDs) Monoclonal Antibody (Clone: TDM-2)	CAC	NM-DND-001
Anti (6-4) photoproducts (6-4PPs) Monoclonal Antibody (Clone: 64M-2)	CAC	NM-DND-002
Anti Dewar photoproducts (DewarPPs) Monoclonal Antibody (Clone: DEM-1)	CAC	NM-DND-003
Anti Acetylaminofluorene(AAF)-DNA adducts Monoclonal Antibody (Clone: AAF-1)	CAC	NM-MA-001
High Sensitivity CPD (Cyclobutane Pyrimidine Dimer) ELISA kit Ver.2	CSR	NM-MA-K003
High Sensitivity 6-4PP ((6-4)Photoproducts) ELISA kit (TMB)	CSR	NM-MA-K004
UVC irradiated DNA samples (0, 2.5, 5, 7.5, 10 J/m <sup>2</sup> )	CSR	NM-MA-R010
PROTAMINE SULFATE COATED ELISA PLATE 96	CSR	NM-MA-P001
PROTAMINE SULFATE COATED ELISA PLATE 96 x 5	CSR	NM-MA-P002
PROTAMINE SULFATE COATED ELISA PLATE 96 x 10	CSR	NM-MA-P003
Anti XPA Monoclonal Antibody (Clone: A-2)	CAC	KUP-TM-M01
Anti XPA Monoclonal Antibody (Clone: 5F12)	BAM	70-032
Anti XPF Monoclonal Antibody (Clone: 19-16)	CAC	KUP-TM-M02
Anti XPG Monoclonal Antibody (Clone: G-26)	CAC	KUP-TM-M03
Anti ERCC1 Monoclonal Antibody (Clone: E1-44)	CAC	KUP-TM-M04
Anti DDB1 Monoclonal Antibody (Clone: 43233-3-1)	CAC	KUP-TM-M05

*For research use only, Not for diagnostic use.*



**COSMO BIO Co., LTD.**

**[JAPAN]**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>



**COSMO BIO USA**

**[Outside Japan]**

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiouusa.com](mailto:info@cosmobiouusa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiouusa.com](http://www.cosmobiouusa.com)