




Composition for catenating single-stranded DNA

Product Description	ssDNA catenating solution (SCS)
Specific Activity	<p>Specific Activity: Circular single-stranded DNA (20 μM) is incubated with >2 μl of SCS in 10 μl of a reaction buffer containing 1~2 mM MgCl₂ at 37 °C for 1 hour.</p> <div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 10px auto;"><p>Agarose gel electrophoresis of catemers (SYBR® Gold (invitrogen) staining) Reaction condition : 20 mM HEPES (pH 7.5), 20 mM NaCl, 1.5 mM DTT, 4% glycerol, 1 mM MgCl₂</p></div> 
Initial Volume	50 μ l
Formulation	20 mM HEPES-NaOH (pH7.5), 0.1 M NaCl, 2 mM DTT, 30% glycerol
Storage & Stability	Store below -70 °C. Stable for 1 year at -70 °C from date of shipment. Aliquot to avoid cycles of freeze/thaw. Recommended to divide into small amounts before first time using.
Reference	Motoki Takaku et. Al. Single-stranded DNA catenation mediated by human EVL and a type I topoisomerase. Nucleic Acids Research 2010 (in press)

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619