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Mount Quick

Mount-Quick "Aqueous"

|                       |                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                    |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Purpose of use</b> | Sample mounting medium / mounter for permanent preservation                                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                    |
| <b>Features</b>       | <b>Mount Quick</b> <ul style="list-style-type: none"> <li>• Rapid drying, highly transparent, and does not form bubbles. Highly soluble in xylene, does not degenerate or stain the specimen. Can withstand long-term storage after mounting.</li> <li>• Tube design provides easy control of amounts used and limits air exposure, preventing variability in concentration.</li> </ul> | <b>Specifications:</b><br>Specific gravity: 0.89sp;<br>Viscosity: 33cp;<br>Refractive index: 1.41                                                                  |
|                       | <b>Mount-Quick "Aqueous"</b> <ul style="list-style-type: none"> <li>• Maintains superior transparency when fully hardened.</li> <li>• Low fluorescence.</li> <li>• Mounting is simple with an easy-to-use dropper bottle.</li> </ul>                                                                                                                                                    | <b>Specifications:</b><br>Mobility: pale yellow, transparent, viscous aqueous solution;<br>Specific gravity: 0.88sp;<br>Viscosity: 33cp;<br>Refractive index: 1.41 |
| <b>How to use</b>     | <b>Mount Quick</b><br>Squeeze out an amount approximately the same size as the sample section at one location on the cover glass or sample.                                                                                                                                                                                                                                             |                                                                                                                                                                    |
|                       | <b>Mount-Quick "Aqueous"</b><br>Drop onto the cover glass or section from above.<br>(Prevent bubble formation by inverting the container and touching the tip on the bottom of the drop.)<br>* Before first use, cut open the container tip at a 45° angle with a box cutter or scissors.                                                                                               |                                                                                                                                                                    |
| <b>Precautions</b>    | After use, always re-cap and store in a cool, dark place.                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                    |

## Application Data

Use Mount Quick Aqueous when samples are not dehydrated and cleared (photograph as the color fades during long-term storage) and use Mount Quick Tube when samples are dehydrated and cleared (can be permanently preserved).

## Staining data chart for Mount Quick Tube and Mount Quick Aqueous

| DM01 Mount Quick Tube<br>(non-water soluble mounting medium for mounting after dehydration and clearing)    |   |
|-------------------------------------------------------------------------------------------------------------|---|
| Hematoxylin-eosin staining                                                                                  | ○ |
| Azan stain                                                                                                  | ○ |
| Grimelius stain                                                                                             | ○ |
| Silver stain                                                                                                | ○ |
| Van Gieson's stain                                                                                          | ○ |
| Periodic acid-Schiff (PAS) stain                                                                            | ○ |
| Giemsa stain                                                                                                | ○ |
| Immunostaining (non-aqueous DAB)                                                                            | ○ |
| DM02 Mount Quick Aqueous<br>(water soluble mounting medium for mounting after staining without dehydration) |   |
| Lipid stains                                                                                                | ○ |
| Immunohistochemical staining (AEC)                                                                          | ○ |
| Blood smear staining                                                                                        | ○ |
| PAS, FITC, DAB staining                                                                                     | ○ |

Both Mount Quick Tube and Mount Quick Aqueous can be used with methyl green and BCIP/NBT staining, but are not suitable for long-term storage. After staining, we recommend photographing as soon as possible, as the color can fade in 24 hours after drying, dehydrating, clearing, and mounting. BCIP/NBT can stain red and black with immunostaining.

**COSMO BIO Co., LTD.**  
Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN  
[International/Outside Japan] [国内連絡先]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com

e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp

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DAIDO SANGYO CO., LTD.