2D-Silver Stain Reagent II

-For Electrophoresis-

Cat. No. DCB-423413

Last Updated: 2016/1/7

www.cosmobio.com

【 I 】 Background

Detection of proteins and nucleic acids migration by polyacrylamide gel electrophoresis and silver staining are now being highlighted as one of the most sensitive staining methods. However the original silver staining method is not necessarily easy to work and takes a long-time to complete. 2D-SILVER STAIN-II is an improved reagent kit developed to provide simple and fast staining.

【 II 】 Features

- 1) Rapid staining results in a short time after electrophoresis
- 2) Staining with higher sensitivity
For Proteins: 50 to 100 times more sensitivity than CBB staining
For Nucleic Acids: 50 to 100 times more sensitivity than EB staining
- 3) Simple preparation and simple operation
- 4) No use of materials such as heavy metals which need to be controlled by regulations.
- 5) Staining process can be stopped any time during development for desired chromatic figures.

【 III 】 Applications

2D-SILVER STAIN II is applicable to the detection of proteins and nucleic acids subjected to polyacrylamide and SDS-polyacrylamide slab gel electrophoresis.

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Reagents	Main Ingredients	Volume
Fixing Reagent	Thiourea	100 ml x 1
Pretreatment reagent	Dithiothreitol, Glutar-aldehyde, Thiourea	100 ml x 1
Staining Solution A	Silver nitrate	100 ml x 1
Staining Solution B	Ammonium hydroxide, Sodium hydroxide	100 ml x 1
Concentrated developer	Citric acid, Formaldehyde, Sodium thiosulfate	100 ml x 1
Stopper	Citric Acid	100 ml x 1

【V】Procedures

For Protein Staining, in addition to the reagents supplied with the kit, the following materials are required.

- 1) Methanol
- 2) Acetic Acid
- 3) Deionized water (less than 10^{-6} mho conductivity)
- 4) Measured cylinder
- 5) Flat-bottomed tray for developing gels
- 6) Measuring pipette

【V-1】Preparation of Reagents

Prepare the following solutions just before staining the gel. The volume of reagents use for preparation shown below is based on the size of 140 mm x 140 mm x 1.0 mm slab gel. For any different sizes of slab gels, the volume of reagents should be prepared in proportion to the size of the gel to be used.

1 Fixing solution I

Mix 100 ml of methanol and 20 ml of acetic acid with 80 ml of deionized water. Stir and use the mixture as fixing solution I.

2 Fixing solution II

Mix 60 ml of methanol, 20 ml of acetic acid and 10 ml of 1 Fixing Reagent with 110 ml of deionized water. Stir and use the mixture as fixing solution II.

3 Pretreatment solution

Mix 100 ml of methanol, 10 ml of 2 Pretreatment Reagent and 90 ml of deionized water. Stir and use the mixture as pretreatment solution.

4 Silver staining solution

Mix 10 ml of 3.) Staining Solution A and 10 ml of 4.) Staining Solution B and stir. Then add 180 ml of deionized water and stir. Use the mixture as silver staining solution.

5 Developer

Add 10 ml of 5.) Concentrated Developer to 190 ml of deionized water and stir. Use the mixture as the developer.

【V-2】Staining Procedure



- 1 Soak the gel in fixing solution I in a clean flat-bottomed tray and shake it for 10 minutes.

Sample: protein gel < 1 mm


Fixing solution I: 200 ml

10 min

Fixing solution II: 200 ml

Aspirator

15 min

- 2 Remove the solution, and then pour 200 ml of fixing solution II to the tray and shake it for 15 minutes.

- 3 Remove the solution, and then pour 200 ml of pretreatment solution to the tray and shake it for 10 minutes.

- 4 Remove the solution, and then pour 200 ml of the deionized water to the tray and shake it for 5 minutes.


- 5 Remove the solution, and then pour 200 ml of silver staining solution to the tray and shake it for 15 minutes

- 6 Remove the silver staining solution into a container. Immediately add 2-3 ml of concentrated hydrochloric acid to the waste solution to convert the silver into silver chloride (This procedure must be followed accurately). Pour 200 ml of deionized water to the tray and shake it for 2 minutes and then remove the water. Repeat this step 3 times.

- 7 Pour 200 ml of developer to the tray and shake it well

- 8 When appropriate chromatic figures are obtained, pour 10 ml of 6.) Stopper solution to the same tray and shake it for a moment

- 9 Leave it for about 10 minutes and then wash the gel with deionized water two or three time and keep it.

【V-3】Notes

- ◆ As the silver staining solution produces explosive silver amide when left to stand, the solution must be treated with hydrochloric acid or sodium chloride to convert the silver into silver chloride (See Procedure 6.)
- ◆ Since Staining Solutions B among the components of 2D-Silver Stain-II, contains ammonium hydroxide and sodium hydroxide as an alkali substance respectively, take care to prevent them from contacting the eye and the skin. In case of eye contact rapidly wash the eye with flowing water, and consult a doctor. In case of contact with the skin and clothes, rapidly wash away with water, and consult a doctor if inflammation occurs.
- ◆ Staining Solution A and other components may also cause a stimulus to the mucous membrane and the skin. In case of contact with the eye, skin, and clothes, rapidly wash away with water, and consult a doctor if inflammation should be caused.
- ◆ Addition of the silver staining solution to the sample gel may color or make turbid the solution. Such coloration and turbidity do not affect the stain. (See Procedure 5.)
- ◆ Deionized water or distilled water of less than 10⁻⁸ mho conductivity is recommended.
- ◆ The tray to be used for the staining should be clean and have a smooth surface.
- ◆ When a sample gel stained by the CBB staining is to be subjected to the silver staining, the sample gel should be washed with deionized water thoroughly. (3 hours or longer while changing deionized water several times prior to the silver staining.)
- ◆ Care should be taken not to damage the gel during operation.
- ◆ Care should be taken to see that the gel may not come out of the surface of solution.
- ◆ Use a pair of disposable gloves in staining operation not to touch the gel.
- ◆ If a proper shaker is not available, shake the tray by hand.
- ◆ The sensitivity of a gel thicker than 1mm declines somewhat as compared with a thinner gel.
- ◆ Keep the reagents in a tightly closed container and refrigerated, and use within 6 months after unsealing.
- ◆ Use the reagents in a tightly closed container and refrigerated, and use within 6 months after unsealing.
- ◆ Use the reagents (methanol, acetic acid) of superior quality

【V-4】Package Size and Storage

Enough to stain 10 slabs of gels (140 mm x 140 mm x 1.0 mm)

The kit should be stored in a dark place at 2-10° C

The term of validity: Printed on the package.

【V-5】References

1. K. Ohsawa et. Al.: Silver Stain for Detecting 10-Femogram Quantities of protein after Polyacrylamide Gel Electrophoresis, Anal. Biochem., 135, 409(1983)
2. S. Irie: A highly sensitive silver Staining for detection of proteins in polyacrylamide gels, biochemistry(Japan), 52, 411(1980)
3. B.R. Oakley et. al.: Visualization of proteins with silver "stain" for detecting proteins in polyacrylamide gels, Anal. Biochem., 105, 361(1980)
4. H.M. Poehling et. al.: Visualization of proteins with silver "stain", a critical analysis, electrophoresis, 2, 141 (1981)
5. U.K. Laemmli: Cleavage of structural proteins during the assembly of the head of bacteriophage T4., Nature, 227, 680(1970)

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