



Viale Piceno 137/f
61032 Fano PU (IT)
Telephone + 39 (0) 721830605
FAX +39 (0) 721837154
e-mail: info@diatheva.com
www.diatheva.com

DPYD IVS14 + 1G>A HRM kit

MBK0048

50 Reactions (42 tests)

Store at -20°C

Intended use

For identification of DPYD IVS14 + 1G>A polymorphism by High Resolution Melt (HRM).

Introduction

DPYD gene (located in 1p22) encodes dihydropyrimidine dehydrogenase protein, a pyrimidine catabolic enzyme representing the initial and rate-limiting factor in the pathway of uracil and thymidine catabolism. It is also the most important enzyme for degradation of 5-FU. Genetic variants causing DPYD deficiency leads to decreased 5-FU catabolism resulting in a significant increase in the effective 5-FU dose in the body. Several variants have been found in DPYD gene and among them IVS14 + 1G>A polymorphism represents one of the most predictive factor for 5-fluorouracil toxicity in cancer treated patients.

Product description

The DPYD IVS14 + 1G>A HRM kit is an easy-to-use master mix dedicated for the use with the Rotor-gene 6000 instrument (Corbett Research) and Rotor-Gene Q (Qiagen). The kit contains reagent, enzyme and genotype controls for the detection of DPYD IVS14 + 1G>A SNP located in the DPYD gene by HRM.

Kit contents

1 x DPYD IVS14 + 1G>A 2X master mix (650 µl)
1 x DNA polymerase (10 µl)
2 X DPYD IVS14 + 1G>A G/G CONTROL (10 µl)
2 X DPYD IVS14 + 1G>A G/A CONTROL (10 µl)
2 X DPYD IVS14 + 1G>A A/A CONTROL (10 µl)

Storage

The product should be stored at -20°C immediately upon arrival. When stored under the recommended conditions and handled correctly, the kit is stable up to the expiry date indicated in the attached label.

Note: The product is guaranteed for two thawing steps. For an intermittently use we suggest to subdivide DPYD IVS14 + 1G>A 2X master mix in small aliquots.

Precautions

The user should always pay attention to:

- use pipette tips with aerosol-preventive filters, deionized DNA-free water and gloves;
- store positive material (plasmidic DNA controls) separately from all other reagents and, if possible, add it to the reaction mix in a separated space;
- do not use the same precision pipettes for reaction mix and DNA;
- thaw all components samples at room temperature before starting an assay;
- when thawed, mix the components and centrifuge briefly.
- Protect DPYD IVS14 + 1G>A 2X master mix from light.

Procedure

PREPARING HRM GENOTYPIC CONTROLS

Genotype controls must be prepared by 1:10 dilution of stock solution in PCR -grade water (not provided).

Genotype controls are provided to characterize the HRM curves of unknown amplicons.

PCR SETUP

Total volume per reaction is 25 µl.

Before each use, thaw all reagents completely, mix and centrifuge.

Pipet mastermix and the Taq polymerase in the quantity needed for the planned reactions into a 1.5 ml reaction tube and mix as indicated in table 1.

Moreover include provided positive controls for each assay as references for genotype assignment.

Table 1

	1 reaction	50 reaction
PCR Mastermix 2X	12.5 µl	625 µl
Taq Diatheva (5U/µl)	0.125 µl	6,25 µl
H ₂ O	10,375	518,75
Total volume	23 µl	1150 µl

Aliquot 23 µl of master mix into each PCR reaction tube before adding 2 µl sample DNA*, or 2 µl positive controls (DPYD IVS14 + 1G>A G/G -A/G-A/A) for each SNP genotype of interest.

*Employ 10 ng DNA/reaction (2 µl of a 5ng/µl concentrated DNA).

DNA extracted employing **QIAamp DNA Blood Mini Kit** and **DNeasy Blood & Tissue Kit** –QIAGEN can be used without quantitation step, simply diluting the sample 1:2 in water.

(DNA quality is a fundamental factor for HRM genotyping: kit has been optimized with DNA extracted by **DNeasy Blood & Tissue Kit** - **QIAamp DNA Blood Mini Kit** -QIAGEN and **Blood genomic DNA isolation mini Kit** – NORGEN).

!! Attention !! Final elution step must be carried out using distilled water instead of the elution buffer supplied with the respective kits.

Alternative isolation methods have to be validated from the user. In any case we suggest to elute purified DNA in water.

THERMAL PROFILE

HRM analysis setting is a prerequisite for accurate results. For details, please refer to the manual provided with your HRM real-time PCR instrument.

Program the HRM instrument according to the operator's manual and table 2

Table 2: optimized cycling protocol for rotor-gene 6000

STEP	TEMPERATURE	TIME	CICLES
Initial denaturation	95°C	10 min	1X
Denaturation	95°C	10 sec	40 X
Annealing / Extension	59°C	40 sec	
Acquire on the GREEN Channel in the Annealing / Extension step			
High Resolution melting	Ramp 82°-93°C	Rise temp 0.2°C/step; wait for 90 sec pre-melt conditions; wait for 2 sec for each step afterwards; For optimal acquisition of fluorescence data, set the gain to 70% of saturation in the highest fluorescent signal.	
Acquire on the HRM channel			

ANALYSIS SETUP

For HRM analysis set confidence level at 90%.

Due to the comparison between controls and unknown samples HRM profiles, the software instrument performs an automatic genotyping showing a confidence level for each sample. Results with a confidence level lower than 90% are defined "variation" and have to be reject.

Troubleshooting

No signal, poor Rn value (PCR or HRM) or signal detected late in PCR	Pipetting error or missing reagent	Check the storage conditions of the reagents, Repeat the assay.
	Problems with starting template DNA	Check the concentration, storage conditions, and quality of the template and control DNA ($\lambda_{260}/\lambda_{280}$ ratio of DNA samples must be over 1,7.) Efficient removal of PCR inhibitors is essential for optimal results; purify nucleic acids from your sample using an appropriate purification method. Insufficient or degraded template DNA, increase the amount of template DNA if possible.

References

Wei X, McLoed HL, McMurrough J, et al. **Molecular basis of the human dihydropyrimidine dehydrogenase deficiency and 5-fluorouracil toxicity.** *J Clin Invest* 98: 610-615 (1996)

Van Kuilenburg ABP, Vreken P, Beex LVAM, et al. **Heterozygosity for a point mutation in an invariant splice site donor of dihydropyrimidine dehydrogenase and severe 5-fluorouracil related toxicity.** *Eur J Cancer* 33: 2258-2264 (1997)

Etienne M.C, Lagrange JL, Dassonville O, et al. **Population study of dihydropyrimidine dehydrogenase in cancer patients.** *J Clin Oncol* 12: 2248- 2253(1994)

Example images for HRM analysis Images





Fig. 1: Graphical visualizations for the identification of the DPYD IVS14+1 G>A polymorphisms.
A) Melting curves show heteroduplexes formation and temperature shift between homozygotes different genotypes.
B) Normalized melting curves of the three genotypes.

Example of experiment set-up for 6 unknown samples:


- 1) Prepare each positive controls and samples as described in procedure paragraph and keep them on ice.
- 2) prepare Master Mix 1X for 10* reaction (the 3 positive controls + 6 unknown samples + 1):

	1 reaction	10 reaction
DPYD IVS14 + 1G>A 2X master mix	12.5 µl	1255 µl
DNA Polymerase (5U/µl)	0.125 µl	1,25 µl
H ₂ O	10.375	103,75 µl
Total volume	23 µl	230 µl

- 3) aliquot 23 µl in each rotor-gene 0,2 ml tubes.



- 4) Dispense in each PCR tubes 2 µl of positive controls or samples prepared in step 1.



- 5) Start thermal cycling as described in procedure paragraph.

*In order to eliminate pipetting errors we suggest to consider a reaction excess for every Master Mix constitution.