



本プロトコールは参考用の資料になります。商品ご購入の際は必ず商品に添付されている資料をご参照ください。



Enabling Discovery in Life Science®

Calcineurin Cellular Activity Assay Kit

A complete colorimetric assay for measuring cellular calcineurin (PP2B) phosphatase activity.

Instruction Manual
BML-AK816

For research use only



assay designs®



Stressgen®

www.enzolifesciences.com

✦ CALCINEURIN CELLULAR ACTIVITY ASSAY KIT - BML-AK816 ✦

BACKGROUND

Calcineurin (CaN) is the neuronal form of the widely distributed Ca^{2+} /calmodulin-dependent Ser/Thr protein phosphatase 2B (PP-2B). CaN is a heterodimer consisting of a catalytic A subunit (57-61 kDa) and a regulatory B subunit (19 kDa). The catalytic A subunit is composed of four functional domains: the catalytic core with sequence homology to PP-1 and PP-2A, binding sites for both calmodulin and CaN B-regulatory subunit, and a C-terminal autoinhibitory domain.

The *Calcineurin Cellular Activity Assay Kit* is a complete colorimetric assay kit for measuring cellular calcineurin (PP2B) phosphatase activity. It employs a convenient 96-well microtiter-plate format with all reagents necessary for measuring calcineurin (PP2B) phosphatase activity in tissue/cellular extracts PLUS, human recombinant calcineurin⁵ is included as a positive control! The RII phosphopeptide substrate, supplied with this kit, is the most efficient and outstanding peptide substrate known for calcineurin^{1,2}. The detection of free-phosphate released is based on the classic Malachite green assay^{3,4} and offers the following advantages: Non-radioactive; convenient 1-step detection; excellent sensitivity.

- Refs: 1. A. Enz *et al. Anal. Biochem.* 1994 **216** 147
 2. A. Donella-Deana *et al. Eur. J. Biochem.* 1994 **219** 109
 3. B. Martin *et al. J. Biol. Chem.* 1985 **260** 14932
 4. K.W. Harder *et al. Biochem. J.* 1994 **298** 395
 5. A. Mondragon *et al. Biochemistry* 1997 **36** 4934

ALSO AVAILABLE SEPARATELY...

PRODUCT	CATALOG #
Calcineurin Phosphatase Assay Kit	BML-AK804
BIOMOL GREEN™ Reagent 250 ml bottle	BML-AK111
Calcineurin (human, recombinant)	BML-SE163
RII phosphopeptide substrate (0.5 mg)	BML-P160
Cyclophilin A (human, recombinant)	BML-SE105
Cyclosporin A	BML-A195
CaN Autoinhibitory Peptide CN412	BML-PR104
Additional Desalting Column and Resin	BML-KI100

PLEASE READ ENTIRE BOOKLET BEFORE PROCEEDING WITH THE ASSAY. CAREFULLY NOTE THE HANDLING AND STORAGE CONDITIONS OF EACH KIT COMPONENT. PLEASE CONTACT ENZO LIFE SCIENCES TECHNICAL SERVICES FOR ASSISTANCE IF NECESSARY.

COMPONENTS

NOTE ON STORAGE: Store all components except the microtiter plate at -70°C for the highest stability. The calcineurin enzyme component BML-SE163-9090 must be handled particularly carefully in order to retain maximal enzymatic activity. Thaw it quickly in a RT water bath or by rubbing between fingers, then immediately store on an ice bath. The remaining unused enzyme should be quickly refrozen by placing at -70°C. To minimize the number of freeze/thaw cycles, aliquot the calcineurin into separate tubes and store at -70°C.

BML-SE163-9090 CALCINEURIN ENZYME (human, recombinant)
 FORM: 50 U/μl in 50 mM Tris, pH 7.5, 100 mM NaCl, 6 mM MgCl_2 , 5 mM DTT, 0.025% NP-40, 0.5 mM CaCl_2 . 1 U=1 pmol/min @30°C.
 STORAGE: -70°C; AVOID FREEZE/THAW CYCLES!
 QUANTITY: 500 U

BML-SE325-9090 CALMODULIN (human, recombinant)
 FORM: 25 μM in dH_2O
 STORAGE: -70°C
 QUANTITY: 100 μl

BML-P160-9090 SUBSTRATE (RII phosphopeptide, sequence Asp-Leu-Asp-Val-Pro-Ile-Pro-Gly-Arg-Phe-Asp-Arg-Val-pSer-Val-Ala-Ala-Glu; MW=2192.0)
 FORM: 1.5 mg net peptide/vial
 STORAGE: -20 or -70°C
 QUANTITY: 1 x 1.5 mg

BML-KI128 2X ASSAY BUFFER
 (100 mM Tris, pH 7.5, 200 mM NaCl, 12 mM MgCl_2 , 1 mM DTT, 0.05% NP-40, 1 mM CaCl_2)
 FORM: Liquid in screw-cap plastic bottle.
 STORAGE: -70°C
 QUANTITY: 20 ml

BML-KI136 2X EGTA BUFFER
 (100 mM Tris, pH 7.5, 20 mM EGTA, 200 mM NaCl, 12 mM MgCl_2 , 1 mM DTT, 0.05% NP-40)
 FORM: Screw-cap microfuge tube.
 STORAGE: -70°C
 QUANTITY: 1 ml

BML-KI135 LYSIS BUFFER
 (50 mM Tris, pH 7.5, 0.1 mM EDTA, 0.1 mM EGTA, 1 mM DTT, 0.2% NP-40)
 FORM: Liquid in screw-cap plastic bottles. (white cap)
 STORAGE: -70°C
 QUANTITY: 2 X 40 ml

BML-KI103 PROTEASE INHIBITORS COCKTAIL
 FORM: Tablet
 STORAGE: 4°C or colder
 QUANTITY: 2 tablets (1 tablet per 10 ml Lysis Buffer BML-KI135)

BML-AK111-9090 BIOMOL GREEN™ REAGENT
 FORM: Liquid in screw-cap plastic bottle.
 STORAGE: 4°C. Long-term at -70°C.
 QUANTITY: 20 ml

BML-KI132 PHOSPHATE STANDARD
 FORM: 80 μM in dH_2O
 STORAGE: -20 or -70°C
 QUANTITY: 0.5 ml

BML-EI181-9090 OKADAIC ACID
 HARMFUL! AVOID ALL FORMS OF CONTACT.
 FORM: 5 μM in 2X Assay Buffer (BML-KI128)
 STORAGE: -70°C.
 QUANTITY: 325 μl

BML-KI100 DESALTING COLUMN AND RESIN

FORM: 5 ml polypropylene disposable column and P6 DG desalting resin
STORAGE: RT. After rehydration store at 4°C.
QUANTITY: 1 column and 1 g of resin

BML-KI101 1/2 VOLUME MICROPLATES

STORAGE: RT
QUANTITY: 2 plates

OTHER MATERIALS REQUIRED

Microplate reader capable of measuring A_{620} to ≥ 3 -decimal accuracy.
Centrifuge capable of 100k x g RCF.
Swing bucket centrifuge capable of 800 x g RCF.
Pipetman capable of pipetting 5-100 μ l accurately.
Multi-channel Pipetman capable of pipetting 100 μ l (optional).
Ice bucket to keep reagents cold until use.
16g needle/syringe
TBS buffer, 100 ml (20 mM Tris, pH 7.2, 150 mM NaCl)
15 and 50 ml conical centrifuge tubes
Biological test material (e.g.: tissue, cells)

EXPERIMENTAL METHODS

NOTE ON HANDLING: Hold all samples on ice bath until use, unless otherwise noted.

PRECAUTIONS: The BIOMOL GREEN™ reagent is a highly sensitive phosphate detection solution. Free phosphate present on labware and in reagent solutions will greatly increase the background absorbance of the assay. This is detected visually as a change in color from yellow to green. Detergents used to clean labware may contain high levels of phosphate. Use caution by either rinsing labware with dH₂O or employ unused plasticware. Do not use phosphate buffered saline (PBS) for any tissue/cell rinses -use TBS (Tris buffered saline)!

PREPARATION OF TISSUE/CELL EXTRACTS

To prepare a tissue/cell extract for calcineurin activity assay:

Calcineurin activity is highly dependent on the experimental conditions and the cell/tissue source. Therefore, the amount of material required for an assay should be determined empirically by the user. Typically, between 0.5-5 μ g of total protein or 5,000 and 50,000 cells per assay will provide sufficient signal for detection.

NOTE ON BIOLOGICAL SAMPLE MATERIAL: The following procedures have been tested for rat and mouse brain tissue. Other tissue or cell culture samples employed may require adjustment to this protocol for satisfactory results.

1. Add protease inhibitor tablets to lysis buffer (BML-KI135) immediately before use (1 tablet/10 ml buffer). Vortex.
2. Obtain tissue, if fresh, excise quickly.
3. Rinse tissue quickly in ice-cold TBS and shake-off/blot excess wetness.
4. Weigh the tissue in centrifuge tube.
5. Add lysis buffer (BML-KI135) with protease inhibitors to tissue. Use 0.33-0.5 ml per gram of tissue.
6. Loosely break up cells by passing them through a 16g needle. Avoid air bubbles.
7. Optional: Sediment at 100-200k x g in centrifuge at 4°C for 45 min. Save supernatant=HSS (high-speed supernatant). Please note that this step will sediment the nucleus and any associated nuclear calcineurin.
8. Freeze immediately at -70°C.

REMOVAL OF FREE PHOSPHATE FROM EXTRACTS

To desalt tissue samples by gel filtration:

NOTE: This procedure is intended to remove excess phosphate and nucleotides (which are slowly hydrolyzed to release free phosphate in the presence of the BIOMOL GREEN™ reagent) in the high speed supernatant (HSS) extract.

1. Rehydrate Desalting Column Resin (BML-KI100) in a 50 ml conical tube by adding 20 ml of phosphate free dH₂O and vortexing briefly. Allow to set for 4 hours at RT or overnight at 4°C.
2. Decant the dH₂O carefully, then add fresh dH₂O at a 1:1 ratio to the rehydrated resin (BML-KI100) (~10 ml).
3. Add rehydrated resin (BML-KI100) to the Chromatography Column (BML-KI100) to obtain a 5 ml settled-bed volume (~5.5 cm bed height). Remove tip from column and allow dH₂O to drain by gravity.
4. Equilibrate column by adding 8 ml of lysis buffer (BML-KI135) (without protease inhibitors) and allow to drain by gravity.
5. Place column in a 15 ml centrifuge tube. Centrifuge at 800 x g for 3 min at 4°C to displace column buffer. Discard flow-through buffer.
6. Place column in a clean 15 ml centrifuge tube.
7. Add up to 350 μ l HSS sample, from above, to column.
8. Centrifuge at 800 x g for 3 min. Save extract flow-through. This is the desalted cell lysate material to be tested for calcineurin activity, below.
9. Freeze sample immediately at -70°C.

TIP: The effective removal of phosphate/nucleotides from the extract should be tested qualitatively by adding 100 μ l BIOMOL GREEN™ reagent to 1 μ l extract, and a separate sample of 1 μ l dH₂O. If no phosphate/nucleotides are present, both samples should remain yellow in color over a time period of 30 min @ RT. The development of a visible green color indicates phosphate contamination, which must be eliminated from the samples before proceeding further!

BIOMOL GREEN™ PHOSPHATASE ASSAY

To prepare reagents for the assay:

1. Thaw all kit components and hold on ice bath, except BIOMOL GREEN™ reagent at RT.
2. Add calmodulin to the 2x assay buffer: Dilute calmodulin (BML-SE325-9090) 1/50 in 2X assay buffer (BML-KI128) to required quantity (25 μ l are required per assay well). For example, add 20 μ l to 980 μ l 2X assay buffer.
3. Reconstitute substrate (R11 phosphopeptide, BML-P160) with dH₂O to 0.75 mM (1.64 mg/ml): Add 915 μ l dH₂O per 1.5 mg vial (10 μ l are needed per assay well).

To prepare phosphate standard curve sample wells:

1. Prepare 1 ml of 1X assay buffer (dilute 500 μ l of 2X assay buffer with 500 μ l of dH₂O)
2. Perform 1:1 serial dilutions of phosphate standard and an assay buffer blank. Concentrations of 40, 20, 10, 5, 2.5, 1.25 and 0.625 μ M correspond to 2, 1, 0.5, 0.25, 0.125, 0.063 and 0.031 nmol PO₄ (see Table 1):
 - a) Add 50 μ l of assay buffer (BML-KI128) to wells A1, and A2 (2nmol PO₄ standards)
 - b) Add 50 μ l of 1X assay buffer (prepared in step 1 above) to wells B1-H1 and wells B2-H2. (remaining standard concentrations)
 - c) Add 50 μ l of 80 μ M phosphate standard to wells A1 and A2 of assay plate. Mix thoroughly by pipetting up and down several times.
 - d) Remove 50 μ l from well A1 and add it to well B1. Mix thoroughly by pipetting up and down several times.
 - e) Remove 50 μ l from well B1 and add it to well C1.

- f) Mix thoroughly and repeat for wells D1-G1. At well G1, remove 50 μ l and discard. DO NOT PROCEED TO WELL H1 (assay buffer blank). Final volume=50 μ l.
- g) Repeat serial dilution for the wells in column 2 (standard curve duplicates)

To prepare calcineurin activity assay sample wells:

(See Tables 1 & 2):

Background (no substrate):

{Control for background phosphate/interfering substances}

1. Add 20 μ l dH₂O to appropriate wells.
2. Add 25 μ l 2X assay buffer (BML-KI128) w/ calmodulin to each well.

Total phosphatase activity wells:

{Total phosphatase activity in the extract }

1. Add 10 μ l dH₂O to each well.
2. Add 25 μ l 2X assay buffer (BML-KI128) w/ calmodulin to each well.

EGTA buffer (Ca²⁺/CaM free):

{Total activity less PP2B (calcineurin)}

1. Add 10 μ l dH₂O to each well.
2. Add 25 μ l 2X EGTA buffer (BML-KI136) to each well.

OA (okadaic acid):

{Total activity less PP1 & PP2A}

1. Add 5 μ l dH₂O to each well.
2. Add 25 μ l 2X assay buffer (BML-KI128) w/ calmodulin to each well.
- 2+. Add 5 μ l okadaic acid (BML-EI181-9090; 5 μ M)

OA + EGTA:

{Total activity less PP1, PP2A & PP2B}

1. Add 5 μ l dH₂O to each well.
2. Add 25 μ l 2X EGTA buffer (BML-KI136) to each well.
- 2+. Add 5 μ l okadaic acid (BML-EI181-9090; 5 μ M)

Positive control (calcineurin enzyme):

{Purified CaN enzyme positive control}

1. Add 10 μ l dH₂O to each well.
2. Add 25 μ l 2X assay buffer (BML-KI128) w/ calmodulin to each well.

Add phosphopeptide substrate:

3. Add 10 μ l phosphopeptide substrate BML-P160 to each well of the calcineurin samples except the "Background" control. DO NOT ADD SUBSTRATE TO THE PHOSPHATE STANDARD CURVE SAMPLES!
4. Equilibrate microtiter plate to reaction temperature (e.g.: 30°C) for 10 min.

To initiate calcineurin assay:

5. Add 5 μ l extract or diluted calcineurin (dilute to 8 U/ μ l prior to use) to appropriate wells. For sample extract wells, it may be necessary to dilute the HSS tissue extract (e.g.: 1/5-1/10 in lysis buffer BML-KI135). For calcineurin "Positive control" add 5 μ l BML-SE163-9090 (40 U/well).
6. Incubate plate at reaction temperature for desired duration (e.g.: 30 min @ 30°C).

To terminate reactions:

7. After incubating wells for desired duration, terminate reactions by adding 100 μ l BIOMOL GREEN™ reagent (BML-AK111-9090) to ALL samples including the phosphate standard curve.
8. Allow color to develop 20-30 minutes, making sure all wells spend approximately the same time with the reagent before reading on microplate reader.

9. Read OD_{620nm} on microplate reader.
10. Perform data analysis (see below).

NOTE: Retain microplate for future use of unused wells!

TABLE 1. EXAMPLE OF MICROTITER PLATE SAMPLES.

Sample†	Std Curve	Extract/CaN Samples
Well #	1,2	3,4
A	2 nmol PO ₄	Background
B	1	Total
C	0.5	EGTA buffer
D	0.25	OA (okadaic acid)
E	0.125	OA + EGTA
F	0.063	Positive Control
G	0.031	
H	0	

† For highest accuracy, perform all samples in duplicate.

TABLE 2. TYPICAL ASSAY COMPONENTS

	H ₂ O	2X Assay Buffer w/CaM	2X EGTA Buffer	OA	Substrate (0.75 mM)	Extract/ CaN (dilute to 8 U/ μ l)
Background	20 μ l	25 μ l	0	0	0	5 μ l ^a
Total	10 μ l	25 μ l	0	0	10 μ l	5 μ l ^a
EGTA buffer	10 μ l	0	25 μ l	0	10 μ l	5 μ l ^a
OA (okadaic acid)	5 μ l	25 μ l	0	5 μ l	10 μ l	5 μ l ^a
OA + EGTA	5 μ l	0	25 μ l	5 μ l	10 μ l	5 μ l ^a
Positive Control	10 μ l	25 μ l	0	0	10 μ l	5 μ l ^b

a Add cellular extract (HSS)

b Add calcineurin enzyme (SE163-9090)


DATA ANALYSIS

Phosphate (PO₄) Standard Curve

1. Plot standard curve data as OD_{620nm} versus nmol PO₄ (see Figure 1).
2. Obtain a line-fit to the data using an appropriate routine.
3. Use the slope and Y-intercept to calculate amount of phosphate released for the experimental data (see below).

NOTE: For highest accuracy, a standard curve must be performed for each new set of assay data. This will normalize for variations in free phosphate in samples, time of incubation with the BIOMOL GREEN™ reagent, and other experimental factors.

FIGURE 1. BIOMOL GREEN™ PHOSPHATE STANDARD CURVE



Conversion of OD_{620nm} to Amount Phosphate Released

1. Convert OD_{620nm} data into the amount of phosphate released using the standard curve line-fit data, from above:

$$\text{Phosphate released} = (\text{OD}_{620\text{nm}} - Y_{\text{int}}) / \text{slope}$$

EXAMPLE:

Std curve slope=0.3 OD_{620nm}/nmol phosphate

Std curve Yint=0.001 OD_{620nm}

Sample OD=0.4

$$\text{Phosphate released} = (0.4 - 0.001) / 0.3 = 1.33 \text{ nmol}$$

Data Reduction to Determine Calcineurin Phosphatase

PRECAUTION: The procedures for data analysis which follow are intended only as a guideline. The individual user must determine the suitability of this analysis for their particular experimental protocol. Additional controls and other samples may be appropriate for accurate analysis.

ANALYSIS DESCRIPTION: This assay uses the RII phosphopeptide, the best known substrate for CaN (PP2B). Nonetheless, in cellular extracts, the phospho-group is cleaved by other competing phosphatases. Thus, a series of conditions must be employed to discriminate between the contribution of other phosphatases. CaN requires calcium for its activity, thus the "EGTA buffer" sample represents total phosphatase activity less CaN. Okadaic acid (OA) at 100 and 500 nM is known to completely inhibit PP1 and PP2A, while it has no effect on CaN (see Figure 3). Finally, OA + EGTA buffer inhibits PP1, PP2A and PP2B, but not PP2C. Thus, for a given biological system, the analysis of these samples allows the quantification of calcineurin (PP2B) activity in a cellular extract. These methods have been described in the literature. Additional experimental modulation of the cell extracts may be desirable. Inhibitors of calcineurin, calmodulin and calcium ion modulators may be appropriate.

Suggested reading

Minireview:

C.B. Klee *et al.* *J. Biol. Chem.* 1998 **273** 13367

Brain cellular extracts:

P.M. Stemmer *et al.* *FEBS Lett.* 1995 **374** 237

X. Wang *et al.* *Nature* 1996 **383** 434

Human recombinant CaN:

A. Mondragon *et al.* *Biochem.* 1997 **36** 4934

CaN Assays:

D.A. Fruman *et al.* *Methods: A Companion to Meth. Enzymol.* 1996 **9** 146

1. Subtract the "Background" phosphate released from each sample except the "Positive control".
2. Plot a graph analogous to Figure 2, below. Use either OD_{620nm} or phosphate released for the Y-axis.
3. Determine the contribution of calcineurin:

$$\text{eq. 1 CaN (PP2B)} = \text{Total} - \text{EGTA buffer}$$

or

$$\text{eq. 2 CaN (PP2B)} = \text{OA} - (\text{OA} + \text{EGTA})$$

Eq. 1 is a conventional method to report CaN activity. However, the user must determine the most appropriate analysis for their specific experimental goal.

Cellular CaN Assay: Mouse Brain Tissue




FIGURE 2. CELLULAR CaN ASSAY USING MOUSE BRAIN

Phosphatase activity from a freshly prepared mouse brain extract. Prior to gel filtration the extract was diluted 1:1 in lysis buffer containing protease inhibitors. After gel filtration and prior to the phosphatase assay, the extract was diluted 1/25 in lysis buffer. Well contents were as in Table 2. The reaction was incubated for 30 min @ 30°C.

CaN Positive Control: Effect of EGTA & OA




FIGURE 3. CaN POSITIVE CONTROL

Purified calcineurin (BML-SE163-9090) positive control. The CaN was incubated 1 hr @ 30°C under various buffer conditions. The results demonstrate that CaN activity is inhibited by the EGTA buffer, but not inhibited by 100 nM or 500 nM concentrations of okadaic acid.

Literature Citations of Calcineurin Assay Kits

B. Mehul *et al. J. Biol. Chem.* 2000 **275** 12841
T. Taigen *et al. Proc. Natl. Acad. Sci.* 2000 **97** 1196
G. Mallert *et al. Cell* 2001 **104** 675
M. Ichida and T. Finkel *J. Biol. Chem.* 2001 **276** 3524
O. Bueno *et al. Proc. Natl. Acad. Sci.* 2002 **99** 4586
K.J. Sellar *et al. Anal. Biochem.* 2006 **358** 104

USE FOR RESEARCH PURPOSES ONLY

Unless otherwise specified expressly on the packaging, all products sold hereunder are intended for and may be used for research purposes only and may not be used for food, drug, cosmetic or household use or for the diagnosis or treatment of human beings. Purchase does not include any right or license to use, develop or otherwise exploit these products commercially. Any commercial use, development or exploitation of these products or development using these products without the express written authorization of Enzo Life Sciences, Inc. is strictly prohibited. Buyer assumes all risk and liability for the use and/or results obtained by the use of the products covered by this invoice whether used singularly or in combination with other products.

LIMITED WARRANTY; DISCLAIMER OF WARRANTIES

These products are offered under a limited warranty. The products are guaranteed to meet all appropriate specifications described in the package insert at the time of shipment. Enzo Life Sciences' sole obligation is to replace the product to the extent of the purchasing price. All claims must be made to Enzo Life Sciences, Inc., within five (5) days of receipt of order. THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTIES OR LIABILITIES, EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND NONINFRINGEMENT OF THE PATENT OR OTHER INTELLECTUAL PROPERTY RIGHTS OF OTHERS, AND ALL SUCH WARRANTIES (AND ANY OTHER WARRANTIES IMPLIED BY LAW) ARE EXPRESSLY DISCLAIMED.





www.enzolifesciences.com

Enabling Discovery in Life Science®

Global Headquarters

Enzo Life Sciences Inc.

10 Executive Blvd

Farmingdale, NY 11735

(p) 1-800-942-0430

(f) 1-631-694-7501

(e) info-usa@enzolifesciences.com

Enzo Life Sciences (ELS) AG

Industriestrasse 17, Postfach

CH-4415 Lause / Switzerland

(p) +41/0 61 926 89 89

(f) +41/0 61 926 89 79

(e) info-ch@enzolifesciences.com

Please visit our website at www.enzolifesciences.com for additional contact information.

