

Polink-1 HRP Rat-NM (No cross react with Mouse) DAB Detection System for Rat Primary Antibody

(Polymer-HRP detection system, biotin-free, Anti-rat primary antibody)
Ready-to-use One Step Polymer Detection System
Super clean when using rat antibody on mouse tissue

Storage: 4-8°C

| | | |
|-------------|----------------------------------|---------------------------------------|
| Catalog No. | <input type="checkbox"/> D35-110 | 110 ml (bulk, w/o chromogen) |
| | <input type="checkbox"/> D35-18 | 18 ml (with DAB, good for 180 slides) |
| | <input type="checkbox"/> D35-6 | 6 ml (with DAB, good for 50 slides) |

Intended Use:

Detecting RAT primary antibody on MOUSE tissue is a very difficult task in research field due to background staining issues. Polink-1 HRP Rat-NM (No-Mouse) DAB Detection kit is specially designed to solve the problem. This technology provides excellent specificity to detect rat primary antibody (user supplied) on mouse tissue. Specimen can be frozen tissues, paraffin-embedded tissues, or freshly prepared monolayer cell smears.

Polink-1 HRP Rat-NM DAB Detection kit is a 1-step polymer detection system that uses polymeric HRP-linked anti-rat secondary antibody to directly detect rat primary antibody bound to the mouse tissue. The secondary antibody is adsorbed to mouse, rabbit and human serum proteins. Besides mouse tissue Polink-1 HRP Rat-NM DAB Detection kit also can be used on human tissue and rabbit tissue as well. It is a biotin-free system, therefore, overcomes the non-specific staining caused by endogenous biotin¹. It is a 1-step detection system is a much faster assay compared to traditional two step methods (Biotinylated 2nd antibody, and then streptavidin-HRP). These advantages provide laboratories the benefit of more accurate and quicker result, less trouble shooting and better cost-saving.

If users need a more sensitive polymer detection system for rat primary antibody on mouse tissue, they may choose a two-step polymer detection system, Polink-2 Plus HRP Rat-NM DAB kit (Cat No. D46-110, D46-18, D46-6). For AEC staining please choose Polink-1 HRP Rat-NM for AEC (D36-110, D36-18, and D36-6).

Kit components:

| Catalog No. | Product Name | Reagent 1: Polymer HRP-linked anti-Rat IgG (No cross react with mouse) (Ready-to-use) | Reagent 2: 2A: DAB Substrate 2B: Chromogen concentrate |
|-------------|---------------------------------------|--|--|
| D35-110 | Polink-1 HRP RAT-NM Bulk for DAB kit | 110ml | Not provided |
| D35-18 | Polink-1 HRP RAT-NM with DAB 18ml kit | 18ml | 30 ml of 2A and 2 ml of 2B |
| D35-6 | Polink-1 HRP RAT-NM with DAB 6ml kit | 6ml | 12 ml of 2A and 1.5 ml of 2B |

Recommended Protocol:

1. Fixation: To ensure the quality of the staining and obtain reproducible performance, users need to supply appropriately fixed tissue and well prepared slides.
2. Tissue needs to be adhered to the slide tightly to avoid tissue falling off.
3. Paraffin embedded section must be deparaffinized with xylene and rehydrated with a graded series of ethanol before staining.
4. Cell smear samples should be made as much monolayer as possible to obtain satisfactory results.
5. Investigator needs to optimize dilution and incubation times for primary antibodies.
6. Three control slides will aid the interpretation of the result: positive tissue control, reagent control (slides treated with Isotype control reagent), and negative control.
7. Proceed IHC staining: DO NOT let specimen or tissue dry from this point on.

| Reagent | Staining Procedure | Incubation Time (Min.) |
|--|---|------------------------------|
| 1. Peroxidase Blocking Reagent Supplied by user | a. Incubate slides in peroxidase blocking reagent (Ready-to-use 3% H ₂ O ₂ solution) for 10 minutes. b. Rinse the slide using distilled water. | 10 |
| 2. HIER Pretreatment: Refer to antibody data sheet. | a. Heat Induced Epitope Retrieval (HIER) may be required for primary antibody suggested by vendor. b. Wash with PBS for 2 minutes each time for 3 times. | Refer to vendor's data sheet |
| 3. Pre-Block (Optional) Not provided | a. Add 2 (100 µL) or more drops of Pre-Block solution to cover the tissue section and Incubate | 10 |

| | | |
|--|---|-----------------|
| | b. Drain or blot off solution. DO NOT RINSE. | |
| 4. Primary antibody: Supplied by user | Notes: Investigator needs to optimize dilution and incubation times a. Apply 2 (100 µL) or more drops of primary antibody to cover the tissue completely. Incubate in moist chamber for 30-60 min. b. Rinse with PBS containing 0.05% Tween-20 for 2 minutes each time for 3 times. | 30-60 |
| 5. Reagent 1: Polymer-HRP anti-rat (Ready-to-use) | a. Apply 2 (100 µL) or more drops of Polymer- HRP anti-Rat 2 nd antibody to cover tissue section and Incubate in moist chamber for 10-15 min. c. Rinse with PBS containing 0.05% Tween-20 for 2 minutes each time for 3 times. | 10-15 |
| 6. Reagents 2A, 2B: 2A: DAB Substrate 2B: DAB Chromogen | a. Adding 1 drop or 2 drops (for higher contrast) of DAB chromogen concentrate (Reagent 2B) in 1ml of DAB substrate buffer (Reagent 2A). Mix well. b. Apply 2 drops (100 µL) or enough volume of pre-mixed DAB Chromogen to completely cover tissue. Incubate for 5 min. use the prepared DAB solution within 5 hours c. When appropriate color is developed, rinse under tap water gently for about 1-2 minutes. | 3-10 |
| 8. Hematoxylin: Supplied by user. | a. Counterstain with 2 (100 µl) or more drops hematoxylin to cover tissue completely and wait about 20 seconds. b. Rinse well with tap water for 1-2 min. c. Put slides in PBS until the color turn blue (about ½ - 1 min.) d. Rinse in distill water, then rinse well with tap water | 20-30 seconds |
| 9. Mounting medium: Supplied by user | Follow the manufacture data sheet procedure for mounting. Recommended product: 1. GB-Mount: Cat. No. E01-18 (18ml), for alcohol soluble substrates (AEC, AP-Red, and AP-blue) 2. O-Mount: Cat. No. E02-18 (18ml), for DAB and BCIP/NBT 3. Simpo-Mount: Cat.No. E03-18 (18ml), or E03-100 (100ml), universal permanent mounting medium. Can be used with or without cover slip | Refer to insert |

Protocol Notes:

1. The fixation, tissue slide thickness, and primary antibody dilution and incubation time affect results significantly. Users need to consider all factors and determine optimal conditions when interpreting the result.
2. Tissue staining is dependent upon the proper handling and processing of tissues prior to staining. Improper tissue preparation may lead to false negative results or inconsistent results.
3. Do not mix reagents from different lot.
4. Do not allow the slides to dry at any time during staining

Related Products:

| Product | Catalog No. | Size | | Product | Catalog No. | Size |
|---------------------------------------|----------------|------------|--|--|-----------------|--------------|
| Polink-1 HRP Broad Bulk kit for DAB | D11-110 | 110ml | | Polink-1 HRP Mouse 18ml, 6ml DAB Kit | D12-18 / D12-6 | 18ml / 6ml |
| Polink-1 HRP Broad 18ml, 6ml DAB Kit | D11-18 / D11-6 | 18ml / 6ml | | *Polink-1 HRP Mouse-NR Bulk kit for DAB | D55-110 | 110ml |
| Polink-1 HRP Rabbit Bulk kit for DAB | D13-110 | 110ml | | *Polink-1 HRP Mouse-NR 18ml, 6ml DAB Kit | D55-18 / D55-6 | 18ml / 6ml |
| Polink-1 HRP Rabbit 18ml, 6ml DAB Kit | D13-18 / D13-6 | 18ml / 6ml | | DAB Kit (2-components) | C09-12 | 12ml +240ml |
| Polink-1 HRP Goat Bulk kit for DAB | D33-110 | 110ml | | O-Mount (Organic) | E02-18 | 18ml |
| Polink-1 HRP Goat 18ml, 6ml DAB Kit | D33-18 / D33-6 | 18ml / 6ml | | Simpo-Mount (Aqueous) | E03-100/ E03-18 | 100ml / 18ml |
| Polink-1 HRP Mouse Bulk kit for DAB | D12-110 | 110ml | | | | |

*Polink -1 HRP Mouse-NR kit does not cross react with Rat primary antibody

Precautions:

AEC may be carcinogenic. Please wear gloves and take other necessary precautions.

Remarks:

For research use only.

References:

1. Bisgaard K, Pluzed KP. *Use of polymer conjugates in immunohistochemistry: A comparative study of a traditional staining method to a staining method utilizing polymer conjugates.* Abstract XXI Intl Cong Intl Acad Pathol and 12th World Cong Acad Environ Pathol. Budapest, Hungary, October 20-25, 1996.
2. Shi ZR, Itzkowitz SH, Kim YS. *A comparison of three immunoperoxidase techniques for antigen detection in colorectal carcinoma tissues.* J Histochem Cytochem 36:317-322,