

## Polink TS-MRR-Ms A Kit for Immunohistochemistry Staining

**Polymer-HRP & AP triple staining kit to detect one mouse and two rabbit primary antibodies on mouse/rat tissue with DAB(Brown), AP-Red(Red), and Emerald(Green)**

Storage: 2-8°C
----------------

Catalog No.:	<input type="checkbox"/>	TS309A-6	*24ml (for 120 slides**)
	<input type="checkbox"/>	TS309A-18	*72ml (for 360 slides**)
	<input type="checkbox"/>	TS309A-60	*240ml (for 1200 slides**)

*\*Volume of polymer conjugate  
\*\* If use 100ul per slide*

### Intended Use:

The **Polink TS-MRR-Ms A Kit** is designed to use with user supplied one mouse and two rabbit primary antibodies to detect three distinct antigens on a single mouse/rat tissue or cell samples. Kit has been tested on tissue specimens that are paraffin embedded. For frozen tissue a lower temperature of 65°C must be used for Antibody Blocker (Reagent 6) to prevent tissue from dissociating from slide. Please read through entire protocol as this protocol requires many step to be done in the defined order.

Triple staining uses traditional methods in immunohistostaining to reveal three distinct antigens and their co-expression on a single tissue<sup>1,2</sup>. **Polink TS-MRR-Ms A Kit** from GBI Labs (Golden Bridge International) supplies polymer enzyme conjugates: polymer-HRP anti-mouse IgG, polymer-HRP anti-rabbit IgG and polymer-AP anti-rabbit IgG with three substrates/chromogens; DAB (brown), Emerald (green), and AP-Red+ (Red). **Polink TS-MRR-Ms A Kit** is a non-biotin system, avoiding non-specific binding caused by endogenous biotin. This kit has been optimized to have no cross detection when detecting more two primary antibodies from the same host species using unique blocking system. Simplified steps allow users to complete triple staining within 5 hours (without antigen retrieval) or 6-7 hours (with antigen retrieval). The well tested protocol provides user a method to permanently mount slides with coverslip.

### Kit Components:

Component No.	Content	TS309A-6	TS309A-18	TS309A-60
Reagent 1	<b>Ms HRP-Primer (RTU)</b>	12ml	36ml	120ml
Reagent 2	<b>Ms HRP-Polymer (RTU)</b>	6ml	18ml	60ml
Reagent 3	<b>Rb AP-Polymer (RTU)</b>	6ml	18ml	60ml
Reagent 4A	<b>DAB substrate (RTU)</b>	12ml	18mlx2	120ml
Reagent 4B	<b>DAB chromogen (20x)</b>	1.5ml	2ml	6ml
Reagent 5A	<b>AP-Red+ Enhancer (40x)</b>	1ml	1ml	2ml
Reagent 5B	<b>AP-Red+ Solution (40x)</b>	1ml	1ml	2ml
Reagent 5C	<b>AP-Red+ Substrate (20x)</b>	4ml	4ml	8ml
Reagent 6	<b>Antibody Blocker (40x)</b>	2 x 15ml	100ml	3 x 100ml
Reagent 7	<b>Rb HRP-Polymer (RTU)</b>	12ml	36ml	120ml
Reagent 8	<b>Emerald Chromogen (RTU)</b>	12ml	36ml	120ml

HRP = Horseradish Peroxidase AP = Alkaline Phosphatase Ms = Mouse Rb = Rabbit

### Protocol Notes:

1. Fixation: To ensure the quality of the staining and obtain reproducible performance, user needs to supply appropriately fixed tissue and well prepared slides.
2. Tissue needs to be adhered to the slide tightly to avoid falling off.
3. Paraffin embedded sections must be deparaffinized with xylene and rehydrated with a graded series of alcohols before staining.
4. Cell smear samples should be prepared as close to a monolayer as possible to obtain satisfactory results.
5. Control slides are recommended for interpretation of results: positive, reagent (slides treated with Isotype control reagent), and negative control.
6. **DO NOT** let specimen or tissue dry during protocol. This will generate false positive and/or false negative signal.
7. **Important:** Incubate 1st primary mouse antibody with rabbit antibody. Never combine two antibodies from the same host species in one incubation step.
8. The fixation, tissue slide thickness, antigen retrieval and primary antibody dilution and incubation time effect results significantly. Investigator needs to consider all factors and determine optimal conditions when interpreting results.

### Precautions:

Please wear gloves, eye protection and take other necessary precautions. If any of the reagent come in contact with skin wash area completely with plenty of water and soap. If irritation develops seek medical attention.

### Equipment or material needed but not provided:

1. Equipment and material for deparaffinization, such as fume absorbing hood, etc.

2. Heat source (microwave or hot plate) for HIER and antigen retrieval buffers
3. Thermometer
4. Timer
5. Beaker
6. Slides warmer (or hair dryer, hybridizer, vacuum oven) that can keep temperature to 85°C during slides heat drying process
7. Wash buffer: 0.01 M PBS with 0.5% Tween20, pH7.4
8. Peroxidase and alkaline phosphatase blocking buffer
9. 100% ethanol
10. 100% Xylene
11. Hematoxylin
12. 100% methanol at -20°C
13. Mounting medium
14. Coverslip

#### Staining protocol selection and limitation of the kit:

- Most antigens will not be destroyed by heat. However, users need to check if there are proteins on the tissue that are heat sensitive before proceeding with the staining.
- You may encounter conditions that 1st mouse antibody and one rabbit antibody need HIER and the 3rd protein detected by 2<sup>nd</sup> rabbit antibody is heat sensitive. In this situation you may download our triple color staining protocol from our web site.
- Please read the following table carefully before you start the experiment to ensure the result.
- This kit is not suitable for the following condition: 2 proteins are heat sensitive and detected by 2 mouse antibodies and one rabbit antibody requires HIER.

#### Staining Protocol TS309A

Steps / Reagent	Staining Protocol	Incubation Time
1. Peroxidase and alkaline phosphatase Blocking Reagent Supplied by user.	<ol style="list-style-type: none"> <li>a. Incubate slides in peroxidase and alkaline phosphatase blocking reagent for 10 minutes.</li> <li>b. Rinse slides using 2 changes of distilled water.</li> </ol>	10 minutes
2. Antigen retrieval: <b>Refer to primary antibody data sheet.</b> Supplied by user.	<ol style="list-style-type: none"> <li>a. Refer to primary antibody data sheet for antigen retrieval methods.</li> <li>b. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	Up to 1 hour
3. Primary Antibody Mix: <b>Mix one Mouse and one Rabbit primary antibody</b>  Supplied by user.	<p><b>Note:</b> Investigator needs to optimize dilution prior to triple staining. <b>DO NOT</b> combine the same host species primary antibodies together at this step.</p> <ol style="list-style-type: none"> <li>a. Apply 2 drops or enough volume of mouse and rabbit primary antibody mixture to cover the tissue completely. Incubate in moist chamber for 30-60 min. Recommend 30min to shorten total protocol time.</li> <li>b. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	30 minutes
4. <b>Reagent 1: Ms HRP-Primer</b>	<p><b>Note:</b> This step is required for activation of Ms HRP Polymer, <b>DO NOT</b> skip.</p> <ol style="list-style-type: none"> <li>a. Apply 1 to 2 drops (50-100µl) of the mixture to cover the tissue completely.</li> <li>b. Incubate in moist chamber for 10 min.</li> <li>c. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	10 minutes
5. <b>Mix Reagent 2: Ms HRP-Polymer (RTU) with Reagent 3: Rb AP-Polymer (RTU)</b>	<p><b>Note:</b> Make sufficient polymer mixture by adding <b>Reagent 2 (Ms HRP-Polymer)</b> and <b>Reagent 3 (Rb AP-Polymer)</b> at 1:1 ratio, mix well. Do Not Mix More than you need for the experiment because the polymer mixture may not be as stable as non-mixed polymer.</p> <ol style="list-style-type: none"> <li>a. Apply 1 to 2 drops (50-100µl) of the mixture to cover the tissue completely.</li> <li>b. Incubate in moist chamber for 30 min.</li> <li>c. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	30 minutes
6. <b>Reagent 4A: DAB Substrate (RTU)</b> <b>Reagent 4B: DAB Chromogen (20x)</b>	<p><b>Note:</b> Make enough DAB mix by adding 1 drop of <b>Reagent 4B (DAB chromogen)</b> in 1ml of <b>Reagent 4A (DAB substrate)</b>. Mix well. Use within 7 hours. Prepare AP-Red Plus at this time (see step 7).</p> <ol style="list-style-type: none"> <li>a. Apply 1 to 2 drops (50-100µl) of your DAB mixture to cover the tissue completely.</li> <li>b. Incubate for 5min.</li> <li>c. Rinse slides in multiple changes of distilled water 3 times, 2 each time or under running tap water for 1minute.</li> </ol>	5 minutes

<p>7. Reagents for AP-Red+ 5A, 5B, 5C  <b>5A: AP-Red+ Enhancer</b> (40x)  <b>5B: AP-Red+ Solution</b> (40x)  <b>5C: AP-Red+ Substrate</b> (20x)</p>	<ol style="list-style-type: none"> <li>a. Add 1 drop (50µl) of <b>Reagent 5A (Enhancer)</b> and 1 drop of <b>Reagent 5B (Solution)</b> to a test tube. Mix well and set at room temperature for 5 minutes.</li> <li>b. Add 2ml of distilled water to the mixture. Mix well.</li> <li>c. Add 4 drops (200µl) of <b>Reagent 5C (Substrate)</b> to the mixture and mix well.</li> <li>d. Apply 2 drops (100µl) or enough volume of AP-Red Plus mixture to cover the tissue completely. Incubate for 20 minutes observe appropriate color development.</li> <li>e. Rinse slides in multiple changes of distilled water 2 minutes 3 times or running tap water for 1 minute.</li> </ol>	<p>20 minutes</p>
<p>8. <b>Reagent 6: Antibody Blocker</b> (40x)</p>	<p><b>Note:</b> This step will block antibodies of previous step so no cross reaction will occur at end of protocol. HIER can be done immediately after <b>Antibody Blocker</b> step if only one primary antibody requires antigen retrieval.</p> <ol style="list-style-type: none"> <li>a. Use hot plate or water bath to heat diluted <b>Reagent 6</b> to 1x solution (1 part of <b>Antibody Blocker</b> in 39 parts of distilled water) to 80-95°C. Make enough volume to cover the tissue in beaker.</li> <li>b. Put slides in heated Antibody Blocker for 10 minutes at 80-95°C.</li> <li>c. Remove slides from the Antibody blocker; cool slides 5 seconds.</li> <li>d. Rinse slides in multiple changes of distilled water.</li> <li>e. Wash with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	<p>10 minutes</p>
<p>9. Antigen retrieval:  <b>Refer to primary antibody data sheet.</b></p>	<ol style="list-style-type: none"> <li>a. Refer to primary antibody data sheet for antigen retrieval methods.</li> <li>b. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	<p>Up to 1 hour</p>
<p>10. 2<sup>nd</sup> Rabbit primary antibody  Supplied by user.</p>	<p><b>Note:</b> Investigator needs to optimize dilution prior to triple staining.</p> <ol style="list-style-type: none"> <li>a. Apply 2 drops or enough volume of the 2<sup>nd</sup> rabbit primary antibody to cover the tissue completely. Incubate in moist chamber for 30-60 min. Recommend 30 minutes to shorten total protocol time.</li> <li>b. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> </ol>	<p>30 minutes</p>
<p>11. <b>Reagent 7:Rb HRP-Polymer</b> (RTU)</p>	<ol style="list-style-type: none"> <li>a. Apply 1 to 2 drops (50-100µl) of <b>Reagent 7 (Rb HRP-Polymer)</b>, to cover the tissue completely. Incubate slides in moist chamber for 15 min.</li> <li>b. Wash slides with PBS/ 0.05% Tween20 for 2 minutes, 3 times.</li> <li>c. Rinse slides with Tap water for 3 times.</li> <li>d. You need to set up temperature (see step 13) for heat drying process at the end of this step.</li> </ol>	<p>15 minutes</p>
<p>12. Counterstain (<b>Optional</b>)  Not provided</p>	<ol style="list-style-type: none"> <li>a. Counterstain dip in diluted hematoxylin for 5 seconds. (Dilute hematoxylin 1:5 in dH2O). <b>DO NOT</b> over stain with hematoxylin.</li> <li>b. Rinse thoroughly with tap water for 15-30 seconds.</li> <li>c. Put slides in PBS for 5 seconds to blue, <b>DO NOT</b> over blue.</li> <li>d. Rinse well in distilled or tap water for 15-30 seconds.</li> </ol>	<p>5 seconds</p>
<p>13. <b>Reagent 8: Emerald Chromogen</b> (RTU)   <b>Do counterstain first.</b></p>	<ol style="list-style-type: none"> <li>a. Apply 1 to 2 drops (50-100µl) of <b>Reagent 8 Emerald Chromogen</b> to cover the tissue completely.</li> <li>b. Incubate slides in humid chamber for 5 minutes.</li> <li>c. Wash slides in tap water for 10 seconds!</li> </ol> <p><b>Important to READ:</b> Emerald Chromogen is water soluble, do counter stain first. <i>Do not leave slides sitting in water.</i> Always stain Emerald chromogen <b>AFTER</b> AP-Red+ stain because AP-Red+ removes the Emerald and after hematoxylin.</p>	<p>5-10 minutes</p>
<p>14. Heat drying slides:  User may choose any drying method to dry slides.</p>	<p>First, wipe off excessive water then air dry slides for a few seconds. Choose one of the four heat drying processes to stabilize chromogens.</p> <p><b>By hybridizer:</b></p> <ol style="list-style-type: none"> <li>a. Set temperature to 85°C.</li> <li>b. Place slides in the hybridizer and heat for 10 minutes.</li> <li>c. Re-rack slides and allow slides to cool for 1 or 2 minutes before going to the <b>step 15. CAUTION: DO NOT</b> touch hot slides with bare hands!</li> </ol> <p><b>By slide warmer:</b></p> <ol style="list-style-type: none"> <li>a. Set temperature to 85°C.</li> <li>b. Place slide in the slide warmer and heat for 10 minutes.</li> <li>c. Re-rack the slides and allow slides to cool for 1 or 2 minutes before going to <b>step 15. CAUTION: DO NOT</b> touch hot slides with bare hands!</li> </ol> <p><b>By vacuum oven:</b></p> <ol style="list-style-type: none"> <li>a. Preheat Vacuum oven to 85°C.</li> <li>b. Place the slides in rack tissue side up in the oven and apply vacuum to &gt;10in of Hg. Heat for 5 minutes.</li> </ol>	<p>10 minutes</p>

	<p>c. Take out slides and allow slides to cool for 1 or 2 minutes before going to <b>step 15. CAUTION: DO NOT</b> touch hot slides with bare hands!</p> <p><b>By hair dryer:</b>  <b>(CAUTION: To prevent electrocution move the slide tray to a dry area in the lab. Hair Dryer and Water do not mix!)</b></p> <p>a. Turn hair dryer on High.  b. Hold hair dryer about 2 inches above slides.  c. Slowly move the Hair dryer back and forth over the slides! <b>DO NOT</b> let the dryer stay over one tissue more than 3 seconds.  d. Blow air over sections until Emerald turns from blue to green. It takes 2 minutes of waving the hair dryer back and forth over the slides.  e. Check every 15 seconds after the first two minutes for green.  f. Re-rack slides and allow them cool for 1 or 2 minutes before going to <b>step 15. CAUTION: DO NOT</b> touch hot slides with bare hands.</p>	
15. Fix slides -20°C methanol supplied by user.	<p>a. Let the hot slides cool down for 30 seconds in a -20°C freezer.  b. Place slides into -20°C methanol for 2 minutes. Do not leave longer! Do not let dry! Move slides directly to 100% Ethanol, see <b>step 16. CAUTION: Do not air dry slides after methanol fixation! It will erase AP Red stain!</b></p>	2 minutes
16. Cover slip slides with permanent Xylene base mountant.	<p>a. 100% Ethanol for <b>10 seconds</b>.  b. 100% Ethanol for <b>10 seconds</b>. Longer may wash away AP Red!  c. Dip <b>3 to 5 seconds</b> in Xylene. Longer may wash away Emerald!  d. Quickly lay slides down.  e. Add Xylene based mountant to slides and apply glass coverslip. Apply force to coverslip to squeeze out any extra mountant and bubble for optimal clarity.</p>	30 seconds

**Trouble shoot:**

<b>Problem</b>	<b>Tips</b>
Uneven stain on 3 primary antibodies	<ol style="list-style-type: none"> <li>1. Need to adjust the titer of each antibody.</li> <li>2. The amount of each protein expressed on tissue may be different.</li> <li>3. Set slides in water too long so that AP-Red is washed away.</li> <li>4. Set slides in Xylene too long so that Emerald is washed away.</li> </ol>
No stain on 1 or 2 antibodies	Missing steps or step reversed.
Emerald did not turn green	Drying process is too short.
Artifacts on slides	Slides are not completely dried before mount.

**Remarks:**

This kit is for research or investigation use only, not for diagnostic or therapeutic use.

## Work Sheet for TS309A kit

We designed work sheet to help you track each step since triple color staining takes many steps. You may use this sheet for our technical support staff to review if needed.

- Used for tester to check “v “each step during the experiment
- Steps follow after de-paraffinization
- Refer to insert for details of each step

<b>Step/ Protocol</b>	<b>Protocol TS309A</b>	<b>Experiment 1 Date:</b>	<b>Experiment 2 Date:</b>	<b>Experiment 3 Date:</b>	<b>Experiment 4 Date:</b>
<b>Step 1</b>	Peroxidase Block				
<b>Step 2</b>	HIER if needed				
<b>Step 3</b>	Ms 1°Ab & Rb 1°Ab mix (30-60 min.)				
<b>Step 4</b>	Ms HRP-Primer (10 min.)				
<b>Step 5</b>	Ms HRP-Polymer & Rb AP- Polymer mix (30 min.)				
<b>Step 6</b>	DAB (5 min.)				
<b>Step 7</b>	AP-Red+ (20 min)				
<b>Step 8</b>	Antibody Blocker (10 min.)				
<b>Step 9</b>	HIER if needed				
<b>Step 10</b>	Rb 1°Ab (30-60 min.)				
<b>Step 11</b>	Rb HRP-Polymer (15 min.)				
<b>Step 12</b>	Counter stain (5 sec.)				
<b>Step 13</b>	Emerald (5-10 min.)				
<b>Step 14</b>	Heat drying (10 min.)				
<b>Step 15</b>	-20 °C MeOH Fix (2 min.)				
<b>Step 16</b>	Mount & coverslip				

Testing result: