CD9/CD9 Exosome ELISA Kit, Human

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

Exosomes are membrane vesicles, which are secreted from many cells and have a diameter of about 30 nm ~ 200 nm. They are present in many bodily fluids, including saliva and blood ^{1) ~ 5)}. It is also known that exosomes are secreted into the mammalian cell culture medium in vitro. Exosomes are enveloped with a lipid bilayer membrane, so like cells, membrane proteins exist at their surface, and they encase various molecular constituents, including proteins and microRNAs. When the target cells receive the exosomes, the proteins or microRNAs inside will fulfill their function, and in this way exosomes are thought to play a role of cell-to-cell communication ⁶⁾. One of the exosomes' structural characteristics is the tetraspanin family which is located on their surface. CD9 and CD63 are member molecules of this family, and known to be exosome surface markers ⁷⁾.

The conventional way to quantify exosomes is indirect. For example, it measures encapsulated protein abundance or performs Nanoparticle Tracking Analysis to determine the size distribution profile of small particles ⁸⁾. The drawback for these methods is that they require exosome purification utilizing ultracentrifugation, for instance. The direct methods to quantitate exosomes in body fluids or culture supernatant is extremely limited, and there were no common methods available until now. The CD9/CD9 Exosome ELISA Kit is a Sandwich ELISA kit, which utilizes high-performance anti-CD9 antibodies. This product detects exosome markers, CD9 molecules that are located on the exosome surface in the body fluids or cell culture supernatant.

【 I – 2 】 Features

- Directly quantitate exosomes in human blood samples or cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Utilize CD9 immobilized beads of 200 nm diameter(CD9 Standard Beads), instead of unstable/hard to store exosome itself, to implement stability and reproducibility.
- Normalization with CD9 standard beads similar to exosomal structures (CD9 Standard Beads) enable to relative quantitate each samples.
- Capture exosomes with solid phase anti-CD9 antibody, then detect using HRP conjugated anti-CD9 antibody.

【 I – 3 】 Kit Principle

An anti-CD9 antibody was immobilized and placed on the plate. First, samples were added onto the plate to capture exosomes by the anti-CD9 antibody. Next, HRP conjugated anti-CD9 antibody will be added to react exosome surface antigen, CD9. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-CD9 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	CD9 Standard Beads	200μL	1tube* ¹
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X) ^{*2}	25mL	1vial
5	HRP Conjugated Anti-CD9 Antibody (500X) ^{*3}	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		2 sheets

*¹ Sufficient to create 4 standard curves with n=2.

*² Crystals may precipitate in the Washing Buffer (10x) during refrigerated storage. Warm the buffer to dissolve it at 45°C before use.

*³ If the kit is not going to be used immediately, remove the labeled antibody from the kit and store it at -20°C.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 µL)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【 II 】 Preparation of Reagents and Samples

【 II – 1 】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to x10 with purified water.
e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【 II – 2 】 Preparation of CD9 Standard Beads solution

	Concentration (ng/mL)	CD9 Standard Beads	Assay Buffer	Dilution factor
A	180			
B	18	50µL of A	450µL	10
C	9	250µL of B	250µL	2
D	4.5	250µL of C	250µL	2
E	2.25	250µL of D	250µL	2
F	1.125	250µL of E	250µL	2
G	0.563	250µL of F	250µL	2

H	0.281	250μL of G	250μL	2
---	-------	------------	-------	---

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of CD9 Standard Beads (Solution A: Vortex Standard Beads for 30 seconds before use.), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 100μL for measurement, using 2 wells for each solution (n=2).
- Diluted CD9 Standard Beads Solution(0.281~18 ng/mL) should be freshly prepared at each time before use.

【II-3】 Preparation of antibody solution

Dilute HRP conjugated anti-CD9 antibody (500x) to 500 folds using Assay Buffer.

e.g. For 1 plate, add 20μL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.

* Diluted antibody solution should be freshly prepared at each time before use.

【II-4】 Preparation of Samples (For cell culture medium supernatant)

The culture medium containing 10% fetal bovine serum (FBS) can be directly used to quantitate exosomes. Grow cells until it reach to the confluent state. Collect the medium supernatant, centrifuge at 2000 xg for 10min, and then used as a sample. If necessary, dilute with Assay Buffer.

【II-5】 Preparation of Samples (For blood sample)

Dilute the serum sample at least four times with Assay Buffer.

* Out of 100μL reaction solution, serum should not be exceeded more than 25μL.

Dilute the plasma sample at least twenty times with Assay Buffer.

* Out of 100μL reaction solution, plasma should not be exceeded more than 5μL.

【II-6】 Sample Storage

Samples should be stored at 2~8°C after preparation until measurement.

【Ⅲ】 Sample measurement procedure

1. Bring anti-CD9 antibody solid phased plate and the reagents to the room temperature.
2. Prepare CD9 Standard Beads solution by serial dilution. (from step 【Ⅱ – 2】)
3. Add 100μL each of serial diluted CD9 Standard Beads solution (0.281 ~ 18ng/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 2hrs for static reaction.
6. Discard all the reaction solution, and then rinse each well with 300μL of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100μL each of diluted HRP conjugated anti-CD9 antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake in the plate shaker.
9. Incubate at room temperature for 2hrs for static reaction.
10. Discard the antibody solution, and then rinse each well with 300μL of Washing Buffer. Repeat this step for 3 times.
11. Add 100μL of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50μL each of Stop Solution.
13. Place into the plate reader, and read the absorbance of each microwell on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve for the CD9 protein by plotting the mean absorbance (y axis) against the CD9 Standard concentration (x axis) (Fig. 1).


Fig.1 Example of Standard curve for CD9 Standard Beads.

【IV】 Example of Results

【IV- 1】 Cell culture medium supernatant

The culture medium of colon cancer cell line, HCT116, was collected, and exosomes were purified by ultracentrifugation. The purified exosomes, 15.6, 31.3, 62.5, 125, 250, and 500 ng/mL each, were added to the well, and measured using this kit (Fig.2).


Fig.2 Exosomes from the culture supernatant of HCT116

【IV- 2】 Standardization of measured values

The Standard curve was created using the measurement of CD9 Standard Beads (Fig. 3a). From here (Fig.3a), assuming 10ng/mL of CD9 Standard Beads as 1 unit/mL, the OD450 for 1unit will be 1.4. For the exosomes purified from HCT116, the protein concentration for 1 unit of OD450, which will correspond to 1.4, is about 200ng/mL (Fig.3b). Therefore, 200 ng/mL of HCT116 derived exosomes can be considered as 1unit/mL of CD9 positive exosomes.

Presenting the exosome measurement by the units, we could standardize it and/or normalize the measurements, and then will be able to compare exosome measurements directly between different samples or different experiments.


Fig.3 Standardization of measured values

[V] Kit expiry date and storage

Expiry date : 6 months after the manufacturing date.

(The manufacturing date is indicated on the kit box label)

Storage : Refrigeration (2-8°C)

[Reference]

- 1) J. Skog, T. Wurdinger, S. van Rijn, D. H. Meijer, L. Gainche, M. Sena-Esteves, W. T. Jr. Curry, B. S. Carter, A. M. Krichevsky and X. O. Breakefield: *Nat Cell Biol.*, **10**, 1470 (2008).
- 2) T. Pisitkun, R. F. Shen and M. A. Knepper: *Proc Natl Acad Sci USA.*, **101**, 13368 (2004).
- 3) S. Runz, S. Keller, C. Rupp, A. Stoeck, Y. Issa, D. Koensgen, A. Mustea, J. Sehouli, G. Kristiansen and P. Altevogt: *Gynecol Oncol.*, **107**, 563 (2007).
- 4) S. Keller, A. K. Konig, F. Marme, S. Runz, S. Wolterink, D. Koensgen, A. Mustea, J. Sehouli and P. Altevogt: *Cancer Lett.*, **278**, 73 (2009).
- 5) C. Lasser, V. Seyed AlikhS. Gabrielsson, J. Lotvall and H. Valadi: *J Transl Med.*, **9**, 9 (2011).
- 6) Y. Naito, Y. Yoshioka, Y. Yamamoto and T. Ochiya: *Cell Mol Life Sci.*, **74**, 697 (2017).

- 7) A. Zoraida and M. Yáñez-Mó: *Front Immunol.*, **5**, 442 (2014).
- 8) V. Filipe, A. Hawe and W. Jiskoot: *Pharm Res.*, **27**, 796 (2010).

COSMO BIO CO., LTD.
Inspiration for Life Science

Licensed by Cosmo Bio Co., LTD.