PD-L1/CD9 Exosome ELISA Kit, Human

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

Tumor cells evade the immune surveillance by up-regulating surface expression of PD-L1, which interacts with PD-1 on T cells to elicit the immune checkpoint response^{1),2)}. In recent years, anti-PD-1 antibody drugs and anti-PD-L1 antibody drugs that inhibit this signal transduction have been confirmed to have clinical efficacy and safety in various solid cancers by enhancing T cell activation^{3), 4)}.

A well-known predictor of the effects of these immune checkpoint inhibitors is PD-L1 expression on cancer cells. The PD-L1 IHC test is used as a companion diagnostics for measuring PD-L1 expression rate in cancer tissues and cells.

As a general matter, because IHC tests are highly invasive, research on liquid biopsies using exosomes is also in progress.

Exosomal PD-L1 has the same membrane topology as the cell surface PD-L1, with its extracellular domain exposed on the surface of the exosomes⁵⁾. Exosomal PD-L1 binds PD-1 in a concentration-dependent manner. The level of PD-L1 on the circulating exosomes was significantly higher in patients with metastatic melanoma than that in healthy donors⁵⁾.

However, these methods require exosome purification using ultracentrifugation, there were no rapid and convenient direct measurement methods until now.

This product is a two-step sandwich ELISA, which utilizes high-performance anti-PD-L1 and anti-CD9 antibodies to directly detect PD-L1 molecules on the surface of exosomes in human blood samples and cell culture supernatants.

【 I – 2 】 Features

- Directly quantitate exosomal PD-L1 in human blood samples or cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Utilize PD-L1/CD9 fusion protein (Standard Protein), instead of unstable/hard to store exosome itself, to implement stability and reproducibility.
- Normalization with a standard curve using PD-L1/CD9 fusion protein (Standard Protein) enable to relative quantitate each sample.
- Detect human PD-L1 positive exosomes by two-step sandwich method using immobilized anti-CD9 antibody and HRP conjugated anti-PD-L1 antibody.

【 I – 3 】 Kit Principle

This ELISA kit uses two-step sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-CD9 antibody.

First, samples were added onto the plate to capture exosomes by the anti-CD9 antibody. Next, HRP conjugated anti-PD-L1 antibody will be added to react with PD-L1 on the surface of exosomes. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-CD9 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	Standard Protein (PD-L1/CD9 Fusion Protein)	100μL	1tube*
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X)	25mL	1vial
5	HRP Conjugated Anti-PD-L1 Antibody (500X)	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		2sheet

* Sufficient to create 2 standard curves with n=2.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μ L)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【Ⅱ】 Preparation of Reagents and Samples

【Ⅱ – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10 \times) to 10 folds with purified water.
e.g. For 1 plate, add 225mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【Ⅱ – 2】 Preparation of Standard Protein solution

	Concentration (ng/mL)	Standard Protein	Assay Buffer	Dilution factor
A	100			
B	10	50 μ L of A	450 μ L	10
C	5	250 μ L of B	250 μ L	2
D	2.5	250 μ L of C	250 μ L	2
E	1.25	250 μ L of D	250 μ L	2
F	0.625	250 μ L of E	250 μ L	2
G	0.313	250 μ L of F	250 μ L	2
H	0.156	250 μ L of G	250 μ L	2

- To prepare Solution B, add 450 μ L of Assay Buffer into 50 μ L of Standard Protein (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250 μ L of Assay Buffer into 250 μ L of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.

- Use 200 μ L for measurement, using 2 wells for each solution (n=2).
- Diluted Standard Protein Solution(0.156~10 ng/mL) should be freshly prepared at each time before use.

【Ⅱ-3】 Preparation of antibody solution

- Dilute HRP conjugated anti-PD-L1 antibody (500x) to 500 folds using Assay Buffer.
e.g. For 1 plate, add 20 μ L of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.
- * Diluted antibody solution should be freshly prepared at each time before use.

【Ⅱ-4】 Preparation of Samples (For Blood samples)

Serum or plasma is measured as a sample diluted 10-fold with Assay Buffer.
Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【Ⅱ-5】 Preparation of Samples (For cell culture medium supernatant)

Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris.
Collect supernatants and assay.
Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【Ⅱ-6】 Sample Storage

Store un-diluted sample at -20°C or below. Avoid repeated freeze-thaw cycles.

【Ⅲ】 Sample measurement procedure


1. Bring anti-CD9 antibody solid phased plate and the reagents to the room temperature.
2. Prepare Standard Protein solution by serial dilution. (from step 【Ⅱ-2】)
3. Add 100 μ L each of serial diluted Standard Protein solution (0.156 ~ 10 ng/mL) or Sample solution into the well.

4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.
6. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100 μ L each of diluted HRP conjugated anti-PD-L1 antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake it in the plate shaker at 800 rpm for 30 sec.
9. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.
10. Discard the reaction solution, and then rinse each well with 300 μ L of Washing Buffer(from step 【II-1】). Repeat this step for 3 times.
11. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the average blank control subtracted absorbance value for each Standard Protein concentration (y axis) against the Standard Protein concentration (x axis).
15. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Determine the PD-L1-positive exosome concentration (Unit/mL) in the sample solution as the PD-L1/CD9 fusion protein (ng/mL) equivalent. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of PD-L1 positive exosome in the sample.

【IV】 Measurement example

【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against standard concentration is drawn as shown in Figure 1. However, draw a new standard curve for each assay to calculate the concentration in the sample.

(Plotted is the value obtained by subtracting the blank absorbance from the absorbance of each standard protein concentration)

Standard Protein (PD-L1/CD9 Fusion Protein) (ng/mL)	Absorbance (450nm)		Average	Calculated concentration (U/mL)
	1	2		
0	0.085	0.090	0.088	
0.156	0.131	0.131	0.131	
0.313	0.178	0.182	0.180	
0.625	0.275	0.294	0.285	
1.25	0.472	0.468	0.470	
2.5	0.839	0.834	0.837	
5	1.580	1.598	1.589	
10	2.871	2.928	2.900	
Exosome of PC9 cells	0.340	-	0.340	0.84

Fig.1 Standard curve and measured values

【Reference】

- 1) H. Dong, S. E. Strome, D. R. Salomao, H. Tamura, F. Hirano, D. B. Flies, P. C. Roche, J. Lu, G. Zhu, K. Tamada, V. A. Lennon, E. Celis, L. Chen: *Nat Med.*, **8**, 793 (2002).
- 2) L. Chen, X. Han: *J Clin Invest.*, **125**, 3384(2015).
- 3) S. L. Topalian, F. S. Hodi, J. R. Brahmer, et al.: *N Engl J Med.*, **366**, 2443 (2012).
- 4) J. R. Brahmer, S. S. Tykodi, L. Q. Chow, et al.: *N Engl J Med.*, **366**, 2455 (2012).
- 5) G. Chen, A. C. Huang, W. Zhang, et al.: *Nature*, **560**, 382 (2018).