使用説明書
Instruction manual

Equol ELISA Kit

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1. Niwa T., Yokoyama S., Osawa T. (2009). Preparation of soy isoflavonoids for the production of anti-equol monoclonal antibody, Phytochemistry Letters 2(4), 220–222.
2. Taylor JL., Grace PB., Bingham SA. (2005). Optimization of conditions for the enzymatic hydrolysis of phytoestrogen conjugates in urine and plasma, Analytical Biochemistry 341, 220–229.

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TABLE OF CONTENTS

1. INTENDED USE.....10
2. TEST PRINCIPLE.....10
3. SAFETY PRECAUTIONS.....11
4. STORAGE AND STABILITY.....11
5. MATERIALS SUPPLIED.....11
6. MATERIALS REQUIRED BUT NOT SUPPLIED11
7. SAMPLE COLLECTION AND STRAGE.....12
8. REAGENT PREPARATION.....12
9. TEST PROCEDURE.....13
- 1 0. CALCULATION OF RESULTS.....13
- 1 1. PERFORMANCE CHARACTERISTICS.....14
- 1 2. PRODUCT LITERATURE REFERENCES15

1. INTENDED USE

Equol ELISA Kit is an immunoassay specifically designed and validated for quantitative measurement of biological products (urine and other body fluid) or food samples. It is intended only for research. It is not for diagnostic use or treatment.

2. TEST PRINCIPLE

This product is based on the competitive ELISA using monoclonal antibody specific for equol¹⁾ as the basic principle.

1. The microplate provided in this kit has been pre-coated with equol.

● ...Solid phased equol


2. Add Equol Detection Antibody (HRP-labeled anti-equol antibody) and samples. HRP-labeled anti-equol antibody competitively binds to solid phased equol and specimen derived equol.

● ...HRP-labeled anti-equol antibody ● ...Specimen derived equol

3. After incubation, HRP-labeled anti-equol antibody bound to specimen derived equol is washed off.

4. Adding a coloring reagent (TMB Substrate Solution), the reaction solution change blue. Add a Stop Solution to develop a yellow color and measure the absorbance.

☀ ...Coloring reagent

3. SAFETY PRECAUTIONS

1. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit.
2. Please use this product within the expiration date. Even within the expiration date, measurement may not be possible with sufficient accuracy when stored at room temperature or frozen.
3. In an assay, please protect your body using white coat, protective eyeglasses, latex gloves etc. as necessary.

4. STORAGE AND STABILITY

The kit should be stored at 2-8 °C. Keep away from direct sunlight and freezing. After measurement, please seal the surplus strip well and refrigerate it in a zip bag. Please keep other reagents refrigerated as well. Please use up stripwell, reagents within 1 week after the kit is opened.

5. MATERIALS SUPPLIED

	Component	Quantity/Size
1	Pre-Coated Microplate (12×8 well strips)	1plate/ 96wells
2	Equol Standard (2500 μg/mL)	1vial / 50 μL
3	Equol Detection Antibody	1vial / 15 μL
4	Reagent Diluent	1 bottle / 50 mL
5	Wash Buffer Concentrate (10×)	1 bottle / 50 mL
6	TMB Substrate Solution	1 bottle / 12 mL
7	Stop Solution	1 bottle / 6 mL
8	Plate Seal	2 sheets

6. MATERIALS REQUIRED BUT NOT SUPPLIED

1. Micropipette (1-20 μL) and tips
2. Micropipette (50-200 μL) and tips
3. 8-Channel Micropipette (50-200 μL) and tips
4. Microplate reader (measurement wavelength 450 nm, reference wavelength 570 - 620 nm)

7. SAMPLE COLLECTION AND STRAGE

Storing the collected specimens, please store under freezing condition at -20 °C or less

1. Urine ... Urine does not need to be diluted before assaying. If the measured value is close to the upper limit of the calibration curve, the quantitiveness will decrease. We recommend remeasuring with increasing dilution factor. When quantifying the total equol amount in urine, please perform deconjugation referring to the literature ²⁾.


Example: 50 μ L of urine sample and 50 μ L of β -glucuronidase (derived from *Helix pomatia*, 125 U / 100 μ L 0.1 M acetate buffer (pH5.0)) are mixed and treated at 37 °C for 1 hour.

2. Serum ... The concentration of equol in serum is lower than urine, so the measured value may be below the limit of quantification. If necessary, concentrate the specimen with reference to literatures.

3. Foods ... The extraction method will change depending on the type of food. Please extract according to literatures.

8. REAGENT PREPARATION

1. Bring all reagents and samples to room temperature (18–25 °C) before use.
2. Briefly spin or centrifuge the Equol Detection Antibody to collect the inner solution to the bottom of the tube. Dilute the Equol Detection Antibody 400-fold with Reagent Diluent
3. Dilute the whole volume of Wash Buffer Concentrate (10 \times) with 450 mL of distilled water.
4. Mix 10 μ L Equol Standard (2500 μ g / mL) with 990 μ L Reagent Diluent (Std. 1). Mix 400 μ L of Std.1 and 600 μ L of Reagent Diluent, and make it Std.2. Repeat the same operation, please prepare up to Std.7. Std.8 uses Reagent Diluent.

9. TEST PROCEDURE

1. Prepare all reagents, standards, and samples as directed in "7. SAMPLE COLLECTION AND STRAGE" and "8. REAGENT PREPARATION".
2. Dilute samples and serial diluted equol standards, divide each 50 μ L of solution into wells with reference to the layout table below.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std.1	Std.1	Sam-1	Sam-1								
B	Std.2	Std.2	Sam-2	Sam-2								
C	Std.3	Std.3	Sam-3	Sam-3								
D	Std.4	Std.4	Sam-4	Sam-4								
E	Std.5	Std.5	Sam-5	Sam-5								
F	Std.6	Std.6	Sam-6	Sam-6								
G	Std.7	Std.7	Sam-7	Sam-7								
H	Std.8	Std.8	Sam-8	Sam-8								

3. Add 50 μ L of diluted Equol Detection Antibody into each wells.
4. Shake the microplate gently, seal the plate and incubate for 1 hour at room temperature.
5. Aspirate the solution of each wells and wash each wells four times with 300 μ L of diluted Wash Buffer. After the last wash, remove the remaining fluid as much as possible (Invert the plate and blot it against clean paper towels to remove excess liquid).
6. Add 100 μ L of TMB Substrate Solution into all wells, seal the plate and incubate for 15 minutes at room temperature in the dark.
7. Add 50 μ L of Stop Solution to all wells, shake and mix the microplate gently, then measure the absorbance with a microplate reader (measurement wavelength 450 nm, reference wavelength 570 - 620 nm)

10. CALCULATION OF RESULTS

The obtained OD of the standards (y-axis, linear) are plotted against their concentration (x-axis, logarithmic).

We recommend using four-parameter logistic regression models.


Conversion of units

$$\text{Equol } (\mu\text{g} / \text{mL}) \times 4.13 = \mu\text{mol} / \text{L}$$

Typical calibration curve

Example. Do not use for calculation.

Equol (μ g/mL)	OD (450nm/570nm)
25.00	0.187
10.00	0.264
4.00	0.377
1.60	0.610
0.64	1.051
0.26	1.670
0.10	2.305
0.00	2.979

11. PERFORMANCE CHARACTERISTICS


Sensitivity

Detection Limit	0.037 $\mu\text{g/mL}$ (0.151 μM)	RSD>30% RSD: relative standard deviation
Quantitative Limit	0.173 $\mu\text{g/mL}$ (0.716 μM)	RSD>10% RSD: relative standard deviation

Precision

	Mean($\mu\text{g/mL}$)	SD($\mu\text{g/mL}$)	CV(%)	N
Intra-Assay	10.0	1.08	10.8	6
	19.3	1.59	8.2	6
	50.3	3.06	6.1	6
Inter Assay	10.4	1.06	10.2	6
	17.8	1.80	10.1	6
	39.6	5.11	12.9	6

Linearity

Recovery

Urine Sample	Endogenous ($\mu\text{g/mL}$)	Added ($\mu\text{g/mL}$)	Expected ($\mu\text{g/mL}$)	Observed ($\mu\text{g/mL}$)	Recovery (%)
1	9.2	23.2	32.3	31.6	97.8
1	9.2	8.2	17.4	17.3	99.5
1	9.2	3.1	12.3	12.3	100.0
2	27.2	23.2	50.3	48.3	96.0
2	27.2	8.2	35.4	38.5	108.9
2	27.2	3.1	30.3	30.8	101.9

Method Comparison

Regression equation	r	N
ELISA = $0.934 \times \text{HPLC} + 1.055$	0.968	40

Antibody Specificity

Compound	Concentration (μM)	Cross Reactivity (%)
S-equol	10	100.0
R-equol	10	88.2
Daidzein	10	ND
Genistein	10	0.1
Glycitein	10	1.6
Genistin	10	ND

ND : Not Detected

Equol Example Ranges

Urine	Mean($\mu\text{g/mL}$)	SD	N
	64.5	60.8	40

12. PRODUCT LITERATURE REFERENCES

1. Niwa T., Yokoyama S., Osawa T. (2009). Preparation of soy isoflavonoids for the production of anti-equol monoclonal antibody, *Phytochemistry Letters* 2(4), 220–222.
2. Taylor JL, Grace PB., Bingham SA. (2005). Optimization of conditions for the enzymatic hydrolysis of phytoestrogen conjugates in urine and plasma, *Analytical Biochemistry* 341, 220–229.

Revision date: November 12, 2025

Change details:

- Update to "7. SAMPLE COLLECTION AND STRAGE"
- In the table "Recovery" under the section "PERFORMANCE CHARACTERISTICS," the unit for the item "Observed" has been changed as follows:
(pg/mL) → ($\mu\text{g/mL}$)