使用説明書  
Instruction manual

# Equol ELISA Kit

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1. Niwa T., Yokoyama S., Osawa T. (2009). Preparation of soy isoflavonoids for the production of anti-equol monoclonal antibody, *Phytochemistry Letters* 2(4), 220–222.
2. Taylor JL, Grace PB, Bingham SA. (2005). Optimization of conditions for the enzymatic hydrolysis of phytoestrogen conjugates in urine and plasma, *Analytical Biochemistry* 341, 220–229.

## TABLE OF CONTENTS

1. INTENDED USE.....10
2. PRINCIPLE OF THE ASSAY.....10
3. SAFETY PRECAUTIONS.....11
4. STORAGE CONDITIONS AND STABILITY.....11
5. MATERIALS SUPPLIED.....11
6. MATERIALS REQUIRED BUT NOT SUPPLIED .....11
7. SAMPLE COLLECTION AND STORAGE.....12
8. REAGENT PREPARATION.....12
9. ASSAY PROCEDURE.....13
10. CALCULATION OF RESULTS.....14
11. PERFORMANCE CHARACTERISTICS.....14
12. REFERENCES .....15

## 1. INTENDED USE

The Equol ELISA Kit is an immunoassay specifically designed and validated for the quantitative measurement of equol in biological samples, including urine and other body fluids, as well as in food samples.

This kit is intended for research use only and is not intended for diagnostic or therapeutic use.

## 2. PRINCIPLE OF THE ASSAY

This product is based on a competitive ELISA employing a monoclonal antibody specific for equol as its underlying principle<sup>1)</sup>

1. The microplate provided in this kit is pre-coated with equol.

 ...Solid phased equol


2. Samples and the Equol Detection Antibody (HRP-labeled anti-equol monoclonal antibody) are added to the wells. The HRP-labeled anti-equol antibody competitively binds to the equol immobilized on the solid phase and to the equol present in the samples.

 ...HRP-labeled anti-equol antibody       ...Equol present in samples

3. After incubation, unbound HRP-labeled anti-equol antibody, including antibody bound to sample-derived equol, is removed by washing.

4. The TMB Substrate Solution is added, resulting in the development of a blue color. The reaction is then stopped by adding Stop Solution, producing a yellow color, and the absorbance is measured.

 ...Coloring reagent

## 3. SAFETY PRECAUTIONS

1. Before starting the assay, read these instructions carefully and in their entirety. Use only the current, valid version of the package insert provided with the kit.
2. Use this product before the expiration date. Even if used within the expiration date, accurate measurement may not be ensured if the product has been stored at room temperature or frozen.
3. During the assay, wear appropriate personal protective equipment (PPE), such as a laboratory coat, protective eyewear, and latex gloves, as required.

## 4. STORAGE CONDITIONS AND STABILITY

Store the kit at 2–8°C. Protect the kit from direct sunlight and do not freeze.

After use, reseal any unused microplate strips and store them refrigerated in a sealed zip bag. Store all other reagents refrigerated as well. Once the kit has been opened, use all microplate strips and reagents within one week.

## 5. MATERIALS SUPPLIED

	Component	Quantity/Size
1	Pre-Coated Microplate (12×8 well strips)	1plate/ 96wells
2	Equol Standard (2500 μg/mL)	1vial / 50 μL
3	Equol Detection Antibody	1vial / 18 μL
4	Reagent Diluent	1 bottle / 50 mL
5	Wash Buffer Concentrate (10×)	1 bottle / 50 mL
6	TMB Substrate Solution	1 bottle / 12 mL
7	Stop Solution 	1 bottle / 6 mL
8	Plate Seal	2 sheets

## 6. MATERIALS REQUIRED BUT NOT SUPPLIED

1. Micropipette (1–20 μL) with compatible tips
2. Micropipette (50–200 μL) with compatible tips
3. 8-channel micropipette (50–200 μL) with compatible tips
4. Microplate reader capable of measuring absorbance at 450 nm, with a reference wavelength of 570–620 nm

## 7. SAMPLE COLLECTION AND STORAGE

Collected specimens should be stored frozen at  $-20^{\circ}\text{C}$  or below until analysis.


1. Urine•••Urine samples do not require dilution prior to the assay. However, if the measured value is close to the upper limit of the calibration curve, quantitative accuracy may decrease. In such cases, remeasurement with an increased dilution factor is recommended. For quantification of total equol in urine, perform deconjugation according to the literature <sup>2)</sup>.

Example: Mix 50  $\mu\text{L}$  of urine sample with 50  $\mu\text{L}$  of  $\beta$ -glucuronidase solution (derived from *Helix pomatia*, 125 U/100  $\mu\text{L}$  in 0.1 M acetate buffer, pH 5.0), and incubate at  $37^{\circ}\text{C}$  for 1 hour.

2. Serum•••The concentration of equol in serum is lower than that in urine, and measured values may be below the limit of quantification. If necessary, concentrate the sample according to the procedures described in the relevant literature.
3. Food Samples•••Extraction methods vary depending on the type of food. Prepare extracts according to the relevant literature.

## 8. REAGENT PREPARATION

1. Bring all reagents and samples to room temperature ( $18$ – $25^{\circ}\text{C}$ ) before use. <sup>※</sup>
2. Briefly centrifuge the Equol Detection Antibody to collect the contents at the bottom of the tube. Dilute the Equol Detection Antibody 1:400 with Reagent Diluent.
3. Dilute the entire volume of the  $10\times$  Wash Buffer Concentrate with 450 mL of distilled water.
4. Prepare the equol standards as follows:  
Mix 10  $\mu\text{L}$  of Equol Standard (2,500  $\mu\text{g}/\text{mL}$ ) with 990  $\mu\text{L}$  of Reagent Diluent to prepare Standard 1 (Std. 1).  
Mix 400  $\mu\text{L}$  of Std. 1 with 600  $\mu\text{L}$  of Reagent Diluent to prepare Standard 2 (Std. 2).  
Repeat this procedure to prepare Standards 3 through 7 (Std. 3–Std. 7).  
Use Reagent Diluent alone as Standard 8 (Std. 8).

※[Equal Standard] Handling precautions

Equol Standard becomes solid below  $19^{\circ}\text{C}$ .

After taking out from the refrigerator, please keep at room temperature for 1 hour or more.

When the room temperature is low, use a  $25^{\circ}\text{C}$  incubator to dissolve it.

Before use, make sure the solution is completely dissolved.

## 9. ASSAY PROCEDURE

1. Prepare all reagents, standards, and samples according to Section 7, "SAMPLE COLLECTION AND STORAGE," and Section 8, "REAGENT PREPARATION."
2. Dispense 50  $\mu\text{L}$  of each diluted sample and serially diluted equol standard into the appropriate wells, following the plate layout shown below.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std.1	Std.1	Sam-1	Sam-1								
B	Std.2	Std.2	Sam-2	Sam-2								
C	Std.3	Std.3	Sam-3	Sam-3								
D	Std.4	Std.4	Sam-4	Sam-4								
E	Std.5	Std.5	Sam-5	Sam-5								
F	Std.6	Std.6	Sam-6	Sam-6								
G	Std.7	Std.7	Sam-7	Sam-7								
H	Std.8	Std.8	Sam-8	Sam-8								


3. Add 50  $\mu\text{L}$  of diluted Equol Detection Antibody to each well.
4. Gently shake the microplate, seal it, and incubate for 1 hour at room temperature.
5. Aspirate the contents of each well and wash the wells four times with 300  $\mu\text{L}$  of diluted Wash Buffer. After the final wash, remove residual liquid as completely as possible by inverting the plate and blotting it against clean paper towels.
6. Add 100  $\mu\text{L}$  of TMB Substrate Solution to each well, seal the plate, and incubate for 15 minutes at room temperature in the dark.
7. Add 50  $\mu\text{L}$  of Stop Solution to each well. Gently shake the microplate to mix, then measure the absorbance using a microplate reader at 450 nm, with a reference wavelength of 570–620 nm.

## 10. CALCULATION OF RESULTS

Plot the optical density (OD) values of the standards on the y-axis (linear scale) against their concentrations on the x-axis (logarithmic scale).

We recommend fitting the data using a four-parameter logistic (4PL) regression model.

Equlol ( $\mu\text{g/mL}$ )	OD (450nm/570nm)
25.00	0.187
10.00	0.264
4.00	0.377
1.60	0.610
0.64	1.051
0.26	1.670
0.10	2.305
0.00	2.979

Unit Conversion

$$\text{Equlol } (\mu\text{g/mL}) \times 4.13 = \mu\text{mol/L}$$

The calibration curve shown is a typical example and should not be used for calculation.

## 11. PERFORMANCE CHARACTERISTICS


Sensitivity

Detection Limit	0.037 $\mu\text{g/mL}$ (0.151 $\mu\text{M}$ )	RSD>30% RSD: relative standard deviation
Quantitative Limit	0.173 $\mu\text{g/mL}$ (0.716 $\mu\text{M}$ )	RSD>10% RSD: relative standard deviation

Precision

	Mean ( $\mu\text{g/mL}$ )	SD ( $\mu\text{g/mL}$ )	CV (%)	N
Intra-Assay	10.0	1.08	10.8	6
	19.3	1.59	8.2	6
	50.3	3.06	6.1	6
Inter Assay	10.4	1.06	10.2	6
	17.8	1.80	10.1	6
	39.6	5.11	12.9	6

Linearity

Recovery

Urine Sample	Endogenous ( $\mu\text{g/mL}$ )	Added ( $\mu\text{g/mL}$ )	Expected ( $\mu\text{g/mL}$ )	Observed ( $\mu\text{g/mL}$ )	Recovery (%)
1	9.2	23.2	32.3	31.6	97.8
1	9.2	8.2	17.4	17.3	99.5
1	9.2	3.1	12.3	12.3	100.0
2	27.2	23.2	50.3	48.3	96.0
2	27.2	8.2	35.4	38.5	108.9
2	27.2	3.1	30.3	30.8	101.9

Method Comparison

Regression equation	r	N
ELISA = 0.934 × HPLC + 1.055	0.968	40

Antibody Specificity

Compound	Concentration ( $\mu\text{M}$ )	Cross Reactivity (%)
S-equal	10	100.0
R-equal	10	88.2
Daidzein	10	ND
Genistein	10	0.1
Glycitein	10	1.6
Genistin	10	ND

ND : Not Detected

## 12. REFERENCES

- Niwa T., Yokoyama S., Osawa T. (2009). Preparation of soy isoflavonoids for the production of anti-equal monoclonal antibody, *Phytochemistry Letters* 2(4), 220–222.
- Taylor JI., Grace PB., Bingham SA. (2005). Optimization of conditions for the enzymatic hydrolysis of phytoestrogen conjugates in urine and plasma, *Analytical Biochemistry* 341, 220–229.