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#### Technical Information Sheet

#### INTENDED USE

The new **R**apid **S**tain **I**dentification (RSID™) test for semen is designed for fast, easy, and reliable detection of human semen from a variety of samples encountered by forensic laboratories including clothing, bedding, vaginal swabs, prophylactics, and stained surfaces.

The test will detect as little as 1 µl of human semen, and test results are complete within 10 minutes.

The detection protocol can be completely integrated into standard forensic laboratory procedures for DNA analysis, prior to STR analysis (see Provided Protocols).

The test sensitivity has been adjusted so that when semen is detected, sufficient biological material should be present to generate an STR profile (exceptions include low sperm count semen or semen from vasectomized men). Typical extraction protocols are provided for your convenience (see Provided Protocols).

This is the first commercially available confirmatory test for human semen. No other human body fluids tested cross react with this procedure (See Specificity below for fluids tested). The immunochromatographic strip test uses dual monoclonal antibodies specific for human semenogelin.

#### Introduction

Rapid Stain Identification of Human Semen (RSID™-Semen) is a lateral flow immunochromatographic strip test designed to detect the presence of human semenogelin. Semenogelin is a protein produced by the seminal vesicles that causes seminal fluid to form a coagulum subsequent to ejaculation. RSID™-Semen uses two anti-human semenogelin monoclonal antibodies in a lateral flow format, which detects the presence of semenogelin.

RSID™-Semen is specific for human semen and has numerous advantages over other methods for semen detection, including increased sensitivity, specificity, and speed. Current identification methods for semen are presumptive (provide a basis for continued analysis of the tested exhibit but are not specific for semen), and are therefore open to legal and scientific challenge.

#### Principle of the Test

RSID™-Semen is an immunochromatographic assay that uses two mouse monoclonal antibodies specific for human semenogelin. One of these antibodies is conjugated to colloidal gold and is deposited on a conjugate pad beneath the sample window. The other antibody is striped onto the "Test line" on a membrane attached to the conjugate pad. The "Control line" on the membrane consists of anti-mouse IgG antibody and is used as a test control.

Attached to the other end of the membrane is the wick, which absorbs the tested fluid and running buffer at the completion of the test thus preventing back-flow of the sample. Once the tested fluid is added to the sample window, the running buffer and sample diffuse through the conjugate pad, re-dissolving the gold-conjugated antibodies. If human semenogelin is present in the sample, an antigen-antibody conjugated to colloidal gold complex will form. Sample and antibodies (complexed and free) are transported by bulk fluid flow to the membrane phase of the strip test. The immobilized anti-semenogelin antibodies on the test line capture the semenogelin-antibody-gold complexes, producing a red line at the Test position. If no human semenogelin is present in the sample, then gold-conjugated antibody-antigen complexes do not form, and colloidal gold will not be accumulated at the Test line. The anti-mouse IgG on the control line captures any mouse antibodies flowing past the test line, producing a red line at the Control position. This demonstrates that the sample fluid was transported through the length of the test, and that the components of the strip test are working correctly.

#### Reagents and Materials Provided

- i) Test cassettes: 25 cassettes, individually wrapped and sealed in moisture-proof foil (a silica gel desiccant pouch has been added for increased shelf life.)
- ii) 5 ml of RSID™-Semen Running Buffer
- iii) 25 ml of RSID™-Semen Extraction Buffer (formulated for use with RSID™-Semen, certified STR free)
- iv) Suggested laboratory and validation protocols, optional validation summary form.

#### Protocol for Positive Control

Positive controls for RSID™-Semen can be produced from 50 µl of human semen deposited on a sterile cotton swab. The semen swab should be extracted in 1 ml of RSID™-Semen Extraction Buffer for 1-2 hours at room temperature; 5 µl of this extract should be diluted in 95 µl of RSID™-Semen Running Buffer (total volume 100 µl). Load all 100 µl into the sample well; this will give a strong positive signal.

#### Suggested Extraction Protocol for Sample Analysis

Forensic samples obtained on cotton swabs (e.g., vaginal swabs, prophylactics, and stained solid surfaces) should be extracted in 200-300 µl of RSID™-Semen Extraction Buffer for 1-2 hours. A 20 µl volume of this extract is added to 80 µl of RSID™-Semen Running Buffer for a total final volume of 100 µl. Stains on fabric or paper should be sampled by taking a punch or cutting (≈ 20 mm<sup>2</sup>) of the item. The punch or cutting should be extracted in 100 µl of RSID™-Semen Extraction Buffer for 1-2 hours. A 10 µl volume of this extract is then added to 90 µl of RSID™-Semen Running Buffer for a total final volume of 100 µl which is then loaded into the sample well on the cassette. The remainder of the extract can then be processed for STR analysis using one of the recommended protocols (see Provided Protocols).

## Strip Test Assay Procedure


Note: Assays should be performed at room temperature, we recommend that a positive and negative control be included with every assay (see Provided Protocols).

1. Remove cassette from the foil pouch. Discard silica gel desiccant.
2. Bring sample to a final volume of 100 µl with supplied RSID™-Semen Running Buffer, typically in a disposable 0.6 ml microfuge tube.
3. Set a timer for 10 minutes.
4. Add sample in running buffer to sample window. Start timer.
5. At 10 minutes, score and record results as shown in the diagram.
6. Due to the high dose Hook effect, samples giving a weak positive or negative result should be diluted 1:20 and re-tested. For example: If 20 µl from a 200 µl swab extract gives a weak positive or negative result, 1 µl from the original extract should be added to 99 µl RSID™-Semen Running Buffer and analyzed on a new cassette (see High Dose Hook Effect below for details).
7. Discard cassette(s). Each cassette can only be used once.

## Scoring Results

RSID™-Semen should be evaluated exactly 10 minutes after the addition of sample. Fig. 1 illustrates expected results:

- i) A visible red line at the Control (C) position only indicates a negative result.  
**No human semen detected.**
- ii) Visible red lines at both the Control (C) and Test (T) positions indicate a positive result.  
**Human semen detected.**
- iii) A visible red line at the Test (T) position only indicates a failed test.  
**Test failure, no conclusion possible.**

## Stability and Storage

RSID™-Semen cassettes should be stored at room temperature. RSID™-Semen Extraction and Running Buffers should be stored at 4°C. Do not use buffers or cassettes after the printed expiration date.

## Specificity

RSID™-Semen is specific for human semenogelin. No cross-reactivity with human saliva, whole blood, vaginal fluid, menstrual blood, breast milk or urine has been observed.

No cross reactivity with animal semen has been observed. Species tested: dog, cat, mouse, cow, horse, pig, goat, and sheep.

## Test Sensitivity

The detection limit for RSID™-Semen, used as suggested is 1 µl of human semen.

Undiluted semen should not be used with RSID™-Semen, as the viscosity of the sample prevents proper release of the conjugate from the conjugate pad. The tested sample should first be deposited on a sterile cotton swab, extracted in RSID™-Semen Extraction Buffer, and diluted as needed in RSID™-Semen Running Buffer before analysis with RSID™-Semen.

## High Dose Hook Effect

A high dose Hook effect refers to weak positive or false negative results seen with immunochromatographic strip tests when very high levels of target are present in the tested sample. Under these conditions, unbound target antigen can reach the test line *before* the colloidal gold-labeled antibody-bound antigen, resulting in a weak positive or false negative result.

We have observed weak positive and false negative results with RSID™-Semen when samples containing large amounts of human semen (≈ 3 to 50 µl) were analyzed. 20-fold dilution of these samples and re-testing with RSID™-Semen eliminated the weak positive and false negative results in all cases (see Validation Report for details).

**User Note:** Under standard laboratory testing, users of RSID™-Semen may observe weak positive or false negative results due to the high dose Hook effect. Therefore, any weak positive or negative result from RSID™-Semen should be confirmed by diluting the sample 1:20 and re-testing. If re-testing of the diluted sample results in a stronger positive signal, the original result was caused by the high dose Hook effect and a large amount of semen is present in the sample. If re-testing of the diluted sample is once again weakly positive or negative, the original result is confirmed.

## Validation

A typical validation summary page and validation test description are supplied. Laboratories should adapt or edit these to suit their own requirements for technical validation.

*Not for in vitro diagnostic use*

*Manufactured by:*

**Independent Forensics**  
4600 ROOSEVELT RD., STE. 201 - HILLSDALE, IL 60162  
TEL. 866-434-2400 FAX 708-978-5115  
WWW.IFI-TEST.COM/RSID.HTML

**Independent Forensics**  
**Rapid Stain Identification**  
**Of Human Semen (RSID™-Semen)**  
**Provided Protocols**

**GENERAL GUIDELINES**

When possible, stains deposited on fabric or other substrates that can be easily cut, should be dissected to preserve at least half of the questioned stain in order to allow for re-testing. Scissors or scalpels used to cut the underlying substrate of questioned stains should be washed in 95% ethanol and deionized water, and then dried with a fresh Kimwipe before use and between each cutting. Special emphasis should be paid to the hinge region of scissors. We recommend immersing scissors and scalpels sequentially in 95% ethanol and deionized water between uses. Use a clean new cutting surface for each sample.

Stains deposited on substrates that cannot be cut (e.g., glass, metal) may be sampled with a clean new swab moistened with sterile deionized water. Use a 'sponge' technique to transfer the stain to the moistened swab; medium pressure may be required on smooth surfaces, use less pressure on rough surfaces to avoid shredding the swab batting. Swabs should be air dried in a protective environment and stored dry at room temperature, protected from light.

Swab batting can be removed from the shaft of the swab using a scalpel or scissors. The swab head can be removed from the shaft of the swab using scissors: place the swab in a microcentrifuge tube and cut the shaft as close as possible to the batting while leaving the swab head in the tube. Once cut, the shaft can be saved or discarded as per laboratory protocol.

-Supplies required: 0.6 or 1.5 ml disposable microcentrifuge tube (Certified DNase, RNase, Pyrogen, RNA/DNA Free or equivalent), filter/barrier pipette tips (low retention or equivalent), disposable transfer pipettes, Kimwipes

-Reagents required: ddH<sub>2</sub>O and 95% EtOH

---

**Extraction Protocol for RSID™-Semen - Positive control:** A Positive Control for RSID™-Semen can be produced from 50 µl of human semen deposited on a sterile cotton swab.

**Protocol**

- 1.) Deposit 50 µl of human semen onto a sterile cotton swab.
  - 2.) Place swab in a 1.5 ml microcentrifuge tube using laboratory clean technique
  - 3.) Add 1 ml of RSID™-Semen Extraction Buffer and close tube.
  - 4.) Vortex tube vigorously to thoroughly wet swab.
  - 5.) Incubate at room temperature for 1-2 hours.
  - 6.) Mix 5 µl of extract with 95 µl of RSID™-Semen Running Buffer (total volume, 100 µl).
  - 7.) Load entire 100 µl into the sample window of the cassette.
  - 8.) At 10 minutes, score and record results as shown (see diagram, Technical Information Sheet, pg. 3); this will give a strong signal at the Test line.
- 

**Extraction Protocol for RSID™-Semen - Negative control:** A Negative control for RSID™-Semen can be produced from a sterile cotton swab.

**Protocol**

- 1.) Deposit 50 µl of sterile ddH<sub>2</sub>O onto a sterile cotton swab.
  - 2.) Place swab in a 1.5 ml microcentrifuge tube using laboratory clean technique
  - 3.) Add 1 ml of RSID™-Semen Extraction Buffer and close tube.
  - 4.) Vortex tube vigorously to thoroughly wet swab.
  - 5.) Incubate at room temperature for 1-2 hours.
  - 6.) Mix 5 µl of extract with 95 µl of RSID™-Semen Running Buffer (total volume, 100 µl).
  - 7.) Load entire 100 µl into the sample window of the cassette.
  - 8.) At 10 minutes, score and record results as shown (see diagram, Technical Information Sheet, pg. 3); this will give no signal at the Test line.
- 

**Extraction Protocol for RSID™-Semen from Punch or Cutting (separate from STR analysis):**

Many laboratories choose to analyze evidence by testing small punches or cuttings from exhibits (swabs, fabric, etc.) before diverting the remainder of the exhibit for DNA testing. This technique can be easily integrated into RSID™-Semen testing.

**Protocol**

- 1.) Remove a 5 mm diameter punch or ≈20 mm<sup>2</sup> cutting from the exhibit and place in a 0.6 ml microcentrifuge tube using laboratory clean technique.

- 2.) Add 100 µl of RSID™-Semen Extraction Buffer and close tube.
  - 3.) Vortex tube vigorously to thoroughly wet punch or cutting.
  - 4.) Incubate at room temperature for 1-2 hours.
  - 5.) Briefly centrifuge extract to pellet punch or cutting.
  - 6.) Remove 20 µl of extract and bring volume to 100 µl with RSID™-Semen Running Buffer (add 80 µl of running buffer). Set a timer for 10 minutes.
  - 7.) Add extract in running buffer (100 µl) to sample well of an RSID™-Semen cassette. Start timer.
  - 8.) At 10 minutes, score and record results as shown (see diagram, Technical Information Sheet, pg. 3).
  - 9.) Due to high dose Hook effect, weak positive and negative results should be re-tested: Remove 1 µl of the original extract and bring volume to 100 µl with RSID™-Semen Running Buffer (add 99 µl of running buffer). Add extract in running buffer (100 µl) to sample well of an RSID™-Semen cassette and record result at 10 minutes.
- 

#### **Extraction Protocol for RSID™-Semen from whole swab (integrated with STR analysis):**

This protocol is designed to extract an exhibit swab such that *both* stain identification and DNA-STR analysis can be performed from the same sample in a single tube.

#### **Protocol**

- 1.) Place swab in 0.6 or 1.5 ml microcentrifuge tube using laboratory clean technique
- 2.) Add 200 µl of RSID™-Semen Extraction Buffer. Close tube.
- 3.) Vortex tube vigorously to thoroughly wet swab or cutting.
- 4.) Incubate at room temperature for 1-2 hours.
- 5.) Remove a 20 µl aliquot of extract for stain ID testing. Add extract to 80 µl of RSID™-Semen Running Buffer to a final volume of 100 µl in a 0.6 ml microcentrifuge tube. Set timer for 10 minutes.
- 6.) Add extract in running buffer (100 µl) to sample well of an RSID™-Semen cassette. Start timer.
- 7.) At 10 minutes, score and record results as shown (see diagram, Technical Information Sheet, pg. 3).
- 8.) Due to the high dose Hook effect, weak positive and negative results should be re-tested: Remove 1 µl of the original extract and bring volume to 100 µl with RSID™-Semen Running Buffer (add 99 µl of running buffer). Add extract in running buffer (100 µl) to sample well of an RSID™-Semen cassette and record result at 10 minutes.

- 9.) Start DNA extraction with remaining extract –
    - (a) Chelex extraction: add Chelex bead solution directly to swab, proceed as per Chelex extraction protocol.
    - (b) Phenol/Chloroform extraction: remove swab from tube and proceed with DNA extraction as per standard laboratory protocol.
    - (c) Other extraction methods: remove swab from tube and proceed with DNA extraction as per supplied protocol (e.g., Qiagen, Promega or similar).
- 

#### **Validation Protocol -**

This section is provided as a convenience and is meant to be a guide to documenting the validation and testing of RSID™-Semen. All laboratories must follow their own validation and testing protocols.

#### **Validation Summary**

- 1.) Record number of sample types tested. Sample types tested should include those samples commonly encountered in forensic case work.
  - 2.) Record total number of samples tested – a minimum of 5 is usually required for most laboratory audits.
  - 3.) Record test results. Include details regarding precision, sensitivity, accuracy, specificity, and reproducibility. Record the date and the initials of the analyst who performed the validation study.
  - 4.) Record approval of appropriate laboratory supervisory personnel.
  - 5.) Record date of test release. This is the date that the test passes laboratory validation and can be used for samples and casework.
- 

*Not for in vitro diagnostic use*

*Manufactured by:*

 **Independent Forensics**  
4600 ROOSEVELT RD., STE. 201 - HILLSDALE, IL 60162  
TEL. 866-434-2400 FAX 708-978-5115  
[WWW.IFI-TEST.COM/RSID.HTML](http://WWW.IFI-TEST.COM/RSID.HTML)