



## CASPASE FAMILY FLUOROMETRIC SUBSTRATE SET

**Code:** 110-BV-1FSS      **Lot no.:**      **Exp.:** 1 year from date of despatch

**Description:** These fluorometric substrates are for assaying activities of members of the caspase family proteases. These substrates are formulated in a kit in a ready-to-use format.  
YVAD is the recognition sequence for caspase-1/ICE. VDVAD is the recognition sequence for caspase-2/ICH-1. DEVD is the recognition sequence for caspase-3/CPP32. WEHD is the recognition sequence for Caspase-5/ICE<sub>rel</sub>III. VEID is the recognition sequence for caspase-6/Mch2. IETD is the recognition sequence for caspase-8/FLICE, and LEHD is the recognition sequence for Caspase-9/Mch6.

**Quantity:** 125 ul (1 mM) of each of the following substrates in DMSO:  
Caspase-1 Substrate, Ac-YVAD-AFC      Caspase-6 Substrate, Ac-VEID- AFC  
Caspase-2 Substrate, Ac-VDVAD- AFC      Caspase-8 Substrate, Ac-IETD- AFC  
Caspase-3 Substrate, Ac-DEVD- AFC      Caspase-9 Substrate, Ac-LEHD- AFC  
Caspase-5 Substrate, Ac-WEHD- AFC

**Storage:** Store at -20°C.

### Assay Procedure:

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture *without* induction. Count cells and pellet  $1-5 \times 10^6$  cells or use 50-200 ug cell lysates if protein concentration has been measured.
3. Resuspend cells in 50 ul of chilled Cell Lysis Buffer (Code.# 110-BV-1CLBB) and incubate cells on ice for 10 minutes.
4. Centrifuge for 1 min in a microcentrifuge (10,000 x g).
5. Add 50 ul of 2X Reaction Buffer (Code.# 110-BV-1RBA) containing 10 mM DTT (Code # 110-BV-IDTT) to each sample.
6. Add 5 ul of the 1 mM AFC conjugated substrates (50 uM final conc.) into each tube individually and incubate at 37°C for 1-2 hours.
- 7.

Read samples in a fluorometer with a 400nm excitation filter and a 505nm emission filter. For a plate reading set-up, transfer the samples to a 96 well plate. Alternatively you may perform the entire assay directly in a 96-well plate.

Fold-increase in caspase activity can be determined by comparing these results with the level of the uninduced control.

Note: Background reading from cell lysates and buffers should be subtracted from the readings of both induced and the uninduced samples before calculating fold increase in caspase activity

**For Research Use Only. Not For Diagnostic or Therapeutic Use.**

**Conditions:** The information disclosed herein is not be construed as a recommendation to use the above product in violation of any patents. ImmunoKontakt will not be held responsible for patent infringement or other violations that may occur with the use of our products.

Rev. 10/04

**UNITED KINGDOM, INTERNATIONAL:**  
**AMS BIOTECHNOLOGY (EUROPE) LTD.**  
63B MILTON PARK, ABINGDON OXON OX14 4RX  
**TEL: +44 (0) 1235 828200**  
**FAX: +44 (0) 1235 820482**

**SWITZERLAND:**  
**AMS BIOTECHNOLOGY (EUROPE) LTD.**  
CENTRO NORD-SUD 2E, CH-6934 BIOGGIO (LUGANO)  
**TEL: +41 (0) 91 6045522**  
**FAX: +41 (0) 91 6051785**

**DEUTSCHLAND:**  
**AMS BIOTECHNOLOGY (EUROPE) LTD.**  
**TEL: +49 (0) 69 779099**  
**FAX: +49 (0) 69 13376880**