
**Advanced Glycation End Products (AGEs)**  
**Anti Pentosidine Monoclonal Antibody (Clone No. PEN-12)**  
**Peroxidase conjugated**

Reaction of protein amino groups with glucose leads, through the early products such as a Schiff base and Amadori rearrangement products, to the formation of advanced glycation end products (AGEs). Recent immunological studies using anti-AGEs antibody (6D12) demonstrated the presence of AGEs-modified proteins in several human tissues: ( i ) human lens (nondiabetic and noncataractous), ( ii ) renal proximal tubules in patients with diabetic nephropathy and chronic renal failure, ( iii ) diabetic retina, ( iv ) peripheral nerves of diabetic neuropathy, ( v ) atherosclerotic lesions of arterial walls, ( vi )  $\beta$ 2-microglobulin forming amyloid fibrils in patients with hemodialysis-related amyloidosis, ( vii ) senile plaques of patients with Alzheimer's disease, ( viii ) the peritoneum of CAPD patients, ( ix ) skin elastin in actinic elastosis, and ( x ) ceroid / lipofuscin deposits. These results suggest a potential role of AGEs-modification in normal aging as well as age-enhanced disease processes. This antibody named as 6D12 has been used to demonstrate AGEs-modified proteins in these human tissues, indicating potential usefulness of this antibody for histochemical identification and biochemical quantification of AGEs-modified proteins.

Pentosidine is one of the Maillard compounds identified by Monnier *et al* in 1989. It has been proved to cross-link Arginine to Lysine residue and be detected in  $\beta$ 2-microglobulin from patients with hemodialysis-related amyloidosis.

Package Size	50 $\mu$ g (200 $\mu$ L/ vial)
Format	Mouse monoclonal antibody, Peroxidase conjugated 0.25 mg/mL
Buffer	Block Ace as a stabilizer, containing 0.1% Proclin as a bacteriostat
Storage	Store below $-20^{\circ}$ C. Once thawed, store at $4^{\circ}$ C. Repeated freeze-thaw cycles should be avoided.
Clone No.	PEN-12
Subclass	IgG1
Purification method	The splenic lymphocytes from BALB/c mouse, immunized with pentosidine-HSA were fused to myeloma P3U1 cells. The cell line (PEN-12) with positive reaction was grown in ascitic fluid of BALB/c mouse, from which the antibody was purified by Protein G affinity chromatography and conjugated.
Working concentration for immunohistochemistry:	5-10 $\mu$ g/mL ; for ELISA: 0.1-1.0 $\mu$ g/mL
Specificity	Reaction to pentosidine-HSA was suppressed by free-pentosidine in competitive ELISA.

Pentosidine

Advanced Glycation End Products (AGEs)  
**Anti Pentosidine Monoclonal Antibody (Clone No. PEN-12)**  
**Peroxidase conjugated**

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1.	Sell, D.R., <i>et al.</i> (1989): Structure Elucidation of a Senescence Cross-link from Human Extracellular Matrix. <i>J.Biol.Chem.</i> 264 : 21597-21602
2.	Miyata T, <i>et al.</i> (1996): Identification of pentosidine as a native structure for advanced glycation end products in $\beta$ 2-microglobulin-containing amyloid fibrils in patients with dialysis-related amyloidosis. <i>Proc.Natl.Acad.Sci.USA</i> 93 : 2353-2358
3.	Jono T, <i>et al.</i> (2002): Accumulation of imidazolone, pentosidine and N(epsilon)-(carboxymethyl)lysine in hippocampal CA4 pyramidal neurons of aged human brain. <i>Pathol Int.</i> 52:563-71※

\* This product was developed in conjunction with Meiji Institute of Health Science, Meiji Milk Product Co.,Ltd

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
1. Sell, D.R., *et al.* (1989): Structure Elucidation of a Senescence Cross-link from Human Extracellular Matrix. *J.Biol.Chem.*264 : 21597-21602
2. Miyata T, *et al.* (1996): Identification of pentosidine as a native structure for advanced glycation end products in  $\beta$  2-microglobulin-containing amyloid fibrils in patients with dialysis-related amyloidosis. *Proc.Natl.Acad.Sci.USA* 93 : 2353-2358
3. Jono T, *et al.* (2002): Accumulation of imidazolone, pentosidine and N(epsilon)-(carboxymethyl) lysine in hippocampal CA4 pyramidal neurons of aged human brain. *Pathol Int.*52:563-71※

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