Research Use Only

# Extracellular Vesicles derived from Cheese Whey

Cat. No. EXCW1000

November 20, 2025

## [I] Background

Exosomes are membrane vesicles with a diameter of approximately 50–150 nm that are secreted by most cells. They encapsulate miRNA, mRNA, proteins, micropeptides, and other molecules, and are suggested to play a role in transmitting information from the originating cell or tissue to distant target cells or tissues. Microorganisms also secrete extracellular vesicles (EVs) similar to exosomes.

This product is a purified fraction derived from cheese whey, which is generated during the cheese-making process using raw milk. It contains both bovine-derived milk exosomes (2,254.8 µg/mL, measured using the Cosmo Bio EXBMEL Bovine Milk Exosome ELISA Kit) and lactic acid bacteria-derived EVs secreted during fermentation ( $1.378 \times 10^7$  particles/mL, measured using the Cosmo Bio EVEL01 Gram-positive bacteria EV Quantification ELISA Kit). This product is purified from cheese whey with three step purification process, ultrafiltration, and ion-exchange chromatography 1 and 2. Through removal of impurities such as lipids and protein components, a highly purified EV fraction is obtained. The product is prepared at high particle density ( $5 \times 10^{12}$  particles/mL) and can be used for a wide range of applications from in vitro to in vivo studies as a DDS carrier. Bulk supply is also available upon request.

## Milk Exosome + LAB EVs


Fig. 1. Overview of Cheese Whey-Derived EVs Production Diagram

### [II] Instructions for Use

1. Return it to room temperature before use.
2. Add 1 mL of sterile distilled water and mix gently until dissolved. After reconstitution, store at 4°C and use within 3 months.
3. Dilute the reconstituted solution 1:1000 to 1:10,000 in culture medium before use.

### [III] Reference materials


Fig. 2. SDS-PAGE analysis and protein quantification


Fig. 3. Particle size distribution (Nano Sight)

Item No.	Product Name	Package	Particle Count	Example Applications	Storage
<b>EXCW1000</b>	Extracellular Vesicles from Cheese Whey	Freeze-dried (1 mL)	$5 \times 10^{12}$ particles/mL	Promotion of epidermal cell proliferation; enhancement of collagen synthesis in fibroblasts, etc.	Refrigerated

### [III] Reference Data

Human epidermal HaCaT cells were seeded at low density ( $3 \times 10^3$  cells/well) in a 96-well plate. EVs samples were added at a 1:1000 dilution. After 48 hours, cells were fixed with 4% paraformaldehyde and stained with Giemsa. Colony area was quantified using ImageJ.

Compared with the control (no addition), the colony area increased in the group treated with cheese whey-derived EVs. Furthermore, the activity approached that of the positive control group treated with exosomes derived from human dental pulp mesenchymal stem cells. Example shown: 1:1000 dilution.


Fig.4. Epidermal Cell Proliferation Assay

### [IV] References

- [1] Maghraby MK, *et al.*, Extracellular vesicles isolated from milk can improve gut barrier dysfunction induced by malnutrition, *Scientific Reports.*, 2021, **11**:7635.
- [2] Cui Z, *et al.*, Potential therapeutic effects of milk-derived exosomes on intestinal diseases, *J. Nanobiotechnology.*, 2023, **21**:496.
- [3] Lu L, *et al.*, Novel roles of bovine milk-derived exosomes in skin antiaging, *J. Cosmet. Dermatol.* 2024, **23**:1374–1385.

[V] Not for Human Use.

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## Extracellular Vesicles derived from Cheese Whey

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